

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new software version in comparison with the previous version.

- New version: 1.2.50
- Old version: 1.1.14

2 New features

Topic	Description
XTR checkbox	<ul style="list-style-type: none">▪ An XTR checkbox displays XTR data stored in the Baselined property if XTR metadata exists for that sample.▪ Samples taken with previous versions of the device firmware may not be recognized as XTR samples.▪ Baselined and XTR data are selected or deselected separately from each other (displays Calibrated for that type if deselected).
XTR options for operating procedures	Operating procedures have XTR options or can be disabled from settings (once an XTR-enabled device has been connected).
XTR metadata	RMNB, BRMS, and SRMP file formats contain representations of XTR metadata.
Disable XTR features	[Disable XTR Features] option in General Settings : Settings must be saved before UI changes are enacted.
X and Y calibration	Added [XTR X-Calibration] button. Y calibration has been updated (service portal only).

Topic	Description
Auto synchronization	<p>[Auto Synchronize] of operating procedures and libraries is selectable in Settings.</p> <ul style="list-style-type: none"> ▪ Added options for when auto synchronization is off to allow custom synchronization. ▪ Users can choose custom synchronization of operating procedures, licensed libraries, and user libraries when [Auto Synchronize] is disabled. ▪ Samples will always synchronize automatically, regardless of other synchronization settings.
Advanced mode	<p>Advanced mode can be password protected.</p> <ul style="list-style-type: none"> ▪ First time enabling this feature, password change is required. ▪ Add/remove users via File ► Advanced ► Manage Advanced Users, or by toggling option in menu. ▪ Users can change their password once logged into advanced mode.
New language packs	
Spelling corrections	"Mira" is updated to "MIRA".

3 Improvements

Topic	Description
Synchronization changes	Improved handling of synchronization protocols from the device to MIRA Cal DS.
Licensed libraries	<p>All licensed libraries are always "viewable".</p> <ul style="list-style-type: none"> ▪ Licensed libraries that are uploaded as individual LRME + LRMC (as opposed to via a certificate package zip file) do not unlock the library on desktop. In this case libraries will be empty/read-only. However, they may be used to create operating procedures. ▪ The "Illicit Materials (SERS)" library (id=79a2ceb9-3c35-4b61-846d-c6cc9807c89f) is a special case and is hard-coded to be added as a read-only library in certain cases.
Operating procedures	<ul style="list-style-type: none"> ▪ Default operating procedure identification threshold changed from 0.85 to 0.75. ▪ Operating procedures with XTR Mode other than "Never" will no longer synchronize to devices without the XTR capabilities. ▪ Operating procedures created on desktop will default to XTR Mode of "Never" if XTR features are not enabled, or "Automatic" if XTR features are enabled. ▪ Operating procedures built on devices that do not have XTR will have their XTR Mode set to "Never" on desktop when a sync operation runs.

.....

4 Fixed bugs

Multiple bug fixes and backend improvements to synchronization operations (device and desktop).

5 Known issues

One cannot view XTR processed data and baselined data for the same spectrum at the same time. One cannot view baselined data of an XTR processed spectrum.