


Measuring Module Analog



6.02101.010

Product manual

8.0108.8012EN / 2022-01-12



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

Measuring Module Analog

6.02101.010

Product manual

8.0108.8012EN /
2022-01-12

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Overview	1
1.1	Measuring Module Analog – Product description	1
1.2	Measuring Module Analog – Product versions	1
1.3	Symbols and conventions	2
1.4	Further information	2
1.5	Accessories	2
2	Safety	4
2.1	Intended use	4
2.2	Responsibility of the operator	4
2.3	Requirements for operating personnel	5
2.4	Safety instructions	5
2.4.1	Danger from electrical potential	5
2.4.2	Danger from biological and chemical hazardous substances	5
2.4.3	Danger from highly flammable substances	6
2.5	Design of warning messages	6
2.6	Meaning of warning signs	7
3	Functional description	9
3.1	Measuring Module Analog – Overview	9
3.2	Measuring Module Analog – Functional description	10
3.3	System – Signals	10
3.4	Measuring Module Analog – Interfaces	10
4	Delivery and packaging	12
4.1	Delivery	12
4.2	Packaging	12
5	Operation and control	13
5.1	Operation	13
6	Maintenance	14
6.1	Maintenance	14
6.2	Cleaning the product surface	14
7	Troubleshooting	16

1 Overview

1.1 Measuring Module Analog – Product description

The Measuring Module Analog is used as a measuring input for analog electrodes on an OMNIS Titrator or an OMNIS Titration Module.


1.2 Measuring Module Analog – Product versions

The product is available in the following versions:

Table 1 Product versions

Article number	Designation
6.02101.010	Measuring Module Analog

The article number and serial number for identifying the product can be found on the type plate:



1	(01) = Article number in accordance with GS1 standard	2	(21) = Serial number
3	(240) = Metrohm article number	4	Certification
5	Certification	6	Technical specifications

-

2 Safety

2.1 Intended use

Metrohm products are used for the analysis and handling of chemicals.

Usage therefore requires the user to have basic knowledge and experience in handling chemicals. Knowledge regarding the application of fire prevention measures prescribed for laboratories is also mandatory.

Adherence to this technical documentation and compliance with the maintenance specifications make up an important part of intended use.

Any utilization in excess of or deviating from the intended use is regarded as misuse.

Specifications regarding the operating values and limit values of individual products are contained in the "Technical specifications" section, if relevant.

Exceeding and/or not observing the mentioned limit values during operation puts people and components at risk. The manufacturer assumes no liability for damage due to non-observance of these limit values.

The EU declaration of conformity loses its validity as soon as modifications are carried out on the products and/or the components.

2.2 Responsibility of the operator

The operator must ensure that basic regulations on occupational safety and accident prevention in chemical laboratories are observed. The operator has the following responsibilities:

- Instruct personnel in the safe handling of the product.
- Train personnel in the use of the product according to the user documentation (e.g. install, operate, clean, eliminate faults).
- Train staff on basic occupational safety and accident prevention regulations.
- Provide personal protective equipment (e.g. protective glasses, gloves).
- Provide suitable tools and equipment to carry out the work safely.

The product may be used only when it is in perfect condition. The following measures are required to ensure the safe operation of the product:

- Check the condition of the product before use.
- Remedy defects and malfunctions immediately.
- Maintain and clean the product regularly.

2.3 Requirements for operating personnel

Only qualified personnel may operate the product. Qualified personnel are persons who meet the following requirements:

- Basic regulations on occupational safety and accident prevention for chemical laboratories are known and complied with.
- Knowledge of handling hazardous chemicals is present. Personnel have the ability to recognize and avoid potential dangers.
- Knowledge regarding the application of fire prevention measures for laboratories is available.
- Safety-relevant information is communicated and understood. The personnel can operate the product safely.
- The user documentation has been read and understood. The personnel operate the product according to the instructions in the user documentation.

2.4 Safety instructions

2.4.1 Danger from electrical potential

Contact with electrical potential can cause serious injuries or death. To avoid danger from electrical potential, observe the following:

- Operate the product only if it is in perfect condition. The housing must also be intact.
- Only use the product with the covers fitted. If covers are damaged or missing, disconnect the product from the energy supply and contact the regional Metrohm service representative.
- Protect live components (e.g. power supply unit, power cord, connection sockets) against moisture.
- Always have maintenance work and repairs on electrical components carried out by a regional Metrohm service representative.
- Disconnect the product from the energy supply immediately if at least one of the following cases occurs:
 - The housing is damaged or open.
 - Live parts are damaged.
 - Moisture penetrates.

2.4.2 Danger from biological and chemical hazardous substances

Contact with biological hazardous substances may cause poisoning from toxins or infections from microorganisms. Contact with aggressive chemical substances may cause poisoning or chemical burns. To avoid danger from biological or chemical hazardous substances, observe the following:

- Label the product according to regulations if it is used for substances that have a potential for chemical hazards and are generally subject to the Hazardous Substances Ordinance.
- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Dispose of hazardous substances in accordance with regulations.
- Clean and disinfect contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.
- Proceed as follows in case of a return shipment to Metrohm AG or a regional Metrohm representative:
 - Decontaminate the product or product component.
 - Remove the labeling for hazardous substances.
 - Create a declaration of decontamination and enclose it with the product.

2.4.3 Danger from highly flammable substances

Using highly flammable substances or gases may cause fires or explosions. To avoid danger from highly flammable substances, observe the following:

- Avoid ignition sources.
- Use protective grounding.
- Use exhaust equipment.

2.5 Design of warning messages

There are 4 hazard levels for warning messages. The following signal words are used for classifying the hazard levels in warning messages:

- **DANGER** indicates a hazardous situation which, if not avoided, will result in serious injury or death.
- **WARNING** indicates a hazardous situation which, if not avoided, could result in serious injury or death.
- **CAUTION** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
- **NOTICE** indicates a hazardous situation which, if not avoided, could result in property damage.

Warning messages differ in design (color and warning sign) depending on the hazard level:



DANGER

Type and source of danger

Consequences when not observing the notice: An irreversible injury that may result in death is very probable.

- Measures to avoid the danger



WARNING

Type or source of danger

Consequences when not observing the notice: A serious injury that may result in death is probable.

- Measures to avoid the danger



CAUTION

Type or source of danger

Consequences when not observing the notice: A minor to moderate injury is probable.






- Measures to avoid the danger

2.6 Meaning of warning signs

This documentation uses the following warning signs:

Table 2 Warning sign according to ISO 7010

Warning sign	Meaning
	General warning sign
	Warning of electrical voltage
	Warning of hand injuries
	Warning of sharp object
	Warning of hot surface
	Warning of biological hazard

Warning sign	Meaning
	Warning of toxic materials
	Warning of flammable materials
	Warning of corrosive substances
	Warning of optical radiation
	Warning of laser beams

Depending on the intended use of the product, the corresponding warning sign stickers must be placed on the product.

3 Functional description

3.1 Measuring Module Analog – Overview




Figure 1 Measuring Module Analog – Overview

1 INPUT 1

Connection socket for potentiometric sensors (green coding), temperature sensors (red coding) and polarizable sensors (blue coding)

2 INPUT 2

Connection socket for potentiometric sensors (green coding) and temperature sensors (red coding)

3 REF

Connection socket for reference electrodes

4 Fastening screws

Fastening screws, left and right. These fasten the measuring module in the housing and ground the electronics.




Figure 2 Measuring Module Analog – Measuring inputs

1	INPUT 1	2	INPUT 2
3	REF		

Measuring inputs INPUT 1 and INPUT 2

The measuring inputs **INPUT 1** and **INPUT 2** are marked with colored circle segments. The markings indicate that only certain types of electrode cables may be plugged into the connection socket:

Table 3 Meaning of the colors

Red	The connector supports temperature sensors.
Blue	The connector supports polarized sensors.
Green	The connector supports potentiometric sensors.

REF measuring input

Reference electrodes can be plugged into the **REF** measuring input.

5 Operation and control

5.1 Operation

The product can be operated via the OMNIS Software. Further information on the OMNIS Software under [OMNIS Help](#).




Severe injuries with possibly fatal consequences.

- Prerequisite:**

- ### Required accessories:

- 3** Clean the connectors with a dry cloth.

8 Disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.



9 Technical specifications

9.1 Ambient conditions

Nominal function range	+5 to +45 °C	at max. 80% relative humidity, non-condensing
Storage	+5 to +45 °C	

9.2 Measuring module – Energy supply

Power consumption	max. 0.6 W	
Energy transmission		inductive coupling

9.3 Measuring module – Dimensions

Measurements		
Width	105 mm	
Height	31 mm	
Depth	72 mm	
Weight	approx. 420 g	



9.4 Measuring module – Housing

Materials

<i>Lid</i>	AW-5754 H12 / H22	aluminum, coated
<i>Back panel</i>	PBT	poly(butylene terephthalate)
<i>Enclosure</i>	GD-ZnAl4Cu1	zinc die cast, nickel-plated

IP degree of protection IP 40

9.5 Measuring Module Analog – Connectors specifications

Measuring inputs

INPUT 1

Socket		round plug 7-pin, size 0, 45°
Potentiometric	pH, ISE, Redox	measuring input for potentiometric electrodes
Temperature	Temp.	measuring input for temperature sensors of the Pt1000 or NTC type for automatic temperature compensation
Polarizer	Pol.	measuring input for polarizable electrodes

INPUT 2

Socket		round plug 7-pin, size 0, 45°
Potentiometric	pH, ISE, Redox	measuring input for potentiometric electrodes

Temperature	Temp.	measuring input for temperature sensors of the Pt1000 or NTC type for automatic temperature compensation
REF		reference potential
Type	2 mm	
(INPUT 1 - INPUT 2)	pH, ISE, Redox	potentiometric differential measurement, with respect to REF

9.6 Measuring Module Analog – Display specifications

Status display	LED	green-red
-----------------------	-----	-----------

9.7 Measuring Module Analog – Measuring specifications

Potentiometric

Measuring range	−2,400 to +2,400 mV	
Resolution	1.56 μV	
Measuring accuracy	±0.5 mV	in the measuring range −2,000 mV to +2,000 mV
Input resistance	≥ 1*10 ¹² Ω	
Offset current	≤ ±1*10 ^{−12} A	

Temperature

<i>Pt1000</i>		
Measuring range	−150 to +250 °C	
Resolution	approx. 0.002 °C	
Measuring accuracy	±0.4 °C	in the measuring range −20.0 to +150.0 °C

Polarizer

 $I_{pol} DC$



Polarization current	−200.0 to +200.0 μ A	adjustable in 0.5 μ A steps
Measuring range	−2,400 to +2,400 mV	
Measuring resolution	0.1 mV	

Potentiometric differential measurement

<i>Measuring range</i>	−2,400 to +2,400 mV	
<i>Measuring resolution</i>	1.56 μ V	
<i>Measuring accuracy</i>	± 1.0 mV	in the measuring range −2,000 mV to +2,000 mV

Reference conditions

<i>Relative humidity</i>	$\leq 60\%$	
<i>Ambient temperature</i>	+25 °C (± 3 °C)	
<i>Instrument status</i>		min. 30 minutes in operation

Measuring accuracy

applies for all measuring ranges without sensor error, under reference conditions, measuring interval 100 ms