

Cylinder unit OMNIS special



6.03004.210 / 6.03004.220

Product manual

8.0108.8030EN / v5 / 2024-06-06



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

Cylinder unit OMNIS special

Product manual

8.0108.8030EN / v5 /
2024-06-06

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Overview	1
1.1	Cylinder unit OMNIS – Product description	1
1.2	Cylinder unit OMNIS – Product versions	1
1.3	About the documentation	2
1.4	Further information	2
1.5	Displaying accessories	3
2	Safety	4
2.1	Intended use	4
2.2	Responsibility of the operator	4
2.3	Requirements for operating personnel	5
2.4	Safety instructions	5
2.4.1	Danger from electrical potential	5
2.4.2	Danger from biological and chemical hazardous substances	5
2.4.3	Danger from highly flammable substances	6
2.4.4	Danger from leaking liquids	6
2.4.5	Danger during transport of the product	7
2.5	Design of warning messages	7
2.6	Meaning of warning signs	8
3	Functional description	9
3.1	Dosing unit – Overview	9
3.1.1	Cylinder unit OMNIS – Overview	11
3.2	Dosing unit – Function	12
3.2.1	Cylinder unit OMNIS – Function	13
3.3	Cylinder unit OMNIS – Resistance to chemicals	13
3.3.1	Cylinder unit OMNIS – Resistance to chemicals of the cylinder housing	14
4	Delivery and packaging	15
4.1	Delivery	15
4.2	Packaging	15
5	Cylinder unit OMNIS – Operation	16
5.1	Attaching the cylinder unit OMNIS	17
5.2	Removing the cylinder unit OMNIS	19

6	Maintenance	22
6.1	Performing maintenance on the cylinder unit OMNIS	22
6.2	Cleaning the cylinder unit OMNIS	23
6.3	Storing the cylinder unit OMNIS	25
6.4	Disassembling the cylinder unit OMNIS	26
6.5	Cleaning and greasing the cylinder unit OMNIS	29
6.6	Checking and replacing the cylinder unit OMNIS	31
6.7	Assembling the cylinder unit OMNIS	32
7	Troubleshooting	36
7.1	Cylinder unit OMNIS – Malfunctions	36
7.2	Cylinder unit OMNIS – Adjusting the piston position	39
7.3	Cylinder unit OMNIS – Clearing the jamming	42
8	Disposal	44
9	Technical specifications	45
9.1	Ambient conditions	45
9.2	Cylinder unit OMNIS – Dimensions	45
9.3	Cylinder unit OMNIS – Housing	45
9.4	Cylinder unit OMNIS – Connectors specifications	46
9.5	Cylinder unit OMNIS – Liquid handling specifications	46
10	Dosing unit – Explanations regarding dosing accuracy	47

1 Overview

1.1 Cylinder unit OMNIS – Product description

The cylinder unit is part of the dosing unit, which comprises, in addition to the cylinder unit, the dosing drive. It provides the liquid volume required for the analysis and can be equipped with different volumes.

The cylinder unit OMNIS special is used in particular for dosing the following:

- Aqueous alkali solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkali solutions
- Permanganate solutions
- EDTA solutions

1.2 Cylinder unit OMNIS – Product versions

The product is available in the following versions:

Table 1 Product versions

Article number	Designation	Version feature
6.03004.210	Cylinder unit OMNIS 10 mL special	Volume 10 mL
6.01508.210	Cylinder unit OMNIS 10 mL special, with- out accessories	Volume 10 mL
6.03004.220	Cylinder unit OMNIS 20 mL special	Volume 20 mL
6.01508.220	Cylinder unit OMNIS 20 mL special, with- out accessories	Volume 20 mL


An antidiffusion tip (6.1543.200) is available as an accessory. The antidiffusion tip is used whenever the titration tip is immersed into the sample. The antidiffusion valve prevents the sample from diffusing into the tip.

If accurate dosing is important, a dosing tip (6.1543.060) can be ordered as an alternative to the antidiffusion tip.

1.5 Displaying accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website.

1 Searching for a product on the website

- Go to <https://www.metrohm.com>.
- Click on .
- Enter the article number of the product (e.g. **2.1001.0010**) into the search field and press **[Enter]**.

The search result is displayed.


2 Displaying product information

- To display the products matching the search term, click on **Product models**.
- Click on the desired product.

Detailed information regarding the product is displayed.

3 Displaying accessories and downloading the accessories list

- To display the accessories, scroll down to **Accessories and more**.
 - The **scope of delivery** is displayed.
 - Click on **[Optional parts]** for the optional accessories.
- To download the accessories list, click on **[Download accessories PDF]** under **Accessories and more**.

 Metrohm recommends keeping the accessories list for reference purposes.

2 Safety

2.1 Intended use

Metrohm products are used for the analysis and handling of chemicals.

Usage therefore requires the user to have basic knowledge and experience in handling chemicals. Knowledge regarding the application of fire prevention measures prescribed for laboratories is also mandatory.

Adherence to this technical documentation and compliance with the maintenance specifications make up an important part of intended use.

Any utilization in excess of, or deviating from, the intended use is regarded as misuse.

Specifications regarding the operating values and limit values of individual products are contained in the "Technical specifications" section, if relevant.

Exceeding and/or not observing the mentioned limit values during operation puts people and components at risk. The manufacturer assumes no liability for damage due to non-observance of these limit values.

The EU declaration of conformity loses its validity as soon as modifications are carried out on the products and/or the components.

2.2 Responsibility of the operator

The operator must ensure that basic regulations on occupational safety and accident prevention in chemical laboratories are observed. The operator has the following responsibilities:

- Instruct personnel in the safe handling of the product.
- Train personnel in the use of the product according to the user documentation (e.g. install, operate, clean, eliminate faults).
- Train staff on basic occupational safety and accident prevention regulations.
- Provide personal protective equipment (e.g. protective glasses, gloves).
- Provide suitable tools and equipment to carry out the work safely.

The product may be used only when it is in perfect condition. The following measures are required to ensure the safe operation of the product:

- Check the condition of the product before use.
- Remedy defects and malfunctions immediately.
- Maintain and clean the product regularly.

2.3 Requirements for operating personnel

Only qualified personnel may operate the product. Qualified personnel are persons who meet the following requirements:

- Basic regulations on occupational safety and accident prevention for chemical laboratories are known and complied with.
- Knowledge of handling hazardous chemicals is present. Personnel have the ability to recognize and avoid potential dangers.
- Knowledge regarding the application of fire prevention measures for laboratories is available.
- Safety-relevant information is communicated and understood. The personnel can operate the product safely.
- The user documentation has been read and understood. The personnel operate the product according to the instructions in the user documentation.

2.4 Safety instructions

2.4.1 Danger from electrical potential

Contact with electrical potential can cause serious injuries or death. To avoid danger from electrical potential, observe the following:

- Operate the product only if it is in perfect condition. The housing must also be intact.
- Only use the product with the covers fitted. If covers are damaged or missing, disconnect the product from the energy supply and contact the regional Metrohm service representative.
- Protect live components (e.g. power supply unit, power cord, connection sockets) against moisture.
- Always have maintenance work and repairs on electrical components carried out by a regional Metrohm service representative.
- Disconnect the product from the energy supply immediately if at least one of the following cases occurs:
 - The housing is damaged or open.
 - Live parts are damaged.
 - Moisture penetrates.

2.4.2 Danger from biological and chemical hazardous substances

Contact with biological hazardous substances may cause poisoning from toxins or infections from microorganisms. Contact with aggressive chemical substances may cause poisoning or chemical burns. To avoid danger from biological or chemical hazardous substances, observe the following:

- Label the product according to regulations if it is used for substances that have a potential for chemical hazards and are generally subject to the Hazardous Substances Ordinance.
- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Dispose of hazardous substances in accordance with regulations.
- Clean and disinfect contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.
- Proceed as follows in case of a return shipment to Metrohm AG or a regional Metrohm representative:
 - Decontaminate the product or product component.
 - Remove the labeling for hazardous substances.
 - Create a declaration of decontamination and enclose it with the product.

2.4.3 Danger from highly flammable substances

Using highly flammable substances or gases may cause fires or explosions. To avoid danger from highly flammable substances, observe the following:

- Avoid ignition sources.
- Use protective grounding.
- Use exhaust equipment.

2.4.4 Danger from leaking liquids

Leaking liquids may cause injuries and may damage the product. To avoid danger from leaking liquids, observe the following:

- Check the product and its accessories for leakages and loose connections.
- Replace leaking parts and connecting elements without delay.
- Tighten loose connecting elements.
- Do not loosen tubing connections under pressure.
- Do not remove aspiration tubing under pressure.
- Carefully pull the ends of the tubing out of the containers.
- Carefully let liquids from tubing drain into suitable containers.
- Insert the buret tips completely into the containers.
- Remove and dispose of leaked liquids in accordance with regulations.
- If you suspect that liquid has penetrated the instrument, disconnect the instrument from the energy supply. Then have the instrument checked by a regional Metrohm service representative.

2.4.5 Danger during transport of the product

Chemical or biological substances may be spilled during the transport of the product. Parts of the product may fall down or may be damaged. There is a risk of injury from chemical or biological substances and pieces of broken glass. To ensure safe transport, observe the following:

- Remove loose parts (e.g. sample racks, sample vessels, bottles) before transport.
- Remove liquids.
- Lift and transport the product with both hands on the base plate.
- Lift and transport heavy products only according to instructions.

2.5 Design of warning messages

The present documentation uses warning messages as follows.

Structure

1. Severity of the danger (signal word)
2. Type and source of danger
3. Consequences of disregarding the danger
4. Measures for averting the danger

Hazard levels

Signal color and signal word designate the hazard level.



DANGER

Indicates an immediate danger. It will result in serious injuries or death if not avoided.



WARNING

Indicates a potential danger. Failure to avoid the danger may result in death or serious injury.



CAUTION

Indicates a potential danger. If not avoided, it may result in light or minor injuries.

NOTICE












Indicates a potentially damaging situation. If not avoided, the product or something in the surrounding area could be damaged.

2.6 Meaning of warning signs

Warning signs on the product or in the documentation indicate potential dangers or draw attention to certain behaviors in order to avoid accidents or damage.

Depending on the application purpose, the operating company attaches additional warning signs to the product. The corresponding instructions of the operator must be followed.

Table 2 Warning signs according to ISO 7010 (examples)

Warning signs / meaning		Warning signs / meaning	
	General warning sign		Warning of hot surface
	Warning of sharp object (cut/puncture)		Warning of hand injuries (crushing)
	Warning of electrical voltage		Warning of corrosive substances
	Warning of optical radiation		Warning of a laser beam
	Warning of flammable materials		Warning of biological hazard
	Warning of toxic materials		

3 Functional description

3.1 Dosing unit – Overview



Figure 1 Dosing unit – Overview

1 Cylinder unit

2 Dosing drive
Not in scope of delivery

3.1.1 **Cylinder unit OMNIS – Overview**




Figure 3 *Cylinder unit – Overview*

1	Distributor with 4 ports	2	Cylinder top piece
3	Cylinder housing	4	Centering tube
5	Dosing cylinder	6	Data chip
7	Piston stopper	8	Dosing piston
9	Spring clip	10	Unlocking button




Figure 4 Cylinder unit – Overview from above

The following table shows the standard use of the 4 ports. The use of the ports can be changed in the OMNIS Software.

Port	Use	Connecting or sealing with
1	Dosing	Dosing tip
2	Filling the cylinder	Chemical bottle
3	Not used	
4	Not used	

3.2 Dosing unit – Function

Software control of the dosing unit can be used for accurate dosing of liquid volumes.

The dosing unit is comprised of the following units:

- Dosing drive
- Cylinder unit

The dosing drive is permanently installed in the housing of the instrument. The dosing drive is controlled via the OMNIS Software and is responsible for the accurate dosing of the solution.

If the cylinder unit is placed on top of the dosing drive, the dosing drive assumes responsibility for the following functions:

- **Raising and lowering the dosing piston:**
Solution is aspirated when the dosing piston is lowered. The cylinder fills up.
Solution is dosed when the dosing piston is raised. The cylinder empties.

- **Rotating the cylinder:**

The rotation of the cylinder controls which of the 4 ports the solution flows through.

The valve disk with an opening is located in the middle of the cylinder base.

The distributor disk with 4 openings corresponding to the 4 ports of the distributor is located at the bottom in the cylinder top piece.

The dosing drive rotates the cylinder by 90° stages so that the opening of the valve disk fits with an opening on the distributor disk. This results in a passage for the solution to the corresponding port of the distributor.

3.2.1 Cylinder unit OMNIS – Function

The cylinder unit is part of the dosing unit. It provides the liquid volume required for the analysis. A distributor with 4 ports enables the filling and emptying of the cylinder.

The cylinder unit OMNIS special is used in particular for dosing the following:

- Aqueous alkali solutions
- Titrant 5
- Silver nitrate solutions
- Nonaqueous alkali solutions
- Permanganate solutions
- EDTA solutions

3.3 Cylinder unit OMNIS – Resistance to chemicals

Conventional reagents and media can be dosed with the cylinder unit. The materials of the single parts that come into contact with the liquid being dosed have been selected for maximum resistance to chemicals and functionality.

However, not all types of aggressive or high-concentration reagents can be conveyed without difficulty. It is the user's own responsibility to determine the resistance of the various single parts to specific, aggressive media.

Observe the following notes to ensure the functional capability of the cylinder unit:

- When using strong inorganic alkalis or concentrated solutions that could crystallize, adhere to the (*see "Cylinder unit OMNIS – Resistance to chemicals of the cylinder housing", chapter 3.3.1, page 14*).
- The temperature of the media must not exceed 50 °C.

5.1 Attaching the cylinder unit OMNIS

Default settings for the ports 1 and 2


Port 1 is defined as dosing port and port 2 as fill port in the data chip default settings of the cylinder unit. The following instructions describe the default setting.

If the ports should be used differently from the default setting, adjust the ports in the OMNIS Software in **Properties** ► **Specific data**.

Preparing for attaching

- 1 Open the **Manual control** of the dosing unit in the OMNIS Software, see [software help](#).
- 2 Start the **Exchange position** function.

Attaching the cylinder unit

 These instructions describe the default installation as defined in the OMNIS Software.

Prerequisite:

- Dosing drive: Valve coupling and piston rod are in the exchange position (port 2 is set).
- Cylinder unit: The piston stopper is flush with the base of the cylinder housing. The centering tube is in the correct position.


Required accessories:

- Wrench (6.2739.000)
- 2 FEP tubings (6.1805.100)
- Titration tip (6.1543.200)

1 Aligning the cylinder unit


Rotate the cylinder unit until the marking with the label **UNLOCK** is in line with the marking on the dosing drive.

2




Set the cylinder unit down onto the two twistlocks straight from above.

3 Locking the cylinder unit



Rotate the cylinder unit to the left until it stops.
Use the marking with the label **LOCK** as a guide.

4 Mounting the tubings



Tighten an FEP tubing (6.1805.100) into port 1.

This FEP tubing is used as dosing tubing. Tighten the other end to the titration tip (6.1543.200).

5 Tighten the other FEP tubing (6.1805.100) into port 2.

This FEP tubing is used as filling tubing. Tighten the other end to the OMNIS Liquid Adapter.

6 Firmly tighten the tubing with the wrench (6.2739.000).

See also

Cylinder unit OMNIS – Overview (chapter 3.1.1, page 11)

5.2 Removing the cylinder unit OMNIS

Preparing the removal

- 1 Open the **Manual control** of the dosing unit in the OMNIS Software, see [software help](#).
- 2 Start the **Empty** function.
- 3 Start the **Exchange position** function.

Removing the cylinder unit

Prerequisite:


- Dosing drive: Valve coupling and piston rod are in the exchange position (port 2 is set).
- Cylinder unit: The piston stopper is flush with the base of the cylinder housing. The centering tube is in the correct position.

1 Removing the tubing



Unscrew the dosing tubing and the filling tubing.

2 Unlocking the cylinder unit



Rotate the cylinder unit to the right as far as the **UNLOCK** position.

3 Raising the cylinder unit



Raise the cylinder unit straight upwards.

See also

Cylinder unit OMNIS – Overview (chapter 3.1.1, page 11)

6 Maintenance

6.1 Performing maintenance on the cylinder unit OMNIS



CAUTION

Damage by chemicals

Chemicals may escape through leaks. Aggressive chemicals damage the data chip and the dosing drive.

- Regularly check the cylinder unit for leaking liquid (under the dosing piston, on the base of the centering tube or the cylinder unit).
- Regularly check the cylinder and dosing piston for wear. (*see "Checking and replacing the cylinder unit OMNIS", chapter 6.6, page 31*)
- Replace a defective cylinder unit immediately. Do not continue using it.



Depending on the application, cylinder, dosing piston and flat stopcock are subject to different mechanical strain. A cylinder unit, for example, that is often used for alkaline, high-concentration or crystallized reagents will be subject to higher wear. This results in shorter maintenance intervals. The cylinder unit must therefore be replaced more often.

Maintenance work	Maintenance interval
Check cylinder housing for contaminations and clean if necessary. (see "Cleaning the cylinder unit OMNIS", chapter 6.2, page 23)	Daily
Check electrical contacts for contaminations and clean if necessary. (see "Cleaning the cylinder unit OMNIS", chapter 6.2, page 23)	Weekly if using:
Clean cylinder top piece and valve disk. Grease cylinder housing and valve disk. (see "Cleaning and greasing the cylinder unit OMNIS", chapter 6.5, page 29)	<ul style="list-style-type: none"> concentrated solutions that tend to crystallize EDTA solutions, ultrapure solvents and ultrapure water organic solvents alkaline (e.g. KOH or isopropyl alcohol), corrosive or high-concentration reagents
Inspect cylinder and dosing piston. (see "Checking and replacing the cylinder unit OMNIS", chapter 6.6, page 31)	Every 3 months if using unproblematic reagents.

6.2 Cleaning the cylinder unit OMNIS



WARNING

Chemical hazardous substances

Contact with aggressive chemical substances may cause poisoning or chemical burns.

- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Clean contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.



The cylinder unit requires appropriate care. Excess contamination of the cylinder unit results in malfunctions and a reduction in the service life.

Prerequisite:

The cylinder unit has been removed from the dosing drive. (*see "Removing the cylinder unit OMNIS", chapter 5.2, page 19*)

1 Cleaning the cylinder housing



The cylinder housing is not dishwasher proof.

Clean the cylinder housing with lukewarm water and dishwashing detergent.

2

If the cylinder top piece is stuck, place the cylinder unit with the cylinder top piece facing down in warm water (possibly with a little dishwashing detergent) for at least 30 minutes. (*see "Cylinder unit OMNIS – Malfunctions", chapter 7.1, page 36*)

- ### 6.3 Storing the cylinder unit OMNIS

Titrants	Cleaning solution
Aqueous alkali solutions	Deionized water
Titrant 5	Methanol
AgNO ₃ solutions	0.1 mol/L HNO ₃
Nonaqueous alkaline solutions	Deionized water
KMnO ₄ solutions	(NH ₄) ₂ Fe(SO ₄) ₂ (1)
EDTA solutions	Ethanol

 If you are using water-sensitive reagents, rinse the cylinder with solvent and then store it empty.

- 1 Connect the cleaning solution to the cylinder unit.
- 2 Execute the "Best practice" operating procedure. This empties the cylinder unit and executes 6 cleaning cycles with the cleaning solution. Afterwards it is ensured that the cylinder unit is in exchange position and filled with cleaning solution.
- 3 If the cylinder unit is to be stored empty,
 - remove the filling tubing from the bottle with rinsing solution and
 - start the **Empty** function.

- 4 Start the **Exchange position** function.
- 5 Store the cylinder unit at ambient temperature and protect it from direct sunlight.

Rinsing the cylinder unit automatically

To automatically rinse the cylinder unit, download the "Best practice" method for automatically rinsing the cylinder unit as a template or create it yourself.

6.4 Disassembling the cylinder unit OMNIS

Prerequisite:

The cylinder is empty and the cylinder unit has been removed from the dosing drive. (see *"Removing the cylinder unit OMNIS", chapter 5.2, page 19*)

 CAUTION

Damage to cylinder unit

Improper handling when disassembling leads to damage of the cylinder unit and/or the dosing piston.


- Do not use force to rotate the cylinder top piece. Instead, place the cylinder unit in water. (*see "Cylinder unit OMNIS – Clearing the jamming", chapter 7.3, page 42*)
- Follow the instructions for taking the cylinder unit apart.
- Do not disconnect the cylinder from the cylinder base.
- Do not remove the dosing piston from the cylinder.

1 CAUTION

Contact with chemicals

Chemicals in the cylinder unit may cause chemical burns.

- Empty and rinse the cylinder unit before disassembling it.
- Wear protective equipment, especially gloves.



Remove the cylinder top piece.



3



Remove the cylinder element (centering tube incl. cylinder).

4



Hold the black cylinder base. Remove the cylinder together with the dosing piston from the centering tube.

Turn the cylinder base with the cylinder around and place it on a flat surface.

 Do not disconnect the cylinder from the cylinder base.
Do not remove the dosing piston from the cylinder.

6.5 Cleaning and greasing the cylinder unit OMNIS

Cleaning the cylinder unit

Prerequisite:

The cylinder top piece and the cylinder element (centering tube incl. cylinder) have been removed. (see *"Disassembling the cylinder unit OMNIS", chapter 6.4, page 26*)


1



Clean the cylinder top piece with water.

i Do not remove the distributor disk from the cylinder top piece.

2



Clean the contact surface of the valve disk with ethanol.

 Do not disconnect the cylinder from the cylinder base.

- 3 Rinse the cylinder housing and centering tube with water and wipe them with ethanol.


Checking the cylinder unit


- 1** (see "Checking and replacing the cylinder unit OMNIS", chapter 6.6, page 31)

Greasing the centering tube and the valve disk

Required accessories:


Grease (6.2803.010 or 6.2803.000)

- 1  Apply a very thin layer of grease. Wipe off excess grease with a cloth.



Grease the centering tube.

2



Grease the valve disk sparingly.

Assembling the cylinder unit

1 *(see "Assembling the cylinder unit OMNIS", chapter 6.7, page 32)*

6.6 Checking and replacing the cylinder unit OMNIS

Prerequisite:

The cylinder unit has been taken apart. *(see "Disassembling the cylinder unit OMNIS", chapter 6.4, page 26)*

1 Inspecting the cylinder

- Are rough areas or scratches visible on the cylinder?

2 Inspecting the dosing piston


- Are scratches visible on the surface of the dosing piston?
- Is any unevenness visible on the sealing lips of the dosing piston?
- Are the cylinder and the dosing piston leakproof?

3 Replacing the cylinder unit

If any of these defects is visible, replace the entire cylinder unit.




Place the cylinder element (centering tube incl. cylinder) on the cylinder top piece.



Rotate the cylinder element (centering tube incl. cylinder) in such a way that the markings on the centering tube and on the cylinder top piece are positioned above one another.

4




Put the cylinder housing into place.


The markings on the cylinder housing, on the centering tube and on the cylinder top piece must be positioned above one another.

The cylinder housing rests on the cylinder top piece.

5



Hold on to the cylinder top piece and rotate the cylinder housing to the left until the unlocking lever snaps in place.



If the piston stopper protrudes the cylinder housing, push the cylinder unit vertically downwards on a flat surface. (*see "Cylinder unit OMNIS – Adjusting the piston position", chapter 7.2, page 39*)

See also

Attaching the cylinder unit OMNIS (chapter 5.1, page 17)

Cylinder unit OMNIS – Adjusting the piston position (chapter 7.2, page 39)

7 Troubleshooting

Messages on malfunctions and errors are displayed in the control software or in the embedded software (e.g. on the display of an instrument) and contain the following information:

- Descriptions of causes of malfunctions (e.g. jammed drive)
- Descriptions of problems with the control (e.g. missing or invalid parameter)
- Information on how to solve the problem

System components with status display elements also indicate malfunctions and errors with a red flashing LED.

Troubleshooting on the product is often only possible with the control software or the embedded software (e.g. initializing, moving to a defined position).

7.1 Cylinder unit OMNIS – Malfunctions

Problem	Cause	Remedy
The entire cylinder unit rotates during dosing.	The friction points have not been greased.	Grease the centering tube and the valve disk. <i>(see "Cleaning and greasing the cylinder unit OMNIS", chapter 6.5, page 29)</i>
There is liquid under the dosing piston, on the base of the centering tube or the cylinder unit.	The dosing piston is worn out or defective.	Replace the cylinder unit.
	The cylinder does not seal.	Replace the cylinder unit.
	The distributor disk does not seal.	Clean valve disk and distributor disk. <i>(see "Cleaning and greasing the cylinder unit OMNIS", chapter 6.5, page 29)</i>
The cylinder housing does not shut.	The spring clip is inserted incorrectly.	Remove the cylinder housing and insert the spring clip correctly.
The cylinder unit can only be removed from the dosing drive with difficulty.	The friction points have not been greased.	Grease the centering tube and the valve disk. <i>(see "Cleaning and greasing the cylinder unit OMNIS", chapter 6.5, page 29)</i>

Problem	Cause	Remedy
	The coupling is contaminated.	Remove the contaminations on the coupling between cylinder unit and drive.
The cylinder unit cannot be removed from the dosing drive. The cylinder top piece cannot be removed from the cylinder unit easily.	The cylinder unit is not in the exchange position.	Start the Exchange position function.
	The cylinder unit is jammed because the valve disk and the distributor disk stick together.	<i>(see "Cylinder unit OMNIS – Clearing the jamming", chapter 7.3, page 42)</i>
The cylinder unit is recognized either not at all or incorrectly.	The cylinder unit was attached or assembled incorrectly.	<ul style="list-style-type: none"> ▪ <i>(see "Removing the cylinder unit OMNIS", chapter 5.2, page 19)</i> ▪ <i>(see "Attaching the cylinder unit OMNIS", chapter 5.1, page 17)</i> ▪ Check the correct placement of the cylinder unit. ▪ Switch the control instrument off and then back on again. ▪ If the problem persists, contact your regional Metrohm representative.
	The data chip is mechanically damaged or impaired by chemicals.	<ul style="list-style-type: none"> ▪ <i>(see "Cleaning the cylinder unit OMNIS", chapter 6.2, page 23)</i> ▪ If the problem persists, contact your regional Metrohm representative.
Air bubbles are in the cylinder or in the dosing tubing.	Air intrusion through leaky connection.	<ul style="list-style-type: none"> ▪ Check the ends of the tubing, in particular the end of the aspiration tubing. ▪ Tighten the tubing connections at the fill port with the wrench (6.2739.000). ▪ Check the correct placement of the OMNIS Liquid Adapter. ▪ Check the tubing connection of the bottle cap multi-use.
	The reagent degasses excessively; the released air forms bubbles.	<ul style="list-style-type: none"> ▪ Start the Preparing function to rinse the cylinder unit and all tubing. ▪ Reduce the filling rate. ▪ Degas the reagent with ultrasound, nitrogen or in a vacuum.

Problem	Cause	Remedy
		<ul style="list-style-type: none"> Carefully loosen the distributor disk from the valve disk, placing the cylinder unit once again in warm water in the event of problems. Rinse all parts with distilled water, taking care while doing so not to remove the piston or to disconnect the cylinder from the base plate. Dry the parts (e.g. with nitrogen). Grease the parts. <i>(see "Cleaning and greasing the cylinder unit OMNIS", chapter 6.5, page 29)</i> Assemble the cylinder unit and attach it to the titrator. Initialize the dosing unit.


7.2 Cylinder unit OMNIS – Adjusting the piston position

If the piston stopper is not flush with the base of the cylinder housing, the piston rod of the dosing drive cannot reach the dosing piston.

Required accessories:

- Piston tongs (6.1546.030)


1 Inserting the piston tongs



Insert the piston tongs into the opening of the cylinder.

- Hold on to the plunger (blue) of the piston tongs and carefully pull up the dosing piston until it stops.


4 Removing the piston tongs



- Press the plunger of the piston tongs (blue) and hold it down.
- Remove the piston tongs.

5 Check the position of the piston stopper

If the piston stopper protrudes the cylinder housing (see zoomed in view), perform the following steps.



- Place the cylinder unit on a flat surface (e.g. laboratory table).
- Carefully push the cylinder unit vertically downwards onto the support surface.

The piston stopper is positioned flush with the cylinder housing. The cylinder unit can now be attached.


7.3 Cylinder unit OMNIS – Clearing the jamming

If it is difficult or impossible to rotate the cylinder top piece, the valve disk and the distributor disk stick to one another. The software displays an error message.

Clearing the jamming of the attached cylinder unit


- 1 Remove the tubings and the stoppers.
- 2 Open the **Manual control** of the dosing unit in the OMNIS Software .
- 3 Start the **Exchange position** function.
- 4 As soon as the **Exchange position** function has been executed successfully, remove the cylinder unit and soak it (*see "Clearing the jamming of the not attached cylinder unit", page 43*).

If the **Exchange position** function is not executed successfully, proceed as follows.
- 5 Remove the liquid from all ports with a syringe.
- 6 Fill every port with deionized water or another suitable solvent with a syringe (with needle). Make sure that the needle reaches the valve disk (remains stuck in the port).
- 7 Allow the cylinder unit to stand for 2 hours.
- 8 Initialize the dosing unit in the OMNIS Software, if possible, or force the valve to switch using the **Filling** function or the **Exchange position** function.

 Do not force the valve to switch multiple times.
- 9 If the cylinder unit remains jammed, repeat steps 5 to 8.

Clearing the jamming of the not attached cylinder unit

- 1 Place the jammed cylinder unit with the cylinder top piece facing down in warm water (possibly with a little dishwashing detergent) for at least 30 minutes.
- 2 Remove the cylinder unit from the water and dry thoroughly.
- 3 Attach the cylinder unit to the dosing drive and lock it.
- 4 Initialize the dosing unit in the OMNIS Software, if possible, or force the valve to switch using the **Dosing** function or the **Exchange position** function.


 Do not force the valve to switch multiple times.
- 5 If the cylinder unit remains jammed, repeat steps 1 to 4.

If the error persists, call the regional Metrohm service representative or replace the entire cylinder unit.

See also

Cylinder unit OMNIS – Overview (chapter 3.1.1, page 11)

8 Disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.



<i>Distributor</i>	PCTFE	Polychlorotrifluorethy- lene
Degree of protection	IP 40	

9.4 Cylinder unit OMNIS – Connectors specifications

Electrical contacts	4	spring contacts
---------------------	---	-----------------

9.5 Cylinder unit OMNIS – Liquid handling specifications

Cylinder unit

<i>Cylinder volume</i>	10, 20 mL
------------------------	-----------

Tubings

<i>Tubing nipple outer thread</i>	M6	
<i>Inner diameter</i>	2 mm	
<i>Material</i>	FEP	Tetrafluoroethylene- hexafluoropropylene copolymer

10 Dosing unit – Explanations regarding dosing accuracy

The dosing unit is equipped with a resolution of 102,400 steps per stroke. Provided that the cylinder is entirely filled, these 102,400 steps can aspirate and dose the following typical whole-number volumes.

Cylinder volume	Examples for volumes that can be dosed with microliter precision	Theoretically smallest volume step
2 mL	5 µL, 10 µL, 15 µL, ...	19.53125 nL
5 mL	25 µL, 50 µL, 75 µL, ...	48.828125 nL
10 mL	25 µL, 50 µL, 75 µL, ...	97.65625 nL
20 mL	25 µL, 50 µL, 75 µL, ...	195.3125 nL
50 mL	125 µL, 250 µL, 375 µL, ...	488.28125 nL

If a volume is dosed or aspirated that is not a multiple of the theoretically smallest volume step, it is rounded to the previous volume step. The maximum deviation from the required volume thus equals the smallest volume step.

Limit values of the dosing unit

The dosing unit fulfills the *systematic error* and the *random error* according to DIN EN ISO 8655-3 "Piston-operated volumetric apparatus – Part 3: Piston burettes".

Metrohm guarantees that the dosing unit is in compliance with the following limit values at the time of shipment:

Cylinder volume	Maximum permissible systematic measurement deviation		Maximum permissible random measurement deviation	
2 mL	± 0.5%	± 10 µL	± 0.1%	± 2 µL
5 mL	± 0.3%	± 15 µL	± 0.1%	± 5 µL
10 mL	± 0.2%	± 20 µL	± 0.07%	± 7 µL
20 mL	± 0.2%	± 40 µL	± 0.07%	± 14 µL
50 mL	± 0.2%	± 100 µL	± 0.05%	± 25 µL

i The regional Metrohm representatives offer the possibility of on-site inspections and certifications of dosing units with respect to accuracy.