# Metrohm White paper

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# Metrohm White paper

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- [1] D. Bertrand et al., La spectroscopie infrarouge et ses applications analytiques, Tec&Doc, 2000.
- [2] R. A. Shaw and H. H. Mantsch, «Near-IR Spectrometers», in Encyclopedia of Spectroscopy and Spectrometry, Academic Press, 1999, pp. 1451-1461.
- [3] T. Meyer et al., «Suppression of mechanical noise and the selection of optimal resolution in FT-NIR spectroscopy», NIR news, 17 (8), pp. 12-14, 2006.
- [4] O. Kolomiets et al., «The influence of spectral resolution on the quantitative near infrared spectroscopic determination of an active ingredient in a solid drug formulation», Journal of Near Infrared Spectroscopy, 12 (5), pp. 271-277, 2004.
- [5] FOSS, NIR Spectrometer Technology Comparison, White paper from FOSS, 2013.
- [6] E. Ciurczak et al., «Examination of NIR Spectrometers: Dispersive vs. Interferometric Types», Amer. Pharm. Rev., 11 (4), pp. 3-5, 2008.
- [7] E. Voigtman, «The Multiplex Disadvantage and Excess Low-Frequency Noise», Applied Spectroscopy, 41 (7) pp. 1182-1184, 1987.
- [8] F. Grandmont, «Développement d'un spectromètre imageur à transformée de Fourier pour l'astronomie», Doctoral Dissertation, Université Laval, 2006.
- [9] J. B. I. Reeves and C. M. Zapf, «Discriminant analysis of selected food ingredients by near infrared diffuse reflectance spectroscopy», J. Near Infrared Spectroscopy, (5), pp. 209-221, 1997.
- [10] D. Cozzolino et al., «The use of visible and near-infrared reflectance spectroscopy to predict colour on both intact and homogenised pork muscle», LWT - Food Science and Technology, 2003.
- [11] P. R. Armstrong et al., «Comparison of Dispersive and Fourier-Transform NIR Instruments for Measuring Grain and Flour Attributes», Applied Engineering in Agriculture, 22 (3), pp. 453-457, 2006.
- [12] A. Kazeminy et al., «A comparison of near infrared method development approaches using a drug product on different spectrophotometers and chemometric software algorithms», Journal of Near Infrared Spectroscopy, 17 (5), pp. 233-245, 2009.
- [13] A. Mouazen, «Near infrared spectroscopy for agricultural materials an instrument comparison», Journal of Near Infrared Spectroscopy, 13 (2), pp. 87-98, 2005.
- [14] P. Chalus et al., «Comparaison de spectromètres proche infrarouge pour la détermination de principes actifs dans des comprimés», Spectra Analyse, 247, pp. 44-49, 2005.