

# Application Bulletin



Of interest to:  
Laboratories in the cement industry


No. 17/2 e

## Potentiometric routine determination of the $\text{CaCO}_3$ content in raw meal from the cement industry

<b>Summary</b>	An automatic potentiometric titration method is described for determining the $\text{CaCO}_3$ contents of raw meal. The sample is mixed with HCl and heated. The excess of HCl is backtitrated with NaOH.																		
<b>Apparatus</b>	<ul style="list-style-type: none"> <li>▶ Combi-Titrator, consisting of               <ul style="list-style-type: none"> <li>2.632.0010 pH-Meter</li> <li>2.614.0010 Impulsomat</li> <li>2.665.0040 Dosimat Titrating Stand</li> <li>2.586.001X Labograph</li> </ul> </li> <li>▶ 6.3006.220 Exchange Unit with 20 mL burette cylindre</li> <li>▶ 6.0202.000 combined glass electrode</li> <li>▶ 6.2115.010 cable (connection 614/632)</li> <li>▶ 3.980.3200 cable (connection 614/665/586)</li> </ul>																		
<b>Reagents</b>	<ul style="list-style-type: none"> <li>▶ Hydrochloric acid : <math>c(\text{HCl}) = 1 \text{ mol/L}</math></li> <li>▶ Sodium hydroxide : <math>c(\text{NaOH}) = 0.5 \text{ mol/L}</math></li> </ul>																		
<b>Method</b>	<p>Thoroughly dry the raw meal sample, then weigh 1.25 g app. into a 250 mL Erlenmeyer flask. Pipette 25 mL of <math>c(\text{HCl}) = 1 \text{ mol/L}</math> to the sample and bring it to boil. After boiling for 10 s, empty the content of the flask into a beaker and rinse well with distilled water. Titrate with <math>c(\text{NaOH}) = 0.5 \text{ mol/L}</math>, the apparatus settings as follows:</p> <table style="margin-left: 40px; border: none;"> <tbody> <tr> <td style="padding-right: 20px;">Impulsomat</td> <td style="padding-right: 20px;">End-point</td> <td>pH = 7.00</td> </tr> <tr> <td></td> <td>Range</td> <td>coarse</td> </tr> <tr> <td></td> <td>dynamics</td> <td>6</td> </tr> <tr> <td></td> <td><math>t_{\text{clock}}</math></td> <td>1.5 s</td> </tr> <tr> <td></td> <td><math>t_{\text{delay}}</math></td> <td>15 s</td> </tr> <tr> <td colspan="2">Burette titration speed</td> <td>"3"</td> </tr> </tbody> </table>	Impulsomat	End-point	pH = 7.00		Range	coarse		dynamics	6		$t_{\text{clock}}$	1.5 s		$t_{\text{delay}}$	15 s	Burette titration speed		"3"
Impulsomat	End-point	pH = 7.00																	
	Range	coarse																	
	dynamics	6																	
	$t_{\text{clock}}$	1.5 s																	
	$t_{\text{delay}}$	15 s																	
Burette titration speed		"3"																	
<b>Calculation</b>	<p> <math>A = \text{mL sample } c(\text{HCl}) = 1 \text{ mol/L}</math>  <math>B = \text{mL consumption } c(\text{NaOH}) = 0.5 \text{ mol/L}</math>  <math>E = \text{weight of sample in mg}</math> </p> $\% \text{ CaCO}_3 = \frac{(A - 2B) \cdot 50.05 \cdot 100}{E}$																		
<b>Remarks</b>	<ul style="list-style-type: none"> <li>▶ HCl can be standardised with pure <math>\text{CaCO}_3</math> (Merck p.a.); 1 g of substance is treated as mentioned above.</li> <li>▶ The accuracy of the titration is <math>\pm 0.1\%</math>.</li> </ul>																		

**Potentiometric routine determination of the  $\text{CaCO}_3$  content in raw meal from the cement industry**

**Example**



V/mL = 9.76

EP pH = 7

coarse

1.5 s / dynamics 6 / t delay 15 s

Burette "3"

Labograph 200 mm/min