

## Application Bulletin 36/3 e

# Half wave potentials of metal ions for the determination by polarography

### Summary

In the following tables the half-wave potentials or peak potentials of 90 metal ions are listed. The half-wave potentials (listed in Volt) are measured at the dropping mercury electrode (DME) at 25 °C unless indicated otherwise.

### Instruments

VA instrument  
capable of operating a Multi-Mode Electrode and supporting differential pulse (DP) measuring mode

### Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Unsilanized mercury drop capillary	6.1226.030
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel Filled with c(KCl) = 3 mol/L	6.1245.010
AE	Separate Pt rod electrode	6.0343.x00
	or Glassy carbon rod with	6.1247.000
	Electrode holder	6.1241.x20

### Legend

w	Well marked wave or peak
i	Poor wave or evaluation impossible (e.g. in the region of increasing background current of the supporting electrolyte)
cw	Catalytic reduction of hydrogen (in some cases suitable for quantitative analysis)
NR	No signal of the respective metal ion (polarographic or voltammetric analysis not possible)
↓	Formation of precipitate. Attention to solubility product!
–	No information available

## Part 1: Various metals and transition metals

### Electrolytes

Acetate	c(CH <sub>3</sub> COOH) = 2 mol/L c(CH <sub>3</sub> COONH <sub>4</sub> ) = 2 mol/L pH 4.6
NH <sub>3</sub> / NH <sub>4</sub> Cl	c(NH <sub>3</sub> ) = 1 mol/L c(NH <sub>4</sub> Cl) = 1 mol/L pH 9.6
KCl	c(KCl) = 1 mol/L
HCl	c(HCl) = 1 mol/L
NaOH	c(NaOH) = 1 mol/L
KNO <sub>3</sub>	c(KNO <sub>3</sub> ) = 0.1 mol/L
KSCN	c(KSCN) = 1 mol/L
TEA	c(Triethanolamine) = 0.3 mol/L c(KOH) = 0.1 mol/L
Tartrate	c(Na <sub>2</sub> -tartrate) = 0.5 mol/L pH = 9.0
Na <sub>2</sub> EDTA	c(Na <sub>2</sub> EDTA) = 0.1 mol/L
Citrate	c(Na <sub>3</sub> -citrate) = 0.1 mol/L c(NaOH) = 0.1 mol/L
Oxalate pH 2	c(Ammonium oxalate) = 0.1 mol/L, pH setting with H <sub>2</sub> SO <sub>4</sub>
Oxalate pH 4	c(Ammonium oxalate) = 0.1 mol/L, pH setting with H <sub>2</sub> SO <sub>4</sub>
Oxalate pH 6	c(Ammonium oxalate) = 0.1 mol/L, pH setting with H <sub>2</sub> SO <sub>4</sub>
Oxalate pH 8	c(Ammonium oxalate) = 0.1 mol/L, pH setting with NH <sub>3</sub>

Metal ion	Acetate	NH <sub>3</sub> /NH <sub>4</sub> Cl	KCl	HCl	NaOH	KNO <sub>3</sub>	KSCN	TEA	Tartrate	Na <sub>2</sub> EDTA	Citrate	Oxalate pH 2	Oxalate pH 4	Oxalate pH 6	Oxalate pH 8
<b>Al(III)</b>	NR	NR	-1.59 i	NR	NR	NR	NR	NR	NR	NR	NR	-	-	-	-
<b>As(III)</b>	-0.95	-1.62 i	-0.07 i	-0.38 w -0.62	-0.25	-0.11 w	NR	-0.16	NR	-1.23 w	-0.22 w	NR	-1.56 w	-1.70 i	-1.58 i
<b>Bi(III)</b>	-0.21 w	NR	-0.06 w	-0.07 w	-0.68	-1.17 i ↓	NR	-0.87 w	-0.31 w cw	-0.54 w	-0.77 w	-0.13 w	-0.21 w	-0.22 w	-0.26 w ↓
<b>Cd(II)</b>	-0.63 w	-0.77 w	-0.62 w	-0.61 w	-0.82 i	-0.55 w	-0.62 w	-0.81	-0.64	NR cw	-0.68 i	-0.56 w	-0.61 w	-0.62 w	-0.62 w
<b>Ce(IV)</b>	NR	NR	NR	NR	NR	NR ↓	NR	NR	NR	NR	NR	-	-	-	-
<b>Co(II)</b>	NR	-1.26	-1.38	NR	NR ↓	-1.27 w -1.54	-1.11	NR	-1.53	NR	NR	NR	NR	NR	-1.52 i
<b>Cr(III)</b>	-1.20 i	-1.34 w ↓	-1.03	-1.00 w	NR ↓	-0.87 -1.01	-1.01 w	NR	-1.58 i	-1.20 w	NR	NR	NR	-1.73 i	NR
<b>Cr(VI)</b>	NR	-0.28 w -1.45 w -1.68	-0.38 -0.88 w	-1.03 i	-0.85 w	-0.32 -1.10	-0.41 -1.00 w	-1.27 w	-0.33 -0.86 -1.71	-1.20 w	-0.73 w	NR	NR	-0.02 w	-1.14 w -1.72 i
<b>Cu(II)</b>	-0.25 w	-0.20 w -0.46 w	-0.18	-0.20	-0.35 -0.43	NR	-0.64	-0.50	-0.08 w	-0.26 w	-0.41 w	-0.05 w	-0.15 w	-0.18 w	-0.18 w
<b>Fe(II)</b>	NR	-1.44 w	-1.57 i cw	i cw	-1.55 i ↓	-1.31 -1.56 i cw	-1.59 w	-1.02 w -1.68	-1.51 i	-0.09 w	-0.83 w -1.57 w	-0.02 w	-0.17	-0.18 w	-0.18 w -1.61 i
<b>Fe(III)</b>	-0.02 -0.24	NR ↓	-1.42 i	NR	NR ↓	-1.29 i	-1.56 w	-0.99 w -1.68	-0.19 -1.53 w	-0.11 w	-0.85 w -1.59 w	-0.18 w	-0.17 w	-0.19 w	-0.19 w
<b>Ge(IV)</b>	-1.34 i	-1.41	NR	NR	NR	-1.26 i cw	i cw	NR	NR	NR	NR	NR	NR	-1.38 w	-1.38 w
<b>In(III)</b>	-0.67 w	-1.08	-0.57	-0.57	-1.12 w	-0.90 i	-0.59 w	-1.51 w	NR	NR	-1.17	-0.61 w	-0.67 w -0.84 w	-0.68 w -0.87 w	-0.70 w -1.18
<b>Mn(II)</b>	NR	-1.59 w	-1.55 w	NR	-1.68 i ↓	-1.46 i	-1.56 i	-0.44	-1.52 w	NR	NR	NR	NR	-1.61 i	-1.60 i

Metal ion	Acetate	NH <sub>3</sub> /NH <sub>4</sub> Cl	KCl	HCl	NaOH	KNO <sub>3</sub>	KSCN	TEA	Tartrate	Na <sub>2</sub> EDTA	Citrate	Oxalate pH 2	Oxalate pH 4	Oxalate pH 6	Oxalate pH 8
<b>Mo(VI)</b>	-0.63 w -1.15 -1.30	NR	NR	-0.08 i	NR	NR	NR	NR	NR	-0.52 -0.75 w	NR	-0.02 -0.11 -0.33 w -0.74 w	-0.17 w -0.29 -0.59	NR	NR
<b>Ni(II)</b>	-1.09 w	-1.07 w	-1.03	NR	NR ↓	-1.00 w	-0.66 w	-1.36	NR	NR	NR	NR	NR	NR	-1.32
<b>Pb(II)</b>	-0.47 w	-0.47 w	-0.41 w	-0.41 w	-0.73 ↓	-0.36 w	-0.39 w	-0.91	-0.59	-1.03 w	-0.69 w	-0.40 w	-0.45 w	-0.56 w	-0.57 w
<b>Pd(II)</b>	-0.53 -0.76	-0.75 w	i cw	NR	-1.33 i	NR	i cw	-1.01 i ↓	NR	NR	NR	-	-	-	-
<b>Sb(III)</b>	-0.39 -0.50 w	-0.81 w	-0.07 -0.15 w	-0.12 w	-0.42 -1.18	-0.19 w	-0.58 w -0.70	-1.31 w	NR	-0.64 w	-0.33 -0.95 w	-0.33 w	-0.46 w	-0.41 w -0.57 w	-0.76 w
<b>Se(IV)</b>	-0.70 w -1.19	-1.54 w	NR	-0.06 -0.47	NR	NR	NR	NR	-1.21	-0.67 w -1.20 w	NR	-0.03 -0.60 w	-0.68 w -1.11	-1.35	-1.35
<b>Sn(II)</b>	-0.13 w -0.63 w	-0.72	-0.43 w	-0.45 w	-0.83 -1.15 w	-0.35 w	-0.45 w -1.59 w	-0.75 w -1.15 w	-0.54 -0.83 -1.06	-0.13	-0.83 w -1.08 w	-0.12 w -0.51 w	-0.15 w -0.62 w	-0.26 w -0.62 w	-0.51 i -0.64 i
<b>Sn(IV)</b>	NR	NR	NR ↓	-0.45	NR ↓	NR	-1.57	NR	NR	-1.16 i	NR	-0.51 w	NR	NR	NR
<b>Ti(III)</b>	-	-	-	-	-	-	-	-	-	-	-	-0.31 w	-0.41 w	-0.52	-0.70
<b>Ti(IV)</b>	-	-	-	-	-	-	-	-	-	-	-	-0.29	-0.42	NR	NR ↓
<b>Tl(I)</b>	-0.44 w	-0.46 w	-0.47	-0.46 w	-0.45	-0.43 w	-0.50 w	-0.45 w	-0.45 w	-0.44 w	-0.47 w	-0.44	-0.44	-0.44	-0.44
<b>U(VI)</b>	-0.42 w	NR	-0.18 i	-0.17	-0.92 i	NR	NR	-0.87	NR	NR	NR	-0.18 w	-0.32	-0.38 w	-0.42 w
<b>V(V)</b>	i cw	-1.12 -1.29 i	-1.10 i cw	-1.11 i	-0.40 i	-1.03 i	-0.50 w	-0.35	-0.33 i	-1.23	-0.76 -1.11 w	NR	NR	-0.90 -1.30 w	-0.19 w -0.88 -1.29 w
<b>W(VI)</b>	-0.26 i	NR	NR	-1.01 i cw	NR	NR	NR	NR	NR	-1.26 i	NR	-	-	-	-
<b>Zn(II)</b>	-1.05	-1.32 w	-0.99 w	-0.99 i	-1.58	-0.97	-1.01 w	-1.58	-1.27 w	NR	-1.38 w	NR	NR	-1.38	-1.29
<b>Zr(IV)</b>	NR	NR	-1.49 w	NR	NR	-0.95 w	-1.58 w	NR	NR cw	NR	NR	-	-	-	-

## Part 2: Alkaline and alkaline earth metals

### Electrolytes

c(TEAI) = 0.1 mol/L	c(tetraethylammonium iodide) = 0.1 mol/L in water
c(TEAI) = 0.1 mol/L, 20% DMF	c(tetraethylammonium iodide) = 0.1 mol/L in water:dimethylformamide 80:20

### Half-wave potentials

Metal ion	Electrolyte	Half-wave potential
Ba(II)	c(TEAI) = 0.1 mol/L	-1.90
Be(II)	c(TEAI) = 0.1 mol/L, pH 2.3	-1.91
Ca(II)	c(TEAI) = 0.1 mol/L	-2.24 i
Cs(I)	c(TEAI) = 0.1 mol/L	-2.04
K(I)	c(TEAI) = 0.1 mol/L	-2.11
Li(I)	c(TEAI) = 0.1 mol/L	-2.31
Mg(II) <sup>1</sup>	c(TEAI) = 0.1 mol/L	<sup>2</sup>
Na(I)	c(TEAI) = 0.1 mol/L	-2.07
Rb(I)	c(TEAI) = 0.1 mol/L	-2.06
Sr(II)	c(TEAI) = 0.1 mol/L, 20% DMF	-1.85

### Comments

Alkaline and alkaline earth metals usually cannot be determined in real samples because the signals are not well separated from each other.

<sup>1</sup> Magnesium is determined indirectly as an oxime or as a complex with solochrome violet RS.

## Part 3: Rare-earth elements

### Electrolytes

c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	c(tetraethylammonium perchlorate) = 0.1 mol/L in dimethylformamide
---	--

### Half-wave potentials

Metal ion	Electrolyte	Half-wave potential <sup>3</sup>
Ce	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.44
Dy	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.08
Er	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.16
Eu	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-0.80 -2.30
Gd	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.12
Ho	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.20
La	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.16
Lu	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.08
Nd	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.22
Pm	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	- (radioactive)
Pr	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.24
Sm	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-1.98
Sc	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-1.76
Tb	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.16
Tm	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.10
Yb	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-1.58
Y	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-2.06

<sup>2</sup> Maximum, non exploitable.

<sup>3</sup> The reference electrode is filled with supporting electrolyte.

### Comments

- For preparation of the standard solutions the oxides of the rare-earth elements are first dissolved in aqua regia, dried in a vacuum desiccator, and then dissolved in DMF.
- With the exception of Ce, Sc, and Yb, the half-wave potentials of all metals become more positive with increasing concentrations.

## Part 4: Actinides and technetium

### Electrolytes

c(LiClO <sub>4</sub> ) = 0.1 mol/L, pH 1.9 ... 3.1	c(LiClO <sub>4</sub> ) = 0.1 mol/L in water. The pH value is adjusted to 1.9 ... 3.1.
Al <sub>2</sub> (SO <sub>4</sub> ) <sub>3</sub> / citric acid, pH 4.5	c(Al <sub>2</sub> (SO <sub>4</sub> ) <sub>3</sub> ) = 0.2 mol/L, c(citric acid) = 1 mol/L in water. The pH value is adjusted to 4.5.
c(KCl) = 0.5 mol/L, pH 2	c(KCl) = 0.5 mol/L in water. The pH value is adjusted to 2.

### Half-wave potentials

Metal ion	Electrolyte	Half-wave potential	Reference
Ac	c(LiClO <sub>4</sub> ) = 0.1 mol/L, pH 1.9 ... 3.1	-1.61	(1)
Am			(2)
Cf			(2)
Cm			(2)
Fm			(3)
Np			(4)
Pu	Al <sub>2</sub> (SO <sub>4</sub> ) <sub>3</sub> / citric acid, pH 4.5	-0.43	(5)
Tc	c(KCl) = 0.5 mol/L, pH 2	-0.15	(6)

### Comments

Since all these elements are highly radioactive, it was not possible to determine them in our laboratory. The following literature references are useful.

### References

- F. David  
Polarographie von Actinium  
C.R. Acad. Sci. (Paris) Ser. C 271 (1970), 440-442  
Ref.: Z.Anal.Chem. 256 (1971), 307
- F. David  
Polarography of radioisotopes 115Eu, 169Yb, 153Gd, 241Am, 244Cm, 252Cf  
Compt. Rend. 270 C (1970), 2112-2115  
Ref.: Electroanal. Abstr. 9 (1971), 140
- F. David; M. Nussonois  
Reduction of Fm(II) in acetic solution by electrolysis at controlled electric tension  
Radiochem. Radioanal. Lett. 11 (1972), 1-5  
Ref.: Electroanal. Abstr. 11 (1973), 328
- V.F. Peretruhin; D.P. Alekseeva  
Polarographic behaviour of Np(VII) in alkaline aqueous solutions  
Izv. Akad. Nauk. SSSR, Ser. Khim. 23 (1974), 1438  
Ref.: Electroanal. Abstr. 13 (1975), 45
- G.W.C. Milner; A.J. Wood  
The determination of plutonium by A.C. polarography  
J. Electroanal. Chem. 7 (1964), 190-205  
Ref.: Z. Anal. Chem. 210 (1965), 70
- G.B.S. Salaria; C.L. Rulfs; P.J. Elving  
Polarographic and coulometric determination of technetium  
Anal. Chem. 35 (1963), 979-982  
Ref.: Z. Anal. Chem. 219 (1965), 440

## Part 5: Miscellaneous

### Electrolytes

c(TBAI) = 0.1 mol/L in 80% ethanol	c(tetrabutylammonium iodide) = 0.1 mol/L in water:ethanol 20:80
c(HCl) = 0.1 mol/L	c(HCl) = 0.1 mol/L in water.
c(KSCN) = 0.1 mol/L, pH 3	c(KSCN) = 0.1 mol/L in water. The pH value is adjusted to 3
c(NaOH) = 2 mol/L	c(NaOH) = 2 mol/L in water.
HCl conc.	w(HCl) = 32%
Ca(OH) <sub>2</sub> sat.	Saturated solution of Ca(OH) <sub>2</sub> in water.
c(KNO <sub>3</sub> ) = 0.3 mol/L	c(KNO <sub>3</sub> ) = 0.3 mol/L in water.
c(H <sub>2</sub> SO <sub>4</sub> ) = 0.1 mol/L	c(H <sub>2</sub> SO <sub>4</sub> ) = 0.1 mol/L in water.
c(Na <sub>2</sub> -NTA) = 0.1 mol/L	c(nitrilotriacetate disodium salt) = 0.1 mol/L in water.
c(NaOH) = 0.1 mol/L / c(KNO <sub>3</sub> ) = 0.1 mol/L	c(NaOH) = 0.1 mol/L + c(KNO <sub>3</sub> ) = 0.1 mol/L in water.
c(EDTA) = 0.1 mol/L	c(ethylenediaminetetraacetic acid) = 0.1 mol/L in water.
c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	c(tetraethylammonium perchlorate) = 0.1 mol/L in dimethylformamide.

### Half-wave potentials

Metal ion	Supporting electrolyte	Half-wave potential
Al(III)	c(TBAI) = 0.1 mol/L in 80% ethanol	-1.67
Au(III)	c(HCl) = 0.1 mol/L	+0.90 <sup>4</sup>
Ga(III)	c(KSCN) = 0.1 mol/L, pH 3	-0.84
Hf		NR
Ir(IV)	c(NaOH) = 2 mol/L	+0.38 <sup>5</sup>
Nb(V)	HCl conc.	-0.60
Os(VIII)	Ca(OH) <sub>2</sub> sat.	-0.42
Pt(IV)	c(KNO <sub>3</sub> ) = 0.3 mol/L	-0.93 <sup>6</sup>
Re(VII)	c(H <sub>2</sub> SO <sub>4</sub> ) = 0.1 mol/L	-0.66
Rh(III)	c(Na <sub>2</sub> -NTA) = 0.1 mol/L	-0.46
Ru(I)	c(NaOH) = 0.1 mol/L / c(KNO <sub>3</sub> ) = 0.1 mol/L	-0.81 w <sup>7</sup>

4 With glassy carbon working electrode.

5 With platinum working electrode.

6 Negative peaks.

Ta		NR
Ti(IV)	c(EDTA) = 0.1 mol/L	-0.34
Th(III)	c(TEA-CIO <sub>4</sub> ) = 0.1 mol/L in DMF	-1.46 -1.71

7 With respect to the Hg/Hg<sub>2</sub>Cl<sub>2</sub>/KNO<sub>3</sub> sat. reference electrode.