

Karl Fischer water determination in non-explosive gases

Of interest to:
General analytical laboratories

Summary

This Application Bulletin describes the determination of water in non-explosive and non-flammable gaseous samples with the coulometric Karl Fischer method. This method is ideal for very low water contents.

Apparatus and accessories

- 2.756.0010 KF Coulometer, cell with diaphragm, including 728 Magnetic Stirrer
- 2.731.0010 Relay Box
- 6.2148.010 Remote Box
- 2.846.0010 Dosing interface
- 6.1439.020 Gas inlet tube with frit
- 6.1805.010 FEP tubing/ M6/ 13 cm
- 6.2811.000 Molecular sieve 0.3 nm
- 6.2151.000 Cable (Dosing Interface – PC)
- 6.2125.100 Cable (Relay Box)
- 2 x Magnetic valve, Lucifer
- T-distribution connector 6 mm MS, Swagelock
- Copper tubing 2 m, 4 x 6 mm
- 4 x Connector 6 mm x ¼" MS, Swagelock
- Connector ¼" x ⅛", Swagelock
- Flow controller Red-y smart controller GSC (Vögtlin Instruments AG, Aesch (CH))
- Pressure meter Bioblock Scientific MP 340A 0...2000 mbarA
- 6.6056.112 *tiamo* full 1.1

Remark: For gases with a high water content (> 50 ppm), a generator electrode without diaphragm can be used. For gases with a water content below 50 ppm, a generator electrode with diaphragm has to be used, otherwise the results will be too high.

Reagents

- Anolyte: Hydranal-Coulomat AG Oven, Sigma-Aldrich 34739
- Catholyte: Hydranal-Coulomat CG, Sigma-Aldrich 34840

Setup

See Appendix

Parameters

Start drift	:	5 µg/min
Extraction time	:	15 min
Drift correction	:	automatic
Stop criterion	:	rel. drift
Relative stop drift:		5 µg/min

Stirrer

Speed : 6

The stirrer speed should be set high enough to obtain thorough mixing of the gas bubbles with the KF reagent. Too strong stirring of the solution may, however, affect the drift value. For vessel contents between 110 and 200 mL, a stirring speed of 6 is recommended.

Karl Fischer reagent

110...200 mL

Volumes larger than 200 mL KF reagent make it difficult to reach a low start drift. Less than 110 mL is not enough to make all the gas react with the solution.

Gas settings

Gas pressure at the cylinder: at least 2 bar (200 kPa)

Gas flow rate: at least 200 mL/min

Analysis

Prior to the analysis, the whole system should be purged with the gas to be measured. Several liters of gas may be needed for a complete rinsing. The gas tubing and valves should be rinsed until a stable water content is reached.

During the extraction time, valve 1 is closed and valve 2 is opened to let the gas enter the titration cell. At the end of the extraction time, valve 2 is closed and valve 1 is opened. In this way the gas lines are constantly rinsed with gas and no water enters the tubing.

Remark: If measurements are performed over several days, conditioning of the reagent should not be interrupted during the whole period.

No blank measurements were carried out as the blank values are equal to the drift. Therefore, only a drift correction was applied.

Calculation

$$V_1 = t \cdot \Phi \quad (1)$$

$$m_{\text{Gas}} [\text{mg}] = 1000 \cdot \frac{V_1 \cdot M \cdot T_0 \cdot P_1}{V_0 \cdot (T_0 + \theta_1) \cdot P_0} \quad (2)$$

$$\text{H}_2\text{O} [\text{ppm}] = 1000 \cdot \frac{m_1}{m_{\text{Gas}}} \quad (3)$$

m_1 = absolute water content [μg]

Φ = gas flow rate [L/min]

t = extraction time [min]

V_1 = added gas volume [L]

V_0 = gas volume at $0^\circ\text{C} = 22.4$ L/mol

T_0 = 273.15 K

θ_1 = gas temperature [$^\circ\text{C}$]

M = molar mass of gas [g/mol]

P_1 = absolute gas pressure at θ_1 [kPa]

P_0 = gas pressure at $T_0 = 101.325$ kPa

Remark: For a correct calculation of the water content, a pressure correction should be applied, especially for gases with a low water content.

Results

The detection limit of this method was found to be around 2 ppm. The determination limit was approximately 5 ppm.

Relative standard deviations with nitrogen gas were found to lie between 5 and 6%.

Remarks

- As magnetic valves are used, this setup is not suitable for explosive gases.
- A good repeatability can only be obtained with a precise flow meter and pressure meter.
- The setting used is different from the one stipulated in the ISO 10101-3:1993 standard. The flowmeter used in the ISO norm contains water, which means that the gas must be dried before entering the flow meter. Accordingly, the flow rate has to be measured after the cell. However, the measured flow rate is reduced by a factor of 2 compared to a measurement before the titration cell.


Literature

ISO 10101-3:1993


Natural gas – Determination of water by the Karl Fischer method – Part 3: Coulometric procedure

Appendix

Instrument setup



Caution: the valves heat up to 75 °C!





License ID 124049905
 Client name RES-NB
 User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

Method parameters

Method Water N2 15 min cond 10
 Method saving date 2005-05-17 09:49:31 UTC+2
 Method version 1
 Method group Main group
 Method status original
 Method saved by (full name)
 Method saved by (short name) esch

START

Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Determination
 Live display 2 Main track
 Statistics on
 Number of single determinations 3
 Conditioning
 Automatic conditioning off

Application note

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2			off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on Sample size
 Value off
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off



License ID 124049905
Client name RES-NB
User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

Mail to
Subject Meldung von tiamo - Methode 'Neue Methode 1' - Befehl
'Hauptspur'
User
Mail from
SMTP Server
POP3 Server
Acoustic signal off
Action off
Stop determination on
Stop determination and series off

Name **Sample position**
Type Number
Assignment on Sample position
Value off
Check at start on
Comment Sample position number
Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on
Message by e-mail off

Mail to
Subject Meldung von tiamo - Methode 'Neue Methode 1' - Befehl
'Hauptspur'
User
Mail from
SMTP Server
POP3 Server
Acoustic signal off
Action off
Stop determination on
Stop determination and series off

Name **Sample size unit**
Type Text
Assignment on Sample size unit
Value off
Check at start on
Comment Sample size unit

Name **ID1**
Type Text
Assignment on ID1



License ID 124049905
 Client name RES-NB
 User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

Value off
 Check at start on
 Comment Sample identification 1

Name **ID2**
 Type Text
 Assignment on ID2
 Value off
 Check at start off
 Comment

Name **ID3**
 Type Text
 Assignment on ID3
 Value off
 Check at start on
 Comment Sample identification 3

CALL conditioning

Call text	Track name	Condition	Condition
determination track	Determination	off	

RECEIVE conditioning is OK

Wait for event/status
 Wait until all conditions are met on
 Wait until one of the conditions is met off

Command	Event message	Comment
KFC	Condok	

Timeout off
 Message

CALL Call valve

Call text	Track name	Condition	Condition
valve 2 on	Valve	off	

SEND start titration

Send event messages

Receiver	Selection	Event message	Comment
Command	KFC	Start Titration	



License ID 124049905
 Client name RES-NB
 User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

RECEIVE wait for end of titration

Wait for event/status

Wait until all conditions are met on
 Wait until one of the conditions is met off

Command	Event message	Comment
KFC	Finished	

Timeout on
 Time 30 min

Message

CALC water content

Result name	Formula	Unit	Decimal places	Assignment	Statistics
water abs	= 'KFC.EP.QTY'	µg	4	RS01	off
N2 mass	= (15.1 * 'CV.Flowrate' * 273.15 * 28.01) / (22.4 * (273.15 + 25))	mg	2	RS02	off
water %	=0.1 * 'KFC.EP.QTY' / 'RS.N2 mass'	%	4	RS04	off
water ppm	=1000 * 'KFC.EP.QTY' / 'RS.N2 mass'	ppm	1	RS03	off
Startdrift	= 'KFC.DRI'	µg/ min	1	RS05	off
H2O	= 'KFC.LP.QTY'	µg	2	RS06	off
Total duration	= 'KFC.DBL'	s	1	RS07	off

Result name **water abs**
 Formula = 'KFC.EP.QTY'
 Unit µg
 Decimal places 4
 Assignment RS01
 Statistics off
 Description RS.'Resultatname'[.VAL]Resultatwert.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **N2 mass**
 Formula = (15.1 * 'CV.Flowrate' * 273.15 * 28.01) / (22.4 * (273.15 + 25))
 Unit mg
 Decimal places 2
 Assignment RS02
 Statistics off
 Description RS.'Resultatname'[.VAL]Resultatwert.



License ID 124049905
 Client name RES-NB
 User res

Program version tiamo 1.1 - 36
 2006-09-04 15:55:52 UTC+2

```

Result monitoring . . . . . off
Save result as common variable . . . . . off
  Name . . . . .
Save result as titer . . . . . off
  Solution name . . . . .
-----
Result name . . . . . water %
Formula . . . . . =0.1 * 'KFC.EP.QTY' / 'RS.N2 mass'
Unit . . . . . %
Decimal places . . . . . 4
Assignment . . . . . RS04
Statistics . . . . . off
Description . . . . . RS.'Resultatname'[.VAL]Resultatwert.
Result monitoring . . . . . off
Save result as common variable . . . . . off
  Name . . . . .
Save result as titer . . . . . off
  Solution name . . . . .
-----
Result name . . . . . water ppm
Formula . . . . . =1000 * 'KFC.EP.QTY' / 'RS.N2 mass'
Unit . . . . . ppm
Decimal places . . . . . 1
Assignment . . . . . RS03
Statistics . . . . . off
Description . . . . . RS.'Resultatname'[.VAL]Resultatwert.
Result monitoring . . . . . off
Save result as common variable . . . . . off
  Name . . . . .
Save result as titer . . . . . off
  Solution name . . . . .
-----
Result name . . . . . Startdrift
Formula . . . . . = 'KFC.DRI'
Unit . . . . . µg/ min
Decimal places . . . . . 1
Assignment . . . . . RS05
Statistics . . . . . off
Description . . . . . RS.'Resultatname'[.VAL]Resultatwert.
Result monitoring . . . . . off
Save result as common variable . . . . . off
  Name . . . . .
Save result as titer . . . . . off
  Solution name . . . . .
-----
Result name . . . . . H2O
Formula . . . . . ='KFC.LP.QTY'
  
```



License ID 124049905
 Client name RES-NB
 User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

Unit µg
 Decimal places 2
 Assignment RS06
 Statistics off
 Description RS.'Resultatname'[.VAL]Resultatwert.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **Total duration**
 Formula = 'KFC.DBL'
 Unit s
 Decimal places 1
 Assignment RS07
 Statistics off
 Description RS.'Resultatname'[.VAL]Resultatwert.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

DATABASE Water N2 copper

Database

Water N2 copper

TRACK Determination

Return immediately on
 Delete all data on

KFC KFC

General/Hardware

Device
 Device name 831_1
 Indicator electrode
 I(pol) 10 µA
 Electrode test on
 Generator electrode
 Generator type with diaphragm
 Generator current auto mA
 Stirrer
 Switch on/off automatically on

Start conditions



License ID 124049905
 Client name RES-NB
 User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

```

Pause
  Pause ..... 0 s
Control parameters
End point
  EP at ..... 50.0 mV
Control
  Dynamics ..... 70 mV
  Max. rate ..... maximum µg/min
  Min. rate ..... 15.0 µg/min
Stop criterion
  Stop criterion ..... rel. drift
  Relative stop drift ..... 5 µg/min
Titration parameters
  Extraction time ..... 900 s
  Temperature ..... 25.0 °C
  Time interval measuring point ..... 2 s
Stop conditions
  Stop time ..... off s
Conditioning
  Conditioning ..... on
  Start drift ..... 10 µg/min
  Drift correction ..... automatic
  Delay after 'Cond ok' ..... 60 s
  Only start titration by a start command from a SEND command ..... on
  Show measured value during conditioning ..... on
Additional evaluations
  Fix end point evaluation ..... off
Additional measured values
  Additional calculated measured values ..... off
  Additional external measured values ..... off

TRACK Valve
  Return immediately ..... on
  Delete all data ..... off

CTRL Valve 2 on
Device
  Device name ..... 808_1
  Remote box ..... 1
Set lines
  Output signal ..... *****1**

WAIT extraction time
Wait
  Stop track and waiting for [Continue] ..... off
  Stop all tracks and waiting for [Continue] ..... off
  Waiting time ..... on
  Time ..... 906
    
```



License ID 124049905
Client name RES-NB
User res

Program version tiamo 1.1 - 36

2006-09-04 15:55:52 UTC+2

Unit s
Message
Record message on
Message by e-mail off
Acoustic signal off

CTRL **Valve 1 on**
Device
Device name 808_1
Remote box 1
Set lines
Output signal *****0**