

# Application Bulletin



Of interest for:  
Trace analysis, environmental analysis

No. 147/1 e

## ***Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry***

<b>Summary</b>	<p>The metals Cd, Co, Cu, Fe, Ni, Pb and Zn are determined in the sub-ppb range (limit of detection 0.05 µg/L) by means of stripping voltammetry. The DP-ASV method is used for Cd, Cu, Pb and Zn whereas Co, Ni and Fe are determined by means of the DP-CSV method (dimethylglyoxime or catechol complexes).</p> <p>Use of the VA Processor and the sample changer allows automatic determination of the above metal ions in one solution. The method has been specially developed for trace analysis in the manufacture of semiconductor chips based on silicon. It can naturally also be employed successfully in environmental analysis.</p>
<b>Instruments</b>	<p>▶ VA Processor 646 with VA Stand 647 and, depending on the degree of automation required, VA Sample Changer 675 and Multi-Dosimats 665.</p>
<b>Reagents</b>	<p>All reagents used must be of at least purity grade puriss. p.a.. The solution must be prepared with ultrapure water (<math>\leq 0.1 \mu\text{S/cm}</math> 20°C).</p> <p>▶ <b>Acetate buffer</b> pH = 4.6 with <math>c(\text{CH}_3\text{COOH}) = 2 \text{ mol/L}</math> and <math>c(\text{NH}_3) = 1 \text{ mol/L}</math>.</p> <p>▶ <b>Ammonium chloride buffer</b> 26 mL <math>w(\text{HCl}) = 0.32</math> and 44 mL <math>w(\text{NH}_3) = 0.25</math> are made up to 100 mL with ultrapure water.</p> <p>▶ <b>Catechol solution:</b> <math>c(\text{Catechol}) = 1 \text{ mol/L}</math> Catechol (5.5 g) is dissolved in 50 mL oxygen-free ultrapure water under a stream of nitrogen. The solution must be tightly stoppered and stored in the dark. The stability of the solution varies between 1 day and up to 3 weeks, depending on the purity of the catechol. (Allow to stand 1 h before use.)</p> <p>▶ <b>Pipes Buffer:</b> <math>c(\text{Pipes}) = 1 \text{ mol/L}</math> 1,4-bis(2-ethanesulfonic acid) piperazine (6.047 g) is mixed with 1 mL <math>w(\text{NaOH}) = 0.3</math> and a little ultrapure water. After adjustment of the pH to 8.0 with <math>w(\text{NH}_3) = 0.25</math>, the solution is made up to 20 mL with ultrapure water.</p> <p>▶ <b>Dimethylglyoxime in triethanolamine</b> Dimethylglyoxime (0.5 g) is mixed with 100 mL triethanolamine/ultrapure water (1:1) over a period of 1 h using a magnetic stirrer. The insoluble components are then filtered off using a fluted filter and the clear filtrate stored in a dark bottle.</p> <p>▶ <b>Standard solutions</b> Solutions weaker than <math>\rho(\text{Me}) = 100 \text{ mg/L}</math> are stored in plastic bottles and should be prepared every week as follows: The concentrated standard solution is diluted with ultrapure water and acidified with 0.4 mL/100 mL conc. <math>\text{HNO}_3</math>.</p>

## Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

**Sample preparation**

▶ One drop  $w(\text{HF}) = 0.4$  (40%) is placed on the chip and evaporated by warming. A white, powdery residue is formed. This is dissolved in a drop of conc. nitric acid and rinsed into the polarographic vessel with ultrapure water. A blank is prepared using the same reagents.

**Method setup**

- ▶ The first step involves simultaneous determination of the metals Zn, Pb, Cd, and Cu at pH = 5.0 according to the DP-ASV method.
- ▶ In the second step, iron is determined in the pipes buffer in the form of a catechol complex at pH = 7.0...7.2 using the DP-CSV method.
- ▶ Finally, Co and Ni are determined as their dimethylglyoxime complexes in  $\text{NH}_3/\text{NH}_4\text{Cl}$  buffer at pH = 9.1.

**Determination**

- ▶ Acetate buffer (50  $\mu\text{L}$ ) is added to the sample solution and the pH adjusted to 5.0 with  $w(\text{NH}_3) = 0.1$ . Duplicate determinations of Zn, Cd, Pb and Cu each with two standard additions are now performed (PAGE 3, lines 1...10).
- ▶ Catechol solution (100  $\mu\text{L}$ ) and pipes buffer (300  $\mu\text{L}$ ) are now added. Iron is also determined in a duplicate analysis with two standard additions (PAGE 3, lines 11...19). The pH value of the solution should be 7.0...7.2.
- ▶ Finally, the pH value is adjusted to 9.1...9.2 by addition of 500  $\mu\text{L}$  ammonium chloride buffer. After addition of 100  $\mu\text{L}$  triethanolamine/dimethylglyoxime solution, duplicate determinations of Co and Ni each with two standard additions are performed (PAGE 3, lines 20...27). A determination requires approx. 100 min when in this form.

▶ **Fig. 1** shows the operation sequence of the method on PAGE 3.

<p style="text-align: center;">Zn,Cd,Pb,Cu,Fe,III und Co im sub-ppb-Bereich MFL 1 EL.TYPE MME</p> <p style="text-align: center;">OPERATIONS/PARAMETERS</p> <pre> 1 PURGE ;STIR ;CADDL ; 300 s 2 PURGE ;STIR ; 30 s 3 OPURGE ;(REP) ; 10 s 4 HMDE ;STIR ;MEAS ; 120 s 4a M.MODE DPH 50 mV 4b T.STEP 600 ms 4c U.SET -1.100 V 5 OSTIR ; 10 s 6 SMP 0 ; 53 s 6a U.END -750 mV 6b U.STEP 4 mV SW.RATE 6.6 mV/ s 7 SMP 1 ; 45 s 7a U.END -450 mV 7b U.STEP 4 mV SW.RATE 6.6 mV/ s 8 SMP 2 ; 36 s                 </pre>	<p style="text-align: center;">METHOD 4 PAGE 3 OPERATION SEQUENCE</p> <p style="text-align: center;">OPERATIONS/PARAMETERS</p> <pre> 8a U.END -210 mV 8b U.STEP 4 mV SW.RATE 6.6 mV/ s 9 SMP 3 ; 61 s 9a U.END 200 mV 9b U.STEP 4 mV SW.RATE 6.6 mV/ s 10 REP ;(BEEP ;AD01J2) ; 11 BEEP ;PURGE ;STIR ; 30 s 12 HOLD ;CADDL ; 13 PURGE ;STIR ; 30 s 14 OPURGE ;(REP) ; 10 s 15 HMDE ;STIR ;MEAS ; 120 s 15a M.MODE DPH -50 mV 15b T.STEP 600 ms 15c U.SET -250 mV 16 OSTIR ; 15 s                 </pre>
--	---

<p style="text-align: center;">Zn,Cd,Pb,Cu,Fe,III und Co im sub-ppb-Bereich MFL 1 EL.TYPE MME</p> <p style="text-align: center;">OPERATIONS/PARAMETERS</p> <pre> 17 SMP 4 ; 60 s 17a U.END -650 mV 17b U.STEP 4 mV SW.RATE 6.6 mV/ s 18 REP ;(BEEP ;AD02J2) ; 19 BEEP ;PURGE ;STIR ; 30 s 20 HOLD ;CADDL ; 21 PURGE ;STIR ; 30 s 22 OPURGE ;(REP) ; 23 HMDE ;STIR ;MEAS ; 30 s 23a M.MODE DPH -75 mV 23b T.STEP 600 ms 23c U.SET -700 mV 24 OSTIR ;MEAS ; 10 s 24a M.MODE DPH -75 mV 24b T.STEP 600 ms 24c U.SET -600 mV                 </pre>	<p style="text-align: center;">METHOD 4 PAGE 3 OPERATION SEQUENCE</p> <p style="text-align: center;">OPERATIONS/PARAMETERS</p> <pre> 25 SMP 5 ; 60 s 25a U.END -1.000 V 25b U.STEP 2 mV SW.RATE 3.3 mV/ s 26 SMP 6 ; 37 s 26a U.END -1.250 V 26b U.STEP 4 mV SW.RATE 6.6 mV/ s 27 REP ;(BEEP ;AD03J2) ; 28 OMEAS ;BEEP ;END ;                 </pre>
---	--

## Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

### Determination (continued)

► For automation with the sample changer, the following additional instruments are needed:

Seven Dosimats with the solutions:

1. Dosimat for the standard addition of Zn, Cd, Pb and Cu
2. Dosimat for the standard addition of Fe
3. Dosimat for the standard addition of Co and Ni
4. Dosimat for the addition of catechol solution (brown glass)

The Dosimats 1...4 are controlled via PAGE 4 (2.665.0010).


5. Dosimat for the addition of pipes buffer
6. Dosimat for the addition of triethanolamine/dimethylglyoxime
7. Dosimat for the addition of ammonium chloride buffer

The Dosimats 5...7 (as well as the two 683 Pump Units – siphoning and rinsing polarographic vessel/electrodes) are controlled via the 677 Drive Unit. The volume to be dispensed is stored in the burette and the addition triggered by a control pulse. Two special 3.980.3350 cables and one 6.2115.010 cable are needed for the connection to the 677 Drive Unit. The ammonium chloride and triethanolamine buffer solutions are controlled via a single cable and are thus both added at the same time.

► Additionally are needed:

- 3 × 6.182.010 Capillary tubing for the three auxiliary Dosimats
- 1 × A.403.009 Stopper  
(3 holes of  $\varnothing = 1.5$  mm must be drilled by the customer).

Fig. 2 Interconnection of the instruments with the cable connections



## Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

### Illustrations

The following Figures provide an overview of the determination of the 7 elements using the sample changer:

**Fig. 3** PAGE 2 of the 646 VA Processors.  
Note particularly items 7...11.

```

1 Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler          METHOD 5 PAGE 2
2 MPL 1          EL.TYPE MME                        GEN.SPECIFICATIONS

PARAMETERS
3 IR.MODE          N
4 SPEED            5
5 D.SIZE           7
6 N.DROPS         5

RECOGNITION
7 SPIKE THRESH    5
8 H.THRESH        2
9 U.TOL           7
10 U.TOL          9
11 ASYM.TOL       8
  
```

**Fig. 4** PAGE 3, now completely filled for work with the sample changer  
(35 items)

```

Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler          METHOD 5 PAGE 3
MPL 1          EL.TYPE MME                        OPERATION SEQUENCE
  
```

OPERATIONS/PARAMETERS		OPERATIONS/PARAMETERS	
1	PURGE ;STIR ;EADDL ; 30 s	8a	U.END -210 mV
2	PURGE ;STIR ; 60 s	8b	U.STEP 4 mV
3	ØPURGE;(REP ; 10 s		SW.RATE 6.6 mV/ s
4	HMDE ;STIR ;MEAS ; 120 s	9	SWP 3 ; 61 s
4a	M.MODE DPN 50 mV	9a	U.END 200 mV
4b	T.STEP 600 ms	9b	U.STEP 4 mV
4c	U.SET -1.109 V		SW.RATE 6.6 mV/ s
5	ØSTIR ; 10 s	10	REP > 1;BEEP ;ADD1J2;
6	SWP 0 ; 53 s	11	PURGE ;STIR ;DOS+ ; 120 s
6a	U.END -750 mV	12	CTRL 3;ØCTRL3;PURGE ; 60.0 s
6b	U.STEP 4 mV	13	PURGE ;STIR ;EADDL ; 60 s
	SW.RATE 6.6 mV/ s	14	PURGE ;STIR ; 120 s
7	SWP 1 ; 45 s	15	ØPURGE;(REP ; 10 s
7a	U.END -450 mV	16	HMDE ;STIR ;MEAS ; 120 s
7b	U.STEP 4 mV	16a	M.MODE DPN -50 mV
	SW.RATE 6.6 mV/ s	16b	T.STEP 600 ms
8	SWP 2 ; 36 s	16c	U.SET -250 mV

```

Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler          METHOD 5 PAGE 3
MPL 1          EL.TYPE MME                        OPERATION SEQUENCE
  
```

OPERATIONS/PARAMETERS		OPERATIONS/PARAMETERS	
17	ØSTIR ; 15 s	25c	U.SET -800 mV
18	SWP 4 ; 60 s	26	SWP 5 ; 63 s
18a	U.END -650 mV	26a	U.END -1.010 V
18b	U.STEP 4 mV	26b	U.STEP 2 mV
	SW.RATE 6.6 mV/ s		SW.RATE 3.3 mV/ s
19	REP > 1;BEEP ;ADD2J2;	27	SWP 6 ; 36 s
20	CTRL 4;ØCTRL4; 20.0 s	27a	U.END -1.250 V
21	PURGE ;STIR ;EADDL ; 120 s	27b	U.STEP 4 mV
22	PURGE ;STIR ; 120 s		SW.RATE 6.6 mV/ s
23	ØPURGE;(REP ; 10 s	28	REP > 1;BEEP ;ADD3J2;
24	HMDE ;STIR ;MEAS ; 30 s	29	ØMEAS ;STIR ;CTRL 1; 12.0 s
24a	M.MODE DPN -30 mV	30	(REP ;
24b	T.STEP 600 ms	31	ØCTRL1;CTRL 2; 20.0 s
24c	U.SET -700 mV	32	ØCTRL2;CTRL 1;REP > 3; 40.0 s
25	ØSTIR ;MEAS ; 10 s	33	ØCTRL1;ØSTIR ;
25a	M.MODE DPN -30 mV	34	CHANGE;PURGE ;STIR ;
25b	T.STEP 600 ms	35	BEEP ;END ;

**Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry**

Illustrations  
(continued)

**Fig. 5** PAGE 4. Please note that minimum 80 µL are dispensed. Only burette cylinders of 1 mL and 5 mL should be employed.

Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler				METHOD 5 PAGE 4	
MPL 1	EL.TYPE MME			ALLOCATIONS	
a	b	c	d	e	f
SOLUTE	U.VERIF	DOS	V.SOLN	m.CONC	m.BLANK
Subst	Ux	Soln	c, v	rho.x	bx
1 Zn	-998 mV	1	c 80 uL	500.0 ug/L	0.000 g
2 Cd	-560 mV	1	c 80 uL	100.0 ug/L	0.000 g
3 Pb	-340 mV	1	c 80 uL	1.000 mg/L	0.000 g
4 Cu	30 mV	1	c 80 uL	100.0 ug/L	0.000 g
5 Fe	-450 mV	2	c 100 uL	500.0 ug/L	0.000 g
6 Ni	-934 mV	3	c 100 uL	1.000 mg/L	0.000 g
7 Co	-1.100 V	3	c 100 uL	40.00 ug/L	0.000 g
8 Catech		4	c 100 uL	110.0 g/L	0.000 g
9 SUPP.ELEC	c(HNO 3)=0.01 mol/L				
10 V.MEAS	20.000 mL				
11 ALIQUOT	1.000				
12 DATE	87-01-12				
13 TIME	11:09				

**Fig. 6** PAGE 5

Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler				METHOD 5 PAGE 5	
MPL 1	EL.TYPE MME			DATA OUTPUT	
SEGMT	a	b	c	d	e
	Y.AXIS/L	Y.AXIS/R	X.AXIS/DIV		
1 SWP 0	0.00 A	300 nA	100 mV		
2 SWP 1	0.00 A	25.0 nA	100 mV		
3 SWP 2	0.00 A	160 nA	100 mV		
4 SWP 3	0.00 A	160 nA	100 mV		
5 SWP 4	0.00 A	-160 nA	100 mV		
6 SWP 5	0.00 A	-400 nA	50.0 mV		
7 SWP 6	-6.00 nA	-20.0 nA	100 mV		
8					
9 RECORD	FR	SXXX			
10					
11 SEND					
12					

**Fig. 7** PAGE 6 with the calculation principles

Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler				METHOD 5 PAGE 6	
MPL 1	EL.TYPE MME			RES.CALCULATION	
ANALYTE	a	b	c	d	
	EVAL	R.QUANT	R.UNIT	SIGNIF.DIG	
1 Zn	N	rho(Zn)	ug/L	4	
2 Cd	N	rho(Cd)	ug/L	4	
3 Pb	N	rho(Pb)	ug/L	4	
4 Cu	N	rho(Cu)	ug/L	4	
5 Fe	N	rho(Fe)	ug/L	4	
6 Ni	N	rho(Ni)	ug/L	4	
7 Co	N	rho(Co)	ug/L	4	
8					
	(EV.QUANT	+ ADDEND)	* FACTOR	/ DIVISOR	
11 Zn	A	0.00000	1.00000 E+ 6	1.00000	
12 Cd	A	0.00000	1.00000 E+ 6	1.00000	
13 Pb	A	0.00000	1.00000 E+ 6	1.00000	
14 Cu	A	0.00000	1.00000 E+ 6	1.00000	
15 Fe	A	0.00000	1.00000 E+ 6	1.00000	
16 Ni	A	0.00000	1.00000 E+ 6	1.00000	
17 Co	A	0.00000	1.00000 E+ 6	1.00000	
18					

## Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

Illustrations  
(continued)

**Fig. 8** Result block of a blank value determination

```

METROHM 646 VA-PROCESSOR (5.646.5041)
Zn,Cd,Pb,Cu,Fe,Ni,Co mit Probenwechsler          METHOD 5
MPL 1      EL.TYPE  MME

SUPP.ELEC  c(HNO3)=0.01 mol/L
V.MEAS     20.000 mL
ALIQOUT    1.000

REMARK     Halbleiter - Analytik
           7 Elemente in einer Methode; Blindw.HNO3 0.01mol/L
NAME       Wittmann
RUN#       -1
    
```

ANALYTE	L R S	U.SUBST	EV.VALUE	DELTA	m.ANALYTE
Zn	A0 0 0	-962 mV	101.5 nA		
	A0 1 0	-966 mV	104.4 nA		
	A1 0 0	-966 mV	124.0 nA		
	A1 1 0	-966 mV	124.0 nA	21.05 nA	
	A2 0 0	-966 mV	150.3 nA		
	A2 1 0	-966 mV	156.9 nA	29.55 nA	
	m.STD	40.00 ng	SLOPE	1.580 ug/uA	
Cd	A0 0 1	-551 mV	! 792.4 pA		
	A0 1 1	-544 mV	! 816.3 pA		
	A1 0 1	-545 mV	7.575 nA		
	A1 1 1	-544 mV	3.575 nA	2.751 nA	
	A2 0 1	-545 mV	6.649 nA		
	A2 1 1	-545 mV	6.898 nA	3.218 nA	
	m.STD	8.000 ng	SLOPE	2.690 ug/uA	
Pb	A0 0 2	-351 mV	41.50 nA		
	A0 1 2	-352 mV	43.90 nA		
	A1 0 2	-351 mV	60.84 nA		
	A1 1 2	-351 mV	61.67 nA	18.55 nA	
	A2 0 2	-351 mV	79.20 nA		
	A2 1 2	-351 mV	81.18 nA	18.43 nA	
	m.STD	80.00 ng	SLOPE	4.325 ug/uA	
Cu	A0 0 3	15 mV	25.33 nA		
	A0 1 3	15 mV	27.07 nA		
	A1 0 3	7.3 mV	35.17 nA		
	A1 1 3	6.7 mV	34.64 nA	8.705 nA	
	A2 0 3	6.7 mV	43.63 nA		
	A2 1 3	5.9 mV	44.21 nA	9.018 nA	
	m.STD	8.000 ng	SLOPE	902.7 ng/uA	
Fe	A0 0 4	-451 mV	16.42 nA		
	A0 1 4	-454 mV	16.51 nA		
	A1 0 4	-451 mV	24.87 nA		
	A1 1 4	-452 mV	24.71 nA	8.325 nA	
	A2 0 4	-445 mV	32.59 nA		
	A2 1 4	-446 mV	33.53 nA	8.275 nA	
	m.STD	50.00 ng	SLOPE	6.023 ug/uA	
Ni	A0 0 5	-948 mV	105.3 nA		
	A0 1 5	-948 mV	109.5 nA		
	A1 0 5	-948 mV	114.9 nA		
	A1 1 5	-948 mV	137.3 nA	18.67 nA	
	A2 0 5	-947 mV	136.3 nA		
	A2 1 5	-947 mV	159.7 nA	21.94 nA	
	m.STD	100.0 ng	SLOPE	4.923 ug/uA	
Co	A0 0 6	-1.089 V	3.047 nA		
	A0 1 6	-1.089 V	3.164 nA		
	A1 0 6	-1.087 V	6.766 nA		
	A1 1 6	-1.087 V	8.159 nA	4.207 nA	
	A2 0 6	-1.086 V	11.32 nA		
	A2 1 6	-1.085 V	13.43 nA	4.916 nA	
	m.STD	4.000 ng	SLOPE	876.8 ng/uA	


Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

Illustrations  
(continued)

rho(Zn)	=	8.829	ug/L
rho(Ni)	=	-----	ug/L
rho(Co)	=	-----	ug/L
rho(Cd)	=	! 97.36	E- 3 ug/L
rho(Fe)	=	-----	ug/L
rho(Pb)	=	9.242	ug/L
rho(Cu)	=	1.180	ug/L
rho(Cd)	=	-----	ug/L
rho(Pb)	=	-----	ug/L
rho(Fe)	=	4.362	ug/L
rho(Zn)	=	-----	ug/L
rho(Ni)	=	25.31	ug/L
rho(Co)	=	.1375	ug/L


SMPL.V,m 20.0000 mL IDENT Test mit Probenw.  
DATE 87-01-09 TIME 12:37

Fig. 9 Voltammograms for the data in Figure 8



## Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

Illustrations  
(continued)



### Remarks

- ▶ The electrodes and polarographic vessels are stored in  $w(\text{HNO}_3) = 0.08$  (8%) and also rinsed with this solution.
- ▶ The addition of triethanolamine increases the sensitivity of the Co determination. At the same time, it suppresses possible interferences due to relatively large Zn contents. The Ni peak is somewhat attenuated.
- ▶ In the standard additions, the peak potential of the iron-catechol complex is shifted slightly. This is due to the acidic standard addition solution. With this shift, which is dependent on pH, interferences caused by relatively large amounts of Cu can also be eliminated. (Cu also forms catechol complexes with peak potentials around -250 mV. For the iron complexes, the following dependences apply: pH = 6.9 -0.40 V // pH = 7.2: -0.45 V.)
- ▶ Pipes has a limited storage life. We found that the substance had decomposed after two years in a broached reagent bottle and was thus unfit for use. Pipes itself also shows a peak which moreover lies near that of the iron-catechol complex. Since this peak is small and flat, however, it can be subtracted as a blank value on PAGE 4 (e.g. 2 ng/20 mL). If buffer solutions are prepared using decomposed pipes, the base peak is steep and large. The blank value is increased by a factor 10 and makes evaluation of the iron impossible (Fig. 10 and 11).



**Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry**

Remarks  
(continued)

Fig. 10






Fig. 11



► The following method can be very useful to check the quality of the pipes:  
Pipes buffer (100 µL) is added to 20 mL ultrapure water, the solution deaerated with nitrogen and a voltammogram recorded in the same manner as in the iron determination (without catechol addition). Fig. 12 shows a good, Fig. 13 a poor, decomposed pipes buffer.


Fig. 12



## Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry

Remarks  
(continued)

Fig. 13




- The catechol solution should be prepared using pure, undecomposed catechol. If this solution has a brownish colour from the very start, filtration (fluted filter) is usually useful to enable the solution to be kept for 3...4 days.

Recrystallization from hot toluene followed by vacuum drying is better. Solutions prepared from catechol purified in this manner can be stored for up to 20 days.

Catechol solutions undergoing decomposition and which are thus unfit for use show the following behaviour: the peaks are shifted and become smaller. Furthermore, the standard additions are no longer linear and the baseline is highly distorted. Fig. 14 shows an example of such a curve.

Fig. 14



**Simultaneous trace determination of 7 metal ions in "electronic grade" materials with the aid of stripping voltammetry****Literature**

- ▶ C.M.G. Van den Berg, Zi Quiang Huang  
*Determination of iron in sea water using cathodic stripping voltammetry preceeded by adsorptive collection with the hanging mercury drop electrode*  
J. Electroanal. Chem. 177, (1984) 269-280
  
- ▶ P. Ostapczuk, M. Goedde, M. Stöppler, H.W. Nürnberg  
*Kontroll- und Routinebestimmung von Zn, Cd, Pb, Cu, Ni und Co mit differenzieller Pulsvoltammetrie in Materialien der deutschen Umweltbank*  
Fresenius, Z. Anal. Chem. 316, (1984) 252-256
  
- ▶ H. Braun, M. Metzger  
*Umweltanalytische Nickel-Bestimmung durch Adsorptionsvoltammetrie mit der Quecksilberfilmelektrode*  
Fresenius, Z. Anal. Chem. 318, (1984) 321-326
  
- ▶ C.M.G. Van den Berg  
*Determining trace concentration of copper in water by cathodic film stripping voltammetry with adsorptive collection*  
Anal. Lett. 17A, (1984) 2141-2157
  
- ▶ C.M.G. Van den Berg, Zi Quiang Huang  
*Direct electrochemical determination of dissolved vanadium in sea water by cathodic stripping voltammetry with hanging mercury drop electrode*  
Anal. Chem. 56, (1984) 2383-2386
  
- ▶ Axel Meyer, Rolf Neeb  
*Determination of cobalt and nickel by adsorption voltammetry in supporting electrolytes containing triethanolamine and dimethylglyoxime*  
Fresenius, Z. Anal. Chem. 315, (1983) 118-120