

Application Bulletin 192/2 e

Determination of thiourea in the lower mg/L and in the µg/L range by polarography and cathodic stripping voltammetry

Summary

Thiourea forms highly insoluble compounds with mercury. The resulting anodic waves are used for the polarographic determination of thiourea. For the analysis of very small quantities (µg/L), cathodic stripping voltammetry (CSV) is used. Differential Pulse measuring mode is used in both cases.

Instruments

VA instrument
capable of operating a Multi-Mode Electrode and supporting differential pulse (DP) measuring mode

Electrodes

| | | |
|----|--|------------|
| WE | Multi-Mode Electrode pro | 6.1246.120 |
| | Mercury drop capillary | 6.1226.030 |
| RE | Ag/AgCl reference electrode | 6.0728.x20 |
| | Ag/AgCl/KCl (3 mol/L) | |
| | Electrolyte vessel Filled with c(KCl) = 3 mol/L | 6.1245.010 |
| AE | Pt rod electrode | 6.0343.x00 |

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Sodium hydroxide, for trace analysis*, CAS 1310-73-2
- Thiourea, for analysis, CAS 62-56-6
- Ultrapure water, resistivity >18 MΩ·cm (25 °C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

Electrolyte solution c(NaOH) = 2 mol/L

Standard solutions

Thiourea standard stock solution β(Thiourea) = 100 mg/L in ultrapure water.
More dilute solutions must be freshly prepared daily from the stock solution.

Method 1: Polarographic determination of thiourea concentrations 0.2 – 2 mg/L

Analysis

10 mL sample solution and 10 mL electrolyte solution are deaerated well with nitrogen in a polarographic vessel. The concentrations are determined by standard addition. After every standard addition purge well with nitrogen.

Measuring solution

10 mL sample solution
+ 10 mL electrolyte solution

Parameters

Voltammetric

| | |
|--------------------------|-------------------------|
| Electrode operating mode | DME |
| Measuring mode | DP – Differential pulse |
| Stirring rate | 2000 min ⁻¹ |


Potentiostatic pretreatment

| | |
|--------------------|-----|
| Equilibration time | 5 s |
|--------------------|-----|

Sweep

| | |
|---------------------|-----------|
| Start potential | -0.45 V |
| End potential | -0.15 V |
| Potential step | 0.008 V |
| Potential step time | 1.0 s |
| Sweep rate | 0.008 V/s |
| Pulse amplitude | 0.05 V |

| | |
|--------------------------|----------|
| Substance | |
| Name | Thiourea |
| Characteristic potential | -0.26 V |

Example

Result

| | |
|-------------|-----------|
| β(thiourea) | 0.31 mg/L |
|-------------|-----------|

Method 2: Stripping voltammetric determination of thiourea concentrations 5 – 60 µg/L

Analysis

10 mL sample solution and 10 mL electrolyte solution are deaerated well with nitrogen in a polarographic vessel. The concentrations are determined by standard addition. After every standard addition purge well with nitrogen.

Measuring solution


10 mL sample solution
+10 mL electrolyte solution

Parameters


| Voltammetric | |
|-----------------------------|-------------------------|
| Electrode operating mode | HMDE |
| Drop size | 9 |
| Measuring mode | DP – Differential pulse |
| Stirring rate | 2000 min ⁻¹ |
| Potentiostatic pretreatment | |
| Potential 1 | -0.2 V |
| Waiting time 1 | 150 s |
| Equilibration time | 20 s |
| Sweep | |

| | |
|---------------------|-----------|
| Start potential | -0.5 V |
| End potential | -0.9 V |
| Potential step | 0.008 V |
| Potential step time | 1.0 s |
| Sweep rate | 0.008 V/s |
| Pulse amplitude | 0.05 V |

| | |
|--------------------------|----------|
| Substance | |
| Name | Thiourea |
| Characteristic potential | -0.75 V |

Example

Result (using a calibration curve)

| | |
|-------------|-----------|
| β(thiourea) | 16.5 µg/L |
|-------------|-----------|

Calibration curve (stripping voltammetry)

Comments

- The largest possible mercury drops should be used
- Chloride ions up to a content of 10 mg/L do not interfere with the stripping voltammetry.

- If the sample contains more than 2 mg/L thiourea, it must be diluted for the polarographic determination (non-linear region).
- More dilute electrolyte solutions or smaller drop sizes restrict the working range. (Neutralize acidic sample solutions beforehand.)

References

- Artl C., Naumann R.
Vorschlag zur Bestimmung des Selens in Trinkwasser
Z. Anal. Chem 282(1976), 463
- Stara V., Kopanica M.
Adsorptive stripping voltammetric determination of thiourea and its thiourea derivatives. Anal. Chim. Acta 159, (1984), 105-110
- Kirchnerova J., Purdy C.
A new simple voltammetric method for thiourea and thiourea dioxide determinations. Anal. Lett. 13/12, (1980), 1031-1040
- Smyth M. R., Osteryoung J. G.
Determination of some thiourea-containing pesticides by pulse voltammetric methods of analysis. Anal. Chem. 49, (1977), 2310-2314

Appendix

Method print for the polarographic determination of thiourea (method 1)

Thioharnstoff ppm-Bereich
MPL 1 EL.TYPE MME

METHOD 18 PAGE 3
OPERATION SEQUENCE

OPERATIONS/PARAMETERS

```

1 STIR ;PURGE ;          300 s
2 ØPURGE;ØSTIR ;
3 CADDL ;
4 STIR ;PURGE ;          60 s
5 ØSTIR ;ØPURGE;         5 s
6 DME ;MEAS ;            5 s
6a M.MODE      DPN      50 mV
6b T.STEP      1.0 s
6c U.SET       -450 mV
7 SWP Ø ;      37 s
7a U.END       -150 mV
7b U.STEP      8 mV
   SW.RATE     8.0 mV/ s
8 ØMEAS ;ADD1J2;         5 s
9 BEEP ;END ;
  
```

Method print for the stripping voltammetric determination of thiourea (method 2)

Thioharnstoff ppb-Bereich
MPL 1 EL.TYPE MME

METHOD 30 PAGE 3
OPERATION SEQUENCE

OPERATIONS/PARAMETERS

```

1 STIR ;PURGE ;          300 s
2 BEEP ;HOLD ;
3 (REP. ;
4 PURGE ;STIR ;          30 s
5 ØPURGE;              10 s
6 HMDE ;MEAS ;         120 s
6a M.MODE      DPN      -50 mV
6b T.STEP      1.0 s
6c U.SET       -200 mV
7 ØSTIR ;           30 s
8 MEAS ;
8a M.MODE      DPN      -50 mV
8b T.STEP      1.0 s
8c U.SET       -500 mV
9 SWP Ø ;           50 s
9a U.END       -900 mV
9b U.STEP      8 mV
  
```

OPERATIONS/PARAMETERS

```

   SW.RATE     8.0 mV/ s
10 ØMEAS ;REP> 2;
11 BEEP ;END ;
  
```