

## Application Bulletin 223/2 e

# Fully automated determination of uranium

### Branch



Energy, power plant

### Keywords

Titration; uranium; Davies and Gray; automation; branch 16

### Summary

This Bulletin describes the automated determination of uranium according to the Davies and Gray method: Uranium(VI) is reduced with iron(II) to uranium(IV) in a highly concentrated phosphoric acid medium (1). The excess of iron(II) is then oxidized to iron(III) with nitric acid in the presence of a molybdenum catalyst (2). The nitrous acid formed in this reaction is destroyed by sulfamic acid (3) and finally the uranium(IV) is titrated with potassium dichromate in the presence of a vanadium catalyst (4).



### Instruments

- Sample changer
- Titrator with DET mode

### Electrodes

Combined Pt-ring electrode	6.0451.100
----------------------------	------------

### Reagents

- Type 1 water (ASTM) (ultrapure water)
- Iron(II)sulfate heptahydrate,  $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ , analytical reagent grade
- Sulfuric acid,  $w(\text{H}_2\text{SO}_4) = 98\%$ ,  $\Phi(\text{H}_2\text{SO}_4) = 50\%$  (v/v),  $c(\text{H}_2\text{SO}_4) = 1 \text{ mol/L}$ ,  $c(\text{H}_2\text{SO}_4) = 0.1 \text{ mol/L}$
- Sulfamic acid,  $\text{NH}_2\text{HSO}_3$ , laboratory reagent grade
- Potassium dichromate,  $\text{K}_2\text{Cr}_2\text{O}_7$ , analytical reagent grade

- Ortho-phosphoric acid,  $w(\text{H}_3\text{PO}_4) = 85\%$ , analytical reagent grade
- Ammonium molybdate  $[(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}]$
- Nitric acid,  $w(\text{HNO}_3) = 65\%$ ,  $c(\text{HNO}_3) = 10 \text{ mol/L}$
- Vanadyl sulfate monohydrate,  $\text{VOSO}_4 \cdot \text{H}_2\text{O}$ , analytical reagent grade
- Ammonium fluoride,  $\text{NH}_4\text{F}$ , analytical reagent grade

### Solutions

Iron(II)sulfate solution $\beta(\text{FeSO}_4 \cdot 7\text{H}_2\text{O}) = 280 \text{ g/L}$	100 mL $w(\text{H}_2\text{SO}_4) = 98\%$ is added to 600 mL ultrapure water. Then $280 \text{ g} \pm 0.1 \text{ g}$ of $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ is added to the diluted sulfuric acid. After dissolving the $\text{FeSO}_4$ , the solution is left to cool down to room temperature and is made up to 1000 mL. This solution may be used for 7 days.
Sulfamic acid solution $\beta(\text{NH}_2\text{HSO}_3) = 150 \text{ g/L}$	$150 \pm 0.1 \text{ g}$ sulfamic acid is dissolved in 700 mL ultrapure water (this may require gentle heating). After the solution is cooled down to room temperature it is diluted to 1 L. This solution may be used for up to one month.
Titratant $\beta(\text{K}_2\text{Cr}_2\text{O}_7) = 12.5 \text{ g/L}$	$25 \pm 0.1 \text{ g}$ $\text{K}_2\text{Cr}_2\text{O}_7$ is weighed out and dissolved in 1 L ultrapure water then diluted to 2 L. Shake to mix each time before use. A freshly prepared solution should be left standing for at least 24 h hours. Shake to mix each time before use. This solution may be used for up to a month.
$\beta(\text{K}_2\text{Cr}_2\text{O}_7) = 0.6 \text{ g/L}$	$0.30 \pm 0.01 \text{ g}$ $\text{K}_2\text{Cr}_2\text{O}_7$ is weighed out and dissolved in 200 mL ultrapure water then diluted to 500 mL. This solution may be used for up

	to six weeks.
Phosphoric acid solution	50 mL $\beta(\text{K}_2\text{Cr}_2\text{O}_7) = 0.6 \text{ g/L}$ is added to 2.5 L ortho-phosphoric acid and mixed thoroughly. This solution may be used for up to one month.
Oxidizing reagent	4 ± 0.1 g ammonium molybdate is dissolved in 400 mL ultrapure water. 500 mL $c(\text{HNO}_3) = 65\%$ and 100 mL of sulfamic acid solution are added. The solution is mixed thoroughly before cooled down to room temperature and made up to 1 L with ultrapure water. This solution may be used for up to seven days.
Vanadyl sulfate solution	$\beta(\text{VOSO}_4 \cdot \text{H}_2\text{O}) = 2.0 \text{ g/L}$ 2.0 ± 0.01 g vanadyl sulfate is weighed out into a 1 L volumetric flask and made up to the mark with $c(\text{H}_2\text{SO}_4) = 1 \text{ mol/L}$ . This solution may be used for up to 1 month.
Ammonium fluoride $\beta(\text{HN}_4\text{F}) = 400 \text{ g/L}$	400 ± 1 g ammonium fluoride is weighed out and dissolve in ultrapure water. Then diluted to 1 L and mixed thoroughly.
Reducing solution	900 mL iron(II)sulfate solution is mixed with 300 mL sulfamic acid solution and 300 mL $\Phi(\text{H}_2\text{SO}_4) = 50\%$ (v/v). This solution may be used for up to three days.

### Sample preparation

For the appropriate sample preparation, see the referenced standards.

### Analysis

The prepared sample is dissolved in 10 mL of  $c(\text{H}_2\text{SO}_4) = 0.1 \text{ mol/L}$  (the solution may need to be warmed). 2 mL of  $\beta(\text{HN}_4\text{F}) = 400 \text{ g/L}$  is added to the sample solution before the addition of 50 mL phosphoric acid solution. After a mixing time of 10 s, 8 mL of reducing solution is added and the sample is stirred for 150 s. Then, 10 mL oxidizing reagent is added and the solution is stirred for 3 min. After that, 50 mL vanadyl sulfate solution is added and the

solution is stirred for 120 s. The solution is then titrated with  $\beta(\text{K}_2\text{Cr}_2\text{O}_7) = 12.5 \text{ g/L}$  until after the first equivalence point.

### Parameters

Mode	MET U
Start volume	17 mL
Dosing rate	15 mL/min
Pause	10 s
Signal drift	40 mV/min
Max. waiting time	32 s
Volume increment	0.05 mL
Dosing rate	50 mL/min
Filling rate	30 mL/min
EP criterion	30 mV
EP recognition	greatest

### Calculation

For the calculation of the results, see the referenced standards.

### Example determination




Fig. 1: Titration curve of an  $\text{UF}_6$  pellet (blue = titration curve, red = dMV)

### Comments

- The combined Pt-ring electrode is stored in  $c(\text{KCl}) = 3 \text{ mol/L}$  (6.2308.020) when not in use.
- At the beginning of a series, it is recommended to place the electrode in concentrated ammonia for 30 min. This ensures that the ceramic pin diaphragm is clean and the electrode is ready for use.

**References**

- Davies, W., Gray, W., Talanta, 11, (1964), 1203
- ISO 7097-1:2004  
Nuclear fuel technology -- Determination of uranium in solutions, uranium hexafluoride and solids -- Part 1: Iron(II) reduction/potassium dichromate oxidation titrimetric method
- ASTM C1267-11  
Standard Test Method for Uranium by Iron (II) Reduction in Phosphoric Acid Followed by Chromium (VI) Titration in the Presence of Vanadium

**Author**

Competence Center Titration

Metrohm International Headquarters