

# IC Application Note No. C-32

**Title:** Noise reduction in non-suppressed ion chromatography using 732 IC Detector – the benefit of the Metrohm auto zero / full scale approach

**Range:** Absolute Range of conductivity measurement (e.g. 0 ... 1000  $\mu\text{S}/\text{cm}$ ).

**Full Scale:** Relative conductivity measurement window. Zero function activated on the Background conductivity.

Example: Background = 712.5  $\mu\text{S}/\text{cm}$ ; 10  $\mu\text{S}/\text{cm}$  Full Scale (corresponding to a 100 fold signal amplification)


Analog Output	absolute Cond.	relative Cond.
0 mV	712.5 $\mu\text{S}/\text{cm}$	0.0 $\mu\text{S}/\text{cm}$
500 mV	717.5 $\mu\text{S}/\text{cm}$	5.0 $\mu\text{S}/\text{cm}$
1000 mV	722.5 $\mu\text{S}/\text{cm}$	10.0 $\mu\text{S}/\text{cm}$

**Zero:** Electronic background compensation of actual conductivity at the time of activation as well as a reduction of electronic noise and drift.

**Column:** 6.1010.000 Metrosep Cation 1-2

**Eluent:** 4 mmol/L tartaric acid, 1 mmol/L dipicolinic acid

**Flow:** 1.0 mL/min



<b>Description:</b>	Duration of the measurement: 4 min.
<b>trace A)</b>	Baseline measured with Range = 1000 $\mu\text{S}/\text{cm}$ and Full Scale = 1000 $\mu\text{S}/\text{cm}$ . Zero Function not activated. Measured noise level: 45 nS/cm
<b>trace B)</b>	same as A) but Zero function activated. Measured noise level: 22 nS/cm The use of the zero function reduces the noise by a factor of 2.
<b>trace C)</b>	Baseline measured with Range = 1000 $\mu\text{S}/\text{cm}$ . Full Scale is set to 10 $\mu\text{S}/\text{cm}$ . This corresponds to an amplification factor of 100 compared to signals A) and B). Zero function is activated. Measured noise level: 2 nS/cm The use of the Zero function together with a lower Full Scale (amplification) yields in far less baseline noise. In the actual setup the noise is reduced by a factor of 20 compared to the direct measurement (signal A).

<b>Summary:</b>	The use of an appropriate Full Scale together with the Zero Function of the 732 IC Detector minimises the Baseline Noise drastically. Much lower detection limits are achieved.
-----------------	---