


# IC Application Note No. C-99

**Title:** Lanthanides in a rock sample by Ion Chromatography applying non-suppressed conductivity detection

**Summary:** Determination of traces of gadolinium, samarium, neodymium, cerium and lanthanum using cation chromatography with direct conductivity detection applying Metrohm Inline Filtration.

**Sample:** Digested rock sample  
**Sample Preparation:** Nitric acid digestion, direct injection after Inline Filtration

**Column:** 6.1007.000 Nucleosil 100-5-SA  
**Eluent:** 12.0 mmol/L 2-hydroxyisobutyric acid  
 2.0 mmol/L ethylene diamine, pH = 4.6 (ammonia)  
**Flow:** 1.0 mL/min  
**Injection Volume:** 200  $\mu$ L



<b>Lanthanides:</b>	Gadolinium mg/kg	Samarium mg/kg	Neodymium mg/kg	Cerium mg/kg	Lanthanum mg/kg
	<b>78.0</b>	<b>60.6</b>	<b>147.1</b>	<b>368.4</b>	<b>6161</b>