

Thermo. Titr. Application Note No. H-022

Title:	Determination of Moisture in Lubricating Oils with TEOF
---------------	--

Scope:	Determination of moisture in lubricating oils with TEOF (triethyl orthoformate)
---------------	---

Principle:	A sample of oil is dissolved in a solvent mixture containing an acid catalyst, and titrated with a solution of triethyl orthoformate in cyclohexane to an endothermic endpoint.
-------------------	---



Reagents:	<p>1.2 mol/L triethyl orthoformate (TEOF)(CAS 122-51-0, Riedel de Haën cat. No. 60388, M=148.20, d=0.89, purity=99%) in dry (HPLC grade) cyclohexane. Pipette 200mL TEOF into a 1000mL volumetric flask, and make to volume with cyclohexane. (Burette #1)</p> <p>Solvent mixture. Mix 800mL dry (HPLC grade) cyclohexane, 150mL dry (HPLC grade) iso-propanol, 100mL dry chloroform, and 20mL methane sulfonic acid 99.5%+ (CAS 75-75-2, Aldrich cat. no. 471356) and store in a tightly stoppered bottle. (Burette #2) This mixture is hygroscopic, and must be protected from atmospheric moisture</p>
------------------	---

