

Thermo. Titr. Application Note No. H-028

Title: Standardization of 0.1 mol/L KOH
in propan-2-ol

Scope: Standardization of 0.1 mol/L in propan-2-ol for use in applications for the determination of weakly acidic species in non-aqueous media.

Principle: A range of volumes of a standard solution of benzoic acid in propan-2-ol are titrated with ~0.1 mol/L KOH in propan-2-ol to a catalytically determined thermometric endpoint.


Reference:
1. M. J. D. Carneiro, M. A. Feres Júnior, and O. E. S. Godinho. Determination of the acidity of oils using paraformaldehyde as a thermometric end-point indicator. *J. Braz. Chem. Soc.* **13** (5) 692-694 (2002)

Reagents: ~0.1 mol/L KOH in propan-2-ol
Paraformaldehyde
Standard benzoic acid. Weigh accurately ~0.5g benzoic acid (e.g., Sigma-Aldrich cat. No. B-3250, minimum 99.5% purity, F.W. = 122.1), and transfer to a 200mL dry volumetric flask. Dissolve in and make to volume with A.R. propan-2-ol.
A.R. propan-2-ol

Method:	Basic Experimental Parameters:	
	Titrant delivery rate (mL/min.)	2
	No. of exothermic endpoints	1
	Data smoothing factor	40
	Stirring speed (802 stirrer)	5
	Delay before start (secs.)	3
	Procedure: Prepare and titrate solutions of standard benzoic acid according to the following:	
	mL benzoic acid soln.	mL propan-2-ol
	20	10
	15	15
	10	20
	5	25
	Prior to titrating, add ~0.5g of paraformaldehyde to each solution, and titrate to a thermometric endpoint with ~0.1 mol/L KOH in propan-2-ol. Calculate the amount of benzoic acid in each aliquot as mmole, and plot mmole benzoic acid (x-axis) against mL ~0.1 mol/L KOH (y-axis). Determine the gradient of the calibration curve. The strength of the KOH solution is computed from the reciprocal of the gradient.	

Example:		
	Mmole benzoic acid	Titre, mL ~0.1 mol/L KOH
	0.4064	4.212
	0.3048	3.173
	0.2032	2.141
	0.1016	1.117
	Gradient = 10.15476	
	Molarity = $1/10.15476 = 0.09848$ mol/L	

Calibration Plot:



Thermometric Titration Plot:

