

# Thermo. Titr. Application Note No. H-070

**Title:** Determination of Ferric and Cupric Ions in Copper Refining Solutions

**Scope:** Determination of  $\text{Fe}^{3+}$  and  $\text{Cu}^{2+}$  in copper refining solutions by thermometric titration. It was found that the conventional approach of masking  $\text{Fe}^{3+}$  to permit the iodometric determination of  $\text{Cu}^{2+}$  is not possible in some copper refining solutions.

**Principle:**  $\text{Fe}^{3+}$  content is determined by titration with fluoride (refer to AN H-069). The combined  $\text{Fe}^{3+}$  and  $\text{Cu}^{2+}$  content is determined by iodometric titration. The  $\text{Cu}^{2+}$  content is computed by subtraction.

**Reagents:**

1.  $\text{Fe}^{3+}$  determination.

*Titant:* 1mol/L standard NaF solution  
*Combined acetate buffer:* Dissolve 130.9g anhydrous potassium acetate and 54.7g anhydrous sodium acetate in 500mL DI water. Add 115mL glacial acetic acid, and make to 1L with DI water. Alternatively, dissolve 164g anhydrous sodium acetate and 75g potassium chloride in 700mL DI water, add 115mL glacial acetic acid and make to 1L with DI water.

2.  $\text{Fe}^{3+} + \text{Cu}^{2+}$  determination.

*Titant:* 1mol/L standard  $\text{Na}_2\text{S}_2\text{O}_3$  solution

- Glacial acetic acid
- 50% w/v KI solution (store in amber bottle in a cool place).
- 0.04mol/L  $\text{KIO}_3$  solution (for standardizing  $\text{Na}_2\text{S}_2\text{O}_3$  titrant)

**Method:**

*Basic Experimental Parameters:*

1.  $\text{Fe}^{3+}$  determination.

Titant delivery rate (mL/min.)	4
No. of exothermic endpoints	1
Data smoothing factor (DSF)	70
Stirring speed (802 stirrer)	10

Iron must be in  $\text{Fe}^{3+}$  form, and sufficiently acidic to prevent hydrolysis of the  $\text{Fe}(\text{H}_2\text{O})_6^{3+}$  aquo ion. Dispense aliquot into titration vessel. Add 10mL combined acetate buffer and make to approximately 30mL with DI water. Titrate to an exothermic endpoint with 1mol/L NaF solution.

*Standardization of NaF titrant.* This may be standardized against standard Al solution prepared from high purity Al metal.

**2.  $\text{Fe}^{3+} + \text{Cu}^{2+}$  determination.**

*Basic Experimental Parameters:*

Titrant delivery rate (mL/min.)	4
No. of exothermic endpoints	1
Delay start of titration (secs.)	20
Data smoothing factor (DSF)	60
Stirring speed (802 stirrer)	10

Iron must be in  $\text{Fe}^{3+}$  form, and sufficiently acidic to prevent hydrolysis of the  $\text{Fe}(\text{H}_2\text{O})_6^{3+}$  aquo ion. Dispense aliquot into titration vessel. Add 2mL glacial acetic acid. Fit titration vessel to titration head and start the analysis sequence. Add 10mL KI solution through a port in the titration head immediately after clicking the "Start" button.

*Standardization of  $\text{Na}_2\text{S}_2\text{O}_3$  titrant.* Pipette aliquots of 5, 10, 15, 20 and 25mL  $\text{KIO}_3$  solution into titration vessels. Add 2mL glacial acetic acid, and make to ~30mL with DI water. Start the titration, and add 10mL KI solution through a port in the titration head immediately after clicking the "Start" button. Plot mmole of  $\text{KIO}_3$  (x-axis) against mL  $\text{Na}_2\text{S}_2\text{O}_3$  titrant (y-axis) and compute the titrant molarity.

**Examples:**

Solutions from copper refinery operation, containing  $\text{Fe}^{3+}$ ,  $\text{Fe}^{2+}$  and  $\text{Cu}^{2+}$

	Sample no.	$\text{Cu}^{2+}$ g/L	$\text{Fe}^{3+}$ g/L
	1	7.50, 7.48	4.89, 4.91
	2	24.35, 24.29	6.44, 6.45
	3	3.35, 3.31	6.23, 6.22
	4	61.23, 61.68	11.55, 11.50
	5	3.46, 3.54	2.43, 2.37

**Calculation Procedure:**

- |  |   |
|--|---|
| <p>1. From fluoride titration, calculate <math>Fe^{3+}</math> g/L</p>  | <p>(1)</p> $Fe^{3+} \text{ g/L} = \frac{((\text{Titre, mL} - \text{blank, mL}) \times \text{NaF mol/L} \times 55.845)}{(\text{sample vol, mL} \times 6)}$                   |
| <p>2. From iodometric titration, calculate (<math>Cu^{2+} + Fe^{3+}</math>), expressed as <math>Fe^{3+}</math> g/L</p> | <p>(2)</p> $(Cu^{2+} + Fe^{3+}) \text{ g/L} = \frac{((\text{Titre, mL} - \text{blank, mL}) \times \text{Na}_2\text{S}_2\text{O}_3 \times 55.845)}{(\text{sample vol, mL})}$ |
| <p>3. Subtract (2) from (1)</p>  |   |
| <p>4. Convert (3) to <math>Cu^{2+}</math> g/L</p>  | $Cu^{2+} \text{ g/L} = (2) - (1) \times 63.546/55.845$  |