


Determination of the water content in tooth paste with MATi 10



This Application Note describes the automated determination of the water content in tooth paste using volumetric Karl Fischer titration (MATi 10).

Method description

Sample

Tooth paste

Sample preparation

Between 30 and 50 mg of sample are weighed into the titration beaker. After adding a magnetic stirring bar, the titration beaker is tightly closed with Al-foil and foil holder and placed on the rack of the 814 Robotic Sample Processor.

Electrodes

Double Pt electrode	6.0340.000
---------------------	------------

Reagents

HYDRANAL®-Water Standard 10.0	Fluka 34849
-------------------------------	-------------

HYDRANAL®-Composite 5	Fluka 34805
-----------------------	-------------

HYDRANAL®-Methanol dry	Fluka 34741
------------------------	-------------

HYDRANAL®-Formamide dry	Fluka 34724
-------------------------	-------------

Solvents

Methanol/formamide mixture	7 : 3 (v/v)
----------------------------	-------------

Instrumentation

MATi 10 system

Analysis

System preparation

To prepare the system, a blank value is run and the result discarded.

Blank determination

Three blank values are determined using empty titration beakers and 20 mL of the methanol/formamide mixture. The mean value is saved as Common Variable and subtracted from the endpoint volumes of all subsequent determinations (titer and sample).

Sample determination

Due to the high water content small sample sizes between approximately 30 and 50 mg are used for the sample determinations. 20 mL of the methanol/formamide mixture are added and the determination subsequently carried out.

Parameters

Except for the calculation, the parameters for all methods are identical.

The following parameters are changed compared to the default values.

Pause 1	8 s
Extraction time	30 s

Results

Blank determination

Mean / [mL] (n = 3)	RSD / [%]
1.471	2.34

Sample determination

Mean / [%] (n = 6)	RSD / [%]
55.82	0.27

Comments

Depending on the tooth paste sample, it might be necessary to adapt the method.