


Anhydrosugars besides sugar alcohols and sugars collected on air filters



Due to residential wood burning, concentrations of levoglucosan, mannosan, and galactosan are typically elevated during the winter months. By contrast, during the summer months, an increased contribution of primarily biological sugar components can be expected. When analyzing the saccharidic tracers collected on air filters, optimal separation and sensitivity are achieved on the Metrosep Carb 2 - 150/4.0 column applying pulsed amperometric detection.

Results

Compound	Conc. [$\mu\text{g/L}$]	Compound	Conc. [$\mu\text{g/L}$]
1 Inositol	10.0	7 Mannosan	50.0
2 Erythritol	10.0	8 Galactosan	50.0
3 Arabitol	20.0	9 Rhamnose	50.0
4 Sorbitol	20.0	10 Glucose	50.0
5 Mannitol	20.0	11 Xylose	50.0
6 Levoglucosan	50.0	12 Sucrose	50.0

Sample

Standard solution

Sample preparation

none

Columns

Metrosep Carb 2 - 150/4.0	6.1090.420
Metrosep Carb 2 Guard/4.0	6.1090.500
Metrosep CO ₃ Trap 1 - 100/4.0	6.1015.300

Solutions

Eluent	10.0 mmol/L sodium hydroxide
--------	------------------------------

Parameters

Flow rate	1.0 mL/min
Injection volume	100 µL
Recording time	38 min
Column temperature	45 °C

PAD Parameters

Cell	Wall-Jet cell
Working electrode	Gold (3 mm)
Reference electrode	Palladium
Spacer	50 µm
Measuring potential	50 mV
Meas. range	200 µA
Meas. duration	100 ms
Cycle duration	550 ms
Temperature	32 °C
Mode	PAD

Analysis

Pulsed amperometric detection

Instrumentation

930 Compact IC Flex Oven	2.930.2100
IC Amperometric Detector	2.850.9110
889 IC Sample Center	2.889.0010
IC equipment Wall-Jet cell: anion analysis (Au, Pd)	6.5337.010

