

# IC Application Note No. S-277

**Title:** Hexafluorosilicate besides standard anions

**Summary:** Determination of fluoride, chloride, nitrate, phosphate, sulfate, silicate and hexafluorosilicate (calculated) using anion chromatography with conductivity detection after chemical suppression and subsequent UV/VIS detection with post-column reaction (see AN U-48). Hexafluorosilicate is hydrolyzed into fluoride and silicate. Both anion concentrations can be used for the calculation of the  $\text{SiF}_6^{2-}$  concentration.

**Sample:** Standard solution

**Sample Preparation:** –

**Column:** 6.1006.630 Metrosep A Supp 7 – 250


**Temperature:** 45 °C

**Eluent:** 3.6 mmol/L sodium carbonate

**Suppressor:** Metrohm Suppressor Module MSM (50 mmol/L  $\text{H}_2\text{SO}_4$ )

**Flow:** 0.7 mL/min

**Injection Volume:** 20  $\mu\text{L}$



<b>Results:</b>	Fluoride mg/L	Chloride mg/L	Bromide mg/L	Nitrate mg/L	Phosphate mg/L	Sulfate mg/L
Standard anions	5.0	5.0	5.0	5.0	5.0	5.0
Hexafluorosilicate (5 mg/L)	3.98	–	–	–	–	–

See AN U-48 for silicate determination.