


Determination of glycolate and lactate in varnish remover



Glycolate and lactate have to be determined in a dual phase varnish remover. Analyzed is only the upper aqueous phase. The separation is achieved on a Metrosep A Supp 16 - 250/4.0 column. The eluent composition is adapted to get a sufficient separation of glycolate and lactate without interference by formate and acetate. Conductivity detection after sequential suppression is applied.

Results

Anion	Conc. [mg/L]	RT [min]	Recovery [%]	Anion	Conc. [mg/L]	RT [min]	Recovery [%]
1 Fluoride	-	6.4	-	4 Acetate	-	9.3	-
2 Glycolate	n.d.	7.6	101	5 Lactate	17.9	9.8	104
3 Formate	-	8.7	-	6 Chloride	-	12.1	-

Concentrations in the diluted solution. Small picture shows the achieved separation under the applied condition.

Sample

Varnish remover

Sample preparation

Aqueous phase is diluted 1:5000 with ultrapure water.

Columns

Metrosep A Supp 16 - 250/4.0	6.1031.430
Metrosep RP 2 Guard/3.5	6.1011.030

Solutions

Eluent	3.75 mmol/L sodium carbonate 0.375 mmol/L sodium hydroxide
Suppressor regenerant	100 mmol/L sulfuric acid
Rinsing solution	STREAM

Analysis

Conductivity detection after sequential suppression

Instrumentation

930 Compact IC Flex Oven/SeS/PP/Deg	2.930.2560
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0020
MSM Rotor A	6.2832.000
Adapter sleeve for Suppressor Vario	6.2842.020

Parameters

Flow rate	0.8 mL/min
Injection volume	20 µL
P _{max}	20 MPa
Recording time	14 min
Column temperature	45 °C

