


Chloride and sulfate impurities in potassium bicarbonate



Within the scope of the modernization of USP, chloride and sulfate are determined as impurities in potassium hydrogen carbonate (bicarbonate). USP41 monograph for potassium bicarbonate does not check for chloride and sulfate. Applying ion chromatography with conductivity detection after sequential suppression allows quantifying these impurities.

Results

Sample	Conc.** [mg/L]	Corr. Coef. NLT = 0.99	Conc. Calib. Level 1 [mg/L]	RSD [%; N = 6] NMT ¹ = 5 %	Signal/Noise NLT ² = 20
1 Chloride	0.041*	0.9989	0.050	0.9	185
2 Sulfate	0.038*	0.9995	0.050	0.8	69

* Below lowest calibration level ** mg/L in solution corresponds to mg/kg in solid sample

¹ NMT = not more than ² NLT = not less than

Sample

Potassium hydrogen carbonate, USP grade

Sample preparation

0.1 g of sample dissolved in 100 mL of ultrapure and CO₂ free water.

Columns

Metrosep A Supp 16 - 150/4.0	6.1031.420
Metrosep A Supp 16 Guard/4.0	6.1031.500

Solutions

Eluent	7.5 mmol/L sodium carbonate 0.75 mmol/L sodium hydroxide
Suppressor regenerant	500 mmol/L sulfuric acid
Rinsing	STREAM

Analysis


Conductivity detection after sequential suppression

Instrumentation

930 Compact IC Flex Oven/SeS/Deg	2.930.2460
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0020
MSM Rotor A	6.2832.000
Adapter sleeve for Suppressor Vario	6.2842.020
IC equipment: Dosino Regeneration	6.5330.190

Parameters

Flow rate	0.8 mL/min
Injection volume	50 µL
P _{max}	20 MPa
Column temperature	45 °C
Run time	22 min
Dosino regeneration	0.1 mL for 9 min



Calibration

Level	Chloride, sulfate [mg/L]
Level 1	0.05
Level 2	0.10
Level 3	0.50
Level 4	0.75
Level 5	1.00
Level 6	1.50

www.metrohm.com

 **Metrohm**