
## Titration Application Note T-143

# Cobalt analysis by automated photometric titration



Cobalt can be analyzed in aqueous solutions by direct titration with EDTA at pH 9. The indicator is murexide, and the equivalence point is determined with the Optrode at a wavelength of 574 nm.



# Method description

## Sample

Aqueous solution containing cobalt (0.05 mol/L)

#### **Sample preparation**

No sample preparation is required.

## Configuration

	907 Titrando	2.907.0020
	815 Robotic USB Sample Processor XL	2.815.0020
	786 Swing head	2.786.0040
	Swing arm	6.1462.070
	Titration head	6.1458.010
	Sample rack 28 x 200 mL	6.2041.830
	800 Dosino, 3 x	2.800.0010
	802 Stirrer	2.802.0020
	5 mL Dosing unit	6.3032.150
	10 mL Dosing unit	6.3032.210
	50 mL Dosing unit	6.3032.250
	Disposable PP sample beaker, 200 mL	6.1459.310
	Optrode	6.1115.000

#### Solutions

EDTA solution	c(Na <sub>2</sub> EDTA) = 0.1 mol/L If possible this solution should be bought from a supplier.
Murexide	0.2 g murexide (1:100 in NaCl) is dissolved in 50 mL deion. water.

#### Analysis (Optrode 574 nm)

5.0~mL sample solution is pipetted into a 200 mL plastic beaker and 90 mL deion. water is added. After the addition of 5 mL buffer pH 9 and 5 mL murexide indicator solution the cobalt is titrated with c(Na<sub>2</sub>EDTA) = 0.1 mol/L until after the endpoint.

#### **Parameters**

Mode	MET U
Pause	30 s
Stirring rate	8
Signal drift	50 mV/min
Min. waiting time	5 s
Max. waiting time	26 s
Volume increment	0.1 mL
EP criterion	15 mV
EP recognition	Greatest
Stop volume	10 mL

#### Results

Mean result (n = 5)

Cd content / (g/L)	3.205
s(rel) / %	1.51

www.metrohm.com

