


Titration Application Note T-150

Nickel determination using automated photometric titration



Nickel can be determined in alkaline media. To visualize the endpoint, murexide is used as indicator. The equivalence point is detected with the Optrode at a wavelength of 574 nm.

Method description

Sample

Aqueous solution of nickel (0.05 mol/L)

Sample preparation

No sample preparation is required.

Configuration

907 Titrand	2.907.0020
815 Robotic USB Sample Processor XL	2.815.0020
786 Swing head	2.786.0040
Swing arm	6.1462.070
Titration head	6.1458.010
Sample rack 28 x 200 mL	6.2041.830
800 Dosino, 4 x	2.800.0010
802 Stirrer	2.802.0020
5 mL Dosing unit	6.3032.150
10 mL Dosing unit, 2 x	6.3032.210
50 mL Dosing unit	6.3032.250
Disposable PP sample beaker, 200 mL	6.1459.310
Optrode	6.1115.000

Solutions

EDTA solution	$c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ If possible this solution should be bought from a supplier.
Murexide	0.2 g murexide (1:100 in NaCl) is dissolved in 50 mL deion. water
Buffer pH 10	54 g NH_4Cl and 300 mL $w(\text{NH}_3) = 25\%$ is given into a 1000 mL volumetric flask and filled up to the mark with deion. water.

Analysis

5 mL sample solution is pipetted into a 200 mL plastic beaker and 90 mL deion. water is added. After the addition of 5 mL buffer pH 10 and 5 mL murexide indicator solution the nickel is titrated with $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ until after the equivalence point.

Parameters

Mode	MET U
Stirring rate	8
Signal drift	20 mV/min
Min. waiting time	0 s
Max. waiting time	38 s
Volume increment	0.05 mL
EP criterion	15 mV
EP recognition	Greatest
Stop volume	10 mL
Mode	MET U

Results

Mean results (n = 6)

Ni content / (g/L)	2.95
s(rel) / %	0.36