


Determination of the bromine number in petroleum products and aliphatic olefins according to ASTM D1159



The bromine number is an important parameter for the determination of aliphatic C=C double bonds in petroleum products. The bromine number is usually determined using electrochemical titration at 5 °C, where the bromine is generated in situ from a bromide/bromate solution. For the titration, a solvent mixture of glacial acetic acid, methanol, and chloroform is used. In this Application Note, the toxic chloroform was replaced with diethyl carbonate.

Method description

Samples

Cyclohexene

Transformer oil spiked with 6.6 g/L cyclohexene

Sample preparation

Based on the expected bromine number a known amount of sample is dissolved in 50 mL diethyl carbonate.

Analysis

Blank

The blank is determined the same way as the sample, just with diethyl carbonate instead of sample.

Sample

110 mL titration solvent and 5 mL sample solution are pipetted into a titration vessel. While stirring, the solution is cooled down to 0–5 °C. The solution is then titrated with $c(\text{Br}_2) = 0.25 \text{ mol/L}$ until after the equivalence point.

Configuration

| | |
|--|-------------|
| OMNIS Advanced Titrator with stirrer | 2.1001.0220 |
| OMNIS Dosing Module | 2.1003.0010 |
| OMNIS 50 mL cylinder unit, 1x (titration solvent) | 6.03001.250 |
| OMNIS 10 mL cylinder unit, 1x (titrant) | 6.03001.210 |
| Analog measuring module | 6.02101.010 |
| Lauda RE 304 circulation thermostat bath | - |
| Lauda E 300 immersion thermostat | - |
| Electrode cable plug-in head G (pol.) / plug P, 0.55 m | 6.02104.040 |
| Electrode cable plug-in head G (temp.) / plug P, 0.55 m | 6.02104.020 |
| Cable MDL PL/SO 1 m | 6.02102.020 |
| Stirring bar / 30 mm | 6.1903.060 |
| 3-way stopper with antidiffusion valve | 6.1543.210 |
| Titration vessel with thermostat jacket / 50-150 mL | 6.1418.250 |
| Titration vessel lid automation | 6.1414.080 |
| OMNIS Stand-alone license (including one instrument license) | 6.06003.010 |
| Double Pt-wire electrode for coulometry | 6.0341.100 |

| | |
|---------------------------|------------|
| Pt1000 Temperature Sensor | 6.1110.100 |
|---------------------------|------------|

Solutions

| | |
|-------------------|---|
| Titration solvent | <p>$c(\text{Br}_2) = 0.25 \text{ mol/L}$</p> <p>$w(\text{KBr}) = 5.1\%$ and $w(\text{KBrO}_3) = 1.392\%$ in deionized water, if possible this solution should be bought from a supplier</p> <p>$\Phi(\text{diethyl carbonate}) = 13.4\%$, (v/v)</p> <p>$\Phi(\text{methanol}) = 13.4\%$, (v/v)</p> <p>$\Phi(w(\text{H}_2\text{SO}_4) = 20\%) = 1.8\%$, (v/v)</p> <p>$\Phi(\text{glacial acetic acid}) = 71.4\%$, (v/v)</p> <p>134 mL $w(\text{toluene}) \geq 99.3\%$, 134 mL $w(\text{methanol}) \geq 99.8\%$ and 18 mL $w(\text{H}_2\text{SO}_4) = 20\%$ are added into a 1000 mL volumetric flask. The flask is then filled up to the mark with $w(\text{glacial acetic acid}) \geq 99.8\%$.</p> |
|-------------------|---|

Parameters

| | |
|-------------------|-------------------|
| Mode | MET Ipol |
| Pause | 60 s |
| Start volume | 2 mL |
| Stirring rate | 6 |
| Volume increment | 0.05 mL |
| Signal drift | 30 mV/min |
| Min. waiting time | 0 s |
| Max. waiting time | 32 s |
| Dosing rate | Maximum |
| I(pol) | 1.0 μA |
| Stop volume | 10 mL |
| Stop EP | 1 |
| Volume after EP | 0.5 mL |
| EP criterion | 200 mV |
| EP recognition | All |

Results

| Sample | Result / (n = 8) | s(rel) |
|-----------------|-----------------------|--------|
| Cyclohexene | 193.57 g bromine/100g | 0.28% |
| Transformer oil | 12.82 g bromine/100g | 0.32% |

www.metrohm.com

 **Metrohm**