


Iodate and iodide in an electroplating bath applying direct UV/VIS detection



The analysis of iodate and iodide in used electroplating baths is a demanding task due to the high content of other ions. Iodate is used as a stabilizer for the bath and needs to be checked for proper electroplating. The use of a sodium chloride eluent, the Metrosep A Supp 5 - 250/4.0 column and direct UV/VIS allows the analysis of these samples without interference.

Results

Diluted sample	spiked	found	Recovery
Iodate	2.0 mg/L	1.9 mg/L	95%
Iodide	4.0 mg/L	43.2 mg/L	101%

Sample

Fresh resp. used nickel electroplating bath

Sample preparation

Dilution 1:20

Columns

Metrosep A Supp 5 - 250/4.0	6.1006.530
-----------------------------	------------

Metrosep A Supp 4/5 Guard/4.0	6.1006.500
-------------------------------	------------

Solutions

Eluent	10 g/L sodium chloride
--------	------------------------

Parameters

Flow rate	0.7 mL/min
-----------	------------

Injection volume	20 µL
------------------	-------

P _{max}	15 MPa
------------------	--------

Recording time	25 min
----------------	--------

Column temperature	55 °C
--------------------	-------

Wavelength	230 nm
------------	--------

Analysis

Direct UV/VIS detection

Instrumentation

940 Professional IC Vario ONE	2.940.1100
-------------------------------	------------

944 Professional UV/VIS Detector	2.944.0010
----------------------------------	------------

858 Professional Sample Processor	2.858.0020
-----------------------------------	------------

