


# Sulfur speciation in mining leachate by ion chromatography applying a perchlorate eluent and UV/VIS detection



In gold mining there is a tendency to switch from cyanide leaching to the much less toxic thiosulfate leaching process. Thiosulfate leaching is a sensitive process that requires more optimization of the components of the leach reaction to maximize gold recovery and reagent loss. Sulfite, thiosulfate, thiocyanate, and tetrathionate are separated on a Metrosep A Supp 5 - 250/4.0 column. Perchlorate is chosen as an eluent as most of the metal perchlorates are soluble in water. This avoids metal precipitation in the IC System.

## Results

Anions	Concentration [mg/kg]	RSD [%, n = 3]
1 Sulfite	< 1	-
2 Thiosulfate	58.8	2.2
3 Thiocyanate	49.7	1.9
4 Tetrathionate	28.6	1.4

### Sample

Mining leachate

### Sample preparation

Dilution 1:2 in ultrapure water, injection after filtration (0.45 µm).

### Columns

Metrosep A Supp 5 - 250/4.0	6.1006.530
Metrosep A Supp 4/5 - Guard/4.0	6.1006.500

### Solutions

Eluent	10 mmol/L sodium perchlorate 1.0 mmol/L sodium hydroxide
--------	---

### Analysis


UV/VIS detection

### Parameters

Flow rate	0.7 mL/min
Injection volume	20 µL
P <sub>max</sub>	15 MPa
Recording time	32 min
Column temperature	35 °C
Wavelength	215 nm

### Instrumentation

940 Professional IC Vario ONE	2.940.1100
944 Professional UV/VIS Detector Vario	2.944.0010
858 Professional Sample Processor	2.858.0020



[www.metrohm.com](http://www.metrohm.com)

 **Metrohm**