

VA Application Note No. V- 2

Title:	Chromium, manganese and titanium in a polyterephthalic acid solution
---------------	---

Summary:	Determination of Cr, Mn and Ti in a PTA solution containing HCl.
-----------------	--

Sample:	PTA in HCl
Sample Preparation:	none


Chromium:	
Electrolyte:	NaAc, DTPA (diethylenetriaminepentaaceticacid), NaNO ₃ . pH = 6.2 with NaOH.
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPASV (-50 mV), HMDE $U_{\text{meas}} = -1000 \text{ mV (60s)}$, $U_{\text{start}} = -1000 \text{ mV}$, $U_{\text{end}} = -1500 \text{ mV}$, Ep (Cr) = -1250 mV

Manganese:	
Electrolyte:	NH ₄ Cl / NH ₃ buffer, borate buffer, Zn standard
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-75 mV), HMDE $U_{\text{meas}} = -1700 \text{ mV (90s)}$, $U_{\text{start}} = -1620 \text{ mV}$, $U_{\text{end}} = -1250 \text{ mV}$, Ep (Mn) = -1400 mV


Titanium:	
Electrolyte:	Mandelic acid. pH = 3 with NH ₃ 25%.
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-50 mV), HMDE $U_{\text{meas}} = -100 \text{ mV (30s)}$, $U_{\text{start}} = -100 \text{ mV}$, $U_{\text{end}} = -900 \text{ mV}$, Ep (Ti) = -750 mV

Results:	Cr μg/L	Mn μg/L	Ti μg/L
	7.2	324.8	91.8

Determination of chromium



Determination of manganese



Determination of titanium

