

VA Application Note No. V- 15

Title:	Nickel, antimony, cadmium, thallium and copper in a neutral, highly concentrated Zn solution
---------------	---

Summary:	Determination of Ni, Sb, Cd, Tl and Cu in a neutral, highly concentrated zinc solution from the plating industry
-----------------	--

Sample:	Zn electrolyte 130-170g/L, Mn 5-10g/L, Mg 5-10g/L
Sample Preparation:	none

Nickel:	
Electrolyte:	Dimethylglyoxime (DMG), NH_4Cl buffer, ascorbic acid, NH_3 , pH = 9.9
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-50 mV), HMDE $U_{\text{meas}} = -750 \text{ mV}$ (30s), $U_{\text{start}} = -750 \text{ mV}$, $U_{\text{end}} = -1100 \text{ mV}$ $E_p(\text{Ni}) = -1000 \text{ mV}$


Antimony:	
Electrolyte:	HCl
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPASV (+50 mV), HMDE $U_{\text{meas}} = -250 \text{ mV}$ (180s), $U_{\text{start}} = -250 \text{ mV}$, $U_{\text{end}} = -50 \text{ mV}$ $E_p(\text{Sb}) = -150 \text{ mV}$

Cadmium, thallium and copper:	
Electrolyte:	HCl, KMnO_4 , ascorbic acid
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPASV (+50 mV), HMDE $U_{\text{meas}} = -800 \text{ mV}$ (60s), $U_{\text{start}} = -800 \text{ mV}$, $U_{\text{end}} = 0 \text{ mV}$ $E_p(\text{Cd}) = -580 \text{ mV}$, $E_p(\text{Tl}) = -450 \text{ mV}$, $E_p(\text{Cu}) = -120 \text{ mV}$

Results:	Ni µg/L	Sb µg/L	Cd µg/L	Cu µg/L	Tl µg/L
	494	5	7.2	< 1	208.8

Determination of nickel

Determination of antimony



Determination of cadmium, thallium and copper

