

VA Application Note No. V- 59

Title:	1-Methyl-nicotinamide hydrochloride in standard
---------------	--

Summary:	Determination of 1-methyl-nicotinamide hydrochloride in standard using Na_2CO_3 as electrolyte
-----------------	--

Sample:	Standard
Sample Preparation:	none

1-Methyl-nicotinamide hydrochloride:	
Electrolyte:	Na_2CO_3
AE:	GC
RE:	Ag/AgCl/KCl 3M
Parameters:	DPPOL (-40 mV), DME $U_{\text{start}} = -1300 \text{ mV}$, $U_{\text{end}} = -1900 \text{ mV}$ E_p (1-methyl-nicotinamide hydrochloride) = - 1550 mV

Results:	1-Methyl-nicotinamide hydrochloride g/L
	1.29

Determination of 1-methyl-nicotinamide hydrochloride

