

VA Application Note No. V - 103

Title:	Chromium in lime (CaCO_3)
---------------	--

Summary:	Cr(VI) is determined at the HMDE in an electrolyte containing ethylenediamine and acetate. Because Cr(III) is electrochemically inactive, all Cr has to be oxidised prior to analysis.
-----------------	--

Sample:	lime (CaCO_3)
Sample preparation:	<p>2 g lime are dissolved in 7 mL of $w(\text{HNO}_3) = 65\%$, suprapur and filled up to 100 mL with ultrapure water.</p> <p>1 mL lime solution + 10 mL water + 1 drop $c(\text{KMnO}_4) = 0.02 \text{ mol/L}$.</p> <p>The mixture is boiled for 10 min. KMnO_4 solution has to be added to keep the pink color. The volume has also kept constant at approx. 10 mL by adding water.</p> <p>Before cooling to room temperature the pH is adjusted to 7-9 with NaOH.</p>

Analysis of Cr(VI)

Electrolyte

Measuring solution

10 mL oxidised sample solution
+ 10 μL ethylene diamine
+ 150 μL $w(\text{acetic acid}) = 100\%$, suprapur
+ 200 μL $w(\text{NH}_3) = 25\%$, suprapur
adjust pH to 6.8 with acetic acid

Auxiliary electrode (AE) Pt

Reference electrode (RE) Ag/AgCl/KCl (3 mol/L)

Parameters	Working electrode	HMDE
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	600 s
	Deposition potential	no deposition
	Deposition time	0 s
	Equilibration time	10 s
	Pulse amplitude	50 mV
	Start potential	70 mV
	End potential	-170 mV
	Voltage step	6 mV
	Voltage step time	0.2 mV
	Sweep rate	30 mV/s
	Peak potential Cr(VI)	-30 mV

Results:	Cr(VI)
	5.5 µg/g

Determination of Cr(VI)

