

VA Application Note No. V - 112

Title: Thiourea in nickel plating baths

Summary: Thiourea is determined by cathodic stripping voltammetry (CSV) at the HMDE in ammonia buffer at pH 8.9. Chloride in the sample does not interfere with this determination.

Sample: Ni plating bath containing chloride

Sample preparation: none

Analysis of thiourea

Electrolyte

Acetate solution:

c(sodium acetate, suprapur) = 1 mol/L

Perchlorate solution:

c(NaClO₄) = 2.5 mol/L

Ammonia buffer:

c(HCl) = 1 mol/L

+ c(NH₃) = 3 mol/L

Measuring solution

10 mL ultrapure water

+ 10 µL sample

+ 250 µL acetate solution

+ 1 mL perchlorate solution

adjust the pH to exactly 8.9 with ammonia buffer

Auxiliary electrode (AE)

Pt

Reference electrode (RE)

Ag/AgCl/KCl (3 mol/L)

Bridge electrolyte: NaCl (3 mol/L)

Parameters

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	+100 mV
Deposition time	30 s
Equilibration time	5 s
Pulse amplitude	50 mV

Start potential	+100 mV
End potential	-800 mV
Voltage step	4 mV
Voltage step time	0.1 s
Sweep rate	40 mV/s
Peak potential thiourea	-370 mV

Results:	thiourea
	6.8 mg/L

Determination of thiourea

