

VA Application Note No. V - 118

Title: Gold in ammonium thiosulfate solution

Summary: Gold can be determined by anodic stripping voltammetry (ASV) in the $\mu\text{g/l}$ range at the Ultra Trace Graphite electrode. The solution should not contain halide ions.

Sample: ammonium thiosulfate solution

Sample preparation: 2 mL sample and are heated in a porcelain cup to approx. 120°C and evaporated to dryness. Then it was ashed in a muffle furnace at 650°C for 2 hours. The residue was mixed with 2 mL $w(\text{HNO}_3) = 65\%$ suprapur and 1 mL $w(\text{H}_2\text{SO}_4) = 96\%$ suprapur and heated until SO_3 vapors are evolved. After cooling once more 2 mL of $w(\text{HNO}_3) = 65\%$ suprapur were added and evaporated almost to dryness. After cooling 25 mL ultrapure water were added and the solution was boiled for 1 min. 1 mL of $c(\text{EDTA-Na}_2) = 0.1\text{ mol/L}$ was added to the cold solution. The solution was filled up with water to 50 mL.

Analysis of Au		
Electrolyte	$c(\text{KNO}_3) = 0.2\text{ mol/L}$ + $c(\text{Na}_2\text{-EDTA}) = 0.04\text{ mol/L}$	
Measuring solution	7.5 mL ultrapure water + 1 mL digestion solution + 2.5 mL electrolyte The pH is adjusted to 3.8 with $c(\text{NaOH}) = 2\text{ mol/L}$	
Auxiliary electrode (AE)	Glassy Carbon (GC)	
Reference electrode (RE)	Ag/AgCl/KCl (3 mol/L) Bridge electrolyte : KNO_3 sat.	
Parameters	Working electrode	Ultra Trace RDE
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Deposition potential	-400 mV
	Deposition time	120 s

Deposition potential (without stirring)	-600 mV
Deposition time without stirring = Equilibration time	25 s
Pulse amplitude	50 mV
Start potential	-400 mV
End potential	+500 mV
Voltage step	4 mV
Voltage step time	0.1 s
Sweep rate	40 mV/s
Cleaning potential	+500 mV
Cleaning time	60 s
Peak potential Au	+150 mV

Results:	Au
	702 µg/L

Determination of Au

