

VA Application Note No. V - 154

Title:	Nitrilotriacetic acid (NTA) in a cyanidic gold bath
---------------	------------------------------------------------------------


Summary:	NTA in a cyanidic gold bath is determined as Bi-NTA complex by polarography. For standard addition a Bi-NTA standard solution is used.
-----------------	----------------------------------------------------------------------------------------------------------------------------------------

Sample:	Cyanidic Au bath
Sample preparation:	Dilution 1:2000

Analysis of NTA		
Supporting electrolyte	β (ascorbic acid) = 40 g/L	
Bi nitrate solution	β (Bi ³⁺) = 2 g/L	
Measuring solution	10 mL supporting electrolyte + 50 μ L Bi nitrate solution + 100 μ L diluted Au bath	
Working electrode (WE)	MME (Multi Mode Electrode) 6.1246.020	
Auxiliary electrode (AE)	Pt 6.0343.000	
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020 Intermediate electrolyte: c(KCl) = 3 mol/L 6.1245.010	
Parameters	Working electrode	DME
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Equilibration time	10 s
	Pulse amplitude	0.05 V
	Start potential	0.1 V
	End potential	-0.4 V
	Voltage step	0.006 V
	Voltage step time	0.6 s
	Sweep rate	0.01 V/s
	Peak potential Bi-NTA	-0.2 V

Results:	NTA
	410.4 g/L

Determination of NTA



NTA
 $c = 410.406 \text{ g/L}$
 $\pm 1.395 \text{ g/L (0.34\%)}$

