

VA Application Note No. V - 158

Title: Indium in a tin bath

Summary: The concentration of In in a Sn bath is determined in a HCl / Urotropin[®] containing electrolyte by anodic stripping voltammetry (ASV). The determination is linear up to approx. 0.5 mg/L with respect to the concentration of In in the measuring vessel. The standard addition solution is also prepared with HCl and Urotropin[®].

Sample: Acid Sn bath
Sample preparation: None

Analysis of In

HCl solution w(HCl) = 30%

Urotropin[®] solution c(Urotropin[®]) = 2 mol/L
Urotropin[®]: Hexamethylenetetramine

Measuring solution 10 mL H₂O
 + 0.05 mL Sn bath
 + 3 mL HCl solution
 + 3 mL Urotropin[®] solution

Working electrode (WE) MME (Multi Mode Electrode) 6.1246.020

Auxiliary electrode (AE) Pt 6.0343.000


Reference electrode (RE) Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020
 Intermediate electrolyte: c(KCl) = 3 mol/L 6.1245.010

Parameters

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	-0.68 V
Deposition time	5 s
Equilibration time	5 s
Pulse amplitude	0.05 V
Start potential	-0.68 V
End potential	-0.53
Voltage step	0.002 V
Voltage step time	0.3 s
Sweep rate	0.0067 V/s
Peak potential In	-0.60 V

Results:	In
	62.0 mg/L

Determination of In



In
 c = 62.021 mg/L
 +/- 1.056 mg/L (1.70%)

