

VA Application Note No. V - 160

Title: Palladium in an activator bath

Summary: The concentration of Pd in activator bath is determined by polarography in ammonium chloride electrolyte.


Sample: Activator solution

Sample preparation: None

Analysis of Pd		
Supporting electrolyte	c(NH ₄ Cl) = 1 mol/L c(NH ₃) = 1 mol/L c(KCl) = 1.5 mol/L	
Measuring solution	15 mL H ₂ O + 5 mL supporting electrolyte + 0.5 mL activator solution	
Working electrode (WE)	MME (Multi Mode Electrode)	6.1246.020
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: c(KCl) = 3 mol/L	6.1245.010
Parameters	Working electrode	DME
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Equilibration time	5 s
	Pulse amplitude	0.05 V
	Start potential	-0.4 V
	End potential	-0.9 V
	Voltage step	0.005 V
	Voltage step time	0.6 s
	Sweep rate	0.0083 V/s
	Peak potential Pd	-0.71 V

Results:	Pd
	64.1 mg/L

Determination of



Pd
c = 64.067 mg/L
+/- 0.287 mg/L (0.45%)

