

VA Application Note No. V - 165

Title: Zinc in a phosphatation bath

Summary: The concentration of Zn in a zinc phosphatation bath is determined by polarography in ammonia buffer pH 9.3.


Sample: Zinc phosphatation bath

Sample preparation: Dilution 1:100

Analysis of Zn		
Ammonia buffer pH 9.3	c(NH ₄ Cl) = 1 mol/L c(NH ₃) = 2 mol/L pH 9.3 adjusted with NH ₃ or HCl	
Measuring solution	10 mL H ₂ O + 5 mL ammonia buffer pH 9.3 + 50 µL diluted phosphatation bath	
Working electrode (WE)	MME (Multi Mode Electrode)	6.1246.020
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: c(KCl) = 3 mol/L	6.1245.010
Parameters	Working electrode	DME
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Equilibration time	5 s
	Pulse amplitude	0.05 V
	Start potential	-1.05 V
	End potential	-1.4 V
	Voltage step	0.006 V
	Voltage step time	0.6 s
	Sweep rate	0.01 V/s
	Peak potential Zn	-1.23 V

Results:	Zn
	18.1 g/L

Determination of Zn



Zn
c = 18.061 g/L
+/- 0.019 g/L (0.10%)

