

# VA Application Note No. V - 181

**Title:** Determination of chromium(VI) in cement

**Summary:** The concentration of Cr(VI) in cement is determined in tartrate electrolyte after acid extraction of the sample.

**Sample:** cement

**Sample preparation:** 10 g cement sample is weighed into a beaker and elutriated with 50 mL nitric acid ( $c(\text{HNO}_3) = 1 \text{ mmol/L}$ ,  $\text{pH} = 3$ ). The mixture is stirred for 1 h and then filtrated through a filter paper. The residue in the filter is rinsed twice with nitric acid. The filtrate is transferred into a 100 mL volumetric flask and made up to the mark with ultrapure water.

## Analysis of Cr(VI)


<b>Ammonium tartrate solution</b>	$c((\text{NH}_4)_2\text{C}_4\text{H}_4\text{O}_6) = 0.1 \text{ mol/L}$ pH 9 adjusted with $\text{NH}_3$	
<b>Measuring solution</b>	10 mL ammonium tartrate solution + 2 mL extraction solution	
<b>Working electrode (WE)</b>	<b>MME</b> (Multi Mode Electrode)	6.1246.020
<b>Auxiliary electrode (AE)</b>	<b>Pt</b>	6.0343.000
<b>Reference electrode (RE)</b>	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: $c(\text{KCl}) = 3 \text{ mol/L}$	6.1245.010

## Parameters

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	0
Deposition time	0
Equilibration time	5 s
Pulse amplitude	0.05 V
Start potential	0 V
End potential	-0.4 V
Voltage step	0.003 V
Voltage step time	0.3 s
Sweep rate	0.01 V/s
Peak potential Cr(VI)	-0.2 V

<b>Results:</b>	Cr(VI)
	17.9 mg/kg

### Determination of Cr(VI)



Cr VI  
 c = 17.903 mg/kg  
 +/- 0.296 mg/kg (1.66%)

