


Determination of dimethylamine in meropenem applying sequential suppression



Meropenem is used to treat a wide variety of bacterial infections. It is known as a carbapenem-type antibiotic, and, as a β -lactam antibiotic, it works by inhibiting the cell wall synthesis of bacteria. Dimethylamine is a precursor in meropenem production and thus occurs as an impurity in the final product. It is determined on a Metrosep C Supp 1 - 250/4.0 column with conductivity detection after sequential suppression.

Results

Cation	Concentration [mg/kg]
Dimethylamine	116

Sample

Meropenem sample

Sample preparation

50 mg sample dissolved in 10 mL diluent and injected after filtration (0.2 µm)

Columns

Metrosep C Supp 1 - 250/4.0	6.1052.430
-----------------------------	------------

Metrosep C Supp 1 Guard/4.0	6.1052.500
-----------------------------	------------

Solutions

Eluent	3.0 mmol/L nitric acid 50 µg/L rubidium 15% acetonitrile
--------	--

Suppressor regenerant	70 mmol/L sodium carbonate 70 mmol/L sodium hydrogen carbonate 20% acetonitrile
-----------------------	---

Rinsing solution	20% acetonitrile in ultrapure water
------------------	-------------------------------------

Diluent	10 mmol/L nitric acid
---------	-----------------------

Analysis

Conductivity detection after sequential suppression

Instrumentation

940 Professional IC Vario ONE/SeS/PP	2.940.1500
--------------------------------------	------------

IC Conductivity Detector	2.850.9010
--------------------------	------------

858 Professional Sample Processor	2.858.0020
-----------------------------------	------------

MSM-HC Rotor C	6.2842.200
----------------	------------

Parameters

Flow rate	1.0 mL/min
-----------	------------

Injection volume	20 µL
------------------	-------

P _{max}	15 MPa
------------------	--------

Recording time	28 min
----------------	--------

Column temperature	40 °C
--------------------	-------

