

Determination of sodium with the ion-selective electrode

Of interest to:

General analytical laboratories; Water analysis; Pharmaceutical industry; Food analysis

A 1, 2, 4, 7, 12, 16

Summary

The determination of sodium with the sodium ISE represents a selective, rapid, accurate and favorably priced method, which is described in this Bulletin. Examples are used to show how determinations can be carried out with the 692 pH/Ion Meter using either direct measurement or the standard addition technique. The sodium concentration has been determined in standard solutions, water samples (tap water, mineral water, waste water), foodstuffs (spinach, baby food) and urine.

The construction, working principles and areas of application of the two Metrohm ion-selective sodium electrodes – the 6.0501.100 glass membrane ISE and the 6.0508.100 polymer membrane ISE – are explained in detail.

Theory

The activity of ions can be determined by means of ion-selective electrodes. In very dilute solutions the activity corresponds approximately to the concentration of the respective ion.

The fact that in an ideal case ISEs only respond to a particular type of ion means that it is possible to specify the present form of bonding of an ion and to determine it selectively (examples: selective determination of a particular oxidation state of an element; determination of free metal ions in the presence of complexly bound metal ions).

However, ion-selective electrodes that have no cross-sensitivity towards ions that are chemically similar do not exist. ISEs can therefore not be said to be specific.

The use of the two ion-selective sodium electrodes

The measuring range of the glass membrane ISE lies between 1×10^{-5} and 1 mol/L Na^+ (corresponds to $0.23 \text{ mg/L} \dots 22.99 \text{ g/L Na}^+$). Even lower Na^+ concentrations can be determined with the polymer membrane ISE: down to $5 \times 10^{-6} \text{ mol/L}$ (0.11 mg/L).

The 6.0501.100 glass membrane ISE is primarily used for samples with a difficult matrix (fruit and vegetable juice; samples containing chlorate, perchlorate and/or chloroacetate ions).

The 6.0508.100 polymer membrane ISE is recommended for the analysis of water samples (ultrapure, rain, ground, surface and drinking water) and for acidic

samples containing fluoride, fluoroborate or fluorosilicate ions.

6.0501.100 glass membrane ISE

Construction of the electrode

The electrode shaft and membrane are made of glass. The selectivity for Na^+ ions is determined by the composition of the glass membrane.

Composition and conditioning of the membrane

The glass membrane (composition: $\text{Li}_2\text{O-Al}_2\text{O}_3\text{-SiO}_2$) is immersed in water before use; this generates a gel-like hydrated silicic acid layer.

Working principle

In aqueous solution, cations are removed from the outer hydrated layer and replaced by Na^+ ions from the sample. This creates a potential at the boundary layer between the sample solution and the membrane that depends on the activity of sodium.

ISA solutions

ISA is the abbreviation for Ionic Strength Adjustor. This is a solution with a high ionic strength that is used to dilute the sample solutions while simultaneously «fixing» their ionic strength. As a result of the almost constant ionic strength, the activity coefficient $\gamma_{\text{meas. ion}}$ remains virtually unchanged, even when measurements are made in very different sample solutions.

Apart from fixing the ionic strength, ISA solutions can also adjust the pH value at the same time. Such solutions are known as **TISAB: Total Ionic Strength Adjustment Buffer**.

In order to remove interfering ions, these sample-conditioning solutions can additionally contain complexing agents. It is thus possible to carry out highly selective measurements even with an electrode that has a relatively low selectivity.

The choice of one of these solutions depends on the particular application. For sodium determinations with the glass membrane ISE it is necessary to use a TISAB solution that also buffers pH values above 9.

6.0508.100 polymer membrane ISE

Construction of the electrode

The electrode consists of a plastic shaft (PVC), a polymer membrane and an inner Ag/AgCl system. It is filled with an internal buffer containing NaCl.

Composition and conditioning of the membrane

The membrane is made of PVC, plasticizer and ionophore + additive.

It is not necessary to immerse the PVC membrane in water before use, i.e. the polymer membrane ISE is immediately ready for use.

Working principle

The ionophore within the polymer membrane is responsible for the selectivity, i.e. it is able to take up Na⁺ ions selectively.

The molecular structure of the ionophore is adapted to suit the coordination number and size of the sodium ion. During the measurement, Na⁺ ions diffuse from the sample solution into the polymer membrane. This results in a concentration gradient and therefore a membrane potential.

The polymer has lipophilic properties, i.e. water or ions with a hydration sphere cannot penetrate into the membrane. In aqueous solution, sodium ions are also surrounded by a hydration sphere. However, as a result of the gain in energy this is removed upon inclusion in the ionophore.

ISA solution

The addition of an ISA solution for «fixing» the ionic strength or adjusting the pH value is not necessary if the polymer membrane ISE is used in a concentration range below 50 ppm Na⁺. At higher concentrations it is the sample matrix that determines whether an ISA solution is required. pH adjustment with a TISAB solution is not necessary. However, we recommend that the polymer membrane ISE is chiefly used in a pH range between 3 and 12, as longer periods of use outside this range can shorten the service life of the electrode.

Direct measurement with calibration

The ionic concentration of the sample is interpolated using a calibration curve. This curve is obtained by measuring standard solutions. The expected ionic concentration should be approximately in the middle of the concentration range of the standard solutions.

Standard addition

For the standard addition technique, a known sample volume is spiked with a defined quantity of the ion to be determined (possibly in several steps). The ionic concentration of the sample is then calculated automatically from the resulting potential differences between the sample solution alone and the spiked solutions.

The volume of the added standard solution should not exceed 25% of the sample volume. The concentration of the standard solution should be as high as possible (so that dilution effects can be neglected).

There are three types of standard addition: *manual*, *auto dos* and *auto*. The quickest method is the automatic standard addition *auto*. However, please note that in order to obtain accurate results, the chosen potential difference ΔU must be at least 12 mV per addition and that at least three standard additions must be performed (ΔU_{total} at least 30 mV). If exactly defined volumes of the standard additions are required but maximum ease of operation is also desirable, choose *auto dos* and enter the individual volumes (see Instructions for use).

Instruments and accessories

2.692.0010	pH/Ion Meter
2.765.0030	Dosimat Titration Stand
6.3014.213	10 mL Exchange Unit
6.2138.010	Connection cable 692–765
6.0501.100	Sodium-selective glass membrane electrode
6.0508.100	Sodium-selective polymer membrane electrode
6.2104.020	Cable for Na ISE
6.0726.100	Double-junction Ag/AgCl reference electrode (for measurements with the glass membrane ISE) Inner electrolyte: c(KCl) = 3 mol/L Bridge electrolyte: TISAB solution (see reagents)
6.0733.100	LL Ag/AgCl reference electrode (for measurements with the polymer membrane ISE) Reference electrolyte: c(KCl) = 3 mol/L
6.2106.020	Cable for reference electrode
6.1110.100	Pt 1000 temperature sensor
6.2104.080	Cable for temperature sensor
6.1415.220	Titration vessel
6.1414.010	Titration vessel lid
6.2125.060	Connection cable 692–PC
6.2125.010	Adapter cable
6.6008.500	VESUV 3.0 Metrodata software

Reagents

a) Measurements with the glass membrane ISE

Standard solution	$\beta(\text{Na}^+) = 2000 \text{ mg/L}$
TISAB solution	$c(\text{tris}(\text{hydroxymethyl})\text{-aminomethane}) = 1 \text{ mol/L}$; pH value adjusted to 8 ... 10 with $c(\text{HNO}_3) = 2 \text{ mol/L}$.
Conditioning solution	$c(\text{NaCl}) = 1 \text{ mol/L}$

b) Measurements with the polymer membrane ISE

Standard solution 1	$\beta(\text{Na}^+) = 2106.9 \text{ mg/L}$
Standard solution 2	$\beta(\text{Na}^+) = 207.42 \text{ mg/L}$
ISA solution (depending on application)	$c(\text{CaCl}_2) = 1 \text{ mol/L}$

Electrode handling

The following points must be observed in order to achieve optimal results.

a) Glass membrane ISE

- We recommend that the electrode is shaken like a clinical thermometer prior to the first measurement.
- The glass membrane ISE should be stored dry if it is out of use for longer periods of time. For short-term storage $c(\text{NaCl}) = 0.1 \text{ mol/L}$ is suitable.
- Before being used for the first time, and after longer periods out of use, the glass membrane must be activated/conditioned by overnight immersion in $c(\text{NaCl}) = 1 \text{ mol/L}$.
- With the reference electrode care should be taken that the internal parts are well wetted. Any air bubbles formed when the electrolyte is filled in must be removed.
- After each measurement the electrode must be thoroughly rinsed with dist. water and wiped off with a slightly moist tissue.
- Potassium interferes if present in 500-fold excess compared to sodium.
- As under acidic conditions this Na ISE responds to H_3O^+ ions like a pH glass electrode, measurements must be carried out under alkaline conditions (pH = 8 ... 10).
- Other ions that interfere are listed in the Instructions for use of the electrode.

b) Polymer membrane ISE

- It is not necessary to shake the polymer membrane ISE before it is used for the first time; the electrode can be used immediately.
- If the electrode is not to be used for some time it should be stored dry.
- A long conditioning period for the membrane before use is not necessary with this Na ISE. Briefly rinsing the electrode with dist. water and drying it with

a soft paper tissue prior to the measurement is completely adequate.

- Avoid touching the polymer membrane with your fingers.
- In contrast to the glass membrane ISE it is not necessary to work with a double-junction reference electrode. A bridge electrolyte (TISAB solution for the glass membrane ISE) is not required. $c(\text{KCl}) = 3 \text{ mol/L}$ as reference electrolyte yields excellent results as regards accuracy and precision.
- There must be no air bubbles on the membrane surface before and during the measurements.
- The polymer membrane ISE should neither be cleaned with nor stored in organic solvents (acetone, ethanol).
- Use at increased temperatures ($>30 \text{ }^\circ\text{C}$) should be avoided if possible.
- The electrode should not be used for longer periods of time at extreme pH values, as solutions whose pH value is not between 3 and 12 shorten the service life of the polymer membrane ISE.
- If a series of determinations is to be carried out then the electrode should be stored in dist. water between the measurements instead of just rinsing and drying it. In this way any Na^+ ions that may have diffused into the membrane during the previous determination can be rinsed out of it again. The working procedure described yields better results with regard to accuracy and precision.

General information

- All solutions must be stored in plastic vessels.
- All standard additions must be performed with stirring (stirrer parameter 692: *ON* or *control.*). Take care that the stirrer is switched on and that there is a stirrer bar in the titration vessel. Additions without stirring result in incorrect measurements!
- Method selection:
Standard addition is a rapid and reliable method, particularly for an undefined sample matrix. A direct measurement is advisable whenever many determinations have to be carried out one after another in similar solutions with a defined and unproblematic sample matrix.
- The addition of TISAB solution before measurements with the glass membrane ISE has the effect that the ionic strength and the pH value of the sample solutions remain constant.
- If the automatic standard addition *auto* is used for solutions with low ionic strength, the obtained potential steps can slightly differ from the defined parameter *delta U*. However, this does not affect the accuracy of the analytical results.
- Choice of the standard concentration:
In order to ensure a correct evaluation of the standard addition by the 692 pH/Ion Meter, the standard concentrations c_{std} for the different buret volumes V_{buret} must be chosen as a function of the

sample concentration c_{smpl} according to the following table:

$V_{\text{buret}}/\text{mL}$	$C_{\text{std}} : C_{\text{smpl}}$
5	40 : 1
10	20 : 1
20	10 : 1
50	5 : 1

If the sample has been greatly diluted with ISA or TISAB solution this ratio must be taken into account.

Example:

Sample concentration c_{smpl}	5 mg/L
Buret volume V_{buret}	10 mL
Sample volume	10 mL
Volume ISA/TISAB	10 mL
Total volume V_{total}	20 mL
Factor $C_{\text{std}} : C_{\text{smpl}}$ from the table	20

This yields a sample concentration of 2.5 mg/L in the titration vessel. The optimal concentration of the standard is therefore $2.5 \text{ mg/L} \times 20 = 50 \text{ mg/L}$.

Please note that this only provides a rough guideline for the standard concentration. Accurate measurements are also possible if the above recommendations are not strictly adhered to.

- The equations for calculating the results are explained in the Instructions for use of the 692 pH/Ion Meter in paragraphs 3.1 and 3.2.
- Low sodium concentrations:
 In this case the response time is distinctly longer. The blank values of the reagents used must be taken into account.

Parameters

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-13      time 08:52:54
parameter
meas.type:          direct
>measuring parameters
ion type:           Na(+1)
conc.unit:         mg/l
meas.input:        1
electr.id:         0501.100
drift              0.5 mV/min
temperature        25.0 C
method id         NA
delta measurement: OFF
stirrer:          control
prestir pause     10 s
stir time         90 s
poststir pause    20 s
>calculation parameters
smpl size         20.0 ml
V total          40.0 ml
factor           1.0
smpl size unit:   ml
>calibration parameters
temperature       25.0 C
drift            0.5 mV/min
report:          full
no.of standards  6
addition:        auto
min.conc.        50.0 mg/l
max.conc.        300.0 mg/l
V init           40.0 ml
no.of Exchange Units 1
conc.1           2000.0 mg/l
V Exchange Unit 1: 10 ml
>analog output
select:          Conc
state:           OFF
>limits Conc
state:           OFF
>limits T
state:           OFF
>plot parameters
left mar.        0.0E+00 mg/l
right mar.       1.0E+30 mg/l
left marg.T      20 C
right marg.T     30 C
>preselections
req.ident:       OFF
req.smpl size:   value
-----
```

Configuration

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-13      time 08:52:24
config
>printer
id1              ISE Na
id2              AB Nr. 83/4
print header:    always
character set:   Epson
>print meas.value
print crit:      drift
date&time:      ON
>auxiliaries
last digit:     ON
dialog:         english
date            96-03-13
time            08:52:25
temp.unit:      C
run number      0
Dosimat:        665
device label
program         692.0020
>RS232 settings
baud rate:      9600
data bit:       8
stop bit:       1
parity:         none
handshake:     HWS
RS control:     ON
-----
```

a) Measurements with the glass membrane ISE

Direct measurement

Please observe the information given in the Instructions for use of the 692 pH/Ion Meter.


The following calibration is an example of a sample with an unproblematic matrix.

Sample:	NaCl solution with 230 mg/L Na^+
Sample volume:	20 mL
Volume of TISAB:	20 mL
Standard:	2000 mg/L Na^+ in 10 mL Exchange Unit
Concentration range:	50 ... 300 mg/L Na^+
No. of calibration points:	6

Calibration

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-12 time 16:37:52 31
Conc NA electr.id 0501.100
conc.calibration
meas.input:          1
ion type             Na(+1)
temperature man.    25.0 C
cal.date 96-03-12 16:37

  conc/mg/l      U/mV      dconc/%
std. 1  5.00E+01  109.5      0.0
std. 2  7.15E+01  118.7     -0.1
std. 3  1.02E+02  127.8      0.0
std. 4  1.47E+02  137.1      0.0
std. 5  2.10E+02  146.2      0.1
std. 6  3.00E+02  155.5     -0.1
```



```
variance      0.001
slope         59.1 mV
E(0)         9.0 mV
=====
```

Results

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-12 time 16:50:36
Conc NA electr.id 0501.100
id1 ISE Na
id2 AB Nr. 83/4
smp1 size 5.0 ml
#35 231 mg/l Na(+1)
96-03-12 16:50:36
=====
```

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-12 time 16:50:50
Conc NA electr.id 0501.100
id1 ISE Na
id2 AB Nr. 83/4
smp1 size 5.0 ml
#36 231 mg/l Na(+1)
96-03-12 16:50:50
=====
```

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-12 time 16:51:06
Conc NA electr.id 0501.100
id1 ISE Na
id2 AB Nr. 83/4
smp1 size 5.0 ml
#37 231 mg/l Na(+1)
96-03-12 16:51:06
=====
```

Result of a threefold determination:

```
Mean value: 231 mg/L Na+
Sabs:      1 mg/L Na+
Srel:      0.4%
```

Standard addition

Please observe the information given in the Instructions for use of the 692 pH/Ion Meter.

Example of determining the concentration by standard addition:

```
Sample:          NaCl solution with 57.5 mg/L Na+
Sample volume:  20 mL
Volume of TISAB: 20 mL
Standard:        2000 mg/L Na+ in 10 mL Exchange Unit
No. of additions: 5
```

Parameters


```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-11 time 15:18:24
parameter
meas.type:          std.add
>measuring parameters
ion type:           Na(+1)
conc.unit:          mg/l
meas.input:        1
electr.id:          0501.100
drift:              0.5 mV/min
temperature:        25.0 C
method id:          NA
stirrer:            control
prestir pause:     10 s
stir time:          90 s
poststir pause:    20 s
>calculation parameters
smp1 size:          20.0 ml
V total:            40.0 ml
factor:             1.0
smp1 size unit:    ml
>standard addition
type:               add
conc.std:           2000.0 mg/l
report:             full
addition:           auto
delta U:            12 mV
dos.rate:           medium
no. of additions:  5
V Exchange Unit:   10
activate pulse:    OFF
stop V:             10.0 ml
>preselections
req.ident:          OFF
req.smp1 size:     value
=====
```

Result

```

692 pH/Ion Meter      ON1/109  692.0020
date 96-03-11 time 14:17:56 27
Conc NA electr.id 0501.100
id1 ISE Na
id2 AB Nr. 83/4
addition/subtraction method
meas.type: std.add
meas.input: 1
temperature man. 25.0 C
conc.std 2000.0 mg/l
V total 40.0 ml
initial voltage 93.5 mV
smpl size 20.0 ml

          dV/ml  U/mV  dU/mV
std.incr. 1  0.349  105.4  11.9
std.incr. 2  0.569  117.3  11.9
std.incr. 3  0.935  129.2  11.9
std.incr. 4  1.463  140.5  11.3
std.incr. 5  2.552  152.2  11.7
  
```



```

variance      0.001
slope         60.0 mV
E(0)         5.1 mV
Na(+1)       59.4 mg/l
  
```

Result:

Measured value: 59.4 mg/L Na⁺
 Theoretical value: 57.5 mg/L Na⁺
 Absolute deviation: 1.9 mg/L Na⁺
 Relative deviation: 3.3%

b) Measurements with the polymer membrane ISE

Direct measurement

The following calibration is an example of a sample with an unproblematic matrix.

Sample: solution with 207.42 mg/L Na⁺
 Sample volume: 40 mL
 Standard: 2106.9 mg/L Na⁺ in 10 mL Exchange Unit
 Concentration range: 50 ... 300 mg/L Na⁺
 No. of calibration points: 6

Parameters

```

692 pH/Ion Meter      692.0021
date 2001-05-22 time 08:58:12
parameter
meas.type: direct
>measuring parameters
ion type: Na(+1)
conc.unit: PPM
meas.input: 1
electr.id: 0508.100
drift 0.1 mV/min
method id *****
delta measurement: OFF
stirrer: control
prestir pause 10 s
stir time 90 s
poststir pause 20 s
>calculation parameters
smpl size 40.0 ml
V total 40.0 ml
factor 1.0
smpl size unit: ml
>calibration parameters
drift 0.5 mV/min
report: full
cal.interval OFF h
no.of standards 6
addition: auto
min.conc. 50.0 PPM
max.conc. 300.0 PPM
V init 40.0 ml
no.of Exchange Units 1
conc.i 2106.9 PPM
V Exchange Unit 1: 10 ml
>analog output
select: Conc
state: ON
0 mV at 1.00E-03 PPM
1 V range 100.00 PPM
>limits Conc
state: OFF
>limits T
state: OFF
>plot parameters
left mar. 0.0E+00 PPM
right mar. 1.0E+00 PPM
left marg.T 20 °C
right marg.T 30 °C
>preselections
req.ident: OFF
req.smpl size: OFF
  
```


Configuration

```


692 pH/Ion Meter          692.0021
date 2001-05-11 time 14:47:09
config
)printer
  id1      Na ISE/ Polymer
  id2      AB Nr. 83/5
  print header:  always
  calibration report:  ON
  character set:  Seiko
)print meas.value
  print crit:  drift
  date&time:  ON
)auxiliaries
  last digit:  ON
  dialog:  english
  date 2001-05-11
  time 14:47:10
  temp.unit:  C
  run number 9
  Dosimat: 665
  device label 692-1
  program 692.0021
)RS232 settings
  baud rate: 9600
  data bit: 7
  stop bit: 1
  parity: even
  handshake: HWS
  RS control: ON
  
```

Calibration

```

692 pH/Ion Meter          692.0021
date 2001-05-22 time 08:59:17 1
Conc ***** electr.id 0508.100
conc.calibration
  meas.input: 1
  ion type Na(+1)
  temperature 23.6 °C
  c.dat. 2001-05-22 08:35

  std. 1 5.00E+01 -16.4 0.8
  std. 2 7.15E+01 -7.1 0.0
  std. 3 1.02E+02 2.2 -1.0
  std. 4 1.46E+02 11.2 -0.6
  std. 5 2.10E+02 20.0 0.1
  std. 6 3.00E+02 29.0 0.7
  
```



Results

```

692 pH/Ion Meter          692.0021
date 2001-05-22 time 09:06:56
Conc ***** electr.id 0508.100
id1 Na ISE/ Polymer
id2 AB Nr. 83/5
smpl size 40.0 ml
#2 206 PPM 23.5 °C Na(+1)
  
```

```

692 pH/Ion Meter          692.0021
date 2001-05-22 time 09:10:50
Conc ***** electr.id 0508.100
id1 Na ISE/ Polymer
id2 AB Nr. 83/5
smpl size 40.0 ml
#3 209 PPM 23.3 °C Na(+1)
2001-05-22 09:10:50
  
```

```

692 pH/Ion Meter          692.0021
date 2001-05-22 time 09:15:45
Conc ***** electr.id 0508.100
id1 Na ISE/ Polymer
id2 AB Nr. 83/5
smpl size 40.0 ml
#4 210 PPM 23.0 °C Na(+1)
2001-05-22 09:15:45
  
```

Result of a sevenfold determination:

Mean value: 210.43 mg/L Na⁺
 S_{abs}: 2.37 mg/L Na⁺
 S_{rel}: 1.13%

Standard addition

Example of determining the concentration by standard addition:

Sample: solution with 46 mg/L Na⁺
 Sample volume: 40 mL
 Standard: 2106.9 mg/L Na⁺ in 10 mL Exchange Unit
 No. of additions: 5

Parameters

```
692 pH/Ion Meter          692.0021
date 2001-05-11 time 11:57:26
parameter
  meas.type:      std.add
>measuring parameters
  ion type:       Na(+1)
  conc.unit:      ppm
  meas.input:     1
  electr.id:      0508.100
  drift:          0.1 mV/min
  method id:      Na3
  stirrer:        control
  prestir pause: 10 s
  stir time:      90 s
  poststir pause: 20 s
>calculation parameters
  smpl size:      40.0 ml
  V total:        40.0 ml
  factor:         1.0
  smpl size unit: ml
>standard addition
  type:           add
  conc.std.:      2106.9 ppm
  report:         full
  addition:       auto
  delta U:        12 mV
  dos.rate:       medium
  no.of additions: 5
  V Exchange Unit: 10
  activate pulse: OFF
  stop V:         10.0 ml
>preselections
  req.ident:      OFF
  req.smpl size:  value
```


Configuration

```
692 pH/Ion Meter          692.0021
date 2001-05-11 time 11:57:24
confis
>printer
  id1:            Na ISE/ Polymer
  id2:            AB Nr. 83/5
  print header:   always
  calibration report: OFF
  character set:   Seiko
>print meas.value
  print crit:     immediate
  date&time:      ON
>auxiliaries
  last digit:     ON
  dialog:         english
  date:           2001-05-11
  time:           11:57:24
  temp.unit:      C
  run number:     6
  Dosimat:        665
  device label:   692-1
  program:        692.0021
>RS232 settings
  baud rate:      9600
  data bit:       7
  stop bit:       1
  parity:         even
  handshake:      HwS
  RS control:     ON
```

Result

```
692 pH/Ion Meter          692.0021
date 2001-05-11 time 09:03:29 2
Conc Na3 electr.id 0508.100
id1 Na ISE/ Polymer
id2 AB Nr. 83/5
addition/subtraction method
  meas.type:      std.add
  meas.input:     1
  temperature:    24.0 C
  conc.std:       2106.9 ppm
  V total:        40.0 ml
  initial voltage: -33.0 mV
  smpl size:      40.0 ml
```

	dV/ml	U/mV	dU/mV
std.incr. 1	0.400	-23.9	9.1
std.incr. 2	0.582	-14.8	9.1
std.incr. 3	0.915	-5.3	9.5
std.incr. 4	1.472	4.4	9.7
std.incr. 5	2.393	14.2	9.8



variance 0.002
 slope 57.3 mV
 E(0) -128.1 mV
 Na(+1) 45.7 ppm

Result of a tenfold determination:

Mean value: 45.66 mg/L Na⁺
 S_{abs}: 0.66 mg/L Na⁺
 S_{rel}: 1.45%

Practical examples

1. Sodium in the inlet of a waste water treatment plant using automatic standard addition

Electrode: glass membrane ISE
 Sample preparation: Sample is diluted 1 : 4 with dist. water (→ factor = 4).
 Analysis: 20 mL sample solution and 20 mL TISAB solution are pipetted into the titration vessel.
 Standard: β(Na⁺) = 2000 mg/L
 No. of additions: 5


Parameters

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-13      time 08:19:20
parameter
  meas.type:          std.add
>measuring parameters
  ion type:           Na(+1)
  conc.unit:          mg/l
  meas.input:         1
  electr.id:          0501.100
  drift:              0.5 mV/min
  temperature:        25.0 C
  method id:          NA
  stirrer:            control
  prestir pause:      10 s
  stir time:          90 s
  poststir pause:    20 s
>calculation parameters
  smpl size:          20.0 ml
  V total:            40.0 ml
  factor:             4.0
  smpl size unit:     ml
>standard addition
  type:               add
  conc.std.:          2000.0 mg/l
  report:             full
  addition:           auto
  delta U:            12 mV
  dos.rate:           medium
  no.of additions:    5
  V Exchange Unit:   10
  activate pulse:     OFF
  stop V:             10.0 ml
>preselections
  req.ident:          OFF
  req.smpl size:      value
```

Result

```
692 pH/Ion Meter      ON1/109 692.0020
date 96-03-11      time 17:08:31 34
Conc      NA      electr.id 0501.100
id1       ISE Na
id2       ARA in
addition/subtraction method
  meas.type:          std.add
  meas.input:         1
  temperature man.:   25.0 C
  conc.std.:          2000.0 mg/l
  V total:            40.0 ml
  initial voltage:    96.8 mV
  smpl size:          5.0 ml
```

	dV/ml	U/mV	dU/mV
std.incr. 1	0.415	108.8	12.0
std.incr. 2	0.675	120.8	12.0
std.incr. 3	1.035	132.2	11.4
std.incr. 4	1.750	143.7	11.5
std.incr. 5	2.957	155.1	11.4



variance 0.001
 slope 60.5 mV
 E(0) 3.4 mV
 Na(+1) 280 mg/l

Result of a fourfold determination:

Mean value: 285 mg/L Na⁺
S_{abs}: 5 mg/L Na⁺
S_{rel}: 1.7%

2. Sodium in the outlet of a waste water treatment plant using automatic standard addition

Electrode: glass membrane ISE
Sample preparation: not required
Analysis: 15 mL sample and 25 mL TISAB solution are pipetted into the titration vessel.
Standard: β(Na⁺) = 2000 mg/L
No. of additions: 5

Parameters

```


692 pH/Ion Meter      ON1/109  692.0020
date 96-03-12      time 13:16:22
parameter
  meas.type:          std.add
>measuring parameters
  ion type:           Na(+1)
  conc.unit:          mg/l
  meas.input:         1
  electr.id:          0501.100
  drift:              0.5 mV/min
  temperature:        25.0 C
  method id:          NA
  stirrer:            control
  prestir pause:      10 s
  stir time:          90 s
  poststir pause:     20 s
>calculation parameters
  smpl size:          15.0 ml
  V total:            40.0 ml
  factor:             1.0
  smpl size unit:     ml
>standard addition
  type:               add
  conc.std:           2000.0 mg/l
  report:             full
  addition:           auto
  delta U:            12 mV
  dos.rate:           medium
  no.of additions:    5
  V Exchange Unit:    10
  activate pulse:     OFF
  stop V:             10.0 ml
>preselections
  req.ident:          OFF
  req.smpl size:      value
  
```

Result

```

692 pH/Ion Meter      ON1/109  692.0020
date 96-03-12      time 10:50:50  8
Conc      NA      electr.id  0501.100
id1       ISE Na
id2       ARA out
addition/subtraction method
  meas.type:          std.add
  meas.input:         1
  temperature man.:   25.0 C
  conc.std:           2000.0 mg/l
  V total:            40.0 ml
  initial voltage:    98.8 mV
  smpl size:          15.0 ml

          dV/ml    U/mV    dU/mV
std.incr. 1    0.406    110.5    11.9
std.incr. 2    0.660    122.4    11.9
std.incr. 3    1.092    134.1    11.7
std.incr. 4    1.855    145.8    11.7
std.incr. 5    3.298    157.4    11.6
  
```



```

variance      0.000
slope         57.9 mV
E(0)          11.0 mV
Na(+1)        86.3 mg/l
  
```

Result of fivefold determination:

Mean value: 87.3 mg/L Na⁺
S_{abs}: 0.7 mg/L Na⁺
S_{rel}: 0.8%

3. Sodium in spinach (salted cream spinach) using automatic standard addition

Electrode: glass membrane ISE
Sample preparation: Approx. 10 g sample (package information: 1% NaCl) is stirred for 15 min with 100 mL dist. water. The sample is filtered and the filtrate made up to 200 mL.
Analysis: 10 mL sample solution (→ factor = 20), 10 mL dist. water and 20 mL TISAB solution are pipetted into the titration vessel.
Standard: β(Na⁺) = 2000 mg/L; this corresponds to w(NaCl) = 0.508%
No. of additions: 4

Parameters

```


692 pH/Ion Meter      ON1/109 692.0020
date 96-03-13 time 08:34:44
parameter
  meas.type:          std.add
>measuring parameters
  ion type:           Na(+1)
  conc.unit:          %
  meas.input:         1
  electr.id:          0501.100
  drift:              0.5 mV/min
  temperature:        25.0 C
  method id:          NA
  stirrer:            control
  prestir pause:      10 s
  stir time:          90 s
  poststir pause:     20 s
>calculation parameters
  smpl size:          9.35 g
  V total:            40.0 ml
  factor:             20.0
  smpl size unit:     g
>standard addition
  type:               add
  conc.std.:          0.508 %
  report:             full
  addition:           auto
  delta U:            12 mV
  dos.rate:           medium
  no.of additions:    4
  V Exchange Unit:    10
  activate pulse:     OFF
  stop V:             10.0 ml
>preselections
  req.ident:          OFF
  req.smpl size:      value
  -----
  
```

Result

```

692 pH/Ion Meter      ON1/109 692.0020
date 96-03-12 time 18:12:58 44
Conc NA electr.id 0501.100
id1 ISE Na
id2 Spinat NaCl
addition/subtraction method
  meas.type:          std.add
  meas.input:         1
  temperature man.:   25.0 C
  conc.std.:          0.508 %
  V total:            40.0 ml
  initial voltage:    104.9 mV
  factor:             20.0
  smpl size:          9.35 g

  std.incr. 1    dV/ml    U/mV    dU/mV
  std.incr. 2    0.532    116.9   12.0
  std.incr. 3    0.868    128.9   12.0
  std.incr. 4    1.416    140.6   11.7
  std.incr. 4    2.398    152.2   11.6
  
```



```

variance          0.001
slope              60.1 mV
E(0)              222.3 mV
Na(+1)            0.958 %
  -----
  
```

Result of a fourfold determination:

```

Mean value: 0.952% NaCl
Sabs:       0.049% NaCl
Srel:       4.9%
  
```

4. Sodium in urine using automatic standard addition

```

Electrode:        glass membrane ISE
Sample preparation: Sample is diluted 1 : 40 with dist.
                  water (→ factor = 40).
Analysis:         20 mL sample solution and
                  20 mL TISAB solution are pipet-
                  ted into the titration vessel.
Standard:         β(Na+) = 2000 mg/L
No. of additions: 4
  
```

Parameters


```

692 pH/Ion Meter      ON1/109 692.0020
date 96-03-14 time 13:39:56
parameter
  meas.type:          std.add
>measuring parameters
  ion type:           Na(+1)
  conc.unit:          g/l
  meas.input:         1
  electr.id:          0501.100
  drift:              0.5 mV/min
  temperature:        25.0 C
  method id:          NA
  stirrer:            control
  prestir pause:      10 s
  stir time:          90 s
  poststir pause:     20 s
>calculation parameters
  smpl size:          20.0 ml
  V total:            40.0 ml
  factor:             40.0
  smpl size unit:     ml
>standard addition
  type:               add
  conc.std.:          2.0 g/l
  report:             full
  addition:           auto
  delta U:            12 mV
  dos.rate:           medium
  no.of additions:    4
  V Exchange Unit:    10
  activate pulse:     OFF
  stop V:             10.0 ml
>preselections
  req.ident:          OFF
  req.smpl size:      value
  -----
  
```

Result

```
692 pH/Ion Meter      ON1/109  692.0020
date 96-03-14      time 13:05:41      8
Conc      NA      electr.id  0501.100
id1      ISE Na
id2      Urin
addition/subtraction method
meas.type:      std.add
meas.input:      1
temperature man. 25.0 C
conc.std      2.0 g/l
V total      40.0 ml
initial voltage 107.1 mV
factor      40.0
smpl size    20.0 ml
```

	dV/ml	U/mV	dU/mV
std.incr. 1	0.545	119.3	12.2
std.incr. 2	0.910	131.5	12.2
std.incr. 3	1.507	143.4	11.9
std.incr. 4	2.602	155.1	11.7



```
variance      0.001
slope      58.9 mV
E(0)      187.6 mV
Na(+1)      3.45 g/l
```

Result of a fivefold determination:

Mean value: 3.44 g/L Na⁺
 Sabs: 0.03 g/L Na⁺
 Srel: 1.0%

5. Sodium in tap water using automatic standard addition

Electrode: polymer membrane ISE
 Sample preparation: not required
 Analysis: 40 mL tap water (taken at 25 °C after 1 min flow time) is pipetted into the titration vessel.
 Standard: β(Na⁺) = 207.42 mg/L
 No. of additions: 5

Parameters

```
692 pH/Ion Meter      692.0021
date 2001-05-09      time 15:48:59
```


```
parameter
meas.type:      std.add
>measuring parameters
ion type:      Na(+1)
conc.unit:      ppm
meas.input:      1
electr.id:      0508.100
drift      0.1 mV/min
method id      Natrium
stirrer:      control
prestir pause  10 s
stir time      90 s
poststir pause 20 s
>calculation parameters
smpl size      40.0 ml
V total      40.0 ml
factor      1.0
smpl size unit: ml
>standard addition
type:      add
conc.std.      207.42 ppm
report:      full
addition:      auto
delta U      12 mV
dos.rate:      medium
no.of additions 5
V Exchange Unit: 10
activate pulse: OFF
stop V      10.0 ml
>preselections
req.ident:      OFF
req.smpl size:  value
```

Configuration

```

692 pH/Ion Meter          692.0021
date 2001-05-09 time 15:49:50
config
>printer
  id1      Na ISE/ Polymer
  id2      AB Nr. 83/5
  print header:    always
  calibration report:  OFF
  character set:    Seiko
>print meas.value
  print crit:      immediate
  date&time:       ON
>auxiliaries
  last digit:      ON
  dialos:          english
  date             2001-05-09
  time            15:49:51
  temp.unit:       C
  run number       1
  Dosimat:        665
  device label     692-1
  program          692.0021
>RS232 settings
  baud rate:       9600
  data bit:        7
  stop bit:        1
  parity:          even
  handshake:       HwS
  RS control:      ON
  
```

	dV/ml	U/mV	dU/mV
std.incr. 1	0.415	-97.9	9.4
std.incr. 2	0.640	-88.1	9.8
std.incr. 3	0.934	-78.6	9.5
std.incr. 4	1.378	-69.4	9.2
std.incr. 5	2.116	-60.3	9.1



```

variance      0.000
slope         59.6 mV
E(0)         -147.6 mV
Na(+1)       4.74 ppm
  
```

Result

```

692 pH/Ion Meter          692.0021
date 2001-05-09 time 14:58:57   4
Conc Natrium electr.id 0508.100
id1   Na ISE/ Polymer
id2   AB Nr. 83/5
addition/subtraction method
  meas.type:      std.add
  meas.input:     1
  temperature     24.8 C
  conc.std        207.42 ppm
  V total         40.0 ml
  initial voltage -107.3 mV
  smpl size       40.0 ml
  
```

Result of an eightfold determination:

```

Mean value:  4.74 mg/L Na+
Sabs:      0.04 mg/L Na+
Srel:      0.86%
  
```

6. Sodium in mineral water using automatic standard addition

Electrode: polymer membrane ISE
 Sample preparation: not required
 Analysis: 40 mL mineral water (without carbon dioxide) is pipetted into the titration vessel.
 Standard: $\beta(\text{Na}^+) = 207.42 \text{ mg/L}$
 No. of additions: 4

Parameters

692 pH/Ion Meter 692.0021
 date 2001-05-09 time 15:48:59

```
parameter
meas.type:      std.add
>measuring parameters
ion type:       Na(+1)
conc.unit:      ppm
meas.input:     1
electr.id:      0508.100
drift           0.1 mV/min
method id      Natrium
stirrer:        control
prestir pause  10 s
stir time       90 s
poststir pause 20 s
```

```
>calculation parameters
sml size        40.0 ml
V total         40.0 ml
factor          1.0
sml size unit:  ml
```

```
>standard addition
type:           add
conc.std.       207.42 ppm
report:         full
addition:       auto
delta U         12 mV
dos.rate:       medium
no.of additions 5
V Exchange Unit 10
activate pulse: OFF
stop V          10.0 ml
```

```
>preselections
req.ident:      OFF
req.sml size:   value
```

Configuration

692 pH/Ion Meter 692.0021
 date 2001-05-09 time 15:49:50

```
config
>printer
id1      Na ISE/ Polymer
id2      AB Nr. 83/5
print header:  always
calibration report: OFF
character set:  Seiko
```

```
>print meas.value
print crit:    immediate
date&time:     ON
```

```
>auxiliaries
last digit:    ON
dialog:        english
date          2001-05-09
time          15:49:51
temp.unit:     C
run number    1
Dosimat:       665
device label   692-1
program       692.0021
```


```
>RS232 settings
baud rate:     9600
data bit:      7
stop bit:      1
parity:        even
handshake:     HWs
RS control:    ON
```

Result

692 pH/Ion Meter 692.0021
 date 2001-05-21 time 15:14:28 11
 Conc Na4 electr.id 0508.100

```
id1      Na ISE/ Polymer
id2      AB Nr. 83/5
addition/subtraction method
meas.type:      std.add
meas.input:     1
temperature     24.4 C
conc.std        207.42 ppm
V total         40.0 ml
initial voltage -64.1 mV
sml size        40.0 ml
```

	dV/ml	U/mV	dU/mV
std.incr. 1	0.790	-53.8	10.3
std.incr. 2	1.215	-43.6	10.2
std.incr. 3	1.964	-33.4	10.2
std.incr. 4	3.216	-23.3	10.1



variance 0.000
 slope 59.2 mV
 E(0) -116.9 mV
 Na(+1) 7.82 ppm

Result of a ninefold determination:

Mean value: 7.83 mg/L Na⁺
 S_{abs}: 0.15 mg/L Na⁺
 S_{rel}: 1.95%

7. Sodium in baby food (milk powder) using automatic standard addition

Electrode: polymer membrane ISE
 Sample preparation: 75.407 g baby food (package information: 0.04 g Na⁺ / 15 g milk powder; this corresponds to 0.27% Na⁺) is dissolved in 1 L warm dist. water.
 Analysis: 10 mL sample solution and 30 mL dist. water are pipetted into the titration vessel.
 Standard: β(Na⁺) = 2106.9 mg/L
 No. of additions: 5

Parameters

```

692 pH/Ion Meter          692.0021
date 2001-05-10 time 15:11:00
parameter
  meas.type:          std.add
>measuring parameters
  ion type:          Na(+1)
  conc.unit:         ppm
  meas.input:        1
  electr.id:         0508.100
  drift:             0.1 mV/min
  method id:         Na2
  stirrer:           control
  prestir pause:     10 s
  stir time:         90 s
  poststir pause:    20 s
>calculation parameters
  smpl size:         10.0 ml
  V total:           40.0 ml
  factor:            1.0
  smpl size unit:    ml
>standard addition
  type:              add
  conc.std:          2106.9 ppm
  report:            full
  addition:          auto
  delta U:           12 mV
  dos.rate:          medium
  no.of additions:   5
  V Exchange Unit:   10
  activate pulse:    OFF
  stop V:            10.0 ml
>preselections
  req.ident:         OFF
  req.smpl size:     value
  
```

Configuration

```

692 pH/Ion Meter          692.0021
date 2001-05-10 time 15:10:23
config
>printer
  id1      Na ISE/ Polymer
  id2      AB Nr. 83/5
  print header:  always
  calibration report: OFF
  character set: Seiko
>print meas.value
  print crit:   immediate
  date&time:    ON
  
```

```


>auxiliaries
last digit:      ON
dialog:         english
date           2001-05-10
time          15:10:24
temp.unit:      C
run number     11
Dosimat:       665
device label   692-1
program       692.0021
>RS232 settings
baud rate:     9600
data bit:      7
stop bit:      1
parity:        even
handshake:     HWs
RS control:    ON
  
```

Result

```

692 pH/Ion Meter      692.0021
date 2001-05-10 time 10:04:32 4
Conc Na2 electr.id 0508.100
idl Na ISE/ Polymer
id2 AB Nr. 83/5
addition/subtraction method
meas.type:           std.add
meas.input:          1
temperature          25.0 C
conc.std             2106.9 ppm
V total              40.0 ml
initial voltage      -48.0 mV
smpl size            10.0 ml
  
```

	dV/ml	U/mV	dU/mV
std.incr. 1	0.371	-37.6	10.4
std.incr. 2	0.562	-27.2	10.4
std.incr. 3	0.863	-16.9	10.3
std.incr. 4	1.374	-6.6	10.3
std.incr. 5	2.292	3.8	10.4



```

variance          0.005
slope             58.6 mV
E(0)              -148.2 mV
Na(+1)            149 ppm
  
```

Result of a tenfold determination:

Sodium concentration in the sample solution:

```

Mean value: 149.20 mg/L Na+
Sabs: 1.99 mg/L Na+
Srel: 1.33%
  
```

Sodium concentration in the original sample (milk powder):

```

Mean value: 0.20% Na+
  
```

Literature

- AOAC 976.25 (1990)
Nutritionally related components. Sodium in foods for special dietary use. Ion-selective electrode method
- L. Ciba
Direct potentiometric determination of sodium in silicic acid and water glass
Chem. Analit. Warsaw 36 (1991) 163–166 (in Polish)
Ref.: Fresenius J. Anal. Chem. 343 (1992) 446
- E. Florence
Determination of sodium in salted foods using an ion-selective electrode
Analyst 111 (1986) 571–573
Ref.: Fresenius Z. Anal. Chem. 326 (1987) 480
- S. A. H. Khalil, G. J. Moody, J. D. R. Thomas
Ion-selective electrode determination of sodium and potassium in blood and urine
Anal. Lett. 19 (1986) 1809–1830
Ref.: Fresenius Z. Anal. Chem. 327 (1987) 636
- E. Rabe
Zur Natrium- und Kaliumbestimmung mit ionensensitiven Elektroden
Z. Lebensm. Unters. Forsch. 176 (1983) 270–274
- S. K. Srivastava, S. Kumar, C. K. Jain
Determination of sodium by ion-selective potentiometry
Analyst 109 (1984) 667–668
- G. Schwedt
Analytische Chemie: Grundlagen, Methoden und Praxis
Georg Thieme Verlag, Stuttgart, New York, 1995
- M. Otto
Analytische Chemie
VCH Verlagsgesellschaft mbH, Weinheim, New York, Basel, Cambridge, Tokyo, 1995
- Instructions for use for ion-selective electrodes
Metrohm Ltd., Herisau, 2000