Application Bulletin 91/2 e



Potentiometric analysis of brass and bronze plating baths

Of interest for: Metals, Electroplating industries A 10

Summary

Methods are described for the potentiometric analysis of the following bath components:

Brass plating bath: copper, zinc, free cyanide, ammonium, carbonate and sulfite

Bronze plating bath: copper, tin and free cyanide

Instruments and accessories

- Titrino or Titrando with Dosino or Dosimat
- Magnetic Stirrer
- Exchange Unit(s)
- 6.0502.140 ion-selective copper electrode (Cu ISE) with 6.2104.020 electrode cable
- 6.0726.107 double-junction Ag/AgCl reference electrode [filled with c(KCl) = 3 mol/L] with 6.2106.020 electrode cable
- 6.0430.100 Ag Titrode with Ag₂S coating
- 6.0255.100 combined LL double-junction pH glass electrode
- 6.0431.100 Pt Titrode

1. Brass plating baths

1.1. Copper and zinc

Reagents

- Acids for digestion: While cooling, carefully add 80 mL H₂SO₄ conc. to 150 mL HNO₃ conc.
- Titrant 1: c(Na₂EDTA) = 0.1 mol/L
- Titrant 2: c(Na₂S₂O₃) = 0.1 mol/L
- Buffer solution pH = 10: Dissolve 54 g NH₄Cl and 350 mL w(NH₃) = 25% in dist. water and fill up to 1 L.
- Sulfuric acid, w(H₂SO₄) ≈ 25%
- Sodium hydroxide solution, c(NaOH) = 2 mol/L
- · Potassium iodide, p.a.

Sample preparation

Work in fume cupboard. Toxic HCN and acid fumes are released!

Add approx. 20 mL dist. water to 10.0 mL bath sample in an Erlenmeyer flask. Carefully add 3 mL digestive acid, warm up the mixture, heating until white sulfuric acid fumes appear. After cooling, rinse with dist. water into a 100 mL graduated flask, fill up to the mark and mix

Iodometric determination of copper

Pipet 10.0 mL of the prepared sample solution (corresponding to 1 mL of original sample) into a glass beaker and dilute with dist. water to approx. 50 mL. After adding 2 mL $w(H_2SO_4) \approx 25\%$ and approx. 1 g KI, titrate the freed iodine with $c(Na_2S_2O_3) = 0.1$ mol/L (Pt Titrode). Consumption of the titrant up to the equivalence point is stored in the titrator as common variable C31.

Calculation

1 mL $c(Na_2S_2O_3)$ = 0.1 mol/L corresponds to 6.3546 mg Cu^{2+}

 $g/L Cu^{2+} = EP1 * C01 / C00$

EP1 = titrant consumption in mL

C00 = 1 (sample volume used in mL original sample)

C01 = 6.3546

Chelatometric (complexometric) determination of the sum of copper and zinc

Dilute 10.0 mL of the prepared sample solution (corresponding to 1 mL of original sample) in a glass beaker with dist. water to approx. 50 mL and pre-neutralize to pH \approx 4 with c(NaOH) = 2 mol/L. After addition of 5 mL buffer solution pH = 10, titrate in MET mode with $c(Na_2EDTA) = 0.1$ mol/L (Cu ISE).



Calculation

1 mL $c(Na_2EDTA) = 0.1$ mol/L corresponds to 6.538 mg Zn^{2+}

 $g/L Zn^{2+} = (EP1 - C31) * C01 / C00$

EP1 = titrant consumption in mL

C00 = 1 (sample volume used in mL original sample)

C01 = 6.538

C31 = titrant consumption in mL during the iodometric titration of Cu²⁺ (common variable)

1.2. Free cyanide

Reagents

Titrant: c(AgNO₃) = 0.1 mol/L

Potassium iodide solution, w(KI) = 10%

• Sodium hydroxide solution, c(NaOH) = 2 mol/L

Analysis

Add 2 mL c(NaOH) = 2 mol/L to approx. 50 mL dist. water in a glass beaker. Add 2.0 mL bath sample as well as 4 mL w(KI) = 10% and titrate with $c(AgNO_3)$ = 0.1 mol/L (Ag Titrode with Ag₂S coating).

Calculation

1 mL c(AgNO₃) = 0.1 mol/L corresponds to 5.204 mg CN⁻ or 9.802 mg NaCN or 13.024 mg KCN

g/L cyanide = EP1 * C01 / C00

EP1 = titrant consumption in mL C00 = 2.0 (sample volume in mL) C01 = 5.204 or 9.802 or 13.024

1.3. Ammonium

Reagents

Titrant: c(HCl) = 0.1 mol/L

Boric acid solution, w(H₃BO₃) = 2%

Sodium hydroxide solution, w(NaOH) = 20%

Iron(II) sulfate solution, w(FeSO₄ x 7 H₂O) = 25%

Analysis

A Kjeldahl distillation apparatus is used for this analysis. Add 5.0 mL bath sample, 4 mL FeSO₄ solution and 15 mL w(NaOH) = 20% to the distillation flask and begin distillation immediately. The cooling apparatus tube submerges in an initial solution of 50 mL $w(\text{H}_3\text{BO}_3) = 2\%$. After 15 min the ammonia is distilled over. Titrate the initial solution with c(HCI) = 0.1 mol/L (combined pH glass electrode). The equivalence point of the titration is at pH ≈ 4.5 .

Calculation

1 mL c(HCI) = 0.1 mol/L corresponds to 1.8038 mg NH₄⁺ or 1.4007 mg N

 $g/L NH_4^+ = EP1 * C01 / C00$

g/L N = EP1 * C02 / C00

EP1 = titrant consumption in mL C00 = 5.0 (sample volume in mL)

C01 = 1.8038C02 = 1.4007

1.4. Carbonate

Reagents

• Titrant: c(HCI) = 1 mol/L

Barium chloride solution, w(BaCl₂) = 25%

Analysis

If the sample still contains free NaOH, two equivalence points are given. Because cyanide is also detected, the titration must be interrupted after this second EP, otherwise *toxic HCN* is released!

To a glass beaker with approx. 50 mL dist. water, add 2.0 bath sample. Add 5 mL $w(BaCl_2)$ = 25% and titrate with c(HCI) = 1 mol/L until shortly after the first (second) equivalence point is reached (combined pH glass electrode).

Calculation

1 mL c(HCI) = 1 mol/L corresponds to 106 mg Na₂CO₃

In presence of free NaOH:

 $g/L Na_2CO_3 = (EP2 - EP1) * C01 / C00$

In absence of free NaOH:

 $g/L Na_2CO_3 = EP1 * C01 / C00$

EP1 = titrant consumption up to the first EP in mL

EP2 = titrant consumption up to the second EP in mL

C00 = 2.0 (sample volume in mL)

C01 = 106

1.5. Sulfite

Reagents

• Titrant: c(I₂) = 0.05 mol/L

Barium chloride solution, w(BaCl₂) = 25%

Acetic acid, w(CH₃COOH) = 96%

• Ammonia, $w(NH_3) = 25\%$



Analysis

Pipet 10.0 mL bath sample into an Erlenmeyer flask and dilute with dist. water to approx. 100 mL. Add 0.25 mL $w(NH_3) = 25\%$, warm up the mixture, then add 30 mL $w(BaCl_2) = 25\%$ and heat to boiling point. Allow to settle, then pass through a folded filter. Rinse precipitation well with dist. water until rinse water reacts neutral. Penetrate the filter and rinse precipitation with dist. water into a glass beaker. While stirring, slowly and carefully add $w(CH_3COOH) = 96\%$ to the mixture to acidify, then titrate immediately with $c(I_2) = 0.05$ mol/L (Pt Titrode).

Calculation

1 mL $c(I_2)$ = 0.05 mol/L corresponds to 6.302 g Na₂SO₃

 $g/L Na_2SO_3 = EP1 * C01 / C00$

EP1 = titrant consumption in mL C00 = 10.0 (sample volume in mL)

C01 = 6.302

2. Bronze plating baths

Sample preparation for copper and tin

Work in fume cupboard. Toxic HCN and acid fumes are released!

In an Erlenmeyer flask, dilute 2.0 mL bath sample with approx. 20 mL dist. water. Carefully add 5 mL conc. $\rm H_2SO_4$ as well as a few drops of $\it w(\rm H_2O_2) = 30\%$ and heat up until white sulfuric acid fumes appear. Cool, then rinse with dist. water into a 100 mL graduated flask, fill up to mark and mix.

2.1. Copper

Reagents

- Titrant: c(Na₂S₂O₃) = 0.1 mol/L
- Sulfuric acid, $w(H_2SO_4) = 25\%$
- Potassium iodide, p.a.

Analysis

Add 5 mL $w(H_2SO_4)$ = 25% and approx. 1 g KI to 50 mL of the prepared sample solution (corresponding to 1 mL original sample) in a glass beaker and then titrate the freed iodine with $c(Na_2S_2O_3)$ = 0.1 mol/L (Pt Titrode).

Calculation

1 mL $c(Na_2S_2O_3)$ = 0.1 mol/L corresponds to 6.3546 mg de Cu²⁺

 $g/L Cu^{2+} = EP1 * C01 / C00$

EP1 = titrant consumption in mL

C00 = 1 (sample volume used in mL original sample)

C01 = 6.3546

2.2. Tin

Reagents

- $c(I_2) = 0.05 \text{ mol/L}$
- Hydrochloric acid, w(HCI) = 36%
- Iron powder, p.a.

Analysis

Add 20 mL w(HCI) = 36% to 50 mL of the prepared sample solution (corresponding to 1 mL original sample) in a glass beaker. While stirring, add approx. 1 g iron powder in small portions and heat shortly to boiling point. Cool the mixture, then pass it through a folded filter to enable the extraction of copper deposits and rinse well with hot dist. water. Add a further 10 mL w(HCI) = 36% as well as 0.5 g iron powder to the filtrate and boil until the iron is completely dissolved. Allow to cool, then titrate immediately with $c(I_2)$ = 0.05 mol/ (Pt Titrode).

Calculation

1 mL $c(I_2)$ = 0.05 mol/L corresponds to 5.9345 mg Sn²⁺

 $g/L Sn^{2+} = EP1 * C01 / C00$

EP1 = titrant consumption in mL

C00 = 1 (sample volume used in mL original sample)

C01 = 5.9345

2.3. Free cyanide

Reagents

- Titrant: c(AgNO₃) = 0.1 mol/L
- Potassium iodide solution, w(KI) = 10%
- Sodium hydroxide solution, c(NaOH) = 2 mol/L

Analysis

Make an initial solution of approx. 50 mL dist. water and 2 mL c(NaOH) = 2 mol/L in a glass beaker. Add 2.0 mL bath sample and 4 mL w(KI) = 10%, then titrate with $c(AgNO_3)$ = 0.1 mol/L (Ag Titrode with Ag_2S coating).



Calculation

1 mL c(AgNO₃) = 0.1 mol/L corresponds to 5.204 mg CN⁻ or 9.802 mg NaCN or 13.024 mg KCN
g/L cyanide = EP1 * C01 / C00
EP1 = titrant consumption in mL

C00 = 2.0 (sample volume in mL)

C01 = 5.204 or 9.802 or 13.024

Literature

- Metrohm Application Bulletin No. 101 Complexometric titrations with the Cu ISE Metrohm Ltd., Herisau
- Metrohm Application Note T-23
 Hydroxide and carbonate in alkaline plating baths for cadmium, copper, lead or zinc
 Metrohm Ltd., Herisau
- Metrohm Application Note T-24
 Metal contents of alkaline plating baths for cadmium, copper, lead or zinc
 Metrohm Ltd., Herisau
- T. W. Jelinek
 Prozessbegleitende Analytik in der Galvanotechnik
 Eugen G. Leuze Verlag, Saulgau, 1999
 ISBN 3-87-480-135-7

Figures




Fig. 1: Parameter settings for the iodometric determination of copper






Fig. 2: Titration curve for the iodometric determination of copper




Fig. 3: Titration curve for the iodometric determination of tin

```
751 GP Titrino 05268 751.0011 date 15-01-03 time 10:59 6
MET U
Parameters
>titration parameters
                       0.05 ml
 V step
 titr.rate
                       max. ml/min
 signal drift
equilibr.time
start V:
                       25 mV/min
34 s
                        OFF
 pause
                           0 s
 dos.element: internal D0
 meas.input:
                        25.0 °C
  temperature
>stop conditions
 stop V:
                        abs.
 stop V
                         10 ml
 stop U
                        OFF mV
 stop EP
 filling rate
                       max. ml/min
>statistics
 status:
                        OFF
>evaluation
                         5.2 mV
 EP recognition:
                        all
                         OFF mV
 fix EP1 at U
 pK/HNP:
                         OFF
>preselections
 req.ident:
                         OFF
 req.smpl size:
                         OFF
  activate pulse:
                         OFF
         =========
```

Fig. 4: Parameter settings for the chelatometric determination of copper/zinc






Fig. 5: Titration curve for the chelatometric determination of copper/zinc




Fig. 6: Titration curve for the determination of free cyanide




Fig. 7: Titration curve for the determination of ammonium






Fig. 8: Titration curve for the determination of carbonate




Fig. 9: Titration curve for the iodometric determination of sulfite