

Application Bulletin

Of interest to: Metals, Electroplating industry A C 10

Titrimetric analysis of cadmium plating baths

Summary

This Bulletin describes titrimetric methods for the determination of cadmium, free sodium hydroxide, sodium carbonate and total cyanide. The free cyanide can be calculated from the total cyanide and the Cd content.

Apparatus and accessories

- Titrino or Titrando with Dosino or Dosimat
- Magnetic swing-out stirrer
- Exchange units
- Photometer (610 nm) or Cu ISE 6.0431.140 with Ag/AgCl reference electrode 6.0726.107 (KCl 3 mol/L) and electrode cable 6.2106.020
- Combined pH glass electrode 6.0255.100 with electrode cable 6.2104.020
- Ag Titrode with Ag₂S coating 6.0430.100

Reagents

These are described under the individual analyses.

1. Determination of cadmium

1.1. By means of photometric titration (Photometer)

Reagents:

- $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- $w(\text{formaldehyde}) = 30\%$
- Buffer solution pH = 10:
Dissolve 114 mL $w(\text{NH}_3) = 25\%$ and 14 g NH₄Cl in dist. H₂O and fill up to 200 mL.
- Colored indicator:
Dissolve 100 mg each eriochrome black T and vitamin C in dist. H₂O and fill up to 100 mL.
- $w(\text{KCN}) = 6.5\%$

Analysis:

Place a 1.0 ... 2.0 mL bath sample, containing approx. 50 mg Cd in a beaker, add 1 mL KCN, approx. 80 mL dist. H₂O, 20 mL buffer solution pH = 10 and 0.25 mL colored indicator. Slowly add 4 mL formaldehyde and, while stirring allow to react for 1 min (to set Cd free from the cyanide complex). Finally titrate with c(Na₂EDTA) = 0.1 mol/L using the Photometer.

Calculations:

$$1 \text{ mL } c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = 11.241 \text{ mg Cd}$$

$$\text{g/L Cd} = \text{EP1} * \text{C01} / \text{C00}$$

C00 = Sample size in mL

C01 = 11.241

Figures:

'pa	>stop conditions	
751 GPD Titrino	05268	751.0011
date 2000-05-30	time 16:01	5
MET U	AB93 Cd	
parameters		
>titration parameters		
V step	0.10 ml	
dos.rate	max. ml/min	
signal drift	30 mV/min	
equilibr.time	32 s	
start V:	OFF	
pause	60 s	
dos.element:	internal D0	
meas.input:	1	
temperature	25.0 °C	
>statistics		
status:		OFF
>evaluation		
EPC		30 mV
EP recognition:	greatest	
fix EP1 at U		OFF mV
pK/HNP:		OFF
>preselections		
req.ident:		OFF
req.smpl.size:		OFF
activate pulse:		OFF
=====		

Fig. 1 Parameter report Titrino for Cd determinations

```
'fr
751 GPD Titrino      05268    751.0011
date 2000-05-30      time 16:01      5
U(init)            325 mV MET U    AB93 Cd
smpl.size          2.0 ml
EP1                7.308 ml      316 mV
Cadmium           41.07 g/l
stop V reached
=====
```




Fig. 2 Titration curve Cd by photometry

1.2. By means of potentiometric titration (Cu ISE)

Reagents:

- $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- $\text{Cu}(\text{NH}_4)_2\text{EDTA}, c = 0.1 \text{ mol/L}$ (Merck No. 105217)
- Buffer solution pH = 10; see under 1.1
- $w(\text{HNO}_3) = 65\%$

Analysis:

Work under a fume cupboard, toxic HCN is set free!!!

Put 5.0 mL bath sample into a Kjeldahl flask and add approx. 10 mL dist. H_2O . Tilting back and forth, carefully add HNO_3 until the solution becomes definitely acidic. Heat up the solution under a fume cupboard and cook until all cyanide is destroyed and entirely removed. After cooling, rinse the solution with dist. H_2O into a 50 mL graduated flask, fill up to the mark and mix.

Pipet 10.0 ... 20.0 mL of the treated sample solution (corresponding to 1 ... 2 mL original bath) into a beaker and complete to approx. 40 mL with dist. H_2O . Add 5 mL buffer solution pH = 10 and 1 mL $\text{Cu}(\text{NH}_4)_2\text{EDTA}$ and, while stirring allow to react for 1 min. Finally titrate with $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ in the MET mode of the titrator, using the Cu ISE.

Calculations:

$$1 \text{ mL } c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = 11.241 \text{ mg Cd}$$
$$\text{g/L Cd} = \text{EP1} * \text{C01} / \text{C00}$$

C00 = Sample size in mL (original sample)

C01 = 11.241

Figures:

```
'fr
751 GPD Titrino          05268   751.0011
date 2000-05-31    time 09:38      2
U(init)      -67 mV MET U      AB93 Cd
smpl size     2.0 ml
EP1          7.438 ml           -110 mV
Cadmium      41.81 g/l
stop V reached
=====
```




Fig. 3 Titration curve Cd by potentiometry

2. Determination of free NaOH and carbonate

Reagents:

- $c(\text{HCl}) = 1 \text{ mol/L}$
- $w(\text{BaCl}_2) = 25\%$

Analysis:

Put approx. 50 mL dist. H₂O and 2.0 mL bath sample into a glass beaker. After addition of 5 mL BaCl₂ solution, titrate with $c(\text{HCl}) = 1 \text{ mol/L}$, using the pH glass electrode, until shortly after the second endpoint.

Calculations:

$$1 \text{ mL } c(\text{HCl}) = 1 \text{ mol/L} = 40.0 \text{ mg NaOH or } 106.0 \text{ mg Na}_2\text{CO}_3$$

$$\text{g/L NaOH} = EP1 * C01 / C00$$

$$\text{g/L Na}_2\text{CO}_3 = (EP2 - EP1) * C02 / C00$$

C00 = Sample size in mL (2)

C01 = 40

C02 = 106

Remarks:

- The titration has to be stopped after reaching the second endpoint, otherwise toxic HCN may be set free. It is best to work under a fume cupboard!

Figures:

'pa	>stop conditions	
751 GPD Titrino	05268	751.0011
date 2000-05-31	time 12:42	3
DET U	AB93 OH-	
parameters		
>titration parameters		
meas.pt.density	4	
min.incr.	10.0 µl	
dos.rate	max. ml/min	
signal drift	25 mV/min	
equilibr.time	34 s	
start V:	OFF	
pause	0 s	
dos.element:	internal D0	
meas.input:	1	
temperature	25.0 °C	
	>statistics	
	status:	OFF
	>evaluation	
	EPC	5
	EP recognition:	2
	fix EP1 at U	OFF mV
	pK/HNP:	OFF
	>preselections	
	req.ident:	OFF
	req.smpl size:	OFF
	activate pulse:	OFF

Fig. 4 Parameter report Titrino for NaOH, Na₂CO₃

'fr	
751 GPD Titrino	05268
date 2000-05-31	time 12:42
U(init)	-309 mV
smpl size	2.0 ml
EP1	1.963 ml
EP2	3.409 ml
NaOH	39.26 g/l
Na ₂ CO ₃	76.64 g/l
manual stop	




Fig. 5 Titration curve NaOH / Na_2CO_3

3. Determination of total cyanide

Reagents:

- $c(\text{AgNO}_3) = 0.1 \text{ mol/L}$
- $c(\text{NaOH}) = 2 \text{ mol/L}$
- $w(\text{potassium iodide}) = 10\%$

Analysis:

Place approx. 50 mL dist. H_2O and 2 mL NaOH into a glass beaker. Then add 1.0 mL bath sample and 2 mL KI solution, and titrate with $c(\text{AgNO}_3) = 0.1 \text{ mol/L}$ using the Ag Titrode.

Calculations:

$$1 \text{ mL } c(\text{AgNO}_3) = 0.1 \text{ mol/L} = 5.204 \text{ mg CN or } 9.802 \text{ mg NaCN or } 13.024 \text{ mg KCN}$$

$$\text{g/L „cyanide“} = EP1 * C01 / C00$$

C00 = Sample size in mL (1)

C01 = 5.204 or 9.802 or 13.024

Figures:

'pa	>stop conditions	
751 GPD Titrino	05268	751.0011
date 2000-05-31	time 13:15	5
DET U	AB93 CN-	
parameters		
>titration parameters		
meas.pt.density	4	
min.incr.	10.0 μl	
dos.rate	max. mL/min	
signal drift	25 mV/min	
equilibr.time	34 s	
start V:	OFF	
pause	0 s	
dos.element:	internal D0	
meas.input:	1	
temperature	25.0 $^{\circ}\text{C}$	
	>statistics	
	status: OFF	
	>evaluation	
	EPC	5
	EP recognition:	greatest
	fix EP1 at U	OFF mV
	pK/HNP:	OFF
	>preselections	
	req.ident:	OFF
	req.smpl size:	OFF
	activate pulse:	OFF
	=====	

Fig. 6 Parameter report Titrino for the cyanide determination

'fr
751 GPD Titrino 05268 751.0011
date 2000-05-31 time 13:15 5
U(init) 204 mV DET U AB93 CN-
smpl size 2.0 ml
EP1 10.019 ml 20 mV
Cyanid 13.03 g/l
stop V reached
=====




Fig. 7 Titration curve total cyanide

4. Calculation of free cyanide

To calculate the free cyanide, the content of cyanide combined with Cd must be subtracted from the total cyanide content. According to the formula $K_2Cd(CN)_4$ this corresponds to 0.926 g CN per g Cd.

$$\text{Free cyanide; g/L CN} = (\text{g/L total-CN}) - (\text{g/L Cd} * 0.926)$$

Literature

- Metrohm Application Bulletin No. 101
- Metrohm Ti Application Note No. T-22, T-23, T-24
- Wild,P.W.
Moderne Analysen für die Galvanik
Eugen G. Leuze Verlag, D-88348 Saulgau/Württ. 1972
- Jelinek,T.W.
Prozessbegleitende Analytik in der Galvanotechnik
Eugen G. Leuze Verlag, D-88348 Saulgau/Württ. 1999
ISBN 3-87-480-135-7