

Determination of mercury at the rotating gold electrode by anodic stripping voltammetry

Of interest to:

General analytical laboratories; water analysis; food analysis; trace analysis

B 1, 2, 7, 9

Summary

This Application Bulletin describes the determination of mercury by anodic stripping voltammetry (ASV) at the rotating gold electrode. With a deposition time of 90 s the calibration curve is linear from 0.4 µg/L to 15 µg/L; the limit of quantification is 0.4 µg/L.

The method has primarily been drawn up for investigating water samples. After appropriate digestion the determination of mercury is possible even in samples with a high load of organic substances (wastewater, food and semi-luxuries, biological fluids, pharmaceuticals).

Instruments and accessories

- 797 VA Computrace
or
746 VA Trace Analyzer + 747 VA Stand
- Electrodes and measuring vessel:
For this application the measuring cell must be completely free from platinum. This means that neither a measuring vessel nor electrodes can be used that have already been used once with a platinum working or auxiliary electrode. This is why it is rgently recommended that electrodes and measuring vessel are only used for the determination of mercury.

WE	Driving axle	6.1204.210
	Au electrode tip	6.1204.140
AE	Electrode holder	6.1241.020
	Glassy carbon rod	6.1247.000
RE	Ag/AgCl reference electrode c(KCl) = 3 mol/L	6.0728.020
	Electrolyte vessel c(NaCl) = 3 mol/L	6.1245.010
Stopper		6.2709.040
Measuring vessel		6.1415.150

- 705 UV Digester

Sample preparation

Voltammetric determinations are sensitive to interference by organic substances. This means that samples that contain organic substances must first be digested.


- Groundwater, drinking water, seawater and mineral water can usually be analyzed directly.
- Water that is lightly to moderately polluted with organic substances is digested using the 705 UV Digester:
10 mL acidified water sample (pH = 2) is treated with 10 µL w(HNO₃) = 65% and 50 µL w(H₂O₂) = 30% and irradiated for 60 min at 90 °C.
- Samples with an organic matrix (food, drugs, etc.) must undergo a wet digestion:
 - High-pressure ashing
 - Microwave digestion
 In both methods the sample is oxidized in a closed vessel using a mixture of concentrated mineral acids.
 - An open wet digestion cannot be recommended because of the volatility of some mercury compounds.

Conditioning the gold electrode

In order to obtain reproducible curves, the gold electrode must be conditioned electrochemically. This should be carried out daily before the start of measurements, and also when the background current alters greatly from measurement to measurement.

The background current should normally be as low as possible. With a well-functioning gold electrode it has an approximate minimum value between 100 nA and 300 nA.

Example of the background current in the conditioning solution of a well-conditioned electrode



Reagents

- Perchloric acid, suprapur, w(HClO₄) = 70%
- Ultrapure water type 1
(electrical resistivity > 18.2 MΩ·cm, TOC < 10 ppb)

Conditioning

Conditioning solution

20 mL ultrapure water
 + 0.2 mL HClO₄

Voltammetric conditioning parameters

Working electrode	RDE
Stirrer speed (rpm)	2000
Mode	DP
Purge time	30 s


Pretreatment	
Hydrodynamic measurement	no
Cleaning potential	0 V
Cleaning time	10 s
Deposition potential	+1.5 V
Deposition time	30 s
Sweep	
Equilibration time	5 s
Start potential	0 V
End potential	+0.8 V
Voltage step	0.006 V
Voltage step time	0.1 s
Sweep rate	0.06 V/s
Pulse amplitude	0.05 V
Pulse time	0.04 s
Cell off after measurement	yes
Determination	
No. of additions	0
No. of replications	10

Potentiostat

Highest current range	10 mA
Lowest current range	1 μA

Example

Changes in the background current of a gold electrode during the first 10 replications of conditioning.



Remarks


- $c(\text{NaCl}) = 3 \text{ mol/L}$ must be used as the electrolyte in the electrolyte vessel (6.1245.010) of the reference electrode. If $c(\text{KCl}) = 3 \text{ mol/L}$ is used then sparingly soluble potassium perchlorate is formed which can then block the reference electrode diaphragm.
- Conditioning with the parameters given above should only be carried out in the electrolyte described. In electrolytes that contain chloride the gold electrode could be irreversibly damaged.
- Conditioning the electrode must be repeated so many times until the baseline no longer changes appreciably. Experience has shown that the background current will be below 400 nA at +0.6 V. If the electrode is in regular use then it is usually sufficient to condition it for about one hour directly before the start of work.
- If conditioning the electrode has to be repeated several times then it is advisable to replace the conditioning solution frequently, particularly at the start of conditioning.
- The increase in background current at approx. +0.35 V becomes smaller as conditioning proceeds. With intensively used electrodes it can disappear completely.
- For cleaning purposes the gold electrode can also be conditioned in $c(\text{NaOH}) = 0.1 \text{ mol/L}$. In this case the following parameters should be used:

Voltammetric parameters for cleaning in $c(\text{NaOH}) = 0.1 \text{ mol/L}$

Working electrode	RDE
Stirrer speed (rpm)	2000
Mode	DP
Purge time	30 s
<i>Pretreatment</i>	
Hydrodynamic measurement	no
Cleaning potential	-1.5 V
Cleaning time	30 s
Deposition potential	+0.7 V
Deposition time	30 s
<i>Sweep</i>	
Equilibration time	5 s
Start potential	0 V
End potential	+0.8 V
Voltage step	0.006 V
Voltage step time	0.1 s
Sweep rate	0.06 V/s
Pulse amplitude	0.05 V
Pulse time	0.04 s

Cell off after measurement	yes
<i>Determination</i>	
No. of additions	0
No. of replications	10

Example of the background current in NaOH



After cleaning in NaOH the electrode must always be thoroughly conditioned in the conditioning solution as described above.

- If the gold electrode is not to be used for more than one day then it is recommended that it is rinsed thoroughly with ultrapure water and stored in a dry condition.
- If after cleaning in NaOH and conditioning in conditioning solution the gold electrode can no longer be regenerated then it can be polished with aluminum oxide (Metrohm Order No. 6.2802.000). The electrode must then again be conditioned thoroughly as described above (time required approx. 0.5 to 1 day). However, after polished the performance of a new electrode can no longer be achieved (background current, detection limits).

Method 1: Determination of mercury in samples with a low chloride concentration

Theory

Mercury is determined by anodic stripping voltammetry (ASV) at the rotating gold electrode. Mercury ions are deposited on the gold electrode by reduction to Hg^0 . In the stripping step the deposited Hg^0 is again oxidized to mercury ions.

This method is only suitable for samples with a low chloride concentration ($c(\text{Cl}^-) < 0.01 \text{ mol/L}$). For samples with a higher chloride concentration or in the presence of other halides or pseudohalides the positive potential range of the method must be limited, as otherwise the gold electrode will be irreversibly damaged by oxidation.

Reagents

- Sulfuric acid, suprapur, $w(\text{H}_2\text{SO}_4) = 96\%$, CAS 7664-93-9
- Ethylenediaminetetraacetic acid disodium salt dihydrate, for analysis, $\text{Na}_2\text{C}_{10}\text{H}_{14}\text{N}_2\text{O}_8 \cdot 2\text{H}_2\text{O}$, CAS 6381-92-6
- Sodium chloride, suprapur, NaCl , CAS 7647-14-5
- Hg standard stock solution, $\beta(\text{Hg}^{2+}) = 1 \text{ g/L}$ (available commercially)
- Nitric acid, suprapur, $w(\text{HNO}_3) = 65\%$
- Ultrapure water type 1 (electrical resistivity $> 18.2 \text{ M}\Omega \cdot \text{cm}$, $\text{TOC} < 10 \text{ ppb}$)

Ready-to-use solutions:

Supporting electrolyte	$c(\text{H}_2\text{SO}_4) = 2 \text{ mol/L}$ $c(\text{Na}_2\text{EDTA}) = 0.02 \text{ mol/L}$ $c(\text{NaCl}) = 0.05 \text{ mol/L}$
	<i>0.372 g Na_2EDTA and 0.146 g NaCl are dissolved in approx. 40 mL ultrapure water and 5.56 mL H_2SO_4. After cooling down to room temperature they are made up to 50 mL with ultrapure water.</i>
Hg standard solution	$\beta(\text{Hg}^{2+}) = 1 \text{ mg/L}$
	<i>0.1 mL Hg standard (1 g/L) and 0.1 mL HNO_3 (65%) are made up to 100 mL with ultrapure water.</i>

Determination

Measuring solution

10 mL sample
+ 1 mL supporting electrolyte

Voltammetric parameters

Working electrode	RDE
Stirrer speed (rpm)	2000
Mode	DP
Purge time	30 s
Pretreatment	
Hydrodynamic measurement	no
Cleaning potential	+0.9 V
Cleaning time	10 s
Deposition potential	+0.37 V
Deposition time	90 s
Sweep	
Equilibration time	5 s
Start potential	+0.4 V
End potential	+0.7 V
Voltage step	0.002 V
Voltage step time	0.1 s
Sweep rate	0.02 V/s
Pulse amplitude	0.05 V
Pulse time	0.04 s
Cell off after measurement	yes
Substances	
Peak position (Hg)	+0.6 V



Potentiostat

Highest current range	10 mA
Lowest current range	1 μA

The concentration in the sample is determined by standard addition.

Example

Determination of mercury in tap water (spiked sample)



Sample volume 10 mL
 Deposition time 90 s
 Result 0.7 $\mu\text{g/L}$ Hg

Remarks

- In order to obtain reproducible curves it is important that the gold electrode is first conditioned. See «Conditioning the gold electrode». Depending on the sample matrix, and if longer deposition times are used, it may be necessary to repeat this conditioning at regular intervals between the determinations.
- $c(\text{NaCl}) = 3 \text{ mol/L}$ must be used as the electrolyte in the electrolyte vessel (6.1245.010) of the reference electrode. If $c(\text{KCl}) = 3 \text{ mol/L}$ is used then sparingly soluble potassium perchlorate is formed which can then block the reference electrode diaphragm.
- If the gold electrode is not to be used for more than one day then it is recommended that it is rinsed thoroughly with ultrapure water and stored in a dry condition.
- The presence of oxygen does not interfere with the determination of mercury, as oxygen is only reduced at a potential of -0.1 V . The measuring solution does not need to be degassed for this reason.
- The pH of the measuring solution should be $\text{pH} < 1$. At $\text{pH} > 1$ the background current increases and from pH 2 completely hides the Hg peak.
- Determination limit and linear range apply under the assumption that a new, well-conditioned gold electrode and 10 mL sample are used, and that the sample matrix does not cause any interference:

	Deposition time		
	30 s	90 s	180 s
Limit of quantification	2 $\mu\text{g/L}$	0.4 $\mu\text{g/L}$	0.2 $\mu\text{g/L}$
Linear range up to	40 $\mu\text{g/L}$	15 $\mu\text{g/L}$	6 $\mu\text{g/L}$

- The supporting electrolyte described above can also be prepared using perchloric acid as an alternative to sulfuric acid. All other parameters remain unchanged.

Supporting electrolyte HClO_4	$c(\text{HClO}_4) = 1 \text{ mol/L}$ $c(\text{Na}_2\text{EDTA}) = 0.02 \text{ mol/L}$ $c(\text{NaCl}) = 0.05 \text{ mol/L}$
	<i>0.372 g Na_2EDTA and 0.146 g NaCl are dissolved in approx. 40 mL ultrapure water and 4.2 mL HClO_4^[1] and made up to 50 ml with ultrapure water.</i>

^[1] Perchloric acid, suprapur, $w(\text{HClO}_4) = 70\%$

- If no interference or matrix effects occur during the determination (e.g. in a standard solution) then the background current at the start of the curve (potential range +0.4 V to +0.5 V) typically lies within the range 150 nA ... 400 nA (with a deposition time of 60 s).
- Possible interferences to the peak of $\beta(\text{Hg}^{2+}) = 10 \mu\text{g/L}$ have been investigated under the conditions given above for the following ions. If nothing else is mentioned then the tests have been carried out up to an excess of 1000:1.

Interfering metals

Element	Remarks
Al^{3+}	Does not interfere.
Bi^{3+}	Does not interfere. Peaks at +0.17 V and +0.29 V
Ca^{2+}	Does not interfere, tested up to 10000:1.
Cd^{2+}	Does not interfere. Peak at -0.2 V
Co^{+2}	Does not interfere.
Cr^{3+}	Does not interfere.
Cr^{6+}	From a ratio of approx. 100:1 shoulder on the Hg peak. Peak at +0.49 V
Cu^{2+}	Does not interfere. Peak at +0.24 V
Fe^{3+}	Does not interfere. Peak at +0.44 V
Mg^{2+}	Does not interfere. From approx. 100 mg/L an increase in the background current becomes apparent in the range +0.4 V to +0.5 V.
Mn^{2+}	Does not interfere.
Ni^{2+}	Does not interfere.
Pb^{2+}	Does not interfere. Peak at -0.1 V
Sb^{3+}	Up to a ratio of approx. 100:1 the interfering peak at +0.56 V can be ignored. Peaks at +0.12 V (main signal), +0.4 V and +0.56 V.
Sb^{5+}	Does not interfere. Peak at +0.12 V
Se^{4+}	Does not interfere, but the Hg peak is slightly displaced to a more negative potential and the background current increases.
Se^{6+}	Does not interfere.
$\text{Sn}^{2+}/\text{Sn}^{4+}$	Does not interfere. Peak at -0.16 V
Zn^{2+}	Does not interfere.

Interfering non-metals

Element	Remarks
Halides, Pseudo-halides	Make oxidation easier and therefore destroy the gold electrode at positive potentials by the formation of gold halides. It may be necessary to adapt the voltammetric parameters. It is important that the positive potential range of the method is limited. This is done by displacing the parameters for «Cleaning potential» and «End potential» to more negative potentials. The size of this displacement depends on the type and concentration of the halide. If the Hg peak is greatly displaced then the «Deposition potential» must also be adapted. The rule is that the deposition potential must be at least 0.2 V more negative than the peak potential of mercury.
Cl^-	Displaces the Hg peak to more negative potentials. With $c(\text{Cl}^-) > 0.01 \text{ mol/L}$ the voltammetric parameters must be adapted (see halides). In cases of doubt use the parameters given under «Determination of mercury in seawater».
I^-	Displaces the Hg peak to more negative potentials. At higher concentrations the voltammetric parameters must be adapted (see halides).
NO_2^-	Peak at +0.6 V. The peak generated by $c(\text{NO}_2^-) = 0.001 \text{ mol/L}$ corresponds to a concentration of about $\beta(\text{Hg}^{2+}) = 40 \mu\text{g/L}$. At higher concentrations nitrous gases are formed.
NO_3^-	Does not interfere, tested up to 2 mol/L.
O_2	Does not interfere. Peak at -0.1 V.
PO_4^{3-}	Does not interfere up to a concentration of about 0.1 mol/L. At higher concentrations the increase in background current produced low-bias results.
SO_4^{2-}	Does not interfere, tested up to 2 mol/L.

Method 2: Determination of mercury in seawater or other matrices containing chloride

Theory

The concentration of chloride in seawater is approx. 0.5 mol/L. This concentration is sufficient to oxidize and therefore destroy the gold electrode under the measuring conditions given in Method 1. This is why the voltammetric parameters must be adapted. The parameters given here can also be used for the determination of other samples with a similar chloride concentration.

Reagents + Determination

See «Method 1: Determination of mercury in samples with a low chloride concentration».

Voltammetric parameters

Working electrode	RDE
Stirrer speed (rpm)	2000
Mode	DP
Purge time	30 s
Pretreatment	
Hydrodynamic measurement	no
Cleaning potential	+0.7 V
Cleaning time	10 s
Deposition potential	0 V
Deposition time	60 s
Sweep	
Equilibration time	5 s
Start potential	+0.15 V
End potential	+0.6 V
Voltage step	0.002 V
Voltage step time	0.1 s
Sweep rate	0.02 V/s
Pulse amplitude	0.05 V
Pulse time	0.04 s
Stand-by potential	no
Substances	
Peak position (Hg)	+0.4 V



Potentiostat

Highest current range	10 mA
Lowest current range	1 μ A

The concentration in the sample is determined by standard addition.

Example

Determination of mercury in low-nutrient seawater (spiked sample)



Sample volume 10 mL
 Deposition time 60 s
 Result 1.1 μ g/L Hg

Remarks

- Under certain circumstances even higher chloride concentrations or the presence of other halides or pseudohalides may make a further adaptation of the parameters necessary (see also the remarks about under «interfering halides» under «Method 1: Determination of mercury in samples with a low chloride concentration»).

Literature

- Sipos L., Golimowski J., Valenta P., Nürnberg H.W.
 New voltammetric procedure for the simultaneous determination of copper and mercury in environmental samples
 Fresenius Z Anal Chem 298(1979) 1-8
- Kapel M., Komaitis M.E.
 Polarographic determination of trace elements in food from a single digest.
 Analyst 104(1979) 124-135
- Nürnberg H.W., Valenta P. Sipos L., Branica M.
 The reliable determination of mercury traces in sea water by subtractive differential pulse voltammetry at the twin gold electrode
 Anal Chim Acta 115(1980) 25-42

- Ahmed R., Valenta P. Nürnberg H.W.
 Voltammetric determination of mercury levels in tuna fish
 Mikrochim Acta (1981) 171-184
- Ireland-Ripert J., Bermond A., Ducauze C.
 Determination of methylmercury in the presence of inorganic mercury by anodic stripping voltammetry
 Anal Chim Acta 143(1982) 249-254
- Golimowski J., Gustavsson I.
 Determination of mercury in fish using differential pulse anodic stripping voltammetry
 Fresenius Z Anal Chem 317(1984) 481
- Leu M., Seiler H.
 AC2 Inversvoltammetrische Bestimmung von Quecksilber in Harn
 Fresenius Z Anal Chem 321(1985) 479-482
- Seritti A., Morelli E., Orsini F., Nannicini L.
 Heavy metals in seawater in front of chlor-alkali plant
 Marine Pollution Bulletin 18/8(1987) 461-463

Annex

Report for example on «Determination of Hg in tap water»

```

===== METROHM 797 VA COMPUTRACE (Version 1.3.0.77) (Serial No. 5195) =====
Determination : 0808141410_Tap water spiked.dth
Sample ID      : Tap water spiked
Creator method: zu           Date : 2008-05-09           Time: 14:20:13
Creator determ.: zu         Date : 2008-08-14           Time: 14:10:51
Modified by   : ---         Date :                   Time:
-----
Method        : AB096_1_Det of Hg with Au RDE.mth
Title         : Determination of Mercury with Gold RDE
Remark1      : 10 mL sample+ 1 mL electrolyte
Remark2      : Electrolyte: c(H2SO4) = 2 M, c(EDTA) = 0.02 M, c(NaCl) = 0.05 M
-----
Sample amount : 10.000 mL
Cell volume   : 11.000 mL
-----
Substance     : Hg
Conc.         : 656.190 ng/L
Conc.dev.     : 12.338 ng/L      ( 1.88%)
Amount       : 7.218 ng
Add.amount    : 10.000 ng
-----

```

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	0.600	44.2	45.0	1.087	0.0	
1 - 2	0.604	45.8				
2 - 1	0.602	106.5	106.8	0.437	61.8	
2 - 2	0.600	107.1				
3 - 1	0.602	170.2	169.7	0.672	62.9	
3 - 2	0.600	169.2				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Hg	std.add.	4.497e-008	6.853e-002	1.259e-009	0.99993
Final results			+/- Res. dev.	%	Comments
Hg: Mercury	=	0.722 ug/L	0.014	1.880	

Report for example on «Determination of Hg in low-nutrient seawater»

```

===== METROHM 797 VA COMPUTRACE (Version 1.3.1.84) (Serial No. 5195) =====
Determination : 0808130959_Low nutrient seawater spiked.dth
Sample ID      : Low nutrient seawater spiked
Creator method :                               Date : 2008-08-12           Time: 16:44:54
Creator determ.: zu                           Date : 2008-08-13           Time: 09:59:52
Modified by    : ---                           Date :                       Time:
  
```

```

Method          : AB096_2 Det of Hg in sea water with Au RDE.mth
Title           : Determination of Mercury in sea water with Gold RDE
Remark1         : 10 mL seawater + 1 mL electrolyte
Remark2         : Electrolyte: c(H2SO4) = 2 M, c(EDTA) = 0.02 M, c(NaCl) = 0.05 M
  
```

```

Sample amount   : 10.000 mL
Cell volume     : 11.000 mL
  
```

```

Substance       : Hg
Conc.           : 1.028 ug/L
Conc.dev.       : 0.003 ug/L      ( 0.30%)
Amount          : 11.311 ng
Add.amount      : 10.000 ng
  
```

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	0.370	44.7	44.7	0.027	0.0	
1 - 2	0.372	44.7				
2 - 1	0.372	83.3	83.3	0.041	38.6	
2 - 2	0.372	83.3				
3 - 1	0.372	121.4	121.6	0.259	38.3	
3 - 2	0.372	121.8				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Hg	std.add.	4.469e-008	4.346e-002	1.625e-010	0.99999
Final results			+/- Res. dev.	%	Comments
Hg: Mercury	=	1.131 ug/L	0.003	0.304	

Method printout for «Conditioning the gold electrode» at the 797 VA Computrace

Method parameters			
Method	:	AB096_Conditioning of Au electrode for Hg determination.mth	
Title	:	Conditioning of Au electrode for Hg determination	
Remark1	:	20 mL H2O + 0.2 mL HClO4 (70%)	
Remark2	:		
Calibration	:	Standard addition	
Technique	:	Batch	
Addition	:	Manual	
Sample ID	:	Cond Au electrode	
Sample amount (mL)	:	20.000	
Cell volume (mL)	:	20.200	
Voltammetric parameters			
Mode	:	DP - Differential Pulse	
Highest current range	:	10 mA	
Lowest current range	:	1 uA	
Electrode	:	SSE/RDE	
Stirrer speed (rpm)	:	2000	
Initial electr. conditioning	:	No	
No. of additions	:	0	
No. of replications	:	10	
Measure blank	:	No	
Addition purge time (s)	:	10	
Initial purge time (s)	:	30	
Conditioning cycles	:		
Start potential (V)	:	-1.200	
End potential (V)	:	-0.100	
No. of cycles	:	0	
Hydrodynamic (measurement)	:	No	
Cleaning potential (V)	:	0.000	
Cleaning time (s)	:	10.000	
Deposition potential (V)	:	1.500	
Deposition time (s)	:	30.000	
Sweep	:		
Equilibration time (s)	:	5.000	
Start potential (V)	:	0.000	
End potential (V)	:	0.800	
Voltage step (V)	:	0.006	
Voltage step time (s)	:	0.100	
Sweep rate (V/s)	:	0.060	
Pulse amplitude (V)	:	0.050	
Pulse time (s)	:	0.040	
Cell off after measurement	:	Yes	
Peak evaluation			
Regression technique	:	Linear Regression	
Peak evaluation	:	Height	
Minimum peak width (V.steps)	:	5	
Minimum peak height (A)	:	1.000e-010	
Reverse peaks	:	No	
Smooth factor	:	4	
Eliminate spikes	:	Yes	
Substances			
Baseline			
Substance Addition	:	automatic start (V) end (V) type	scope

Method printout for «Cleaning the gold electrode in NaOH» at the 797 VA Computrace

Method parameters			
Method	:	AB096_Cleaning of Au electrode in NaOH.mth	
Title	:	Cleaning of Au electrode in NaOH	
Remark1	:	20 mL NaOH	
Remark2	:	c(NaOH) = 0.1 mol/L	
Calibration	:	Standard addition	
Technique	:	Batch	
Addition	:	Manual	
Sample ID	:	Cleaning Au-RDE	
Sample amount (mL)	:	20.000	
Cell volume (mL)	:	20.000	
Voltammetric parameters			
Mode	:	DP - Differential Pulse	
Highest current range	:	10 mA	
Lowest current range	:	1 uA	
Electrode	:	SSE/RDE	
Stirrer speed (rpm)	:	2000	
Initial electr. conditioning	:	No	
No. of additions	:	0	
No. of replications	:	10	
Measure blank	:	No	
Addition purge time (s)	:	10	
Initial purge time (s)	:	30	
Conditioning cycles	:		
Start potential (V)	:	-1.200	
End potential (V)	:	-0.100	
No. of cycles	:	0	
Hydrodynamic (measurement)	:	No	
Cleaning potential (V)	:	-1.500	
Cleaning time (s)	:	30.000	
Deposition potential (V)	:	0.700	
Deposition time (s)	:	30.000	
Sweep	:		
Equilibration time (s)	:	5.000	
Start potential (V)	:	0.000	
End potential (V)	:	0.800	
Voltage step (V)	:	0.006	
Voltage step time (s)	:	0.100	
Sweep rate (V/s)	:	0.060	
Pulse amplitude (V)	:	0.050	
Pulse time (s)	:	0.040	
Cell off after measurement	:	Yes	
Peak evaluation			
Regression technique	:	Linear Regression	
Peak evaluation	:	Height	
Minimum peak width (V.steps)	:	5	
Minimum peak height (A)	:	1.000e-010	
Reverse peaks	:	No	
Smooth factor	:	4	
Eliminate spikes	:	Yes	
Substances			
Baseline			
Substance Addition	:	automatic start (V) end (V) type	scope

Method printout for «Determination of mercury in samples with a low chloride concentration » at the 797 VA Computrace

```

-----
Method parameters
-----
Method       : AB096_1 Det of Hg with Au RDE.mth
Title        : Determination of Mercury with Gold RDE
Remark1      : 10 mL sample + 1 mL electrolyte
Remark2      : Electrolyte: c(H2SO4) = 2 M, c(EDTA) = 0.02 M, c(NaCl) = 0.05 M

Calibration  : Standard addition
Technique    : Batch
Addition     : Manual

Sample ID    : sample
Sample amount (mL): 10.000
Cell volume (mL): 11.000

-----
Voltammetric parameters
-----
Mode          : DP - Differential Pulse
Highest current range : 10 mA
Lowest current range  : 1 uA
Electrode      : SSE/RDE
Stirrer speed (rpm)  : 2000
Initial electr. conditioning : No

No. of additions      : 2
No. of replications   : 2

Measure blank        : No
Addition purge time (s) : 20
Initial purge time (s) : 30

Conditioning cycles
Start potential (V)   : 0.300
End potential (V)     : 0.700
No. of cycles         : 0

Hydrodynamic (measurement) : No
Cleaning potential (V)      : 0.900
Cleaning time (s)           : 10.000
Deposition potential (V)    : 0.370
Deposition time (s)         : 90.000
Sweep
Equilibration time (s)      : 10.000
Start potential (V)         : 0.400
End potential (V)           : 0.700
Voltage step (V)            : 0.002
Voltage step time (s)       : 0.100
Sweep rate (V/s)           : 0.020
Pulse amplitude (V)        : 0.050
Pulse time (s)             : 0.040
Cell off after measurement  : Yes

-----
Peak evaluation
-----
Regression technique      : Linear Regression
Peak evaluation           : Height
Minimum peak width (V.steps) : 10
Minimum peak height (A)   : 1.000e-010
Reverse peaks             : No
Smooth factor             : 1
Eliminate spikes          : Yes

-----
Substances
-----
Hg                        : 0.600 V +/- 0.050 V

Standard solution         : 1 1.000 mg/L
Addition volume (mL)     : 0.010

Mercury                   : Final result (Hg) =
                          Conc * (11 / 10) * (1e+006 / 1) + 0 - 0

-----
Baseline
-----
Substance Addition      automatic start (V) end (V) type      scope
-----
Hg      Sample          yes      ---      ---      linear      wholePeak
        Addition 1      yes      ---      ---      linear      wholePeak
        Addition 2      yes      ---      ---      linear      wholePeak
-----

```

Method printout for «Conditioning the gold electrode» at the 746 VA Trace Analyzer

```

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB96_CON.mth          OPERATION SEQUENCE
Title : Conditioning of Au electrode for Hg determination
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  DOS>M
2  DOS>M
3  STIR          30.0  Rot.speed       2000 /min
4  (REP
5  SEGMENT
6  SEGMENT
7  REP)8
8  OMEAS
9  END

Method: AB96_CON          SEGMENT
                          COND_Au
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  RDE          3.0  Rot.speed       2000 /min
2  DPMODE
   U.ampl       50 mV
   t.step       0.10 s
3  MEAS         10.0  U.meas          0 mV
4  MEAS         30.0  U.meas          1500 mV
5  OSTIR        5.0
6  SWEEP        13.6  U.start         0 mV
   U.end        800 mV
   U.standby    mV
7  OMEAS
8  END

Method: AB96_CON          DOCUMENTATION
-----
COPY  Reports, Curves  TO Destination
-----
Curve Smth Sgmt:COND_Au VR:**  Printer
  
```

Method printout for «Cleaning the gold electrode in NaOH» at the 746 VA Trace Analyzer

```

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB96_cle.mth         OPERATION SEQUENCE
Title : Cleaning of Au electrode in NaOH
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  DOS>M
2  STIR          30.0  Rot.speed       2000 /min
3  (REP
4  SEGMENT
5  SEGMENT
6  REP)8
7  OMEAS
8  END

Method: AB96_cle          SEGMENT
                          Clean_Au
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  RDE          3.0  Rot.speed       2000 /min
2  DPMODE
   U.ampl       50 mV
   t.step       0.10 s
3  MEAS         30.0  U.meas          -1500 mV
4  MEAS         30.0  U.meas          700 mV
5  OSTIR        5.0
6  SWEEP        13.6  U.start         0 mV
   U.end        800 mV
   U.standby    mV
7  OMEAS
8  END

Method: AB96_cle          DOCUMENTATION
-----
COPY  Reports, Curves  TO Destination
-----
Curve Smth Sgmt:Clean_Au VR:**  Printer
  
```

Method printout for «Determination of mercury in samples with a low chloride concentration » at the 746 VA Trace Analyzer

```

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB96_1 .mth OPERATION SEQUENCE
Title : Determination of Mercury with Gold RDE
-----
Instructions t/s Main parameters Auxiliary parameters
-----
1 SMPL>M V.fraction mL V.total L
2 DOS>M Soln.name electrol V.add 1.000 mL
3 STIR 20.0 Rot.speed 2000 /min
4 (ADD
5 STIR 10.0 Rot.speed 2000 /min
6 (REP
7 SEGMENT Segm.name ASV_Hg
8 REP)1
9 ADD>M Soln.name Hg_std V.add 0.010 mL
10 ADD)2
11 OMEAS U.standby mV
12 END

Method: AB96_1 SEGMENT
ASV_Hg
-----
Instructions t/s Main parameters Auxiliary parameters
-----
1 RDE 3.0 Rot.speed 2000 /min
2 DPMODE U.ampl 50 mV t.meas 20.0 ms
t.step 0.10 s t.pulse 40.0 ms
3 MEAS 10.0 U.meas 900 mV
4 MEAS 90.0 U.meas 370 mV
5 OSTIR 10.0
6 SWEEP 15.3 U.start 400 mV U.step 2 mV
U.end 700 mV Sweep rate 20 mV/s
7 OMEAS U.standby mV
8 END

Method: AB96_1 DOCUMENTATION
-----
COPY Reports, Curves TO Destination
-----
Report Full Printer
Curve Smth Sbst:Mercury VR:** Printer
Curve Cal Sbst:Mercury Printer

Method: AB96_1 SUBSTANCES
Mercury - ASV_Hg
-----
Recognition Display / Plot
-----
U.verify 600 mV I.scale auto
U.tol (+/-) 50 mV U.div 50.00 mV/cm
U.width min 10 mV U.begin mV
U.width max 200 mV U.end mV
I.threshold 200 pA

Baseline Evaluation
-----
Type linear Mode VA
Scope whole Quantity I.peak
dU.front auto Sign. digits 4
S.front auto
dU.rear auto
S.rear auto

Calibration 2008-09-30 09:57:01 Coefficients
-----
Technique std.add. Y.reg 5.976e-08
Curve type linear Slope 0.08036
Nonlin.
Mean dev. 8.48e-10

Additions
-----
Soln.name Hg_std
-----
Mass conc. 1 mg/L g/L g/L g/L
Range min g/L g/L g/L g/L
Range max g/L g/L g/L g/L
M.conc./cm g/L g/L g/L g/L

Method: AB96_1 CALCULATION
max. 15 lines
-----
Quantity Formula (R##, C##, A##) Res.unit Sig.dig.
-----
Mercury R1000=MC:Mercury #g/L 5
  
```