

Determination of manganese in water samples by anodic stripping voltammetry

Summary

A sensitive methods to determine manganese is described. It is primarily suitable for the investigation of ground, drinking and surface waters, in which the concentration of manganese is important. The method can naturally also be used for trace analysis in other matrices.

Manganese is determined in an alkaline borate buffer by the anodic stripping voltammetry (ASV). Interference by intermetallic compounds is prevented by the addition of zinc ions in the sample. The limit of determination lies at $\beta(\text{Mn}) = 2 \mu\text{g/L}$.

Instruments

VA instrument capable of operating a Multi-Mode Electrode and supporting differential pulse (DP) measuring mode	
909 UV Digester	2.909.0014

Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Mercury drop capillary	6.1226.030
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel Filled with $c(\text{KCl}) = 3 \text{ mol/L}$	6.1245.010
AE	Pt rod electrode	6.0343.x00

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Hydrochloric acid, $w(\text{HCl}) = 30 \%$, for trace analysis*, CAS 7647-01-0
- Ammonia solution, $w(\text{NH}_3) = 25 \%$, for trace analysis*, CAS 1336-21-6
- Sodium hydroxide solution, $w(\text{NaOH}) = 30\%$, for trace analysis*, CAS 1310-73-2
- Sodium tetraborate decahydrate, $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10 \text{ H}_2\text{O}$, for analysis*, CAS 1303-96-4

- Mn standard stock solution, $\beta(\text{Mn}^{2+}) = 1 \text{ g/L}$ commercially available
 - Zn standard stock solution, $\beta(\text{Zn}^{2+}) = 1 \text{ g/L}$ commercially available
 - Ultrapure water, resistivity $>18 \text{ M}\Omega \cdot \text{cm}$ (25°C), type I grade (ASTM D1193)
- * e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

Electrolyte	$c(\text{Na}_2\text{B}_4\text{O}_7) = 0.1 \text{ mol/L}$ $c(\text{NaOH}) = 0.3 \text{ mol/L}$ 3.81 g $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10 \text{ H}_2\text{O}$ are dissolved in 50 mL water. Add 3 mL sodium hydroxide solution and fill up to 100 mL with ultrapure water.
Zinc solution	$\beta(\text{Zn}) = 100 \text{ mg/L}$ Add 100 μL conc. hydrochloric acid to 10 mL Zn stock solution and fill up to 100 mL with water.

Standard solution

Mn standard solution	$\beta(\text{Mn}) = 10 \text{ mg/L}$ Add 100 μL conc. hydrochloric acid to 1 mL Mn stock solution and fill up to 100 mL with water.
----------------------	---

Sample preparation

- Ground water, surface waters, mineral waters and drinking waters can usually be analyzed without pretreatment.

Organic matter often interferes with voltammetric determinations and therefore sample solutions usually have to be digested.

- Low polluted waste waters can be digested with the 909 UV Digester.

- o Add 50 μL hydrogen peroxide solution $w(\text{H}_2\text{O}_2) = 30\%$ and 10 μL hydrochloric acid $w(\text{HCl}) = 30\%$ to 10 mL acidified sample ($\text{pH} = 2$) and irradiate for 90 minutes at 90°C .
- Samples with organic matter (foods, pharmaceuticals etc.) must be digested.
 - o High-pressure asher
 - o Microwave digestion
 - o Both techniques oxidize the sample in a closed digestion vessel by means of a mixture of concentrated mineral acids.
 - o Open wet digestion with H_2SO_4 and H_2O_2 according to Application Bulletin 113.

Analysis

Measuring solution

10 mL (diluted) sample

50 μL ammonia solution

2.5 mL supporting electrolyte

10 μL Zn solution



If necessary, the pH is adjusted to 9.5 ... 10 with ammonia solution.

The concentration is determined by the standard addition method.

Parameters

Voltammetric	
Electrode operating mode	HMDE
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
Potentiostatic pretreatment	
Potential 1	-1.7 V
Waiting time 1	90 s
Equilibration time	5 s
Sweep	
Start potential	-1.62 V
End potential	-1.25 V
Potential step	0.004 V
Potential step time	0.5 s
Sweep rate	0.008 V/s
Pulse amplitude	-0.075 V
Substance	
Name	Mn
Characteristic potential	-1.44 V

Example



Results

Sample	Tap water
Sample size	10.0 mL
$\beta(\text{Mn})$	57.8 $\mu\text{g/L}$

Comments

It is essential to use a negative pulse amplitude for this determination.

References

- Monien H., Jacob P.
Voltammetrische Bestimmung kleiner Eisenmengen ohne Abtrennung der Matrix.
Fresenius Z. Anal. Chem. 260, (1972) 195-202
- Davidson W.
Comparison of differential pulse and D.C. sampled polarography for the determination of ferrous and manganese ions in lake water.
J. Electroanal. Chem. 72, (1976) 229-237
- Colombini M.P., Fluoco R.
Determination of manganese in ng/ml levels in natural waters by differential pulse polarography.
Talanta 30, (1983) 901-905; Ref: Metrohm Info 3/84, 14
- Huang Z.Q., Van den Berg C.M.G.
Determination of iron in seawater using cathodic stripping voltammetry preceded by adsorptive collection with the hanging mercury drop electrode.
J. Electroanal. Chem. 177, (1984) 269-280; Ref: Metrohm-Info 1/85, 19
- Weidenauer M., Lieser K.H.
Bestimmung von Spurenelementen in Flusswasser mit Hilfe der Voltammetrie.
Fresenius, Z. Anal. Chem. 320, (1985) 550-555
- Wang J., Mahmoud J.
Chelate adsorption of trace voltammetric measurements of iron(III).
Fresenius, Z. Anal., Chem. 327, (1987) 789-793
- Labuda J., Vanickova M., Beinrohr E.
Determination of dissolved manganese in natural waters by differential pulse cathodic stripping voltammetry.
Mikrochim. Acta 1989/I, 113-120; Ref: Fresenius, Z. Anal. Chem. 335, (1989) 433

Appendix

Report for the example determination of manganese in tap water

```
===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Determ.      : 05180901      User:      Date: 1999-05-18
Modified     : 2000-11-30 18:27:52 Run : 0   Time: 09:01:26
Sample table: -
```

Pos.	Ident.1/S1 water	Ident.2/S2	Ident.3/S3	Method.call	Sample size/S0 10 mL
------	---------------------	------------	------------	-------------	-------------------------

```
-----
Method : AB123
Title  : Determination of Manganese in Waters. AB123
Remark1: Determination of manganese in water
Remark2: 10 ml sample + NH3 buffer + borate buffer + Zn standard
-----
```

Substance	Mn	Mass conc.	57.83 ug/L	Mass	578.3 ng	Comments
MC.dev.	1.28 ug/L (2.22%)	Add.mass	250 ng			
Cal.dev.	-	V0.sample	10 mL			

VR	U/mV	I/nA	I.mean	Std.dev.	I.delta	Comments
00	-1444	-74.42	-74.10	0.4504		
01	-1445	-73.78				
10	-1445	-104.7	-104.4	0.4077	-30.29	
11	-1445	-104.1				
20	-1445	-137.2	-137.2	0.1228	-32.86	
21	-1445	-137.3				

Substance	Techn.	Y.reg/offset	Slope	Nonlin.	Mean deviat.
Mn	std.add.	-7.365e-08	-0.001600		8.389e-10

SOLUTIONS
max. 40

Soln.name	Pos.	Std.subst.	Mass conc.	Remark
PbStd	-	Pb	1 g/L	from: Pb

C# Workg.com.var Remark

Final results	+/-	Res.dev.	%	Comments
Mn = 57.825 ug/L	1.28	2.22		

Method print for the determination of manganese

```
===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB123 .mth      OPERATION SEQUENCE
Title : Determination of Manganese in Waters. AB123
-----
```

Instructions	t/s	Main parameters	Auxiliary parameters
1 SMPL>M		V.fraction mL	V.total L
2 DOS>M		Soln.name ammonia	V.add 0.050 mL
3 DOS>M		Soln.name borate	V.add 2.500 mL
4 DOS>M		Soln.name Zn-std	V.add 0.010 mL
5 PURGE			
6 STIR	300.0	Rot.speed 2000 /min	
7 (ADD			
8 PURGE			
9 STIR	20.0	Rot.speed 2000 /min	
10 ØPURGE			
11 SEGMENT		Segm.name DPAdSV	
12 ADD>M		Soln.name MnStd	V.add 0.025 mL
13 ADD)2			
14 END			

Method: AB123_2 SEGMENT
DPAdSV

	Instructions	t/s	Main parameters		Auxiliary parameters	
1	(REP					
2	STIR		Rot.speed	2000 /min		
3	DME	10.0				
4	HMDE		Drop size	7	Meas.cell	normal
5	HMDE		Drop size	7	Meas.cell	normal
6	HMDE		Drop size	7	Meas.cell	normal
7	HMDE		Drop size	4	Meas.cell	normal
8	DPMODE		U.ampl	-75 mV	t.meas	20.0 ms
			t.step	0.50 s	t.pulse	40.0 ms
9	MEAS	90.0	U.meas	-1700 mV		
10	OSTIR	5.0				
11	SWEEP	48.0	U.start	-1620 mV	U.step	4 mV
			U.end	-1250 mV	Sweep rate	8 mV/s
12	OMEAS		U.standby	mV		
13	REP)1					
14	END					

Method: AB123_2

SUBSTANCES
Mn - DPAdSV

Recognition

U.verify	-1440 mV
U.tol (+/-)	50 mV
U.width min	10 mV
U.width max	200 mV
I.threshold	200 pA

Display / Plot

I.scale	auto
U.div	50.00 mV/cm
U.begin	mV
U.end	mV

Baseline

Type	linear
Scope	whole
dU.front	auto
S.front	auto
dU.rear	auto
S.rear	auto

Evaluation

Mode	VA
Quantity	I.peak
Sign. digits	4

Calibration 2000-11-30 18:31:53

Technique std.add.
Curve type linear

Coefficients

Y.reg	-7.365e-08
Slope	-0.0016
Nonlin.	
Mean dev.	8.389e-10

Additions

Soln.name	MnStd			
Mass conc.	10 mg/L	g/L	g/L	g/L
Range min	g/L	g/L	g/L	g/L
Range max	g/L	g/L	g/L	g/L
M.conc./cm	g/L	g/L	g/L	g/L

Method: AB123_2

CALCULATION
max. 15 lines

Quantity	Formula (R##, C##, A##)	Res.unit	Sig.dig.
Mn	R1000=MC:Mn	#g/L	5