Application Bulletin

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Electrochemical detection in HPLC *

General considerations and application examples of the oxidative determination of organic compounds

Introduction

The Operating Instructions for the METROHM 656 Electrochemical Detector provide the user with all basic information concerning the functioning and operation of the instrument, together with treatment of the electrodes. Hints on separating system requirements and on the causes and cures of detection disturbances are also to be found therein.

The present Application Bulletin is intended on the one hand to give a synoptic view of the most important classes of substances and compounds which can be well detected by oxidative methods, i.e. with detection limits in the pg (picogramme, 10⁻¹² g) range, and on the other hand to illustrate the possible working conditions for separation and electrochemical detection by means of some examples.

HPLC apparatus with electrochemical detection

In principle, the 656 Electrochemical Detector can be used in conjunction with all HPLC systems.

For successful results, particularly in the pg range, all components, such as pump, dampening system, injector, column, capillaries, tube fittings, etc., must be in perfect condition. With electrochemical detection, disturbances due to components of the HPLC system which would pass unnoticed with UV detection may be recorded. These are often due to contamination from the system (column) or to pressure fluctuations due to imperfect fittings.

Particular attention should be paid to the following points:

▶ **Cleanliness**

The usual precautions for trace analysis should be observed. All corrosion of steel parts should be avoided, since ferrous (Fe²⁺) ions can interfere. Eluents should not contain complex-forming substances.

▶ **Pulsation**

"Noise" in the detector signal may be produced by pressure fluctuations in the detector cell. An additional pulse dampener working on the principle of the Bourdon tube has shown itself to be of value in such cases, and functions quite adequately even where a single-piston pump is used.

▶ **Eluent**

Many of the solvents normally used for HPLC can be used in making up the eluent. The extra conductivity required for electrochemical detection is obtained by adding an "electrolyte", generally in concentrations from 1...10 g/L, giving conductivities of about 1...10 mS/cm. Among others, the following can be used: sulfates, nitrates, phosphates, acetic acid, sulfuric acid, perchloric acid and lithium perchlorate.

Chlorides and hydroxyacids should be avoided, while ion-pairing reagents can be used in certain cases.

Gradient elution is possible in the upper sensitivity ranges (current ranges greater than 100 nA).

The eluent is degassed in the usual manner by vacuum or by passing through helium.

▶ **Stationary phase**

In principle, all reversed-phase and ion-exchange materials can be used.

▶ **Injection solution**

Wherever possible, the eluent should be used as a solvent to avoid disturbances due to change of the environment (lowest possible solvent front).

The mass of solute injected should not exceed 100 ng in order to avoid surface filming of the electrode.


* High-Pressure Liquid Chromatography

Electrochemical detection in HPLC *

Important substance-classes

1. Aromatic alcohols

- 1.1 Phenols
- 1.2 Halogenated Phenols
- 1.3 Hydroxyphenyls
- 1.4 Catechols
- 1.5 Methoxyphenols
- 1.6 Hydroxycoumarins
- 1.7 Flavones
- 1.8 Estrogens
- 1.9 Tocopherols
- 1.10 Antioxidants
- 1.11 Others


Voltage U_{pot} *

- + 1200 mV
- + 1200 mV
- + 800 mV
- + 800 mV
- + 800 mV
- + 1000 mV
- + 1000 mV
- + 1000 mV
- + 800 mV
- + 800/1000/1200 mV

2. Aromatic amines


- 2.1 Anilines
- 2.2 Benzidines
- 2.3 Sulfonamides

- + 1000 mV
- + 600 mV
- + 1200 mV


3. Indoles

- 3.1 Indolyl-3-derivatives
- 3.2 5-Hydroxy-indoles


- + 1000 mV
- + 800 mV

4. Phenothiazines

- + 1000 mV

5. Mercaptans

- + 800 mV

6. Others

- 6.1 Ascorbic acid
- 6.2 Vitamin A
- 6.3 Carotenes
- 6.4 Purines

- + 800 mV
- + 1000 mV
- + 800 mV
- + 800 1000 mV


* approximative values

Electrochemical detection in HPLC *

Some examples of compounds belonging to the substance-classes


1. Aromatic alcohols

1.1 Phenols


4-Hydroxybenzoic acid
 4-Hydroximandelic acid
 Hydroxyephedrine
 Phenylephrine
 Salicylic acids
 Synephrin
 Tyrosine
 Tyramine
 Thyroxine
 Thyronine
 2-Phenylphenol

1.2 Halogenated Phenols


Chlorophenols

1.3 Hydroxybiphenyls


Chlorinated Hydroxybiphenyls

1.4 Catechols


Adrenaline/Noradrenaline
 Caffeic acid
 Chlorogenic acid
 Dopamine
 L-Dopa
 Gentisic acids
 Homogentisic acids
 Dihydroxymandelic acid
 Isoprenaline
 Isoproterenol
 Terbutaline

1.5 Methoxyphenols

Homovanillic acid
 Metanephrine
 3-Methoxy-4-hydroxyphenyl-ethyleneglycol
 Normethanephrine
 Vanillin
 Vanillic acid
 Vanillylmandelic acid
 α -Methyldopa, -dopamine


1.6 Hydroxycoumarins

Esculetin
 Esculin
 Scopoletin
 Dicoumarol


Electrochemical detection in HPLC *

1.7 Flavones


Quercetin
Catechins
Cyanidine
Hesperitin
Rutin

1.8 Estrogens

Estrone
Estradiol
Estriol
Ethinylestradiol
Diethylstilbestrol
Stilboestrol
Dienestrol
Hexestrol
Zearalenon
Zeranol

1.9 Tocopherols

α -, β -, γ -, δ -Tocopherols
Tocol

1.10 Antioxidants


tert-Butylhydroxyanisoles (BHA)
Butylphenols
tert-Butylhydroxytoluenes (BHT)
Gallates
Nordihydroguaiaretic acid (NDGA)
tert-Butylhydroxyquinones (TBHQ)

1.11 Others

Morphine
Dihydromorphine
Apomorphine
Tetrahydrocannabinols
Cinchonine
8-Hydroxycarteolol
Salsolinol
 β -Cetotetrine
Folic acid


2. Aromatic amines

2.1 Anilines

Chloroanilines
Bromoanilines
Dichloroanilines
Toluidines
Amino-chlorophenols
p-Phenylenediamine


2.2 Benzidines

Benzidine
substituted Benzidines

Electrochemical detection in HPLC *


2.3 Sulfonamides

Sulfanilamide
Sulfamethoxydiazine
and others


3. **Indoles**

3.1 Indolyl-3-derivatives


Indolyl-3-acetic acid
3-Methylindole
Tryptophan
Tryptamine
Melatonin

3.2 5-Hydroxyindoles


Serotonin
5-Hydroxyindolylacetic acid
5-Hydroxytryptophan

4. **Phenothiazines**

Chlorpromazine
Promethazine
Perphenazine
and others


5. **Mercaptans**


Cysteine
Penicillamine
Glutathione

6. **Others**


6.1 Ascorbic acid

6.2 Vitamin A


6.3 Carotenes

α-Carotene
β-Carotene
Lycopene
Phytoene
Xanthophyll

6.4 Purines

Uric acid
Xanthine
Guanine

Electrochemical detection in HPLC *

Instrumentation

METROHM ELCD System:	641 VA Detector 656 Electrochemical Detector Auxiliary electrode (AE) : Glassy carbon Reference electrode: Silver/silver chloride c(KCl) = 3 mol/L
HPLC System:	Pump: Altex Mod. 110 Pulse dampening system: Bourdon tube principle Lit.: Ventura et al., Anal. Chem. <u>50</u> , 1017 (1978) Loop injection volume: 20 μ L Temperature: 20 °C

Names and symbols

Volume fraction of substance B: $V_B/\sum_i V_i$	ϕ_B
Mass concentration: m_i/V	ρ_i
Volume	V
Glassy carbon	GC

List of the single compounds in the chromatograms

Compound	substance-class	chromatogram/page
Acetaminophen	1.1	10
L-Adrenaline	1.4	13/15
Ascorbic acid	6.1	29/30
3-tert-Butylhydroxyanisole	1.10	21
tert-Butylhydroxytoluene	1.10	22
β -Carotene	6.3	31/32
Catechin	1.7	17
3-Chloroaniline	2.1	23
Chlorpromazine	4.	26
o-Cresol	1.1	11
m-Cresol	1.1	11
p-Cresol	1.1	11
L-Cysteine	5.	27
L-Dopa	1.4	13
Dopamine	1.4	13/14/15
3,4-Dihydroxyphenylacetic acid	1.4	25
2,5-Dimethylphenol	1.1	11/12
3,5-Dimethylphenol	1.1	11

Electrochemical detection in HPLC *

Compound	substance-class	chromatogram/page
Esculin	1.6	16
Estradiol	1.8	19
Estriol	1.8	19
Estrone	1.8	19
L-Glutathione	5.	27
Homovanillic acid	1.5	25
L-Noradrenaline	1.4	13/14/15
Nordihydroguaiaretic acid	1.10	21
Penicillamine	5.	27/28
Pentachlorophenol	1.2	12
Phenol	1.1	11/12
p-Phenylenediamine	2.1	23
4-Phenylphenol	1.1	12
Propyl gallate	1.10	21
Rutin	1.7	17
Quercetin	1.7	18
Sulfanilamide	2.3	24
Sulfamethoxydiazine	2.3	24
2,4,5-Trichlorophenol	1.2	12
2,3,4,5-Tetrachlorophenol	1.2	12
α -Tocopherol	1.9	20/31/32
β -, γ -Tocopherol	1.9	20
δ -Tocopherol	1.9	20
o-Toluidine	2.1	24
D,L-Tryptophan	3.1	25
Uric acid	6.4	29
Vitamine-A-acetate	6.2	31

HPLC/ELCD-Application**Substance**Acetaminophen in serum,
 $\rho = 200 \text{ ng/mL}$ **Phenols****Liquid Chromatographic System**


Column: LiChrosorb RP-18, $5 \mu\text{m}$, $250 \text{ mm} \times 4.6 \text{ mm ID}$
Eluent: methanol-water (140:18), with ammonium acetate and acetic acid
 $\phi(\text{MeOH}) = 0.885$, $\rho(\text{CH}_3\text{COONH}_4) = 5.2 \text{ g/L}$, $\rho(\text{HAc}) = 7 \text{ g/L}$
(pH 5.2)
Flow rate: 0.8 mL/min
Injection volume: $20 \mu\text{L}$

Detector

Working electrode: GC
Voltage: $+800 \text{ mV}$
Range: 5 nA

Sampleserum, $V = 200 \mu\text{L}$, extracted with dichlormethane-2-propanol-ether**Literature**

D.J. Miner et al.
J. Pharm. Sci. 68, 96 (1979)

HPLC/ELCD-Application


Substance

Phenols in wastewater

- 1 Phenol
- 2 p-, m-Cresol
- 3 o-Cresol
- 4 3,5-Dimethylphenol
- 5 2,5-Dimethylphenol

Liquid Chromatographic System

Column: Nucleosil 5 C₁₈, 5 μm, 250 mm × 4.6 mm ID
 Eluent: methanol-water (1:1), with potassium nitrate and sulfuric acid
 $\phi(\text{MeOH}) = 0.5$, $\rho(\text{KNO}_3) = 2 \text{ g L}$, $\rho(\text{H}_2\text{SO}_4) = 0.05 \text{ g L}$
 Flow rate: 0.7 mL/min
 Mass of solute injected: standard: 1 ng each
 sample: filtrated wastewater, $\rho(\text{Phenol}) = 10 \text{ μg/L}$


HPLC/ELCD-Application**Substance****Halogenated Phenols**

- 1 Phenol
- 2 2,5-Dimethylphenol
- 3 4-Phenylphenol
- 4 2,4,5-Trichlorophenol
- 5 2,3,4,5-Tetrachlorophenol
- 6 Pentachlorophenol

Liquid Chromatographic System

Column: Nucleosil 5 C₁₈, 5 μ m, 250 mm \times 4.6 mm ID
Eluent: methanol-water (7:3), with potassium nitrate and sulfuric acid
 $\phi(\text{MeOH}) = 0.7$, $\rho(\text{KNO}_3) = 2 \text{ g/L}$, $\rho(\text{H}_2\text{SO}_4) = 0.05 \text{ g/L}$
Flow rate: 0.9 mL min
Mass of solute injected: 20 ng each

HPLC/ELCD-Application**Substance****Catecholamine**


- 1 L-Noradrenaline
- 2 L-Adrenaline
- 3 L-Dopa
- 4 Dopamine-HCl

Liquid Chromatographic System

Column: LiChrosorb RP-18, 5 μm , 250 mm \times 4.6 mm ID
Eluent: water, with ammonium sulfate and acetic acid (pH 3.2)
 $\rho(\text{NH}_4)_2\text{SO}_4 = 5 \text{ g/L}$, $\rho(\text{HAc}) = 3 \text{ g/L}$
Flow rate: 1 mL/min
Mass of solute injected: 200 μg each

Detector

Working electrode: GC
Voltage: + 800 mV
Range: 0.5 nA

HPLC/ELCD-Application**Substance****Catecholamines from tissue**

- | | |
|-------------------|---------------------|
| 1 L-Noradrenaline | $m = 15 \text{ pg}$ |
| 2 Dopamine | $m = 40 \text{ pg}$ |

Liquid Chromatographic System

Column: Nucleosil 5 SA, $5 \mu\text{m}$, $150 \text{ mm} \times 4.6 \text{ mm ID}$
Eluent: acetate-citrate-buffer, pH 5.2
Flow rate: 1 mL min
Injection volume: $20 \mu\text{L}$

Detector


Working electrode: GC
Voltage: $+800 \text{ mV}$
Range: 0.5 nA

Sample preparation

Isolation of Catecholamines
from tissue by alumina (Al_2O_3)
extraction

Literature

- R. Keller et al.
Life Sci. 19, 995 (1976)
S. Allemark et al.
J. Liq. Chrom. 2, 277 (1979)

HPLC/ELCD-Application

Substance

Catecholamines from urine

- 1 L-Noradrenaline (NA)
- 2 L-Adrenaline (A)
- 3 Dopamine (DA)

Liquid Chromatographic System


Column: Nucleosil 10 SA, 10 μ m, 250 mm \times 4.6 mm ID
 Eluent: water, with sodium acetate and sodium citrate, pH 4.2
 Flow rate: 2 mL/min
 Injection volume: 20 μ L

Detector

Working electrode: GC
 Voltage: + 800 mV
 Range: 50 nA

Sample preparation

Isolation of catecholamines from urines ($V = 5$ mL) by alumina (Al_2O_3) extraction

HPLC/ELCD-Application**Substance**


Esculin

Hydroxycoumarins**Liquid Chromatographic System**

Column: LiChrosorb RP-18, 5 μm , 250 mm \times 4.6 mm ID
Eluent: methanol-water (2:3), with acetic and lithium perchlorate
 $\phi(\text{MeOH}) = 0.4$, $\rho(\text{HAc}) = 5 \text{ g/L}$, $\rho(\text{LiClO}_4) = 5 \text{ g/L}$
Flow rate: 0.5 mL/min
Mass of solute injected: 20 ng

Detector

Working electrode: GC
Voltage: +1000 mV
Range: 50 nA

HPLC/ELCD-Application

Substance

- 1 Catechin
- 2 Rutin



Flavones

Liquid Chromatographic System

Column: LiChrosorb RP-18, 5 μ m, 250 mm \times 4.6 mm ID
 Eluent: methanol-water, with lithium perchlorate and acetic acid
 $\rho(\text{LiClO}_4) = 5 \text{ g L}$, $\rho(\text{HAc}) = 5 \text{ g L}$
 Flow rate: 0.5 mL/min
 Mass of solute injected: 20 ng each

Detector

Working electrode: GC
 Voltage: + 1000 mV
 Range: 50 nA or 10 nA

HPLC/ELCD-Application**Substance**

Quercetin
from hawthorn (*crataegus oxyacanthae*)
tincture

Flavones**Liquid Chromatographic System**


Column: LiChrosorb RP-2, 5 μm , 250 mm \times 4.6 mm ID
Eluent: methanol-water (1:1), with lithium perchlorate and acetic acid
 $\phi(\text{MeOH}) = 0.5$, $\rho(\text{LiClO}_4) = 2 \text{ g/L}$, $\rho(\text{HAc}) = 10 \text{ g/L}$
Flow rate: 0.9 mL min
Injection volume: 20 μL

Detector

Working electrode: GC
Voltage: + 1000 mV
Range: 50 nA

Sample preparation

tincture diluted with eluent
 $\phi(\text{tincture}) = 0.1$

HPLC/ELCD-Application

Substance

Estrogens


- 1 Estriol
- 2 Estradiol
- 3 Estrone

Liquid Chromatographic System

Column: LiChrosorb RP-18, 5 μ m. 250 mm \times 4.6 mm ID
 Eluent: acetone-water (1:1), with lithium perchlorate and acetic acid
 $\phi(\text{acetone}) = 0.5$, $\rho(\text{LiClO}_4) = 5 \text{ g/L}$, $\rho(\text{HAc}) = 5 \text{ g/L}$
 Flow rate: 0.7 mL/min
 Mass of solute injected: 20 ng each

Detector

Working electrode: GC
 Voltage: + 1000 mV
 Range: 10 nA

HPLC/ELCD-Application**Substance**

- 1 δ -Tocopherol
- 2 β, γ -Tocopherol
- 3 α -Tocopherol

Tocopherols from soya oil**Liquid Chromatographic System**


Column: LiChrosorb RP-18, 5 μm , 250 mm \times 4.6 mm ID
Eluent: methanol, with lithium perchlorate and acetic acid
 $\rho(\text{LiClO}_4) = 5 \text{ g/L}$, $\rho(\text{HAc}) = 1 \text{ g/L}$
Flow rate: 1.5 mL/min
Injection volume: 20 μL

Detector

Working electrode: GC
Voltage: + 800 mV
Range: 100 nA

Sample preparation

soya oil $m = 182 \text{ mg}$
extracted with eluent $V = 20 \text{ mL}$

HPLC/ELCD-Application**Substance**


- 1 Propyl gallate
- 2 Nordihydroguaiaretic acid
- 3 3-tert-Butyl-4-hydroxyanisole

Antioxidants**Liquid Chromatographic System**

Column: LiChrosorb RP-18, 5 μ m, 250 mm \times 4.6 mm ID
Eluent: methanol-water (7:3), with lithium perchlorate and acetic acid
 ϕ (MeOH) = 0.7, ρ (LiClO₄) = 5 g/L, ρ (HAc) = 5 g/L
Flow rate: 0.6 mL/min
Mass of solute injected: 20 ng each

Detector

Working electrode: GC
Voltage: +800 mV
Range: 50 nA

HPLC/ELCD-Application**Substance**


tert-Butylhydroxytoluene

Antioxidants**Liquid Chromatographic System**

Column: LiChrosorb RP-18, 5 μm , 250 mm \times 4.6 mm ID
Eluent: methanol-water (9:1), with lithium perchlorate and acetic acid
 $\phi(\text{MeOH}) = 0.9$, $\rho(\text{LiClO}_4) = 5 \text{ g/L}$, $\rho(\text{HAc}) = 5 \text{ g/L}$
Flow rate: 1 mL/min
Mass of solute injected: 20 ng

Detector

Working electrode: GC
Voltage: + 1200 mV
Range: 50 nA

HPLC/ELCD-Application

Substance

- 1 p-Phenylenediamine
- 2 3-Chloroaniline



Anilines

Liquid Chromatographic System

Column: Nucleosil 5 C₁₈, 5 μm, 250 mm × 4.6 mm ID
 Eluent: methanol-water (1:1), with potassium nitrate and sulfuric acid
 $\phi(\text{MeOH}) = 0.5$, $\rho(\text{KNO}_3) = 2 \text{ g/L}$, $\rho(\text{H}_2\text{SO}_4) = 0.05 \text{ g/L}$
 Flow rate: 0.7 mL/min
 Mass of solute injected: 1: 20 ng, 2: 25 ng

Detector

Working electrode: GC
 Voltage: + 1000 mV
 Range: 100 nA

HPLC/ELCD-Application**Substance****Anilines
Sulfonamides**



- 1 o-Toluidine
- 2 Sulfanilamide
- 3 Sulfamethoxydiazine

Liquid Chromatographic System

Column: Nucleosil 5 C₁₈, 5 μ m, 250 mm \times 4.6 mm ID
Eluent: methanol-water (1:1), with potassium nitrate and sulfuric acid
 $\phi(\text{MeOH}) = 0.5$, $\rho(\text{KNO}_3) = 2 \text{ g/L}$, $\rho(\text{H}_2\text{SO}_4) = 0.05 \text{ g/L}$
Flow rate: 0.7 mL/min
Mass of solute injected: 20 ng each

Detector

Working electrode: GC
Voltage: +1000 mV / 1200 mV
Range: 100 nA

HPLC/ELCD-Application

Substance

*Indoles
Aromatic alcohols*



- 1 3,4-Dihydroxyphenylacetic acid
- 2 D,L-Tryptophan
- 3 Homovanillic acid

Liquid Chromatographic System

Column: Nucleosil 5 C₁₈, 5 μm, 250 mm × 4.6 mm ID
 Eluent: methanol-water (1:1), with potassium nitrate and sulfuric acid
 $\phi(\text{MeOH}) = 0.3$, $\rho(\text{CH}_3\text{COONH}_4) = 6 \text{ g/L}$, $\rho(\text{H}_2\text{SO}_4) = 1.5 \text{ g/L}$
 Flow rate: 0.7 mL/min
 Mass of solute injected: 40 ng each

Detector

Working electrode: GC
 Voltage: + 1000 mV / 800 mV
 Range: 100 nA

HPLC/ELCD-Application**Substance**


Chlorpromazine-HCl

Phenothiazines**Liquid Chromatographic System**

Column: LiChrosorb RP-2, 5 μ m, 250 mm \times 4.6 mm ID
Eluent: methanol-water (65:35), with lithium perchlorate and acetic acid
 $\phi(\text{MeOH}) = 0.65$, $\rho(\text{LiClO}_4) = 2 \text{ g/L}$, $\rho(\text{HAc}) = 1 \text{ g/L}$
Flow rate: 1 mL/min
Mass of solute injected: 20 ng

Detector

Working electrode: GC
Voltage: + 1000 mV
Range: 10 nA

HPLC/ELCD-Application

Substance

- 1 L-Cysteine
- 2 L-Glutathione
- 3 Penicillamine

Mercaptans

Liquid Chromatographic System


Column: Nucleosil 5 SA, 5 μm , 150 mm \times 4.6 mm ID
Eluent: water, with ammonium citrate and phosphoric acid
 $\rho(\text{C}_6\text{H}_6\text{O}_7(\text{NH}_4)_2 \cdot \text{H}_2\text{O}) = 4.5 \text{ g/L}$, $\rho(\text{H}_3\text{PO}_4) = 6 \text{ g/L}$
Flow rate: 1 mL/min
Mass of solute injected: 20 ng each

Detector

Working electrode: Au
Voltage: + 800 mV
Range: 50 nA

Literature

D.L. Rabenstein and R. Saetre
Anal. Chem. 49, 1036 (1977)

HPLC/ELCD-Application**Substance**Penicillamine from urine,
 $\rho = 20 \mu\text{g/mL}$ **Mercaptans****Liquid Chromatographic System**

Column: Nucleosil 5 SA, $5 \mu\text{m}$, $150 \text{ mm} \times 4.6 \text{ mm ID}$
Eluent: water, with ammonium citrate and phosphoric acid, pH 2.2
 $\rho(\text{C}_6\text{H}_6\text{O}_7(\text{NH}_4)_2 \cdot \text{H}_2\text{O}) = 4.5 \text{ g/L}$, $\rho(\text{H}_3\text{PO}_4) = 6 \text{ g/L}$
Flow rate: 1 mL/min
Injection volume: $20 \mu\text{L}$

Detector

Working electrode: Au
Voltage: $+800 \text{ mV}$
Range: 100 nA


Sample

urine diluted with eluent
 $\phi(\text{urine}) = 0.1$

Literature

D.L. Rabenstein and R. Saetre
Anal. Chem. 49, 1036 (1977)

F. Kreuzig
private communication

HPLC/ELCD-Application

Substance

- 1 Ascorbic acid
- 2 Uric acid from urine or serum

Vitamins
Purines

Liquid Chromatographic System

Column: LiChrosorb RP-18, 5 μm , 250 mm \times 4.6 mm ID
Eluent: water, with meta-phosphoric acid
 $\rho(\text{HPO}_3) = 8 \text{ g L}^{-1}$
Flow rate: 1 mL/min
Injection volume: 20 μL

Detector


Working electrode: GC
Voltage: + 800 mV
Range: 100 nA

Sample preparation

urine: diluted with eluent
 $\phi(\text{urine}) = 0.0025$
serum: deproteinised serum
diluted with eluent
 $\phi(\text{serum}) = 0.1$

Literature

E.S. Wagner et al.
J. Chrom. 163, 225 (1979)

HPLC/ELCD-Application**Substance**

Ascorbic acid
from fruit drinks, fruit juices,
pharmaceutical preparations

Vitamins**Liquid Chromatographic System**

Column: LiChrosorb RP-18, 5 μm , 250 mm \times 4.6 mm ID

Eluent: water, with meta-phosphoric acid
 $\rho(\text{HPO}_3) = 8 \text{ g/L}$

Flow rate: 1 mL/min

Injection volume: 20 μL

Detector

Working electrode: GC


Voltage: + 800 mV

Range: 100 nA

Literature

E.S. Wagner et al.
J. Chrom. 163, 225 (1979)

L.A. Pachla et al.
Methods in Enzymology 62, 15 (1979)

HPLC/ELCD-Application**Substance**


- 1 Vitamin-A-acetate
- 2 α -Tocopherol
- 3 β -Carotene

Vitamins**Liquid Chromatographic System**

Column: Nucleosil 5 C₁₈, 5 μ m, 250 mm \times 4.6 mm ID
Eluent: methanol, with lithium perchlorate and acetic acid
 $\rho(\text{LiClO}_4) = 2 \text{ g/L}$, $\rho(\text{HAc}) = 1 \text{ g/L}$
Flow rate: 1.5 mL min
Mass of solute injected: 20 ng each

Detector

Working electrode: GC
Voltage: + 1000 mV
Range: 100 nA

HPLC/ELCD-Application

Substance

Carotene
Tocopherol

- 1 α -Tocopherol
- 2 β -Carotene
in carrot juice

Liquid Chromatographic System

Column: Nucleosil 5 C₁₈, 5 μ m, 250 mm \times 4.6 mm ID
 Eluent: methanol, with lithium perchlorate and acetic acid
 $\rho(\text{LiClO}_4) = 2 \text{ g L}$, $\rho(\text{HAc}) = 1 \text{ g L}$
 Flow rate: 1.5 mL/min
 Mass of solute injected: standard: 1: 4 ng, 2: 8 ng

Sample preparation extraction with methanol

Detector

Working electrode: GC
 Voltage: + 800 mV
 Range: 10 nA