

## Application Bulletin 191/3 e

# Determination of cysteine and cystine simultaneously by polarography

### Summary

After the degradation of biological samples (e.g. milk, wool, etc.), it is often important to know the cystine/cysteine ratio. This Application Bulletin describes the simultaneous, polarographic determination of the two amino acids. The determination is performed in perchloric acid solution at the DME. Samples with a high protein content require that the determination is performed in an alkaline solution.

### Instruments

VA instrument  
capable of operating a Multi-Mode Electrode and supporting differential pulse (DP) measuring mode

### Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Mercury drop capillary	6.1226.030
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel	6.1245.010
	Filled with $c(\text{CH}_3\text{COOLi}) = 1 \text{ mol/L}$ (method 1) or $c(\text{KCl}) = 3 \text{ mol/L}$ (method 2)	
AE	Pt rod electrode	6.0343.x00

### Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis\*).

- Cystine, for analysis, CAS 56-89-3
- Cysteine, for analysis, CAS 52-90-4
- Ultrapure water, resistivity  $>18 \text{ M}\Omega\cdot\text{cm}$  (25 °C), type I grade (ASTM D1193)

### For method 1

- Perchloric acid,  $w(\text{HClO}_4) = 70\%$ , for analysis, CAS 7601-90-3
- Lithium acetate dihydrate,  $w(\text{LiOOCCH}_3\cdot 2 \text{ H}_2\text{O}) \geq 99.0\%$ , for analysis, CAS 6108-17-4

### For method 2

- Ammonium hydroxide solution,  $w(\text{NH}_3) = 25\%$ , for analysis, CAS 1336-21-6
- Hydrochloric acid,  $w(\text{HCl}) = 30\%$ , for analysis, CAS 7647-01-0
- Sodium hydroxide solution,  $c(\text{NaOH}) = 2.0 \text{ mol/L}$ , CAS 1310-73-2

## Method 1: Samples containing few proteins

### Summary

Samples without proteins can be determined in perchloric acid solution. The determination limit is heavily dependent on the matrix and lies at approx. 1 mg/L for both substances.

### Solutions

Electrolyte	$c(\text{HClO}_4) = 0.1 \text{ mol/L}$ 4.16 mL perchloric acid are diluted to 500 mL with ultrapure water.
Li-acetate solution	$c(\text{CH}_3\text{COOLi}) = 1 \text{ mol/L}$ for reference electrode

### Standard solutions

Cysteine standard solution	$c(\text{cysteine}) = 1 \text{ g/L}$ in perchloric acid Mix 0.1 g cysteine to a slurry in 20 mL dist. water and dissolve by adding 0.86 mL perchloric acid, then fill up to 100 mL with dist. water.
Cystine standard solution	$c(\text{cystine}) = 1 \text{ g/L}$ in perchloric acid Mix 0.1 g cystine to a slurry in 20 mL dist. water and dissolve by adding 0.86 mL perchloric acid, then fill up to 100 mL with dist. water.

## Analysis

### Measuring solution

10 mL diluted sample





10 mL supporting electrolyte

The concentrations are determined by standard addition.

### Parameters

Voltammetric	
Electrode operating mode	DME
Measuring mode	DP – Differential pulse
Stirring rate	2000 min <sup>-1</sup>
Equilibration time	30 s
Sweep	
Start potential	0.2 V
End potential	-0.8 V
Potential step	0.006 V
Potential step time	0.6 s
Sweep rate	0.01 V/s
Pulse amplitude	0.05 V
Substance	
Name	Cysteine
Characteristic potential	-0.08 V
Name	Cystine
Characteristic potential	-0.49 V

### Example



## Result

Sample size	10.0 mL
$\beta$ (cysteine)	2.9 mg/L
$\beta$ (cystine)	5.8 mg/L

## Comments

- The peak potentials of both substances shift to slightly more negative values with increasing concentrations. The value must perhaps be corrected.
- If the concentrations of cystine and cysteine differ strongly, it is possibly necessary to perform two analyses in various dilutions for the respective substance.

## Method 2: Samples with a high protein content

### Summary

Samples with a high protein content cannot be determined using perchloric acid as electrolyte, since proteins are precipitated in acidic solutions. Therefore the determination is carried out in ammonium buffer at pH 9.6.

The limit of detection for cysteine is approx. 0.05 mg/L, for cystine approx. 1 mg/L. The determination runs linear for cysteine up to 180 mg/L, for cystine up to 300 mg/L.

### Solutions

Ammonia buffer pH 9.6	$c(\text{NH}_4\text{Cl}) = 1 \text{ mol/L}$ $c(\text{NH}_3) = 2 \text{ mol/L (pH 9.6)}$ 112.5 mL $\text{NH}_3$ + 53 mL HCl, filled up to 500 mL with ultrapure water.
--------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Standard solutions

Cysteine standard solution	$c(\text{cysteine}) = 1 \text{ g/L}$ Dissolve 0.100 g cysteine in 20 mL ammonium buffer and fill up to 100 mL with dist. water.
Cystine standard solution	$c(\text{cystine}) = 1 \text{ g/L}$ Dissolve 0.100 g cystine in 5 mL sodium hydroxide and fill up to 100 mL with dist. water.

## Analysis

### Measuring solution

10 mL diluted sample

1 mL ammonia buffer pH 9.6


The concentrations are determined by standard addition.

The peak potential of cystine is strongly dependent on the pH value and shifts in negative direction with increasing pH value.

## Parameters

Voltammetric	
Electrode operating mode	DME
Measuring mode	DP – Differential pulse
Stirring rate	2000 $\text{min}^{-1}$
Equilibration time	30 s
Sweep	
Start potential	-0.25 V
End potential	-1.75 V
Potential step	0.006 V
Potential step time	0.8 s
Sweep rate	0.0075 V/s
Pulse amplitude	0.05 V
Substance	
Name	Cysteine
Characteristic potential	-1.0 V
Name	Cystine
Characteristic potential	-0.5 V


## Example



**Cysteine**


c = 0.031 g/l

+/- 0.001 g/l (1.92%)


**Cystine**

c = 0.456 g/l

+/- 0.006 g/l (1.22%)


**Result**

Sample	Sodium caseinate
Sample size	0.1 mL
$\beta$ (Cysteine)	0.031 g/L
$\beta$ (Cystine)	0.456 g/L

**Comments**

- If the pH value is over 10, the cystine peak becomes very flat and insensitive.
- If the concentrations of cystine and cysteine differ greatly, it will possibly be necessary to perform two analyses in various dilutions for the respective substance.

## Appendix

### Report for the example determination of cysteine and cystine according to method 1

```

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Determ.      : 06111056          User:  mj          Date: 1999-06-11
Modified     : 1999-06-11 11:12:06  Run :  18          Time: 10:56:21
Sample table: -
    
```

Pos.	Ident.1/S1	Ident.2/S2	Ident.3/S3	Method.call	Sample size/S0
	cysHClO4				10.0 mL

```

Method      : AB191/1
Title       : Determination of Cystine and Cysteine
Remark1     : 10 mL supporting electrolyte + 10 mL sample
Remark2     : supporting electrolyte: HClO4 c = 0.1 mol/L
    
```

Substance	Mass conc.:	MC.dev.	Cal.dev.	Mass	Add.mass	V0.sample	Comments
Cysteine	2.873 mg/L	0.121 mg/L (4.22%)	-	28.73 µg	25 µg	10 mL	-----

VR	U/mV	I/nA	I.mean	Std.dev.	I.delta	Comments
00	-80	-92.07	-92.09	0.0219		-----
01	-79	-92.11				
10	-88	-175.9	-176.9	1.432	-84.84	
11	-87	-177.9				
20	-92	-254.5	-252.3	3.016	-75.41	
21	-91	-250.2				

Substance	Mass conc.:	MC.dev.	Cal.dev.	Mass	Add.mass	V0.sample	Comments
Cystine	5.798 mg/L	0.090 mg/L (1.55%)	-	57.98 µg	25 µg	10 mL	-----

VR	U/mV	I/nA	I.mean	Std.dev.	I.delta	Comments
00	-483	-48.28	-48.25	0.0341		-----
01	-482	-48.23				
10	-487	-68.13	-68.30	0.2406	-20.05	
11	-487	-68.47				
20	-489	-89.29	-88.88	0.5862	-20.58	
21	-492	-88.46				

Substance	Techn.	Y.reg/offset	Slope	Nonlin.	Mean deviat.
Cysteine	std.add.	-9.352e-08	-3.284e-05		3.016e-09
Cystine	std.add.	-4.813e-08	-8.373e-06		3.830e-10

Final results	+/-	Res.dev.	%	Comments
Cysteine =	2.8726 mg/L	0.121	4.22	
Cystine =	5.7979 mg/L	0.090	1.55	

### Method print for the determination cysteine and cystine according to method 1

```

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB191_1 .mth          OPERATION SEQUENCE
Title : Det of cystine and cystine (low protein content)
    
```

	Instructions	t/s	Main parameters	Auxiliary parameters
1	SMPL>M		V.fraction mL	V.total L
2	DOS>M		Soln.name	V.add 10 mL
3	PURGE			
4	STIR	300.0	Rot.speed 2000 /min	
5	(ADD			
6	PURGE			
7	STIR	30.0	Rot.speed 2000 /min	
8	0PURGE			
9	0STIR	20.0		
10	(REP			
11	SEGMENT		Segm.name POL	
12	REP)1			
13	PURGE			
14	ADD>M		Soln.name Cystine	V.add 0.025 mL
15	ADD>M		Soln.name Cysteine	V.add 0.025 mL
16	ADD)2			
17	END			

Method: AB191\_1 SEGMENT  
POL

	Instructions	t/s	Main parameters	Auxiliary parameters
1	DME			
2	DPMODE		U.ampl -50 mV	t.meas 20.0 ms
			t.step 0.60 s	t.pulse 40.0 ms
3	SWEEP	101.4	U.start 200 mV	U.step 6 mV
			U.end -800 mV	Sweep rate 10 mV/s
4	OMEAS		U.standby mV	
5	END			

**Report for the example determination of cysteine and cystine according to method 2**

```

===== METROHM 797 VA COMPUTRACE (Version 1.0.0.1) (Serial No. 0) =====
Determination : 06150859_Na-caseinat.dth
Sample ID      : Na-caseinat
Creator method :                               Date :                               Time:
Creator determ.:                               Date : 1999-06-15           Time: 08:59:33
Modified by    :                               Date : 2017-07-20           Time: 16:46:51
  
```

```

Method      : AB191_Det of Cystin and Cystein in alkaline solution.mth
Title       : Bestimmung von Cystin und Cystein
Remark1     : background electrolyte: NH3/NH4Cl pH 9.6
Remark2     :
  
```

```

Sample amount : 0.100 mL
Cell volume   : 10.100 mL
  
```

```

Substance    : Cystine
Conc.        : 4.513 mg/L
Conc.dev.    : 0.055 mg/L ( 1.22%)
Amount       : 45.583 ug
Add.amount   : 50.000 ug
  
```

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-1.018	-45.3	-45.5	0.251	0.0	
1 - 2	-1.006	-45.7				
2 - 1	-1.006	-95.2	-95.3	0.166	-49.8	
2 - 2	-1.006	-95.4				
3 - 1	-1.000	-141.9	-143.0	1.527	-47.7	
3 - 2	-1.000	-144.1				

```

Substance    : Cysteine
Conc.        : 305.007 ug/L
Conc.dev.    : 5.841 ug/L ( 1.92%)
Amount       : 3.081 ug
Add.amount   : 2.500 ug
  
```

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.482	-49.7	-49.5	0.249	0.0	
1 - 2	-0.482	-49.4				
2 - 1	-0.488	-90.1	-90.3	0.346	-40.8	
2 - 2	-0.488	-90.6				
3 - 1	-0.488	-128.1	-127.5	0.746	-37.2	
3 - 2	-0.488	-127.0				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Cystine	std.add.	-4.553e-008	-1.009e-005	7.436e-010	0.99983
Cysteine	std.add.	-4.968e-008	-1.629e-004	1.092e-009	0.99968

Final results		+/-	Res. dev.	%	Comments
Cystine: default	=	0.456	g/l	0.006	1.215
Cysteine: default	=	0.031	g/l	0.001	1.915

**Method print for the determination cysteine and cystine according to method 2**

Method parameters

Method : AB191\_2 Det of Cystine and Cysteine (high protein content).mth  
 Title : Determination of Cystine and Cysteine (high protein content)  
 Remark1 : 10 mL sample + 1 mL ammonia buffer (pH 9.6)  
 Remark2 :

Calibration : Standard addition  
 Technique : Batch  
 Addition : Manual

Sample ID : Na-caseinat  
 Sample amount (mL): 10.000  
 Cell volume (mL): 11.000

## Voltammetric parameters

-----  
 Mode : DP - Differential Pulse  
 Highest current range : 1 mA  
 Lowest current range : 100 nA  
 Electrode : DME  
 Stirrer speed (rpm) : 2000  
 Initial electr. conditioning : No  
 No. of additions : 2  
 No. of replications : 2  
 Measure blank : No  
 Addition purge time (s) : 30  
 Initial purge time (s) : 300  
 Sweep  
 Equilibration time (s) : 5.000  
 Start potential (V) : -0.250  
 End potential (V) : -1.750  
 Voltage step (V) : 0.006  
 Voltage step time (s) : 0.800  
 Sweep rate (V/s) : 0.007  
 Pulse amplitude (V) : 0.050  
 Pulse time (s) : 0.040  
 Cell off after measurement : Yes

## Peak evaluation

-----  
 Regression technique : Linear Regression  
 Peak evaluation : Height  
 Minimum peak width (V.steps) : 5  
 Minimum peak height (A) : 1.000e-010  
 Reverse peaks : No  
 Smooth factor : 4  
 Eliminate spikes : Yes

## Substances

-----  
 Cystin : -1.000 V +/- 0.050 V  
 Standard solution : 1 1.000 g/L  
 Addition volume (mL) : 0.050  
 Cystine : Final result (Cystin) =  
           Conc \* (11 / 10) \* (1000 / 1) + 0 - 0  
 Cystei : -0.480 V +/- 0.050 V  
 Standard solution : 2 0.100 g/L  
 Addition volume (mL) : 0.025  
 Cysteine : Final result (Cystei) =  
           Conc \* (11 / 10) \* (1000 / 1) + 0 - 0

## Baseline

Substance Addition		automatic	start (V)	end (V)	type	scope
Cystin	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak
Cystei	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak