


Determination of folic acid by polarography

Summary

This Application Bulletin describes the polarographic determination of folic acid, a vitamin of the B series, also known as vitamin B₉ or vitamin Bc. Instructions for the determination in solutions (e.g. fruit juice), vitamin capsules and multivitamin tablets are given. The linear range of the determination is also specified. The limit of detection is approx. 75 µg/L folic acid.

Folic acid is reduced in alkaline solution (pH~ 9) at the DME to tetrahydrofolic acid following the below mechanism:



The reaction has only limited reversibility.

Instruments

VA instrument
capable of operating a Multi-Mode
Electrode and supporting differential
pulse (DP) measuring mode

Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Mercury drop capillary	6.1226.030
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel Filled with c(KCl) = 3 mol/L	6.1245.010
AE	Pt rod electrode	6.0343.x00

Reagents

All of the used reagents must be of purest quality possible (for analysis).

- Boric acid, for analysis, CAS 10043-35-3
- Sodium hydroxide, for analysis, CAS 1310-73-2
- Folic acid, for analysis, CAS 59-30-3
- Ultrapure water, resistivity >18 MΩ·cm (25 °C), type I grade (ASTM D1193)

Solutions

Electrolyte	c(boric acid) = 0.1 mol/L c(NaOH) = 0.05 mol/L 100 mL ultrapure water are added to 6.2 g boric acid in a beaker. After addition of 2 g NaOH, the mixture is stirred until everything is dissolved. After dilution with dist. water to 950 mL, the pH is adjusted to 11.1 - 11.2 with c(NaOH) = 2 mol/L. The solution is then made up to 1000 mL with ultrapure water. Use only CO ₂ free ultrapure water!
2 mol/L NaOH solution	c(NaOH) = 2 mol/L 80 g/L NaOH in CO ₂ free ultrapure water.
0.1 mol/L NaOH solution	c(NaOH) = 0.1 mol/L 4 g/L NaOH in CO ₂ free ultrapure water.

Standard solutions

Folic acid standard stock solution	$\beta(\text{folic acid}) = 0.5 \text{ g/L}$ The stock solution of the standard is prepared with a concentration of 0.5 g/L. The actual content of the folic acid used must be taken into account here. Example: 127.5 mg 98% folic acid = 125 mg 100% folic acid are weighed into a beaker and 80 mL ultrapure water added. $c(\text{NaOH}) = 0.1 \text{ mol/L}$ is added with stirring until the solution becomes clear and a pH value of 8.0 has been reached. An additional 120 mL ultrapure water are added and the pH value readjusted to 8.0. The solution is added to a 250 mL volumetric flask, which is then filled to the mark with ultrapure water and the contents mixed. Store solution in a brown glass bottle in a refrigerator. It is best to prepare fresh solutions daily and used CO_2 free ultrapure water.
Folic acid standard solution	This is prepared when needed from the stock solution by dilution with ultrapure water pH = 8.

If need be, a pH value of 8.0 is set by addition of diluted NaOH.

Sample preparation

Injection solutions

These are adjusted to pH = 8 with NaOH and can then be used directly.

Vitamin tablets (monovitamin tablets)

10 tablets are weighed out to obtain the average weight and then ground to a powder (grinder, mortar). 200 mg of the resulting powder are weighed into a beaker and 30 ml water added. The pH value is adjusted to 8.0 by addition of $c(\text{NaOH}) = 0.1 \text{ mol/L}$ with stirring. The solution is stirred for a further 15 min and the pH value again adjusted to 8.0 if need be. After filtration of the solution through a paper filter into a dark glass bottle, the filter is washed four times with 5 mL aliquots of dist. water pH = 8.0. The filtrate is transferred to a 100 mL volumetric flask, which is filled to the mark with

ultrapure water (pH = 8.0) and the contents mixed. Use solution for the analysis immediately.

Analysis

0.50 mL of the solution obtained from the tablets are mixed in a polarographic vessel with 19.5 mL electrolyte.

Measuring solution

0.5 mL sample solution


19.5 mL electrolyte

The concentration is determined by standard addition.

Parameters

Voltammetric	
Electrode operating mode	DME
Measuring mode	DP – Differential pulse
Stirring rate	2000 min ⁻¹
Equilibration time	5 s
Sweep	
Start potential	-0.7 V
End potential	-1.2 V
Potential step	0.004 V
Potential step time	0.8 s
Sweep rate	0.005 V/s
Pulse amplitude	0.05 V
Substance	
Name	Folic acid
Characteristic potential	-0.97 V

Example



Result

Sample	Folic acid vitamin tablet
$\beta(\text{Folic acid})$	4.7 mg/tablet

Comments

- The linearity range lies between 75 µg/L and 8.75 mg/L folic acid.
- Ensure that the concentration of the folic acid in the polarographic vessel including the standard additions does not exceed the linear range of 8.75 mg/L.
- The presence of ascorbic acid, ascorbates and iron(II) compounds (e.g. iron(II) fumarate) does not disturb the determination of folic acid.
- Folic acid has to be extracted from foodstuffs and fodder.

References

- Ben-Bassat, A. H. I., Frydam-Kupfer, G. et al., Polarographic rapid microdetermination of folic acid with hanging drop electrode and rapid sweep, *Polarography* 1964, (1966) 993
- Chen, H. Y., Han, J. L. et al., Study of alternating current adsorptive stripping voltammetry in flowing systems for determination of ultratrace amounts of folic acid, *Fresenius, Z. Anal. Chem.* 334/7, (1989) 621
- El Maali, N. A., Vire, J. C. et al., Cathodic adsorptive stripping square wave voltammetry of folic acid (Vitamin B-9), *Analisis* 17, (1989) 213
- Göbbler, K. H., Breinlich, J., Quantitative wechselstrompolarographische Simultanbestimmung von Vitaminen der B-Reihe, *Pharm. Ztg.* 48, (1972) 1859
- Jacobsen, E., Björnsen, M. W., Polarographic determination of folic acid in pharmaceutical preparations, *Anal. Chim. Acta* 96, (1978) 345
- Kretzschmar, K., Jaenicke, W., The redox mechanism of the system folic acid - dihydrofolic acid – tetrahydrofolic acid(II), *Z. Naturforsch.* 26/10, (1971) 999
- Kruze, J. E., Polarographic determination of ascorbic acid, folic acid and riboflavin. *Farmatsiya* 18/4, (1969) 59
- Luo, D. B., Determination of folic acid by adsorptive stripping voltammetry at the static mercury drop electrode, *Anal. Chim. Acta* 189, (1986) 277
- Rozanski, L. Polarographic determination of folic acid in tablets containing iron(II) sulphate, *Analyst* 103, (1978) 950
- Solunina, L. A., Dewjatinin, V.A. et al., Polarographic determination of folic acid in pharmaceuticals, *Farmatsiya* 17/2, (1968) 45

Appendix

Report for the example determination of folic acid in a vitamin tablet

METROHM 646 VA-PROCESSOR (5.846.6041)

Detn. of Folic Acid (Vit. Bc) in Monovitamin Tab. METHOD

MPL 1 EL.TYPE MME

SUPP.ELEC Borate buff.11.2pH

V.MEAS 20.000 mL

ALIQOT 1.000

REMARK Detn.folic acid in monovitamin tablet (5mg/tab)

Ag/AgCl (3M KCl) reference electrode

NAME Prof.J.G.Dick

RUN# 1

ANALYTE	L R S	U.SUBST	EV.VALUE	DELTA	m.ANALYTE
Niamde	A0 0 0	-968 mV	36.79 nA		
	A0 1 0	-970 mV	36.97 nA		
	A1 0 0	-970 mV	59.18 nA		
	A1 1 0	-971 mV	58.54 nA	21.97 nA	
	A2 0 0	-971 mV	79.78 nA		
	A2 1 0	-972 mV	79.04 nA	20.54 nA	
m.STD	25.00 ug	SLOPE	1.175 mg/uA		43.64 ug

rho(fol.ac = 4.735 mg/g

SMPL.V,m 1.00000 mg

IDENT APO folic ac.5mg/tab

DATE 91-05-25 TIME 14:31

Method print for the determination of folic acid

Detn. of Folic Acid (Vit. Bc) in Monovitamin Tab.

METHOD 4 PAGE 3

MPL 1 EL.TYPE MME

OPERATION SEQUENCE

OPERATIONS/PARAMETERS

OPERATIONS/PARAMETERS

```

1 PURGE ;STIR ; 5 s
2 [ADDL ;OPURGE;OSTIR ; 5 s
3 (REP ;
4 DME ;MEAS ; 5 s
4a M.MODE DPN -50 mV
4b T.STEP 800 ms
4c U.SET -700 mV
5 SWP 0 ; 100 s
5a U.END -1.200 V
5b U.STEP 4 mV
    SW.RATE 5.0 mV/ s
6 OMEAS ;
7 REP) 1;
8 OMEAS ;PURGE ;STIR ;
9 BEEP ;ADD1]2; 30 s
10 OMEAS ;OPURGE;OSTIR ;
11 BEEP ;END ;

```