

Application Bulletin 231/3 e

Determination of zinc, cadmium, lead, copper, thallium, nickel and cobalt in water samples by anodic and adsorptive stripping voltammetry according to DIN 38406-16

Summary

The standard method in DIN 38406-16, describes the determination of Zn, Cd, Pb, Cu, Tl, Ni and Co in drinking, ground, surface and precipitating (e.g. rain) waters. Because the presence of organic substances in the water samples can strongly interfere with the voltammetric determination, a pretreatment with UV digestion using hydrogen peroxide is necessary. This digestion ensures the elimination of all organic substances without introduction of blank values.

These methods can, of course, also be applied for trace analysis in other materials, e. g. trace analysis in the production of semiconductor chips on silicon basis.

Zn, Cd, Pb, Cu and Tl are determined on the HMDE by means of anodic stripping voltammetry (ASV), Ni and Co by means of adsorptive stripping voltammetry (AdSV).

These metals can be determined in the following ranges:

Element	Determination limit	Max. concentration
Zn	1.0 µg/L	50 mg/L
Cd	0.1 µg/L	50 mg/L
Pb	0.1 µg/L	50 mg/L
Cu	1.0 µg/L	50 mg/L
Tl	0.1 µg/L	50 mg/L
Ni	0.1 µg/L	10 mg/L
Co	0.1 µg/L	10 mg/L

Instruments

VA instrument capable of operating a Multi-Mode Electrode and supporting differential pulse (DP) measuring mode	
909 UV Digester	2.909.0014

Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Mercury drop capillary	6.1226.030
	or	6.1226.050
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel Filled with c(KCl) = 3 mol/L	6.1245.010
AE	Pt rod electrode	6.0343.x00

Sample preparation

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Nitric acid, w(HNO₃) = 65%, for trace analysis*, CAS 7697-37-2

or

- Hydrochloric acid, w(HCl) = 30%, for trace analysis*, CAS 7647-01-0

If UV digestion is required:

- Hydrogen peroxide, w(H₂O₂) = 30%, for trace analysis*, CAS 7722-84-1

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Sample preparation

- Immediately upon taking the water sample, filter it through a microfilter (0.45 µm). The filtrates and filters with residues are treated separately. Acidify the filtrates, per liter, with 1 mL conc. HNO₃ or HCl. The pH value of the acidified filtrates should lie between 1.7 and 2.0. If this is not the case, more acid must be added.

- Tap water, surface waters, mineral waters and drinking waters can usually be analyzed without pretreatment.

Organic matter often interferes with voltammetric determinations and therefore sample solutions usually have to be digested.

- Low polluted waste waters can be digested with the 909 UV Digester.
 - Add 50 - 100 μL hydrogen peroxide $w(\text{H}_2\text{O}_2) = 30\%$ to 10 mL acidified samples ($\text{pH} = 2$). The quartz tubes are irradiated for 90 min at 90 °C. After cooling to room temperature, the digested samples can be transferred directly to the polarographic vessel. The blank value of this digestion is relatively small.
- Filters, filter residues and samples with organic matter (foods, pharmaceuticals etc.) must be digested.
 - High-pressure asher
 - Microwave digestion

Both techniques oxidize the samples in a closed digestion vessel by means of a mixture of concentrated mineral acids.

 - According to Application Bulletin 113, open wet digestion with H_2SO_4 and H_2O_2 .
- Preparation of sample for the analysis of semiconductor chips: Place a drop of hydrofluoric acid $w(\text{HF}) = 40\%$ on chip and fume off by heating. A white, powdery residue forms. Dissolve this with a drop of nitric acid $w(\text{HNO}_3) = 65\%$ and rinse into a polarographic vessel. Prepare a blank value with the same reagents

Comments

- The analysis of the digested filters and that of the filtrates is identical.
- The blank values of the digestions and the chemicals must be taken into consideration.
- If samples contain higher contents of humic acids, these must be digested immediately to avoid formation of precipitation.
- If Fe must also be determined, e.g. for trace analysis for the production of semiconductor chips, please refer to Application Bulletin 74.

Method 1: Determination of Zn, Cd, Pb and Cu

Zinc, cadmium, lead and copper are determined in KCl sodium acetate solution at the HMDE by anodic stripping voltammetry (ASV).

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Sodium hydroxide solution, for trace analysis*, $w(\text{NaOH}) = 30\%$, CAS 1310-73-2
- Acetic acid, for trace analysis*, $w(\text{CH}_3\text{COOH}) = 100\%$, CAS 64-19-7
- Potassium chloride, KCl, for trace analysis*, CAS 7447-40-7
- Zn standard stock solution, $\beta(\text{Zn}^{2+}) = 1 \text{ g/L}$ (commercially available)
- Cd standard stock solution, $\beta(\text{Cd}^{2+}) = 1 \text{ g/L}$ (commercially available)
- Pb standard stock solution, $\beta(\text{Pb}^{2+}) = 1 \text{ g/L}$ (commercially available)
- Cu standard stock solution, $\beta(\text{Cu}^{2+}) = 1 \text{ g/L}$ (commercially available)
- Ultrapure water, resistivity $>18 \text{ M}\Omega \cdot \text{cm}$ (25 °C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

KCl sodium acetate solution	$c(\text{KCl}) = 1.5 \text{ mol/L}$ $c(\text{CH}_3\text{COONa}) = 0.5 \text{ mol/L}$ 55.9 g KCl + 25 mL NaOH + 14.2 mL CH_3COOH are filled up to 500 mL with ultrapure water.
-----------------------------	---

Standard solutions

Zn standard solution	$\beta(\text{Zn}^{2+}) = 10 \text{ mg/L}$ The solution is prepared using $c(\text{HNO}_3) = 0.014 \text{ mol/L}$.
Cd standard solution	$\beta(\text{Cd}^{2+}) = 0.1 \text{ mg/L}$ The solution is prepared using $c(\text{HNO}_3) = 0.014 \text{ mol/L}$.
Pb standard solution	$\beta(\text{Pb}^{2+}) = 0.5 \text{ mg/L}$

	The solution is prepared using $c(\text{HNO}_3) = 0.014 \text{ mol/L}$.
Cu standard solution	$\beta(\text{Cu}^{2+}) = 2.5 \text{ mg/L}$ The solution is prepared using $c(\text{HNO}_3) = 0.014 \text{ mol/L}$.

Analysis

Measuring solution

10 mL (diluted) sample

1 mL KCl sodium acetate solution


If necessary, adjust the pH of the solution to 4.6 ± 0.2 .

The concentration is determined by standard addition.


Parameters

Voltammetric	
Electrode operating mode	HMDE
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
Potentiostatic pretreatment	
Potential 1	-1.15 V
Waiting time 1	90 s
Equilibration time	10 s
Sweep	
Start potential	-1.15 V
End potential	0.05 V
Potential step	0.006 V
Potential step time	0.1 s
Sweep rate	0.06 V/s
Pulse amplitude	0.05 V
Substance	
Name	Zn
Characteristic potential	-0.98 V
Name	Cd
Characteristic potential	-0.56 V
Name	Pb
Characteristic potential	-0.38 V
Name	Cu
Characteristic potential	-0.1 V


Example




Zn
 $c = 250.125 \text{ ug/L}$
 $+/- 9.603 \text{ ug/L (3.84\%)}$




Cd
 $c = \text{---}$



Pb
 c = 7.233 µg/l
 +/- 0.321 µg/l (4.44%)



Cu
 c = 7.035 µg/l
 +/- 0.254 µg/l (3.61%)



Results

Sample	Tap water
Sample size	10.0 mL
$\beta(\text{Zn})$	250.1 µg/L
$\beta(\text{Cd})$	n/a
$\beta(\text{Pb})$	7.2 µg/L
$\beta(\text{Cu})$	7.0 µg/L

Comments

- Ammonium acetate buffer pH 4.6 can be used instead of KCl sodium acetate solution. If the sample contains no chloride, it is preferable to use this buffer instead of the KCl sodium acetate solution. Cu achieves peaks with more positive potential (Peak potential Cu: +0.05 V), which are easier to evaluate.
 - NH₄ acetate buffer pH 4.6:
 - c(CH₃COOH) = 2 mol/L
 - c(NH₃) = 1 mol/L
 - 55.5 mL w(CH₃COOH) = 100 % + 37 mL w(NH₃) = 25 % filled up to 500 mL.

- For samples with higher concentrations the time for deposition (Waiting time 1) can be reduced to 60 s or even 30 s.
- Samples with higher concentrations can also be diluted.
- Metals with concentrations >50 mg/L can also be determined on the DME or SMDE without deposition. In this case start and end potential have to be swapped reversing the sweep from positive to negative potentials.
- Zn is often present in tap water in higher concentrations as the other metals. In this case, Zn should be determined separately with shorter deposition time or after being diluted to a greater extent.
- Seeing that the difference of the half wave potentials of Tl and Pb is only approx. 50 mV, Tl and Pb are determined together in acetate solutions. For this reason, following the analysis of Zn, Cd, Pb and Cu, the Tl content is then determined according to method 2. If the concentration of Tl is less than 50% of the concentration of Pb the error for the Pb determination is less than 10%.
- If the Tl concentration is more than 50% of the Pb concentration it is not recommended to use method 1 for determination of Pb. An alternative is a post-electrolysis where a large quantity of the Tl deposited on the electrode is removed before the determination of Pb.
- Another possibility is the determination of Pb and Tl at a pH>12 (NH₄Cl buffer according to method 3 + NaOH). With this electrolyte it is possible to separate the peaks of Pb and Tl. But in this case Cd interferes the determination of Pb. The resulting peak potentials are:

Peak potential Tl	-0.43 V
-------------------	---------

Peak potential Pb	-0.61 V
-------------------	---------

Peak potential Cd	-0.68 V
-------------------	---------

Method 2: Determination of Tl

Thallium is determined in KCl sodium acetate solution at the HMDE by anodic stripping voltammetry (ASV). Cd and Pb, which interfere, are masked by addition of EDTA.

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Sodium hydroxide solution, for trace analysis*, $w(\text{NaOH}) = 30\%$, CAS 1310-73-2
- Acetic acid, for trace analysis*, $w(\text{CH}_3\text{COOH}) = 100\%$, CAS 64-19-7
- Potassium chloride, KCl, for trace analysis*, CAS 7447-40-7
- Ethylenediaminetetraacetic acid disodium salt dihydrate, $\text{Na}_2\text{-EDTA}\cdot 2\text{H}_2\text{O}$, for analysis, CAS 6381-92-6
- Tl standard stock solution, $\beta(\text{Tl}^+) = 1\text{ g/L}$ (commercially available)
- Ultrapure water, resistivity $>18\text{ M}\Omega\cdot\text{cm}$ (25 °C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

KCl sodium acetate solution	$c(\text{KCl}) = 1.5\text{ mol/L}$ $c(\text{CH}_3\text{COONa}) = 0.5\text{ mol/L}$ 55.9 g KCl + 25 mL NaOH + 14.2 mL CH_3COOH are filled up to 500 mL with high purity water.
EDTA solution	$c(\text{EDTA}) = 0.1\text{ mol/L}$ Dissolve 0.93 g $\text{Na}_2\text{-EDTA}\cdot 2\text{H}_2\text{O}$ in ultrapure water and fill up to 25 mL.

Standard solutions

Tl standard solution	$\beta(\text{Tl}^+) = 0.5\text{ mg/L}$ The solution is prepared using $c(\text{HNO}_3) = 0.014\text{ mol/L}$.
----------------------	---

Analysis

Measuring solution

10 mL (diluted) sample

1 mL KCl sodium acetate solution

(or measuring solution from method 1)


0.2 mL EDTA solution

If necessary, adjust the pH of the solution to 4.6 ± 0.2 .


The concentration is determined by standard addition.

Parameters

Voltammetric	
Electrode operating mode	HMDE
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
Potentiostatic pretreatment	
Potential 1	-0.8 V
Waiting time 1	180 s
Equilibration time	10 s
Sweep	
Start potential	-0.8 V
End potential	-0.2 V
Potential step	0.006 V
Potential step time	0.1 s
Sweep rate	0.06 V/s
Pulse amplitude	0.05 V
Substance	
Name	Tl
Characteristic potential	-0.45 V

Example


Tl
 c = 0.305 µg/l
 +/- 0.028 µg/l (9.04%)


Result

Sample	Tap water
Sample size	10.0 mL
β (Tl)	0.3 µg/L

Comments

- Ammonium acetate buffer pH 4.6 can be used instead of KCl sodium acetate solution.
 - NH₄ acetate buffer pH 4.6:
 c(CH₃COOH) = 2 mol/L
 c(NH₃) = 1 mol/L
 55.5 mL w(CH₃COOH) = 100 % + 37 mL w(NH₃) = 25 % filled up to 500 mL.

- For samples with higher concentrations the time for deposition (Waiting time 1) can be reduced to 60 s or even 30 s.
- Samples with higher concentrations can also be diluted.
- Metals with concentrations >50 mg/L can also be determined on the DME or SMDE without deposition. In this case start and end potential have to be swapped reversing the sweep from positive to negative potentials.
- Large excess of Cu makes the determination of thallium impossible. The Cu-EDTA peak is located close to the thallium signal. With high copper concentrations the Tl peak is vanishing in the slope of the Cu-EDTA peak.

Method 3: Determination of Ni and Co

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Hydrochloric acid, w(HCl) = 30%, for trace analysis*, CAS 7647-01-0
- Ammonia solution, for trace analysis*, w(NH₃) = 25%, CAS 1336-21-6
- Dimethylglyoxime disodium salt octahydrate, DMG, for analysis, CAS 75006-64-3
- Ni standard stock solution, β (Ni²⁺) = 1 g/L (commercially available)
- Co standard stock solution, β (Co²⁺) = 1 g/L (commercially available)
- Ultrapure water, resistivity >18 MΩ·cm (25 °C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

NH ₄ Cl buffer pH 9.5	c(NH ₄ Cl) = 1 mol/L c(NH ₃) = 2 mol/L 112.5 mL NH ₃ + 53 mL HCl, filled up to 500 mL with ultrapure water.
----------------------------------	---

DMG solution	c(Na ₂ DMG) = 0.1 mol/L Dissolve 0.304 g dimethylglyoxime disodium salt in 10 mL ultrapure water.
--------------	---

Standard solutions

Ni standard solution	$\beta(\text{Ni}^{2+}) = 1 \text{ mg/L}$ The solution is prepared using $c(\text{HNO}_3) = 0.014 \text{ mol/L}$.
Co standard solution	$\beta(\text{Co}^{2+}) = 0.1 \text{ mg/L}$ The solution is prepared using $c(\text{HNO}_3) = 0.014 \text{ mol/L}$.

Analysis
Measuring solution


10 mL (diluted) sample
0.5 mL NH_4Cl buffer
0.1 mL DMG solution

If necessary, adjust the pH of the solution to 9.5 ± 0.2 .


The concentration is determined by standard addition.

Parameters


Voltammetric	
Electrode operating mode	HMDE
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
Potentiostatic pretreatment	
Potential 1	-0.7 V
Waiting time 1	90 s
Equilibration time	10 s
Sweep	
Start potential	-0.8 V
End potential	-1.25 V
Potential step	0.004 V
Potential step time	0.3 s
Sweep rate	0.013 V/s
Pulse amplitude	0.05 V
Substance	
Name	Ni
Characteristic potential	-0.97 V
Name	Co
Characteristic potential	-1.13 V

Example


Ni
c = 0.611 $\mu\text{g/l}$
+/- 0.027 $\mu\text{g/l}$ (4.41%)



Co
c = 1.860 ng/l
+/- 3.947 ng/l (212.22%)



Result

Sample	Sea water
Sample size	2.0 mL
$\beta(\text{Ni})$	0.61 $\mu\text{g/L}$
$\beta(\text{Co})$	0.002 $\mu\text{g/L}$

Comments

- For samples with higher concentrations the time for deposition (Waiting time 1) can be reduced to 60 s or even 30 s.
- Samples with higher concentrations can also be diluted.
- If the disodium salt of the dimethylglyoxime should not be used, DMG can be dissolved in ethanol or triethanolamine.
 - DMG in ethanol: $c(\text{DMG}) = 0.1 \text{ mol/L}$
Dissolve 0.29 g dimethylglyoxime in ethanol and fill up to 25 mL.
 - or
 - DMG in triethanolamine: $c(\text{DMG}) = \text{saturated}$
Dissolve 500 mg dimethylglyoxime in 100 mL TEA:H₂O (1:1). This solution is used when a big excess of Zn disturbs the Co determination.
- Higher Zn concentrations (low mg/L-range) interfere with the determination of Co.
 - In the DIN standard $c(\text{Na}_2\text{-Iminodiacetate}) = 0.2 \text{ mol/L}$ is used to mask Zn in the Co determination. Our tests have shown that the determination of Ni besides Co is then not possible anymore because Ni is also complexed.
 - To determine Co in the presence of excess Zn the complexing agent α -furildioxime (CAS 6035-71-8)

can be used, see Application Note V-172. With this complexing agent a simultaneous determination of Ni is not possible.

References

- DIN 38406-16
German standard methods for the examination of water, waste water and sludge; cations; determination of 7 metals (zinc, cadmium, lead, copper, thallium, nickel, cobalt) by voltammetry
- Ostapczuk, M. Goedde, M. Stöppler, H. W. Nürnberg Kontroll- und Routinebestimmung von Zn, Cd, Pb, Cu, Ni und Co mit differentieller Pulsvoltammetrie in Materialien der deutschen Umweltbank, Fresenius Z. Anal. Chem. 318 (1984) 321-326
- Braun, M. Metzger, Umweltanalytische Nickel-Bestimmung durch Adsorptionsvoltammetrie mit der Quecksilberfilmelektrode, Fresenius Z. Anal. Chem. 318 (1984) 321-326
- G. Van den Berg, Determining trace concentration of copper in water by cathodic film stripping voltammetry with adsorptive collection, Anal. Lett. 17 A, (1984) 2141-2157 und Application Bulletin 74
- Axel Meyer, Rolf Neeb, Determination of cobalt and nickel by adsorption voltammetry in supporting electrolytes containing triethanolamine and dimethylglyoxime, Fresenius Z. Anal. Chem. 315 (1983) 118-120
- Metrohm Application Bulletin 113

Appendix

Report for the example determination of Zn, Cd, Pb and Cu in tap water according to method 1

===== METROHM 797 VA COMPUTRACE (Version 1.0.0.1) (Serial No. 0) =====

Determination : 06051547_tap water_ZnCdPbCu.dth

Sample ID :

Creator method :

Date :

Time:

Creator determ.:

Date : 2001-06-05

Time: 15:47:51

Modified by :

Date : 2017-07-11

Time: 14:36:09

Method : AB231_1_Det of Zn Cd Pb Cu with HMDE.mth

Title :

Remark1 : 10 ml sample + 1 ml acetate buffer (1.5mol/l KCl + 0.5mol/l NaAc)

Remark2 :

Sample amount : 10.000 mL

Cell volume : 11.000 mL

Substance : Zn

Conc. : 227.386 ug/L

Conc.dev. : 8.730 ug/L (3.84%)

Amount : 2.501 ug

Add.amount : 1.000 ug

VR	V	uA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.989	1.501	1.483	0.024	0.000	
1 - 2	-0.989	1.466				
2 - 1	-0.989	2.024	2.042	0.026	0.559	Ovl. in sweep
2 - 2	-0.989	2.061				
3 - 1	-0.989	2.589	2.626	0.052	0.583	Ovl. in sweep
3 - 2	-0.989	2.663				Ovl. in sweep

Substance : Cd

Conc. : ---

Conc.dev. : ---

Add.amount : 10.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	---	---	---	---		No peak found
1 - 2	---	---				No peak found
2 - 1	-0.596	4.175	4.218	0.061		Ovl. in sweep
2 - 2	-0.596	4.261				
3 - 1	-0.596	7.922	8.077	0.219	3.859	Ovl. in sweep
3 - 2	-0.596	8.232				Ovl. in sweep

Substance : Pb

Conc. : 6.575 ug/L

Conc.dev. : 0.292 ug/L (4.44%)

Amount : 72.328 ng

Add.amount : 50.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.400	15.96	15.71	0.348	0.00	
1 - 2	-0.400	15.46				
2 - 1	-0.400	26.55	26.85	0.427	11.14	Ovl. in sweep
2 - 2	-0.400	27.16				
3 - 1	-0.400	35.85	36.54	0.978	9.69	Ovl. in sweep
3 - 2	-0.400	37.23				Ovl. in sweep

Substance : Cu

Conc. : 6.396 ug/L

Conc.dev. : 0.231 ug/L (3.61%)

Amount : 70.355 ng

Add.amount : 250.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.132	16.8	16.5	0.336	0.0	
1 - 2	-0.132	16.3				
2 - 1	-0.120	76.6	77.0	0.662	60.5	Ovl. in sweep
2 - 2	-0.120	77.5				
3 - 1	-0.114	127.2	128.8	2.307	51.8	Ovl. in sweep
3 - 2	-0.114	130.4				Ovl. in sweep

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Zn	std.add.	1.480e-006	6.511e-003	4.062e-008	0.99842
Cd	std.add.	2.884e-010	4.362e-003	2.276e-010	0.99827

Pb	std.add.	1.578e-008	2.400e-003	7.292e-010	0.99801
Cu	std.add.	1.658e-008	2.592e-003	1.201e-009	0.99900

Final results			+/-	Res. dev.	%	Comments
Zn:						
default	=	250.125 ug/l	9.603		3.839	
Cd:						
default	=	--- ug/l				No result found
Pb:						
default	=	7.233 ug/l	0.321		4.439	
Cu:						
default	=	7.035 ug/l	0.254		3.612	

Method print for the determination of Zn, Cd, Pb, Cu according to method 1

Method parameters

Method : AB231_1_Det of Zn Cd Pb Cu.mth
 Title : Determination of Zinc, Cadmium, Lead, and Copper acc. to DIN38406/16
 Remark1 : 10 mL sample + 1 mL acetate buffer (1.5mol/L KCl + 0.5mol/L NaAc)
 Remark2 :

Calibration : Standard addition
 Technique : Batch
 Addition : Manual

Sample ID : tap water
 Sample amount (mL): 10.000
 Cell volume (mL): 11.000

Voltammetric parameters

Mode : DP - Differential Pulse

Highest current range : 10 mA
 Lowest current range : 100 nA

Electrode : HMDE
 Drop size (1..9) : 4
 Stirrer speed (rpm) : 2000

Initial electr. conditioning : No

No. of additions : 2
 No. of replications : 2

Measure blank : No
 Addition purge time (s) : 10

Initial purge time (s) : 300

Conditioning cycles :
 Start potential (V) : 0.000
 End potential (V) : 0.000
 No. of cycles : 0

Hydrodynamic (measurement) : No
 Cleaning potential (V) : 0.000
 Cleaning time (s) : 0.000
 Deposition potential (V) : -1.150
 Deposition time (s) : 90.000

Sweep :
 Equilibration time (s) : 10.000
 Start potential (V) : -1.150
 End potential (V) : 0.050
 Voltage step (V) : 0.006
 Voltage step time (s) : 0.100
 Sweep rate (V/s) : 0.060
 Pulse amplitude (V) : 0.050
 Pulse time (s) : 0.040

Cell off after measurement : Yes

Peak evaluation

```

-----
Regression technique      : Linear Regression
Peak evaluation          : Height
Minimum peak width (V.steps) : 10
Minimum peak height (A)   : 1.000e-010
Reverse peaks            : No
Smooth factor            : 4
Eliminate spikes         : Yes
  
```

Substances

```

-----
Zn      : -0.980 V +/- 0.050 V

Standard solution      : 1 10.000 mg/L
Addition volume (mL)  : 0.100

default               : Final result (Zn) =
                      Conc * (11 / 10) * (1e+006 / 1) + 0 - 0

Cd      : -0.560 V +/- 0.050 V

Standard solution      : 1 0.100 mg/L
Addition volume (mL)  : 0.100

default               : Final result (Cd) =
                      Conc * (11 / 10) * (1e+006 / 1) + 0 - 0

Pb      : -0.380 V +/- 0.050 V

Standard solution      : 1 0.500 mg/L
Addition volume (mL)  : 0.100

default               : Final result (Pb) =
                      Conc * (11 / 10) * (1e+006 / 1) + 0 - 0

Cu      : -0.100 V +/- 0.050 V

Standard solution      : 1 2.500 mg/L
Addition volume (mL)  : 0.100

default               : Final result (Cu) =
                      Conc * (11 / 10) * (1e+006 / 1) + 0 - 0
  
```

Baseline

Substance	Addition	automatic	start (V)	end (V)	type	scope
Zn	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak
Cd	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak
Pb	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak
Cu	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak

Report for the example determination of Tl according to method 2

```

===== METROHM 797 VA COMPUTRACE (Version 1.0.0.1) (Serial No. 0) =====
Determination : 06051653_tap water_Tl.dth
Sample ID     : tap water_Tl
Creator method : Date : Time:
Creator determ.: Date : 2001-06-05 Time: 16:53:19
Modified by   : --- Date : Time:

Method       : AB231_2_Det of Tl.mth
Title        : Determination of Thallium acc. to DIN38406/16
Remark1      : 10 ml sample + 1 ml acetate buffer (1.5mol/l KCl + 0.5mol/l Na
Ac)
Remark2      : + 200µl EDTA-Na2 (0.1 mol/l)

Sample amount : 10.000 mL
Cell volume   : 11.200 mL
  
```

Substance : Tl
 Conc. : 272.059 ng/L
 Conc.dev. : 24.589 ng/L (9.04%)
 Amount : 3.047 ng
 Add.amount : 10.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.467	0.549	0.514	0.049	0.000	
1 - 2	-0.467	0.480				
2 - 1	-0.461	2.063	2.112	0.069	1.598	
2 - 2	-0.461	2.161				
3 - 1	-0.455	3.809	3.786	0.033	1.674	
3 - 2	-0.461	3.763				

Final results	+/-	Res. dev.	%	Comments
Tl: default	=	0.305 ug/l	0.028	9.038

Method print for the determination of Tl according to method 2

Method parameters

Method : AB231_2_Det of Tl.mth
 Title : Determination of Thallium acc. to DIN38406/16
 Remark1 : 10 mL sample + 1 mL acetate buffer (1.5mol/L KCl + 0.5mol/L Na Ac)
 Remark2 : + 200µl EDTA-Na2 (0.1 mol/L)
 Calibration : Standard addition
 Technique : Batch
 Addition : Manual

Sample ID : sample
 Sample amount (mL): 10.000
 Cell volume (mL): 11.200

Voltammetric parameters

Mode : DP - Differential Pulse
 Highest current range : 10 mA
 Lowest current range : 100 nA
 Electrode : HMDE
 Drop size (1..9) : 4
 Stirrer speed (rpm) : 2000
 Initial electr. conditioning : No
 No. of additions : 2
 No. of replications : 2
 Measure blank : No
 Addition purge time (s) : 10
 Initial purge time (s) : 300
 Conditioning cycles :
 Start potential (V) : 0.000
 End potential (V) : 0.000
 No. of cycles : 0
 Hydrodynamic (measurement) : No
 Cleaning potential (V) : 0.000
 Cleaning time (s) : 0.000
 Deposition potential (V) : -0.800
 Deposition time (s) : 180.000
 Sweep :
 Equilibration time (s) : 10.000
 Start potential (V) : -0.800
 End potential (V) : -0.200
 Voltage step (V) : 0.006
 Voltage step time (s) : 0.100
 Sweep rate (V/s) : 0.060
 Pulse amplitude (V) : 0.050
 Pulse time (s) : 0.040

Cell off after measurement : Yes

Peak evaluation

 Regression technique : Linear Regression
 Peak evaluation : Height
 Minimum peak width (V.steps) : 10
 Minimum peak height (A) : 1.000e-010
 Reverse peaks : No
 Smooth factor : 4
 Eliminate spikes : Yes

Substances

 Tl : -0.450 V +/- 0.050 V
 Standard solution : 1 0.500 mg/L
 Addition volume (mL) : 0.020
 default : Final result (Tl) =
 Conc * (11.2 / 10) * (1e+006 / 1) + 0 - 0

Baseline

Substance	Addition	automatic	start (V)	end (V)	type	scope
Tl	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak

Report for the example determination of Ni and Co according to method 3

===== METROHM 797 VA COMPUTRACE (Version 1.0.0.1) (Serial No. 0) =====
 Determination : 09151108_NiCoMeerwasser.dth
 Sample ID : Meerwasser
 Creator method : Date : Time:
 Creator determ.: Date : 1999-09-15 Time: 11:08:54
 Modified by : --- Date : Time:

 Method : Ni-Co Meerwasser.mth
 Title : Determination of Nickel and Cobalt acc. to DIN38406/16
 Remark1 : 10 ml sample + 0.5 mL ammonia buffer
 Remark2 : + 0.1 mL DMG solution

 Sample amount : 10.000 mL
 Cell volume : 11.050 mL

 Substance : Ni
 Conc. : 553.296 ng/L
 Conc.dev. : 24.389 ng/L (4.41%)
 Amount : 6.114 ng
 Add.amount : 10.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.927	-29.5	-30.4	1.276	0.0	
1 - 2	-0.927	-31.3				
2 - 1	-0.927	-78.8	-80.8	2.915	-50.4	
2 - 2	-0.927	-82.9				
3 - 1	-0.927	-126.3	-128.6	3.258	-47.8	
3 - 2	-0.923	-130.9				

 Substance : Co
 Conc. : 1.683 ng/L
 Conc.dev. : 3.572 ng/L (212.22%)
 Amount : 18.597 pg
 Add.amount : 1.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-1.089	-0.31	-0.34	0.046	0.00	
1 - 2	-1.086	-0.37				
2 - 1	-1.086	-11.62	-12.26	0.900	-11.92	
2 - 2	-1.086	-12.89				
3 - 1	-1.082	-24.10	-25.02	1.290	-12.76	
3 - 2	-1.082	-25.93				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Ni	std.add.	-3.047e-008	-5.507e-002	2.146e-009	0.99883
Co	std.add.	-2.280e-010	-1.355e-001	1.012e-009	0.99778

Final results		+/-	Res. dev.	%	Comments
Ni: default	=	0.611 ug/l	0.027	4.408	
Co: default	=	1.860 ng/l	3.947	212.215	

Method print for the determination of Ni and Co according to method 3

Method parameters	
Method	: AB231_3_Det of Ni Co.mth
Title	: Determination of Nickel and Cobalt acc. to DIN38406/16
Remark1	: 10 ml sample + 0.5 mL ammonia buffer
Remark2	: + 0.1 mL DMG solution
Calibration	: Standard addition
Technique	: Batch
Addition	: Manual
Sample ID	: sample
Sample amount (mL)	: 10.000
Cell volume (mL)	: 10.600
Voltammetric parameters	
Mode	: DP - Differential Pulse
Highest current range	: 10 mA
Lowest current range	: 100 nA
Electrode	: HMDE
Drop size (1..9)	: 4
Stirrer speed (rpm)	: 2000
Initial electr. conditioning	: No
No. of additions	: 2
No. of replications	: 2
Measure blank	: No
Addition purge time (s)	: 10
Initial purge time (s)	: 300
Conditioning cycles	
Start potential (V)	: 0.000
End potential (V)	: 0.000
No. of cycles	: 0
Hydrodynamic (measurement)	: No
Cleaning potential (V)	: 0.000
Cleaning time (s)	: 0.000
Deposition potential (V)	: -0.700
Deposition time (s)	: 90.000
Sweep	
Equilibration time (s)	: 10.000
Start potential (V)	: -0.800
End potential (V)	: -1.250
Voltage step (V)	: 0.004
Voltage step time (s)	: 0.300
Sweep rate (V/s)	: 0.013
Pulse amplitude (V)	: 0.050
Pulse time (s)	: 0.040
Cell off after measurement	: Yes
Peak evaluation	
Regression technique	: Linear Regression
Peak evaluation	: Height
Minimum peak width (V.steps)	: 10
Minimum peak height (A)	: 1.000e-010
Reverse peaks	: No
Smooth factor	: 4
Eliminate spikes	: Yes

Substances

```

-----
Ni
: -0.970 V +/- 0.050 V

Standard solution : 1 1.000 mg/L
Addition volume (mL) : 0.010

default : Final result (Ni) =
          Conc * (10.6 / 10) * (1e+006 / 1) + 0 - 0

Co
: -1.130 V +/- 0.050 V

Standard solution : 1 0.100 mg/L
Addition volume (mL) : 0.010

default : Final result (Co) =
          Conc * (10.6 / 10) * (1e+006 / 1) + 0 - 0
  
```

Baseline

```

-----
Substance Addition automatic start (V) end (V) type scope
-----
Ni Sample yes --- --- linear wholePeak
  Addition 1 yes --- --- linear wholePeak
  Addition 2 yes --- --- linear wholePeak
-----
Co Sample yes --- --- linear wholePeak
  Addition 1 yes --- --- linear wholePeak
  Addition 2 yes --- --- linear wholePeak
-----
  
```