

Application Bulletin 241/2 e

Determination of cadmium and lead by anodic stripping voltammetry at a mercury film electrode

Summary

This Application Bulletin describes the determination of cadmium and lead at a mercury film electrode (MFE) by anodic stripping voltammetry (ASV). The mercury film is plated *ex situ* on a glassy carbon electrode and can be used for up to one day. With a deposition time of 30 s, the limit of detection is $\beta(\text{Cd}^{2+}) = 0.02 \mu\text{g/L}$ and $\beta(\text{Pb}^{2+}) = 0.05 \mu\text{g/L}$. The linear working range for both elements goes up to approx. 50 $\mu\text{g/L}$ using the same deposition time.

Samples

Surface water, ground water, sea water, waste water

Instruments

797 VA Computrace	2.797.0020
<i>Accessories</i>	
Driving belt	6.1244.020
Measuring vessel 5 mL	6.1415.150
Stopper	6.2709.040
Polishing set for solid-state electrodes	6.2802.000

Electrodes

WE	Glassy carbon electrode tip	6.1204.600
	Driving axle for RDE	6.1204.210
RE	Ag/AgCl reference electrode	6.0728.020
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel	6.1245.010
	Filled with $c(\text{KCl}) = 3 \text{ mol/L}$	
AE	Glassy carbon rod	6.1247.000
	Electrode holder	6.1241.020

Reagents

- Cd standard stock solution, $\beta(\text{Cd}^{2+}) = 1 \text{ g/L}$
- Pb standard stock solution, $\beta(\text{Pb}^{2+}) = 1 \text{ g/L}$
- Hg standard stock solution, $\beta(\text{Hg}^{2+}) = 1 \text{ g/L}$
- Nitric acid, $w(\text{HNO}_3) = 65\%$, for trace analysis*, CAS 7697-37-2

- Hydrochloric acid, $w(\text{HCl}) = 30\%$, for trace analysis*, CAS 7647-01-0
- Acetic acid, $w(\text{CH}_3\text{COOH}) = 100\%$, for trace analysis*, CAS 64-19-7
- Sulfuric acid, $w(\text{H}_2\text{SO}_4) = 96\%$, for analysis, CAS 7664-93-9
- Ammonium hydroxide solution, $w(\text{NH}_3) = 25\%$, for trace analysis*, CAS 1336-21-6
- Potassium chloride, KCl, for trace analysis*, CAS 7447-40-7
- Ultrapure water, resistivity $>18 \text{ M}\Omega \cdot \text{cm}$ (25 °C), type I grade (ASTM D1193)

In addition, if UV digestion is required:

- Hydrogen peroxide solution, $w(\text{H}_2\text{O}_2) = 30\%$, for trace analysis*, CAS 7722-84-1

* e.g., Merck suprapur®, Sigma-Aldrich TraceSelect® or equivalent

Solutions

Hg plating solution	$\beta(\text{Hg}^{2+}) = 20 \text{ mg/L}$ $c(\text{HCl}) = 0.1 \text{ mol/L}$ 0.4 mL mercury standard stock solution and 0.2 mL hydrochloric acid are diluted to 20 mL with ultrapure water.
Electrolyte	$c(\text{CH}_3\text{COOH}) = 2 \text{ mol/L}$ $c(\text{NH}_3) = 1 \text{ mol/L}$ $c(\text{KCl}) = 1.5 \text{ mol/L}$ 55.9 g KCl are dissolved in approx. 400 mL water. 55.5 mL acetic acid and 37 mL ammonium hydroxide solution are added and filled up to 500 mL with ultrapure water.

Standard solutions

Cd standard solution 1 mg/L	$\beta(\text{Cd}^{2+}) = 1 \text{ mg/L}$ Approx. 40 mL ultrapure water is filled into a 50 mL volumetric flask. 0.05 mL $w(\text{HNO}_3) = 65\%$ and 0.05 mL Cd standard stock solution are added. The solution is made up to the mark with ultrapure water.
Cd standard solution 0.1 mg/L	$\beta(\text{Cd}^{2+}) = 0.1 \text{ mg/L}$ Approx. 40 mL ultrapure water is filled into a 50 mL volumetric flask. 0.05 mL $w(\text{HNO}_3) = 65\%$ and 5 mL Cd standard solution 1 mg/L are added. The solution is made up to the mark with ultrapure water.
Pb standard solution 1 mg/L	$\beta(\text{Pb}^{2+}) = 1 \text{ mg/L}$ Approx. 40 mL ultrapure water is filled into a 50 mL volumetric flask. 0.05 mL $w(\text{HNO}_3) = 65\%$ and 0.05 mL Pb standard stock solution are added. The solution is made up to the mark with ultrapure water.
Pb standard solution 0.1 mg/L	$\beta(\text{Pb}^{2+}) = 0.1 \text{ mg/L}$ Approx. 40 mL ultrapure water is filled into a 50 mL volumetric flask. 0.05 mL $w(\text{HNO}_3) = 65\%$ and 5 mL Pb standard solution 1 mg/L are added. The solution is made up to the mark with ultrapure water.

Sample preparation

- Ground water, drinking water, sea water, and mineral water can usually be analyzed directly.
- Water that contains interfering organic substances is digested using the 909 UV Digester:
10 mL acidified water sample ($\text{pH} = 2$) with 10 μL $w(\text{HNO}_3) = 65\%$ and 100 μL $w(\text{H}_2\text{O}_2) = 30\%$ are irradiated for 90 min at 90 °C.

Comments

- No platinum auxiliary electrode should be used in combination with a glassy carbon working electrode.

- Measuring vessel and reference electrode have to be free of platinum traces. Therefore it is recommended to have accessories dedicated to this application.

Mercury film plating

Preparation of the glassy carbon electrode

The glassy carbon electrode has to be clean and free of any insulating material before the mercury film can be plated. For cleaning the following procedure is recommended.

- With a soft tissue deposit, e.g. an old mercury film, is wiped off the electrode surface.
- For polishing the glassy carbon electrode a small amount of aluminum oxide (5 ... 10 mg) is mixed with a few drops of water (0.2 ... 0.5 mL) on the polishing cloth (part of the polishing set 6.2802.000) to form a slurry. The electrode is then polished with small 8-shaped movements on the polishing cloth.
- The slurry is rinsed off the electrode with ultrapure water.
- The electrode is then rinsed with ethanol and afterwards thoroughly with ultrapure water.

This procedure should be carried out each time before plating the mercury film.

Plating the mercury film

20 mL Hg plating solution are transferred into the measuring vessel. The plating is carried out using the parameters given under «Parameters for mercury film plating».

Measuring solution


20 mL Hg plating solution

Parameters for mercury film plating

Determination	
No. of additions	0
No. of replications	2
Voltammetric	
Electrode	SSE/RDE
Measuring mode	DP
Stirring speed	2000 min^{-1}
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-1.2 V

End potential	-0.1 V
No. of cycles	50
Pretreatment	
Cleaning potential	-0.15 V
Cleaning time	2 s
Deposition potential	-1.3 V
Deposition time	180 s
Sweep	
Equilibration time	5 s
Start potential	-0.8 V
End potential	-0.1 V
Pulse amplitude	0.05 V
Pulse time	0.04 s
Potential step	0.006 V
Potential step time	0.1 s
Sweep rate	0.06 V/s

Example curves for mercury film plating



The peaks at -0.5 V and -0.3 V are related to lead and copper contamination of the Hg plating solution. These signals increase the more often the solution is used.

Comments

- The Hg plating solution can be reused several times (approx. 20 times). The solution itself is stable for at least one month. But with the number of plating processes the concentration of Hg in the solution decreases and the amount of contamination increases. Therefore the solution has to be replaced when no proper mercury film can be plated anymore.

- The waste Hg plating solution still contains hazardous amounts of mercury. Therefore, care has to be taken for an appropriate disposal, in accordance with the local legislation.

Assessment of the mercury film

To assess the quality of the mercury film it is recommended to run a blank determination. The measurement is carried out using the parameters given under «Parameters for mercury film assessment».

Measuring solution

10 mL ultrapure water

1 mL electrolyte

Standard addition

0.1 mL $\beta(\text{Cd}^{2+}) = 0.1 \text{ mg/L}$


0.1 mL $\beta(\text{Pb}^{2+}) = 0.1 \text{ mg/L}$

Parameters for mercury film assessment


Determination	
No. of additions	2
No. of replications	2
Voltammetric	
Electrode	SSE/RDE
Measuring mode	DP
Stirring speed	2000 min ⁻¹
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-0.8 V
End potential	-0.4 V
No. of cycles	25
Pretreatment	
Cleaning potential	-0.4 V
Cleaning time	2 s
Deposition potential	-1.0 V
Deposition time	90 s
Sweep	
Equilibration time	5 s
Start potential	-0.85 V

End potential	-0.35 V
Pulse amplitude	0.05 V
Pulse time	0.04 s
Potential step	0.004 V
Potential step time	0.1 s
Sweep rate	0.04 V/s
Substance + calibration	
Calibration	Standard addition
Name	Cd
Peak potential	-0.7 V
Name	Pb
Peak potential	-0.5 V

Example for a blank determination with a good mercury film



Example for a blank determination with a bad mercury film



Assessment criteria

- The curve of the background current should be smooth.
- With a proper mercury film the background current should be smaller than 200 nA.
- When the measuring solution contains equal mass concentrations of Cd and Pb, the Cd peak should have approximately double the height of the Pb peak.

If the background current is not smooth or too high it is recommended to remove the mercury film, carry out the cleaning and polishing procedure as described under «Preparation of the glassy carbon electrode» and plate a fresh mercury film.

If the ratio of Cd and Pb peak is not ok the reason is usually that the mercury film is too thin. Repeat the plating as described under «Plating the mercury film». If that does not improve the situation, remove the mercury film, carry out the cleaning and polishing procedure as described under «Preparation of the glassy carbon electrode» and plate a fresh mercury film.

Determination of Cd and Pb

10 mL (diluted) sample and 1 mL electrolyte are pipetted into the measuring vessel. The determination is carried out using the parameters given under «Parameters for determination of cadmium and lead».

The concentration of cadmium and lead is quantified by two additions of Cd and Pb standard solutions.

Measuring solution

10 mL (diluted) sample

1 mL electrolyte

Parameters for determination of cadmium and lead

Determination	
No. of additions	2
No. of replications	2
Voltammetric	
Electrode	SSE/RDE
Measuring mode	DP
Stirring speed	2000 min ⁻¹
Hydrodynamic measurement	No
Conditioning cycles	
Start potential	-0.8 V

End potential	-0.4 V
No. of cycles	25
Pretreatment	
Cleaning potential	-0.4 V
Cleaning time	2 s
Deposition potential	-1.0 V
Deposition time	90 s
Sweep	
Equilibration time	5 s
Start potential	-0.85 V
End potential	-0.35 V
Pulse amplitude	0.05 V
Pulse time	0.04 s
Potential step	0.004 V
Potential step time	0.1 s
Sweep rate	0.04 V/s
Substance + calibration	
Calibration	Standard addition
Name	Cd
Peak potential	-0.7 V
Name	Pb
Peak potential	-0.5 V

Results


Average of 3 determinations with relative standard deviation.

Sample	$\beta(\text{Cd})$	$\beta(\text{Pb})$
Tap water	n/a	0.35 $\pm 9.1\%$
Ground water reference material BCR-609	0.19 $\mu\text{g/L}$ $\pm 3.4\%$	1.36 $\mu\text{g/L}$ $\pm 4.3\%$
Ground water reference material BCR-610	3.14 $\mu\text{g/L}$ $\pm 1.3\%$	7.02 $\mu\text{g/L}$ $\pm 1.3\%$
Sea water spiked with 1 $\mu\text{g/L}$ Cd and Pb each	1.10 $\mu\text{g/L}$ $\pm 3.4\%$	1.09 $\mu\text{g/L}$ $\pm 1.0\%$


The following table shows the certified concentrations for the ground water reference materials and the measured result relative to the certified value.


Sample	Certified value	Result / certified value
BCR-609	$\beta(\text{Cd}) = 0.164 \mu\text{g/kg}$	115.8%
	$\beta(\text{Pb}) = 1.63 \mu\text{g/kg}$	83.4%
BCR-610	$\beta(\text{Cd}) = 2.94 \mu\text{g/kg}$	106.8%
	$\beta(\text{Pb}) = 7.78 \mu\text{g/kg}$	90.2%

Example for Cd and Pb determination in ground water reference material BCR-610



Cd
 $c = 3.182 \mu\text{g/L}$
 $\pm 0.023 \mu\text{g/L} (0.72\%)$



Pb
c = 7.116 $\mu\text{g/L}$
+/- 0.032 $\mu\text{g/L}$ (0.45%)


Limit of detection and linear working range

The limit of detection was determined using the «regression approach» [1], where the limit of detection is calculated as « $3 \cdot s_y$ », with s_y as the residual standard deviation of a linear regression.

The following table gives an overview over limit of detection depending on the deposition time.

Deposition time	Limit of detection Cd	Limit of detection Pb
30 s	0.02 $\mu\text{g/L}$	0.05 $\mu\text{g/L}$
90 s	0.007 $\mu\text{g/L}$	0.025 $\mu\text{g/L}$

The linear working range was read out from a calibration curve. With a deposition time of 30 s the linear working range for both elements goes up to approx. 50 $\mu\text{g/L}$. For concentrations above 50 $\mu\text{g/L}$ the calibration curve gets slightly convex, why recoveries of < 90% have to be expected.

The limit of detection and linear working range were determined in standard solutions. Depending on the condition of the electrode or the sample matrix, these values can differ.

Interferences

Tl	Broad peak at approx. -0.6 V, mainly overlapping with the cadmium peak. Shows about 25% of the sensitivity of cadmium and 50% of the lead.
Sn(VI)	Peak at approx. -0.6 V. Shows about 5% of the sensitivity of cadmium and 10% of lead.

Pt, Rh

Platinum and other platinum metals reduce the overpotential for the hydrogen reduction at the mercury film. Therefore even concentrations in the $\mu\text{g/L}$ range show an interfering background current caused by hydrogen formation.

Comments

- Deposition time and standard addition (volume and/or concentration of the standard solution) have to be adapted to the concentration of analyte in the measuring vessel. In this context it is important to consider the linear working range.
- The mercury film can be used for half a day up to one day. But effects from the sample matrix can reduce the lifetime of the mercury film. Indications for a fading mercury film are:
 - bad reproducibility of the replications
 - decreasing sensitivity of the additions (slope) by more than one third
 - increasing or noisy background current

The mentioned indications can also be caused by matrix interferences. Therefore, in case of a suspicious mercury film, it is recommended to run a blank determination as described under «Assessment of the mercury film» and to decide based on these curves whether the film needs to be renewed.

- If the mercury film is not used it is recommended to keep either diluted electrolyte or diluted nitric acid ($c(\text{HNO}_3) \sim 0.05 \text{ mol/L}$) in the measuring vessel. Storing the mercury film with ultrapure water is not recommended since it passivates the film.
- From time to time, electrodes and measuring vessel should be cleaned with sulfuric acid. For that, fill 1–2 mL concentrated sulfuric acid into the slightly wet measuring vessel. Attention, the solution gets hot! Carefully turn the measuring vessel to wet the entire surface with H_2SO_4 . Then add 30–40 mL ultrapure water. Be very careful of splashing of the hot solution. Place the measuring vessel under the measuring head and turn on the stirring until the solution cooled down to room temperature. Empty the measuring vessel and rinse everything thoroughly with ultrapure water.
- Traces of metals can be removed with diluted nitric acid ($c(\text{HNO}_3) \sim 0.1 \text{ mol/L}$). For cleaning, place 20–30 mL diluted nitric acid in the measuring vessel and switch on the stirrer for 10–30 min. Empty the measuring vessel and rinse everything thoroughly with ultrapure water.

References

- [1] J. Mocak, A. Bond, S. Mitchell and G. Scollary, "A statistical overview of standard (IUPAC and ACS) and new procedures for determining the limits of detection and quantification: Application to voltammetric and stripping techniques," *Pure and Applied Chemistry*, vol. 69, no. 2, pp. 297-328, 1997.
- [2] A. M. Bond, *Modern Polarographic Methods in Analytical Chemistry*, Marcel Dekker, 1980.
- [3] K. Brainina and E. Neyman, *Electroanalytical Stripping Methods*, John Wiley & Sons, 1993.
- [4] H. Nürnberg, *Electroanalytical Chemistry*, John Wiley & Sons, 1974.
- [5] M. Smyth and J. Vos, *Analytical Voltammetry*, Elsevier Science Publishers, 1992.

Appendix

Method print «Parameters for mercury film plating»

Method parameters						
Method	: AB 241_Plating of mercury film.mth					
Title	: Plating Hg film on GC electrode 6.1204.600					
Remark1	: 20 mL Hg plating solution					
Remark2	: Plating solution: c(Hg) = 20 mg/L, c(HCl) = 0.1 mol/L					
Calibration	: Standard addition					
Technique	: Batch					
Addition	: Manual					
Sample ID	: Hg film plating					
Sample amount (mL)	: 20.000					
Cell volume (mL)	: 20.000					
Voltammetric parameters						
Mode	: DP - Differential Pulse					
Highest current range	: 10 mA					
Lowest current range	: 100 nA					
Electrode	: SSE/RDE					
Stirrer speed (rpm)	: 2000					
Initial electr. conditioning	: No					
No. of additions	: 0					
No. of replications	: 2					
Measure blank	: No					
Addition purge time (s)	: 20					
Initial purge time (s)	: 60					
Conditioning cycles						
Start potential (V)	: -1.200					
End potential (V)	: -0.100					
No. of cycles	: 50					
Hydrodynamic (measurement)	: No					
Cleaning potential (V)	: -0.150					
Cleaning time (s)	: 2.000					
Deposition potential (V)	: -1.300					
Deposition time (s)	: 180.000					
Sweep						
Equilibration time (s)	: 5.000					
Start potential (V)	: -0.800					
End potential (V)	: -0.100					
Voltage step (V)	: 0.006					
Voltage step time (s)	: 0.100					
Sweep rate (V/s)	: 0.060					
Pulse amplitude (V)	: 0.050					
Pulse time (s)	: 0.040					
Cell off after measurement	: Yes					
Peak evaluation						
Regression technique	: Linear Regression					
Peak evaluation	: Height					
Minimum peak width (V.steps)	: 10					
Minimum peak height (A)	: 1.000e-010					
Reverse peaks	: No					
Smooth factor	: 4					
Eliminate spikes	: Yes					
Substances						
Baseline						
Substance	Addition	automatic start	(V)	end (V)	type	scope

Method print «Parameters for determination of cadmium and lead»

These parameters are identical to «Parameters for mercury film assessment»

```

-----
Method parameters
-----
Method      : AB 241_Determination of Cd Pb at MFE.mth
Title       : Determination of Cd Pb at mercury film electrode (MFE)
Remark1     : 10 mL sample + 1 mL electrolyte
Remark2     : Electrolyte: c(HAC) = 2 mol/L, c(NH3) = 1 mol/L, c(KCl) = 1.5 mol/L

Calibration : Standard addition
Technique   : Batch
Addition    : Manual

Sample ID   : Sample
Sample amount (mL): 10.000
Cell volume (mL): 11.000

Voltammetric parameters
-----
Mode                : DP - Differential Pulse

Highest current range : 10 mA
Lowest current range  : 100 nA

Electrode           : SSE/RDE
Stirrer speed (rpm) : 2000

Initial electr.conditioning : No

No. of additions    : 2
No. of replications : 2

Measure blank       : No
Addition purge time (s) : 20

Initial purge time (s) : 300

Conditioning cycles
Start potential (V)   : -0.800
End potential (V)     : -0.400
No. of cycles         : 25

Hydrodynamic (measurement) : No
Cleaning potential (V)   : -0.400
Cleaning time (s)       : 2.000
Deposition potential (V) : -1.000
Deposition time (s)     : 90.000

Sweep
Equilibration time (s) : 5.000
Start potential (V)     : -0.850
End potential (V)       : -0.350
Voltage step (V)        : 0.004
Voltage step time (s)   : 0.100
Sweep rate (V/s)       : 0.040
Pulse amplitude (V)     : 0.050
Pulse time (s)         : 0.040

Cell off after measurement : Yes

Peak evaluation
-----
Regression technique : Linear Regression
Peak evaluation      : Height
Minimum peak width (V.steps) : 5
Minimum peak height (A) : 1.000e-010
Reverse peaks        : No
Smooth factor        : 2
Eliminate spikes     : Yes

Substances
-----
Cd                : -0.700 V +/- 0.050 V

Standard solution : 1 0.100 mg/L
Addition volume (mL) : 0.100

Cadmium           : Final result (Cd) =
                   Conc * (11 / 10) * (1e+006 / 1) + 0 - 0


Pb                : -0.500 V +/- 0.050 V
  
```

Standard solution : 1 0.100 mg/L
 Addition volume (mL) : 0.100

Lead : Final result (Pb) =
 $\text{Conc} * (11 / 10) * (1e+006 / 1) + 0 - 0$

Baseline

Substance	Addition	automatic	start (V)	end (V)	type	scope
Cd	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak
Pb	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak

Example for the determination of cadmium and lead in tap water


===== METROHM 797 VA COMPUTRACE (Version 1.3.2.85) (Serial No. 18154) =====

Determination : 1410170922_Tap water.dth
 Sample ID : Tap water
 Creator method : zu Date : 2014-07-10 Time: 11:30:05
 Creator determ. : zu Date : 2014-10-17 Time: 09:22:23
 Modified by : --- Date : Time:

Method : Determination Cd Pb on MFE auto.mth
 Title : Determination of Cd Pb on mercury film electrode
 Remark1 : 5mL sample
 Remark2 : + 0.5 mL acetate buffer / KCl

Sample amount : 10.000 mL
 Cell volume : 11.000 mL

Substance : Cd
 Conc. : ---
 Conc.dev. : ---
 Add.amount : 2.500 ng


VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	---	---	---	---	---	No peak found
1 - 2	---	---	---	---	---	No peak found
2 - 1	-0.691	36.41	36.63	0.314		
2 - 2	-0.691	36.86				
3 - 1	-0.691	71.46	70.99	0.654	34.36	
3 - 2	-0.691	70.53				

Substance : Pb
 Conc. : 346.173 ng/L
 Conc.dev. : 17.744 ng/L (5.13%)
 Amount : 3.808 ng
 Add.amount : 2.500 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.505	24.90	24.07	1.174	0.00	
1 - 2	-0.505	23.24				
2 - 1	-0.505	39.15	39.28	0.185	15.20	
2 - 2	-0.505	39.41				
3 - 1	-0.505	55.63	55.65	0.031	16.37	
3 - 2	-0.505	55.67				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Cd	std.add.	2.119e-009	1.522e-001	5.793e-010	0.99978
Pb	std.add.	2.400e-008	6.932e-002	1.230e-009	0.99905

Final results	+/-	Res. dev.	%	Comments
Cd: Cadmium	=	---	ug/L	No result found
Pb: Lead	=	0.381	ug/L	0.020 5.126

Example for the determination of cadmium and lead in groundwater reference material BCR-609


===== METROHM 797 VA COMPUTRACE (Version 1.3.2.85) (Serial No. 2173) =====

Determination : 1405200849_BCR-609_direkt.dth
 Sample ID : BCR-609_diRekt
 Creator method : ka Date : 2014-04-07 Time: 12:59:42
 Creator determ.: ka Date : 2014-05-20 Time: 08:49:24
 Modified by : --- Date : Time:

Method : GC Hg-ex situ.mth
 Title : GC 61204600 Electrode
 Remark1 : BCR-609 Direkt
 Remark2 : 10 mL_BCR-609_direkt

Sample amount : 10.000 mL
 Cell volume : 12.500 mL

Substance : Cd
 Conc. : 153.061 ng/L
 Conc.dev. : 6.417 ng/L (4.19%)
 Amount : 1.913 ng
 Add.amount : 10.000 ng


VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.687	24.0	23.3	1.025	0.0	
1 - 2	-0.683	22.6				
2 - 1	-0.683	151.7	149.2	3.540	125.9	
2 - 2	-0.683	146.7				
3 - 1	-0.683	264.6	260.6	5.703	111.4	
3 - 2	-0.683	256.5				

Substance : Pb
 Conc. : 1.123 ug/L
 Conc.dev. : 0.035 ug/L (3.13%)
 Amount : 14.037 ng
 Add.amount : 10.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.509	111.4	111.4	0.071	0.0	
1 - 2	-0.509	111.5				
2 - 1	-0.509	193.3	194.9	2.256	83.5	
2 - 2	-0.509	196.5				
3 - 1	-0.509	266.4	267.9	2.053	73.0	
3 - 2	-0.509	269.3				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Cd	std.add.	2.333e-008	1.524e-001	2.096e-009	0.99899
Pb	std.add.	1.120e-007	9.970e-002	3.719e-009	0.99907

Final results		+/-	Res. dev.	%	Comments
Cd:					
Cadmium	=	0.191 ug/L	0.008	4.192	
Pb:					
Lead	=	1.404 ug/L	0.044	3.131	

Example for the determination of cadmium and lead in groundwater reference material BCR-610


===== METROHM 797 VA COMPUTRACE (Version 1.3.2.85) (Serial No. 2173) =====

Determination : 1405221026_BCR-610_direkt.dth
 Sample ID : BCR-610_diRekt
 Creator method : ka Date : 2014-04-07 Time: 12:59:42
 Creator determ.: zu Date : 2014-10-10 Time: 12:52:10
 Modified by : --- Date : Time:

Method : GC Hg-ex situ.mth
 Title : GC 61204600 Electrode
 Remark1 : BCR-610 direkt
 Remark2 : 10 mL_BCR-610_direkt

Sample amount : 10.000 mL
 Cell volume : 12.500 mL

Substance : Cd
 Conc. : 2.546 ug/L
 Conc.dev. : 0.018 ug/L (0.72%)
 Amount : 31.820 ng
 Add.amount : 50.000 ng


VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.691	126.4	126.6	0.250	0.0	
1 - 2	-0.691	126.7				
2 - 1	-0.691	321.4	322.9	2.127	196.3	
2 - 2	-0.687	324.4				
3 - 1	-0.687	519.7	521.2	2.180	198.4	
3 - 2	-0.687	522.8				

Substance : Pb
 Conc. : 5.693 ug/L
 Conc.dev. : 0.026 ug/L (0.45%)
 Amount : 71.161 ng
 Add.amount : 50.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.513	159.7	159.7	0.126	0.0	
1 - 2	-0.513	159.6				
2 - 1	-0.513	269.5	270.1	0.836	110.5	
2 - 2	-0.513	270.7				
3 - 1	-0.513	381.6	381.3	0.439	111.2	
3 - 2	-0.513	381.0				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Cd	std.add.	1.265e-007	4.969e-002	1.427e-009	0.99996
Pb	std.add.	1.596e-007	2.803e-002	7.512e-010	0.99998

Final results		+/-	Res. dev.	%	Comments
Cd:					
Zinc	=	3.182	ug/L	0.023	0.719
Pb:					
Cadmium	=	7.116	ug/L	0.032	0.450

Example for the determination of cadmium and lead in spiked sea water


===== METROHM 797 VA COMPUTRACE (Version 1.3.2.85) (Serial No. 2173) =====

Determination : 1405271451 Meerwasser_dot lppb_direkt.dth
 Sample ID : Meerwasser_dot lppb_direkt
 Creator method : ka Date : 2014-04-07 Time: 12:59:42
 Creator determ.: ka Date : 2014-05-27 Time: 14:51:48
 Modified by : --- Date : Time:

Method : GC Hg-ex situ.mth
 Title : GC_61204600 Electrode
 Remark1 : Meerwasser_pH 2_dot lppb_direkt
 Remark2 : 10 mL_Meerwasser_dot lppb Ac_KCl_DIREKT

Sample amount : 10.000 mL
 Cell volume : 12.510 mL

Substance : Cd
 Conc. : 878.253 ng/L
 Conc.dev. : 37.375 ng/L (4.26%)
 Amount : 10.987 ng
 Add.amount : 10.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.707	85.8	83.7	2.952	0.0	
1 - 2	-0.703	81.6				
2 - 1	-0.703	160.9	161.0	0.070	77.3	
2 - 2	-0.703	161.0				
3 - 1	-0.707	240.7	235.0	7.948	74.1	
3 - 2	-0.711	229.4				

Substance : Pb
 Conc. : 873.709 ng/L
 Conc.dev. : 20.666 ng/L (2.37%)
 Amount : 10.930 ng
 Add.amount : 10.000 ng

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.520	65.0	64.6	0.516	0.0	
1 - 2	-0.520	64.2				
2 - 1	-0.520	126.2	126.0	0.370	61.4	
2 - 2	-0.520	125.7				
3 - 1	-0.520	181.3	181.3	0.001	55.3	
3 - 2	-0.520	181.3				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
Cd	std.add.	8.384e-008	9.546e-002	4.339e-009	0.99837
Pb	std.add.	6.485e-008	7.422e-002	1.876e-009	0.99955

Final results	+/-	Res. dev.	%	Comments
Cd: Cadmium	=	1.099 ug/L	0.047	4.256
Pb: Lead	=	1.093 ug/L	0.026	2.365