

## Comprehensive tap water analysis with TitrIC 2

Branch: water, waste water, environmental protection

### Keywords

TitrIC 2 / 861 / 855 / 712 / 800 / Aquatrode plus / 6.0257.000 / Dosino / tiamo / 6.1006.430 / Metrosep SUPP 4 – 250 / 6.1010.210 / Metrosep C 2 – 100 / 2.145.0320 / branch 2

### Summary

The determination of anions / cations in tap water was done on the **Metrosep SUPP 4 – 250 / Metrosep C 2 – 100** using carbonate / tartaric acid – dipicolinic acid eluent with chemical suppression. **TitrIC 2** enables a fully automated and parallel analysis of anions, cations, titration (p- and m- value), and direct measurements (temperature, pH, conductivity). In this application only the standard ions are determined. For a more complete analysis including the oxygenated halides please refer to **AB 303 e**.

### Sample

Tap water Herisau

### Reagents

Eluent: **2.4 mM sodium carbonate, 2.0 mM sodium bicarbonate, 5% acetone** in ultra pure water (resistivity >18 MΩ)  
**7 mM tartaric acid, 1.5 mM dipicolinic acid** in ultra pure water

Suppressor solutions:  
50 mmol/L sulfuric acid  
ultra pure water

### Standards [ppm] in ultra pure water

Anions:

Analyte / Level	1	2	3	4	5
fluoride	0.04	0.06	0.08	0.1	0.12
chloride	6	9	12	15	18
nitrate	6	9	12	15	18
sulfate	4	6	8	10	12

Cations:

Analyte / Level	1	2	3	4	5
sodium	1	2	4	6	8
potassium	0.5	1	2	3	4
calcium	12.5	25	50	75	100
magnesium	3.75	7.5	15	22.5	30

### Apparatus and Accessories

#### - delivered with TitrIC 2 package:

861 Adv. Compact IC with seq. supp.	2.861.0020
861 Adv. Compact IC	2.861.0010
855 Robotic Titrosampler	2.855.0020
712 Conductometer	2.712.0010
Cond. measuring cell Pt 1000	6.0912.110
Conductivity standard	6.2301.060
Aquatrode plus	6.0257.000
3 Dosions	2.800.0010
tiamo 1.1 full	6.6056.112
USB Converter Edgeport/4	2.145.0320
802 Rod Stirrer	2.802.0020

#### - used optional accessories for TitrIC 2:

Metrosep A SUPP 4 – 250	6.1006.430
Metrosep 4/5 Guard	6.1006.500
Metrosep C 2 – 100	6.1010.210
Metrosep C 2 Guard	6.1010.200
Sample rack 59 x 120 mL	6.2041.840
Sample beaker 120 mL (x100)	6.1459.300



### Sample Preparation

The water samples were injected directly.

### Analysis

The standards and the sample were injected using two 800 dosinos (10 and 50 mL Dosing unit) with automatic acidification for the cation run prior to injection. Loop volume anions: **20 µl**, cations **10 µl**.

The two (delivered with the TitrIC 2 package) *tiamo* methods *TitrIC 2 – complete run* and *TitrIC 2 – IC calibrator anion & cation* were used for all analysis.

### Calculation (IC)

Automatic integration with IC Net 2.3 software using peak area.

## Parameters

### Anions

#### STARTUP HARDWARE:

```

RECODER   METHOD  Asupp5.mtw
RECODER   DATA    Data acquisition [Cond]
861 Adv. Compact Unit version 2
861 Adv. Compact Polarity +
861 Adv. Compact Supp. autostep yes
861 Adv. Compact Autostep with Fill
861 Adv. Compact Flow 2.00 mL/min
861 Adv. Compact Pmax 15.0 MPa
861 Adv. Compact Pmin 0.0 MPa
861 Adv. Compact FullScale 50 uS/cm
861 Adv. Compact Remote 00001000

```

### Cations

#### STARTUP HARDWARE:

```

RECODER   METHOD  C2_100.mtw
RECODER   DATA    Data acquisition [Cond]
861 Adv. Compact Polarity -
861 Adv. Compact Supp. autostep no
861 Adv. Compact Autostep with Fill
861 Adv. Compact Flow 1.3 mL/min
861 Adv. Compact Pmax 25.0 MPa
861 Adv. Compact Pmin 0.0 MPa
861 Adv. Compact FullScale 1000uS/cm
861 Adv. Compact Remote 00000000

```

#### START WITH DETERMINATION

```

0.20 861 Adv. Compact Valve Inject
6.50 861 Adv. Compact Valve Fill

```

#### START WITH DETERMINATION

```

0.20 861 Adv. Compact Valve Inject
5.70 861 Adv. Compact Valve Fill

```

#### START WITH INJECT:

```
0.0 RECODER START
```

#### START WITH INJECT:

```
0.0 RECODER START
```

## SET pH 4,3(HCl) – Titration (m-value)

### Control parameters

EP1 at pH	4.3
Titration rate	user
Dynamics pH	1.00
Max. rate	10.00 mL/min
Min rate	15.00 µL/min
Stop criterion	drift
Stop drift	30 µL/min
End point 2	off

## MEAS pH

### Measuring parameters:

Measurement with drift control	on
Signal drift	10.0 mV/min
Min. waiting time	15 s
Max. waiting time	52 s
Measuring interval	2.0 s
Stop measured value pH	off
Measurement without drift control	off
Temperature	25.0 °C

## SET pH 8,2 (HCl) – Titration (p-/m-value)

### Control parameters:

EP1 at pH	8.2
Titration rate	user
Dynamics pH	1.00
Max. rate	10.00 mL/min
Min rate	15.00 µL/min
Stop criterion	drift
Stop drift	30 µL/min
End point 2	on
EP2 at pH	4.3
Dynamics pH	1.00
Max. rate	10.00 mL/min
Min rate	15.00 µL/min
Stop criterion	drift
Stop drift	30 µL/min

## MEAS conductivity

### Measuring parameters:

Measurement frequency	auto
Measuring time	35 s
Measuring interval	2.0 s
Stop measured value	off mS/cm
Temperature	20.0 °C

## Results

### TitrIC-Report

Print date	4/12/2006
<b>General data</b>	
ID	Tap water Herisau
Record date	3/09/2006
Record time	19:15:16
User comment	
System comment	recalculated, within limit
<b>Titration data</b>	
pH	7.72
Cond. [ $\mu\text{S}/\text{cm}$ ]	530
Temp. [°C]	21.72
m value [mmol/L]	5.15
p value [mmol/L]	0

<b>Anions</b>	
F [mg/L]	0.061
Cl [mg/L]	14.088
NO <sub>3</sub> [mg/L]	10.606
SO <sub>4</sub> [mg/L]	5.491
Anion balance [mEq/L]	5.6158
<b>Cations</b>	
Na [mg/L]	7.091
K [mg/L]	1.736
Mg [mg/L]	17.537
Ca [mg/L]	80.367
Cation balance [mEq/L]	5.8064
ionic difference [mEq/L]	0.1907
<b>ionic difference [%]</b>	1.67

Tap water Herisau

**Reproducibility:** Measurement of 10 tap water Herisau samples:

IC:

anions	fluoride	chloride	nitrate	sulfate
<b>Average (N=10)</b>	<b>0.062</b>	<b>13.795</b>	<b>7.463</b>	<b>5.809</b>
Standard Deviation	0.002	0.284	0.230	0.163
% RSD (N=10)	<b>2.92</b>	<b>2.06</b>	<b>3.09</b>	<b>2.81</b>

cations	sodium	potassium	calcium	magnesium
<b>Average (N=10)</b>	<b>7.823</b>	<b>1.834</b>	<b>86.420</b>	<b>18.452</b>
Standard Deviation	0.278	0.060	3.104	0.573
% RSD (N=10)	<b>3.55</b>	<b>3.26</b>	<b>3.59</b>	<b>3.11</b>

Titration / ionic balance:

	pH	Cond. [ $\mu\text{S}/\text{cm}$ ]	Temp. [°C]	m value [mmol/L]	p value [mmol/L]	sum anions [meq/L]	sum cations [meq/L]	ionic difference [%]
<b>Average (N=10)</b>	<b>7.74</b>	<b>572.0</b>	<b>22.33</b>	<b>5.51</b>	<b>0.00</b>	<b>6.143</b>	<b>6.217</b>	<b>+0.74</b>
Standard Deviation	0.039	12.1	0.268	0.167	0.00	0.064	0.059	-
% RSD (N=10)	<b>1.21</b>	<b>2.12</b>	<b>1.20</b>	<b>3.03</b>	-	<b>1.04</b>	<b>0.95</b>	-

### Comments

This Application Bulletin is based on AW CH6-0865-042006.

Overview of the determined parameters:

- direct measurement:  
conductivity, pH, temperature
- titrations:  
p- and m-value (= acidic capacity)
- IC:  
fluoride, chloride, nitrate, sulfate  
sodium, potassium, calcium, magnesium

TitrIC is absolutely flexible and can be adapted to determine the desired variables.

This Application work shows how to obtain the most important parameters for tap water analysis (conductivity, pH, m- and p-value, concentrations of F<sup>-</sup>, Cl<sup>-</sup>, NO<sub>3</sub><sup>-</sup>, SO<sub>4</sub><sup>2-</sup>, Na<sup>+</sup>, K<sup>+</sup>, Ca<sup>2+</sup>, Mg<sup>2+</sup>) in 9 minutes – thanks to the execution of the titration while the IC is already running. This parallel working enables a very high throughput of samples per day.

To obtain a fast anion analysis the A SUPP 4 – 250 mm and a special eluent was chosen with which one is able to determine all 7 standard anions in 6 minutes.

No sample preparation is necessary – even the acidification for the cation-IC-run is done automatically by a Dosino.

The calibration results are good which means that the filling of a loop with dosinos is working properly.

The ionic difference shown in the TitrIC-Report is calculated by comparison of total equivalent of anions (the carbonate concentration is calculated from the p- and m-values) and cations from the IC results:

$$\text{Ionic d. [\%]} = \frac{100 * (\text{sum cations [meq/L]} - \text{sum anions [meq/L]})}{(\text{sum cations [meq/L]} + \text{sum anions [meq/L]})}$$

Example: If 5 ppm sulfate is found, this corresponds to 5 mg/L → 0.052 mmol/L (M (SO<sub>4</sub>) = 96 mg/mmol) → -0.104 meq/L (charge SO<sub>4</sub>: -2).

For more detailed information about the system setup of TitrIC 2 please consult the document **Installation Instructions for TitrIC 2**.

The Aquatrode should be stored in the storage solution of Metrohm (6.2323.000).

### Appendix


Calibration curves, reproducibility, chromatograms with peak tables and titration curves

### Date

Thursday, 2006-04-12

### Name

A. Rumi; IC Marketing,  
Metrohm Ltd.; CH-9101 Herisau




**TitriC 2: System setup**


### Calibration

#### Anions:


fluoride: RSD: 1.516 %  
 $Q = -0.933 \cdot A^2 + 1.012 \cdot A - 0.00364$




chloride: RSD: 0.845 %  
 $Q = -0.0123 \cdot A^2 + 1.241 \cdot A + 0.2877$



nitrate: RSD: 0.646 %  
 $Q = -0.0803 \cdot A^2 + 4.1324 \cdot A - 0.1295$




sulfate: RSD: 0.619 %:  
 $Q = -0.178151 \cdot A^2 + 4.727 \cdot A - 0.1549$




#### Cations:


sodium: RSD: 1.94 %  
 $Q = 0.394518 \cdot A - 0.108419$




potassium: RSD: 2.45 %  
 $Q = 0.905112 \cdot A - 0.0977494$




calcium: RSD: 0.83 %  
 $Q = 0.494963 \cdot A - 0.564849$



magnesium: RSD: 0.44 %  
 $Q = 0.208015 \cdot A - 0.32116$




**Chromatograms with peak tables**




Multi-Anion Standard  
File q3091802  
Volume: 20.0  $\mu\text{l}$   
Dilution: 1.00  
Amount: 1.00

No	Retention min	Height $\mu\text{S}/\text{cm}$	Area $\mu\text{S}/\text{cm}^*\text{sec}$	Conc. mg/L	Name
1	1.84	0.09	0.477	0.0840	fluoride
2	2.56	0.15	0.922	0.4749	chloride
3	2.96	0.08	0.570	0.0000	nitrite
4	3.60	0.09	0.732	0.0000	bromide
5	4.01	0.10	0.906	0.2737	nitrate
6	5.11	0.03	0.493	0.0000	phosphate
7	5.86	0.10	1.382	0.2942	sulfate



Multi-Cation Standard  
File q3091821  
Volume: 10.0  $\mu\text{l}$   
Dilution: 1.00  
Amount: 1.00


No	Retention min	Height $\mu\text{S}/\text{cm}$	Area $\mu\text{S}/\text{cm}^*\text{sec}$	Conc. mg/L	Name
1	1.67	0.33	1.605	0.0000	lithium
2	2.04	0.42	2.405	0.8405	sodium
3	2.27	0.36	2.762	0.0000	ammonium
4	2.93	0.20	1.624	1.3726	potassium
5	3.60	0.34	3.874	1.3526	calcium
6	4.86	0.55	7.847	1.3112	magnesium



Standard 3

File q3091119  
Volume: 20.0  $\mu\text{l}$   
Dilution: 1.00  
Amount: 1.00

No	Retention min	Height $\mu\text{S}/\text{cm}$	Area $\mu\text{S}/\text{cm}^*\text{sec}$	Conc. mg/L	Name
1	1.83	0.09	0.482	0.0790	fluoride
2	2.54	10.47	60.282	11.9343	chloride
3	3.92	3.11	30.706	11.9327	nitrate
4	5.84	1.84	26.313	7.9521	sulfate




Tap water Herisau


File: q3091950  
Volume: 20.0  $\mu\text{l}$   
Dilution: 1.00  
Amount: 1.00

No	Retention min	Height $\mu\text{S}/\text{cm}$	Area $\mu\text{S}/\text{cm}^*\text{sec}$	Conc. mg/L	Name
1	1.84	0.07	0.442	0.0611	fluoride
2	2.57	12.00	73.252	13.4092	chloride
3	4.19	1.86	27.263	7.2891	nitrate
4	5.88	1.28	18.017	5.6201	sulfate

**Cations:**




No	Retention min	Height uS/cm	Area uS/cm*sec	Conc. mg/L	Name
1	2.03	1.69	10.590	4.0694	sodium
2	2.90	0.24	2.245	1.9342	potassium
3	3.40	6.96	102.136	49.9889	calcium
4	4.64	3.89	73.842	15.0391	magnesium



No	Retention min	Height uS/cm	Area uS/cm*sec	Conc. mg/L	Name
1	2.04	3.44	20.854	8.1187	sodium
2	2.91	0.20	1.711	1.4505	potassium
3	3.32	10.49	176.793	86.9410	calcium
4	4.60	4.68	95.884	19.6242	magnesium

**Titration curves – Titration of m-value with HCl**



**Fig. 1:** SET pH of tap water Herisau