

Determination of iron, copper and vanadium by adsorptive stripping voltammetry

Summary

This Application Bulletin describes the voltammetric determination of the elements iron, copper and vanadium. Fe as well as Cu and V can be determined as catechol complex at the HMDE by adsorptive stripping voltammetry (AdSV). Fe(II) and Fe(III) are determined as Fe(total) with the same sensitivity for both species in either phosphate buffer or PIPES electrolyte. Cu and V can be determined in PIPES buffer.

The methods are primarily suitable for the investigation of ground, drinking and surface waters, in which the concentration of these metals is important. But the methods can naturally also be used for trace analysis in other matrices.

The limit of detection for all three elements in PIPES buffer is 0.5 ... 1 µg/L, for iron in phosphate buffer it is approx. 5 µg/L.

Instruments

VA instrument
capable of operating a Multi-Mode
Electrode and supporting DP mode

909 UV Digester

Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Mercury drop capillary	6.1226.030
	or	6.1226.050
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel	6.1245.010
	Filled with c(KCl) = 3 mol/L	
AE	Pt rod electrode	6.0343.x00

Sample preparation

- Ground water, surface waters, mineral waters and drinking waters can usually be analyzed without pre-treatment.

Organic matter often interferes with voltammetric determinations and therefore sample solutions usually have to be digested.

- Low polluted waste waters can be digested with the 909 UV Digester:

Add 50 - 100 µL w(H₂O₂) = 30% and 10 µL w(HCl) = 30% to 10 mL acidified sample (pH = 2) and irradiate for 1 h at 90 °C. After cooling to room temperature, the digested sample can be transferred directly to the polarographic vessel.

- Samples with organic matter (foods, pharmaceuticals etc.) must be digested.

- High-pressure asher
- Microwave digestion

Both techniques oxidize the sample in a closed digestion vessel by means of a mixture of concentrated mineral acids.

- Open wet digestion with H₂SO₄ and H₂O₂ according to Application Bulletin 113.

Method 1: Determination of Fe with phosphate buffer

Iron is determined as a complex with catechol. With adsorptive stripping voltammetry at the HMDE concentrations down to 5 µg/L can be determined. For concentrations >100 µg/L the polarographic determination at the SMDE is recommended.

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Sodium dihydrogen phosphate, water free, NaH₂PO₄, for trace analysis*, CAS 7558-80-7
- Disodium hydrogen phosphate, water free, Na₂HPO₄, for trace analysis*, CAS 7558-79-4

- Catechol (Pyrocatechol, Brenzkatechin, 1,2-Dihydroxybenzene), for analysis, CAS 120-80-9
- Iron standard stock solution, $\beta(\text{Fe}^{3+}) = 1 \text{ g/L}$, commercially available
- Hydrochloric acid, suprapur: $w(\text{HCl}) = 30\%$, for trace analysis*, CAS 7647-01-0
- Ultrapure water, resistivity $>18 \text{ M}\Omega\cdot\text{cm}$ (25 °C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

Phosphate buffer pH 7.0	$c(\text{NaH}_2\text{PO}_4) = 0.2 \text{ mol/L}$ $c(\text{Na}_2\text{HPO}_4) = 0.3 \text{ mol/L}$ 23.996 g NaH_2PO_4 and 42.594 g Na_2HPO_4 are dissolved and filled up to 1 L with ultrapure water. The solution has pH 6.6 and is stable for approx. 1 month.
Catechol solution	$c(\text{catechol}) = 1 \text{ mol/L}$ 2.75 g catechol are dissolved in 25 mL oxygen-free ultrapure water under passage of nitrogen. The solution must be stored in the dark and allowed to stand for 1 h before use. The stability of the solution depends on the purity of the substance and ranges from one to several days.

Standard solutions

Fe standard solution	$\beta(\text{Fe}^{3+}) = 10 \text{ mg/L}$ Add 0.1 mL conc. HCl to 1.00 mL Fe standard stock solution and fill up to 100 mL with ultrapure water.
----------------------	---

Analysis

Measuring solution

10 mL (diluted) sample or digestion solution

100 μL catechol solution or some catechol crystals

1 mL phosphate buffer pH 7.0

The measuring solution should have $\text{pH } 7 \pm 0.1$.

The concentration is determined by standard addition.

Parameters for 5 to 100 $\mu\text{g/L}$ Fe

Voltammetric	
Electrode operating mode	HMDE
Drop size	7
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
<i>Potentiostatic pretreatment</i>	
Potential 1	-0.3 V
Waiting time 1	60 s
Equilibration time	5 s
<i>Sweep</i>	
Start potential	-0.2 V
End potential	-0.55 V
Potential step	0.004 V
Potential step time	0.4 s
Sweep rate	0.01 V/s
Pulse amplitude	0.05 V
Substance	
Name	Fe
Characteristic potential	-0.38 V


Parameters for 0.1 to 100 mg/L Fe

Voltammetric	
Electrode operating mode	SMDE
Drop size	4
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
Equilibration time	10 s
<i>Sweep</i>	
Start potential	-0.2 V
End potential	-0.55 V
Potential step	0.004 V
Potential step time	0.4 s
Sweep rate	0.01 V/s
Pulse amplitude	0.05 V
Substance	
Name	Fe
Characteristic potential	-0.38 V

Parameters for $>100 \text{ mg/L}$ Fe

Samples with concentrations over 100 mg/L Fe have to be diluted.

Example



Result

Sample	Water
Sample size	10.0 mL
$\beta(\text{Fe})$	128 µg/L

Comments

- Phosphate buffer has to be added after the addition of catechol. Otherwise insoluble iron phosphates can be formed, which cannot be determined voltammetrically.
- Catechol can be purified by recrystallization from toluene. It is also possible to make a sublimation of the catechol. Catechol Sigma-Aldrich No. 15890 can be used directly.

- Using this method, purity of the applied reagents is of utmost importance. Small concentrations of surface-active substances (detergents, organic complexing agents such as humic acids etc.) are also interfering. These must be destroyed before beginning the analysis (see "Sample Preparation").

References

- Monien H., Jacob P.
Voltammetrische Bestimmung kleiner Eisenmengen ohne Abtrennung der Matrix.
Fresenius Z. Anal. Chem. 260, (1972) 195-202
- Davidson W.
Comparison of differential pulse and D.C. sampled polarography for the determination of ferrous and manganese ions in lake water.
J. Electroanal. Chem. 72, (1976) 229-237
- Huang Z.Q., Van den Berg C.M.G.
Determination of iron in seawater using cathodic stripping voltammetry preceded by adsorptive collection with the hanging mercury drop electrode.
J. Electroanal. Chem. 177, (1984) 269-280; Ref: *Metrohm-Info* 1/85, 19
- Weidenauer M., Lieser K.H.
Bestimmung von Spurenelementen in Flusswasser mit Hilfe der Voltammetrie.
Fresenius, Z. Anal. Chem. 320, (1985) 550-555
- Wang J., Mahmoud J.
Chelate adsorption of trace voltammetric measurements of iron(III).
Fresenius, Z. Anal., Chem. 327, (1987) 789-793

Method 2: Determination of Fe with PIPES buffer

Iron is determined as a complex with catechol. With adsorptive stripping voltammetry at the HMDE concentrations down to 0.5 µg/L can be determined.

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Sodium hydroxide solution, w(NaOH) = 30%, for trace analysis*, CAS 1310-73-2

- Ammonia solution, $w(\text{NH}_3) = 25\%$, for trace analysis*, CAS 1336-21-6
- Catechol (Pyrocatechol, Brenzkatechin, 1,2-Dihydroxybenzene), for analysis, CAS 120-80-9
- PIPES, Piperazine-1,4-bis-2-ethane sulfonic acid, CAS 5625-37-6
- Iron standard stock solution, $\beta(\text{Fe}^{3+}) = 1 \text{ g/L}$, commercially available
- Nitric acid, $w(\text{HNO}_3) = 65\%$, for trace analysis*, CAS 7697-37-2
- Ultrapure water, resistivity $>18 \text{ M}\Omega \cdot \text{cm}$ (25°C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

PIPES buffer	$c(\text{PIPES}) = 1 \text{ mol/L}$ Mix 6.05 g PIPES with 1 mL sodium hydroxide solution and a little high purity water. Adjust the pH value to 8.0 with ammonia solution and fill up to 20 mL with high purity water. The solution has to be prepared fresh daily. This buffer does not contain any Fe.
Catechol solution	$c(\text{Catechol}) = 1 \text{ mol/L}$ Dissolve 2.75 g catechol in 25 mL oxygen-free high purity water, which had been degassed by purging with nitrogen. Store the solution in a tightly closed dark bottle. Stability of the solution depends on the purity of the used catechol and can vary from 1 day to 1 month. The solution must rest for a few hours before use.
Ammonia solution diluted	$w(\text{NH}_3) = 10\%$

Standard solutions

Iron standard solution	$\beta(\text{Fe}^{3+}) = 1 \text{ mg/L}$ The solution is diluted with $c(\text{HNO}_3) = 0.015 \text{ mol/L}$. It is stable for max. 1 week.
------------------------	--

Analysis

Measuring solution

10 mL acidified (diluted) sample

0.05 mL catechol solution

0.250 mL PIPES buffer

If necessary, adjust the pH value to 7.0 ± 0.1 with $w(\text{NH}_3) = 10\%$.


The measuring vessel should contain at least 10 ng but not more than 100 ng Fe.

The concentration is determined by 2 standard additions.

Parameters

Voltammetric	
Electrode operating mode	HMDE
Drop size	7
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
<i>Potentiostatic pretreatment</i>	
Potential 1	-0.25 V
Waiting time 1	60 s
Equilibration time	10 s
<i>Sweep</i>	
Start potential	-0.25 V
End potential	-0.6 V
Potential step	0.006 V
Potential step time	0.1 s
Sweep rate	0.06 V/s
Pulse amplitude	0.05 V
Substance	
Name	Fe
Characteristic potential	-0.4 V

Example



Results

Sample	Sea water
Sample size	10 mL
$\beta(\text{Fe})$	41.9 $\mu\text{g/L}$

Comments

- Catechol should be added to the acidic sample solution before the buffer is added. Otherwise iron may form hydroxo complexes, which cannot be determined any more. This results in low recovery rates.
- The peak potential of Fe is strongly dependent on the exact pH value and the sample matrix and may have to be adjusted.
- For amounts above 100 ng in the polarographic vessel, standard additions are no longer linear. Then, either the amount of the sample must be reduced or a DP polarographic analysis performed on the SMDE.

- If the total ion concentration is very high (e.g. sea water), the iron determination is much less sensitive and only relatively small peaks are achieved. Linear standard addition curves are achieved, however, with proportionally large standard additions and a reduction of the enrichment time.
- Catechol can be purified by recrystallization from toluene. It is also possible to make a sublimation of the catechol. Catechol Sigma-Aldrich No. 15890 can be used directly.

References

- Monien H., Jacob P.,
Voltammetrische Bestimmung kleiner Eisenmengen ohne Abtrennung der Matrix.,
Fresenius, Z. Anal. Chem. 260, (1972) 195-202
- Davidson W.,
Comparison of differential pulse and D.C. sampled polarography for the determination of ferrous and manganese ions in lake water.,
J. Electroanal. Chem. 72, (1976) 229-237
- Huang Z. Q., Van den Berg C. M. G.,
Determination of iron in seawater using cathodic stripping voltammetry preceded by adsorptive collection with the hanging mercury drop electrode.,
J. Electroanal. Chem. 177, (1984) 269-280
Ref: *Metrohm-Info* 1/85, 19
- Weidenauer M., Lieser K.H.,
Bestimmung von Spurenelementen in Flusswasser mit Hilfe der Voltammetrie.,
Fresenius, Z. Anal. Chem. 320, (1985) 550-555
- Wang J., Mahmoud J.,
Chelate adsorption of trace voltammetric measurements of iron(III).,
Fresenius, Z. Anal., Chem. 327, (1987) 789-793

Method 3: Determination of Fe, Cu and V

Iron, copper and vanadium are determined as a complex with catechol. With adsorptive stripping voltammetry at the HMDE concentrations down to 0.5 ... 1 µg/L can be determined.

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Sodium hydroxide solution, $w(\text{NaOH}) = 30\%$, for trace analysis*, CAS 1310-73-2
- Ammonia solution, $w(\text{NH}_3) = 25\%$, for trace analysis*, CAS 1336-21-6
- Catechol (Pyrocatechol, Brenzkatechin, 1,2-Dihydroxybenzene), for analysis, CAS 120-80-9
- PIPES, Piperazine-1,4-bis-2-ethane sulfonic acid, CAS 5625-37-6
- Iron standard stock solution, $\beta(\text{Fe}^{3+}) = 1 \text{ g/L}$, commercially available
- Copper standard stock solution, $\beta(\text{Cu}^{2+}) = 1 \text{ g/L}$, commercially available.
- Vanadium standard stock solution, $\beta(\text{V}^{5+}) = 1 \text{ g/L}$, commercially available.
- Nitric acid, $w(\text{HNO}_3) = 65\%$, for trace analysis*, CAS 7697-37-2
- Ultrapure water, resistivity $>18 \text{ M}\Omega \cdot \text{cm}$ (25 °C), type I grade (ASTM D1193)

* e.g., Merck suprapur®, Honeywell Fluka TraceSelect® or equivalent

Solutions

PIPES buffer	$c(\text{PIPES}) = 1 \text{ mol/L}$ Mix 6.05 g PIPES with 1 mL sodium hydroxide solution and a little high purity water. Adjust the pH value to 8.0 with ammonia solution and fill up to 20 mL with high purity water. The solution has to be prepared fresh daily. This buffer does not contain any Fe.
Catechol solution	$c(\text{Catechol}) = 1 \text{ mol/L}$ Dissolve 2.75 g catechol in 25 mL oxygen-free high purity water, which had been degassed by purging with nitrogen. Store the

solution in a tightly closed dark bottle. Stability of the solution depends on the purity of the used catechol and can vary from 1 day to 1 month. The solution must rest for a few hours before use.

Ammonia solution diluted

$w(\text{NH}_3) = 10\%$

Standard solutions

Iron standard solution	$\beta(\text{Fe}^{3+}) = 1 \text{ mg/L}$
Copper standard solution	$\beta(\text{Cu}^{2+}) = 1 \text{ mg/L}$
Vanadium standard solution	$\beta(\text{V}^{5+}) = 1 \text{ mg/L}$
	The solution are diluted with $c(\text{HNO}_3) = 0.015 \text{ mol/L}$. The are stable for max. 1 week.

Analysis

Measuring solution

10 mL acidified (diluted) sample

0.05 mL catechol solution

0.250 mL PIPES buffer

If necessary, adjust the pH value to 7.0 ± 0.1 with $w(\text{NH}_3) = 10\%$.

Parameters

Voltammetric	
Electrode operating mode	HMDE
Drop size	7
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
<i>Potentiostatic pretreatment</i>	
Potential 1	0 V
Waiting time 1	60 s
Equilibration time	10 s
<i>Sweep</i>	
Start potential	0 V
End potential	-0.8 V
Potential step	0.006 V
Potential step time	0.1 s
Sweep rate	0.06 V/s

Pulse amplitude 0.05 V

Substance

Name Cu

Characteristic potential -0.16 V

Name Fe





Characteristic potential -0.34 V

Name V

Characteristic potential -0.6 V

Example

Determination of Cu, Fe and V in salt solution



Results

Sample	1% salt solution
Sample size	10 mL
$\beta(\text{Cu})$	1.1 µg/L
$\beta(\text{Fe})$	5.7 µg/L
$\beta(\text{V})$	< 1 µg/L

Comments

- Catechol should be added to the acidic sample solution before the buffer is added. Otherwise iron may form hydroxo complexes, which cannot be determined any more. This results in low recovery rates.
- Vanadium(IV) cannot be determined by voltammetry.
- If a large excess of Cu is present, it is recommended not to preconcentrate this, otherwise the Hg drop is overloaded and neither Fe, nor V can be determined accurately.
- The peak potentials are strongly dependent on the exact pH value and the sample matrix and may have to be adjusted.

Appendix

Report of the example determination Fe in a water sample according to method 1

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====							
Determ.	: 06281559	User:		Date:	1999-06-28		
Modified	: 1999-06-28 15:59:32	Run :	0	Time:	15:59:19		
Sample table:	-						
Pos.	Ident.1/S1 water	Ident.2/S2	Ident.3/S3	Method.call	Sample size/S0 10.0 mL		
Method	: AB431_1						
Title	: Determination of Iron in Waters. AB431 method 1						
Remark1	: Determination of iron in water						
Remark2	: 10mL sample + 50uL catechol (1 mol/L) + 0.5mL phosph buff						
Substance	: Fe			Comments			
Mass conc.:	128.1 ug/L	Mass	: 1.281 ug	-----			
MC.dev.	: 6.12 ug/L (4.78%)	Add.mass	: 750 ng				
Cal.dev.	: -	V0.sample:	: 10 mL				
	VR	U/mV	I/nA	I.mean	Std.dev.	I.delta	Comments
	00	-369	-20.48	-20.40	0.1208		
	01	-369	-20.31				
	10	-367	-33.59	-33.27	0.4580	-12.87	
	11	-367	-32.94				
	20	-364	-44.01	-44.02	0.0158	-10.75	
	21	-364	-44.03				
Substance	Techn.	Y.reg/offset	Slope	Nonlin.	Mean deviat.		
Fe	std.add.	-2.072e-08	-1.706e-04		6.153e-10		
Final results				+/- Res.dev.	%	Comments	
	Fe =	128.13 ug/L		6.12	4.78		

Method print for the determination of Fe according to method 1

===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====							
Method: AB431_1 .mth		OPERATION SEQUENCE					
Title : Determination of Iron in Waters. AB431 method 1							

	Instructions	t/s	Main parameters			Auxiliary parameters	
1	DOS/M		V.added	0.550 mL			
2	REM		PIPES buffer				
3	SMPL/M		V.fraction	mL	V.total	L	
4	PURGE						
5	STIR	300.0	Rot.speed	2000 /min			
6	(ADD						
7	PURGE						
8	STIR	60.0	Rot.speed	2000 /min			
9	0PURGE						
10	SEGMENT		Segm.name	DPAdSV			
11	ADD>M		Soln.name	FeStd	V.add	0.050 mL	
12	ADD)2						
13	END						
Method: AB431_1		SEGMENT DPAdSV					

	Instructions	t/s	Main parameters			Auxiliary parameters	
1	(REP						
2	STIR	5.0	Rot.speed	2000 /min			
3	HMDE		Drop size	7	Meas.cell	normal	
4	DPMODE		U.ampl	-50 mV	t.meas	20.0 ms	
			t.step	0.40 s	t.pulse	40.0 ms	
5	MEAS	60.0	U.meas	-300 mV			
6	OSTIR	5.0					
7	FSWEEP	31.2	U.start	-200 mV	U.step	4 mV	
			U.end	-500 mV	Sweep rate	10 mV/s	
8	OMEAS		U.standby	mV			
9	REP)1						
10	END						

Report of the example determination Fe in sea water according to method 2

```
===== METROHM 746 VA TRACE ANALYZER (5.746.0100) =====
Determ.      : 09291112      User:      Date: 99-09-29
Modified     : no           Run : 1      Time: 11:12:16
Sample table: -
```

```
-----
Pos.  Ident.1/S1  Ident.2/S2  Ident.3/S3  Method.call  Sample size/S0
      Sea Water                               10 mL
-----
Method : AB431_1
Title  : Determination of Iron
Remark1: 10 ml sample + 50 µl catechol + 250 µl buffer --> pH 7.0
Remark2 :
```

```
-----
Substance : Fe
Mass conc.: 41.94 ug/L      Mass      : 419.4 ng
MC.dev.   : 0.573 ug/L (1.37%)  Add.mass : 250 ng
Cal.dev.   : -              V0.sample: 10 mL
-----
```

```
-----
VR  U/mV  I/nA  I.mean  Std.dev.  I.delta  Comments
-----
00 -406 -7.561 -7.489  0.1006
01 -407 -7.418
10 -403 -11.76 -11.75  0.0127  -4.258
11 -403 -11.74
20 -398 -15.67 -15.68  0.0139  -3.934
21 -400 -15.69
-----
```

```
-----
Substance  Techn.  Y.reg/offset  Slope  Nonlin.  Mean deviat.
-----
Fe          std.add.  -7.511e-09  -1.845e-04  -  6.494e-11
-----
```

```
-----
Final results      +/-  Res.dev.  %  Comments
-----
Fe = 41.935 ug/L      0.573  1.37
-----
```

Method print for the determination of Fe according to method 2

```
===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB431_1 .mth      OPERATION SEQUENCE
Title : Determination of Fe. AB 431 method 1
-----
```

```
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  SMPL>M
2  DOS>M      V.fraction      mL      V.total      L
3  DOS>M      Soln.name       Catechol V.add        0.050 mL
4  PURGE      Soln.name       Buffer      V.add        0.250 mL
5  STIR        300.0  Rot.speed      2000 /min
6  (ADD
7  PURGE
8  STIR        60.0  Rot.speed      2000 /min
9  OPURGE
10 (REP
11 SEGMENT      Segm.name      Fe_AdSV
12 REP)1
13 PURGE
14 ADD>M      Soln.name      Fe_Std      V.add        0.250 mL
15 ADD)2
16 END
-----
```

```
Method: AB431_1      SEGMENT
                      Fe_AdSV
-----
```

```
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  STIR        5.0  Rot.speed      2000 /min
2  HMDE        Drop size      7
3  DPMODE      U.ampl      -50 mV      Meas.cell     normal
                t.step      0.10 s      t.meas        20.0 ms
                U.meas      -250 mV      t.pulse       40.0 ms
4  MEAS        60.0
5  OSTIR       5.0
6  FSWEPP      6.2  U.start      -250 mV      U.step        6 mV
                U.end      -600 mV      Sweep rate    60 mV/s
7  OMEAS
8  END
-----
```

Report of the example determination Cu, Fe, V in a salt solution according to method 3

```
===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Determ.      : 02021647      User:      Date: 1998-02-02
Modified     : 1998-02-03 07:48:39 Run : 0 Time: 16:47:28
Sample table: -
```

```
-----
Pos.  Ident.1/S1  Ident.2/S2  Ident.3/S3  Method.call  Sample size/S0
      FeCuVrecal                      10 mL
-----
Method : AB431_3
Title  : Determination of Iron and Vanadium AB431 method 3
Remark1: 1% salt solution
Remark2:
-----
```

```
Substance : Fe
Mass conc.: 5.669 ug/L      Mass      : 56.69 ng
MC.dev.   : 0.847 ug/L (14.9%) Add.mass : 50 ng
Cal.dev.  : -              V0.sample: 10 mL
```

```
-----
VR  U/mV  I/nA  I.mean  Std.dev.  I.delta  Comments
-----
00 -341 -12.79 -12.79  0.0041
01 -341 -12.78
10 -329 -26.73 -27.11  0.5436 -14.33 front overlapping
11 -328 -27.50
20 -318 -35.53 -36.04  0.7203 -8.924 front overlapping
21 -318 -36.55 front overlapping
-----
```

```
Substance : Vanadium
Mass conc.: 192.9 ng/L      Mass      : 1.929 ng
MC.dev.   : 216 ng/L (112%) Add.mass : 50 ng
Cal.dev.  : -              V0.sample: 10 mL
```

```
-----
VR  U/mV  I/nA  I.mean  Std.dev.  I.delta  Comments
-----
00 -607 -0.0845 -0.0966  0.0170
01 -609 -0.1086
10 -602 -8.838 -8.819  0.0270 -8.722
11 -601 -8.800
20 -599 -15.67 -16.05  0.5327 -7.232
21 -600 -16.43
-----
```

```
Substance : Copper
Mass conc.: 1.103 ug/L      Mass      : 11.03 ng
MC.dev.   : 0.047 ug/L (4.23%) Add.mass : 10 ng
Cal.dev.  : -              V0.sample: 10 mL
```

```
-----
VR  U/mV  I/nA  I.mean  Std.dev.  I.delta  Comments
-----
00 -167 -16.00 -15.89  0.1642
01 -168 -15.77
10 -159 -30.57 -30.45  0.1661 -14.57 rear overlapping
11 -159 -30.34 rear overlapping
20 -148 -43.32 -43.97  0.9232 -13.51 rear overlapping
21 -148 -44.62 rear overlapping
-----
```

```
Substance  Techn.  Y.reg/offset  Slope  Nonlin.  Mean deviat.
-----
Fe          std.add.  -1.364e-08  -0.002478  1.572e-09
Vanadium    std.add.  -3.145e-10  -0.001680  4.705e-10
Copper      std.add.  -1.601e-08  -0.01495   5.295e-10
```

```
Final results      +/-  Res.dev.  %  Comments
-----
Fe = 5.6686 ug/L      0.847  14.9
Vanadium = 192.87 ng/L 216.  112.
Copper = 1.1031 ug/L  0.047  4.23
```

Method print for the determination of Cu, Fe, V according to method 3

```
===== METROHM 746 VA TRACE ANALYZER (5.746.0101) =====
Method: AB431_3.mth      OPERATION SEQUENCE
Title : Determination of Cu, Fe and V. AB 431 method 3
-----
```

```
-----
Instructions  t/s  Main parameters  Auxiliary parameters
-----
1  SMPL>M      V.fraction      mL  V.total      L
2  DOS>M      Soln.name      Catechol      V.add         0.050 mL
3  DOS>M      Soln.name      Buffer          V.add         0.250 mL
4  PURGE
5  STIR        300.0  Rot.speed      2000 /min
```

```

6      (ADD
7      PURGE
8      STIR          60.0   Rot.speed      2000 /min
9      OPURGE
10     (REP
11     SEGMENT          Segm.name   AdSV
12     REP)1
13     PURGE
14     ADD>M          Soln.name   Cu_Std      V.add      0.010 mL
15     ADD>M          Soln.name   Fe_Std      V.add      0.050 mL
16     ADD>M          Soln.name   V_Std       V.add      0.050 mL
17     ADD)2
18     END

```

Method: AB431_3

SEGMENT
AdSV

Instructions		t/s	Main parameters		Auxiliary parameters	
1	STIR	5.0	Rot.speed	2000 /min		
2	HMDE		Drop size	7	Meas.cell	normal
3	DPMODE		U.ampl	-50 mV	t.meas	20.0 ms
			t.step	0.10 s	t.pulse	40.0 ms
4	MEAS	60.0	U.meas	0 mV		
5	OSTIR	5.0				
6	FSWEEP	13.6	U.start	0 mV	U.step	6 mV
			U.end	-800 mV	Sweep rate	60 mV/s
7	OMEAS		U.standby	mV		
8	END					