

IC Application Note No. C-100

Title:	Choline in infant milk powder using online
	dialysis

Summary: Determination of choline in infant milk powder using cation chromatography with direct conductivity detection applying Metrohm Inline Dialysis.

Sample:Infant milk powderSample Preparation:500 mg dissolved in 100 mL ultrapure water, direct injection after Inline Dialysis


Column: 6.1010.210 Metrosep C 2 – 100

Eluent: 4.0 mmol/L tartaric acid

0.75 mmol/L dipicolinic acid, 2% acetone

Flow: 1.0 mL/min

Injection Volume: 20 μL



Results:	Choline
Concentration (mg/L)	11.04
RSD n = 6 (%)	1.6